REPORT DOCUMENTATION PAGE

Form Approved OMB No. 0704-0188

The public reporting burden for this collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing the burden, to Department of Defense, Washington Headquarters Services, Directorate for Information Operations and Reports (0704-0188), 1215 Jefferson Davis Highway, Suite 1204, Arlington, VA 22202-4302. Respondents should be aware that notwithstanding any other provision of law, no person shall be subject to any penalty for failing to comply with a collection of information if it does not display a currently valid OMB control number.

PLEASE DO NO	T RETURN YOU	R FORM TO TH	E ABOVE ADDRESS.			
	TE-(DD-MM-YY -11-2006	YY) 2. REPO	ORT TYPE THESIS			3. DATES COVERED (From - To)
	TION INTO TH	HE DEPTH OF PRATIVE MAT	CURE OF RESIN-MITERIALS.	ODIFIED		ANT NUMBER
					5c. PRO	GRAM ELEMENT NUMBER
6. AUTHOR(S) COL ROBER	TS HOWARD	W			5d. PRO	DJECT NUMBER
					5e. TAS	K NUMBER
					5f. WOF	RK UNIT NUMBER
	I G ORGANIZATI E UNIVERSIT		ND ADDRESS(ES)			8. PERFORMING ORGANIZATION REPORT NUMBER CI07-0007
THE DEPAR	TMENT OF TH	AGENCY NAM HE AIR FORC	E(S) AND ADDRESS(ES) E			10. SPONSOR/MONITOR'S ACRONYM(S)
AFIT/ENEL, 1 2275 D STRE WPAFB OH 4	ET					11. SPONSOR/MONITOR'S REPORT NUMBER(S)
12. DISTRIBUTION/AVAILABILITY STATEMENT Unlimited distribution In Accordance With AFI 35-205/AFIT Sup 1 DISTRIBUTION STATEMENT A						
13. SUPPLEMENTARY NOTES Approved for Public Release Distribution Unlimited						
14. ABSTRACT						
15. SUBJECT T	ERMS					
16 CECURITY	OI ACCIFICATIO	N OF	17. LIMITATION OF	18. NUMBER	100 NAS	AE OE DECDONCIDI E DEDCON
a. REPORT	b. ABSTRACT	c. THIS PAGE	ABSTRACT	OF PAGES		ME OF RESPONSIBLE PERSON
				219	19b. TELI	EPHONE NUMBER (Include area code)

INVESTIGATION INTO THE DEPTH OF CURE

OF

RESIN-MODIFIED GLASS-IONOMER RESTORATIVE MATERIALS

by

Howard W. Roberts, DMD

A Dissertation submitted to the Faculty of the Graduate School, Marquette University, in Partial Fulfillment of the Requirements for the Degree of Master of Science

Approved for Public Release
Distribution Unlimited

Milwaukee, Wisconsin August 2006

i

ABSTRACT

This investigation involved an attempt at delineating the depth of cure of resin-modified glassionomer restorative dental materials. Samples of different thicknesses using Vitremer Core Material and Restorative (3M/ESPE), Fuji II LC (GC America), and Photac-Fil Quick (3M/ESPE) were evaluated as to solubility, Knoop hardness, and thermal analysis techniques that included specific heat determination as well as differential scanning calorimetry thermal scans. Specimens were evaluated at time periods that included immediately after fabrication, 24 hours, one week, one month, and at three months. Overall results found that the solubility method chosen for this investigation did not provide the necessary sensitivity for depth of cure analysis of resin-modified glass-ionomer restorative materials. Hardness and thermal analysis provided evidence of a continuing, post-photopolymerization reaction that resulted in increased hardness, specific heat, and thermal requirements over the storage times. Furthermore, the resin-modified glass-ionomer restorative materials demonstrated water storage behavior similar to conventional glass-ionomer materials, in that water gained by the materials became more bound as storage time increased. Individual instances were observed in which the physical properties of 3 mm thick specimens were similar to that observed of 2 mm specimens, however these findings were not consistent throughout the investigation. Based on the conditions of this study, it is recommended that resinmodified glass-ionomer restorative materials should not be cured in thicknesses greater than two millimeters.

ABSTRACT

i

This investigation involved an attempt at delineating the depth of cure of resin-modified glassionomer restorative dental materials. Samples of different thicknesses using Vitremer Core Material and Restorative (3M/ESPE), Fuji II LC (GC America), and Photac-Fil Quick (3M/ESPE) were evaluated as to solubility, Knoop hardness, and thermal analysis techniques that included specific heat determination as well as differential scanning calorimetry thermal scans. Specimens were evaluated at time periods that included immediately after fabrication, 24 hours, one week, one month, and at three months. Overall results found that the solubility method chosen for this investigation did not provide the necessary sensitivity for depth of cure analysis of resin-modified glass-ionomer restorative materials. Hardness and thermal analysis provided evidence of a continuing, post-photopolymerization reaction that resulted in increased hardness, specific heat, and thermal requirements over the storage times. Furthermore, the resin-modified glass-ionomer restorative materials demonstrated water storage behavior similar to conventional glass-ionomer materials, in that water gained by the materials became more bound as storage time increased. Individual instances were observed in which the physical properties of 3 mm thick specimens were similar to that observed of 2 mm specimens, however these findings were not consistent throughout the investigation. Based on the conditions of this study, it is recommended that resinmodified glass-ionomer restorative materials should not be cured in thicknesses greater than two millimeters.

Any Opinions Expressed In This Thesis

Are Those of the Author Only

And Does Not Express the Official Opinion

Of the United States Air Force,

The Department of Defense, or

The United States Government

Marquette University

This is to certify that we have examined this copy of the Master's thesis by

Howard W. Roberts, D.M.D.

and have found that it is complete and satisfactory in all respects.

The thesis has been approved by:
Daniel Blown
David W. Berzins, Ph.D., Phesis Director, Department of General Dental
Sciences
Jeffx Ma Gill 3 Ph. P.
Jeffrey M. Toth, Ph.D.
Liaig 1. Vandewalle
Kraig S. Vandewalle, DDS, MS

ACKNOWLEDGEMENTS

This thesis is the result of clinical questions that I have pondered over years of clinically using resin-modified glass-ionomer restorative materials. Unfortunately, as with all research, as this thesis comes to a conclusion many questions still remain.

Enduring, structured, and meaningful work cannot be accomplished by a single entity. I have many people to whom I am indebted for their support and encouragement. First and foremost, I wish to thank my best friend and life partner:

vocation and military life even at times when you do not understand some of my choices. Without your love and support, none of this would have been possible.

give many thanks for understanding the times that dad was "busy" and had to do his "homework".

I wish to sincerely extend my appreciation to the dental manufacturers who graciously donated the restorative materials for this research. I hope that the investment that 3M/ESPE and GC America made will be worth their expense. If it were not for dental manufacturers such as these, dental research as we know today would not exist. I desire to also express my gratitude to the US Navy Institute for Dental and Biomedical Research and the US Air Force Dental Evaluation and Consultation Service for allowing me the additional office space and access to laboratory equipment.

My eternal gratitude will always be with my thesis committee, who provided wise direction and instruction. Finally, I wish to voice admiration for Drs. David Murchison and David Charlton, who both equally instilled the "fire in the belly" regarding the importance of research, dental materials, and how it relates into the "big picture" of dentistry.

The state of the s

AL LANGE CONTRACTOR OF THE PARTY OF THE PART

TABLE OF CONTENTS

ABSTRACT			1
ACKNOWLDGEMENTS			ii
LIST O	LIST OF TABLES		
LIST O	F FI	GURES	X
CHAPT	ER		
	I.	INTRODUCTION AND LITERATURE REVIEW	1
	II.	METHODS AND MATERIALS	10
	III.	RESULTS	16
	IV.	DISCUSSION	57
	V.	CONCLUSIONS	125
REFERENCES			128
APPEN	IDIC	ES.	
	A.	SOLUBILITY DATA	137
	B.	HARDNESS DATA	168
	C.	SPECIFIC HEAT THERMAL ANALYSIS DATA	181
	D.	DIFFERENTIAL SCANNING CALORIMETRY THERMAL SCAN THERMAL ANALYSIS DATA	185
	E.	THERMAL WEIGHT CHANGE THERMAL ANALYSIS DATA	189

LISTS OF TABLES

TABLE	PAGE
1. Vitremer Mean Solubility Results	17
2. Fuji II LC Mean Solubility Results	19
3. Photac-Fil Quick Mean Solubility Results	21
4. Vitremer Mean Hardness Results	22
 Scheffe Multi Comparisons For Vitremer Top Surface Hardness Over Storage Time 	23
 Scheffe Multi Comparisons For Vitremer Bottom Surface Hardness Over Storage Time 	23
7. Vitremer Mean Bottom/Top KHN Ratios	24
8. Fuji II LC Mean Hardness Results	24
 Scheffe Multi Comparisons For Fuji II LC Bottom Surface Hardness Over Storage Time 	25
10. Fuji II LC Mean Bottom/Top KHN Ratio	25
 Scheffe Mutli Comparisons For Fuji II LC Bottom/Top Surface KHN Ratio Over Storage Time 	26
12. Photac-Fil Quick Hardness Results	26
 Scheffe Multi Comparisons For Photac-Fil Quik Top Surface Hardness Over Storage Time 	27
 Scheffe Multi Comparisons For Photac-Fil Quik Bottom Surface Hardness Over Storage Time 	27
15. Photac-Fil Quick Mean Bottom/Top KHN Ratios	28
 Scheffe Multi Comparisons For Photac-Fil Quick Mean Bottom/Top KHN Ratios 	28
17. Vitremer Mean Specific Heat Results	28
 Scheffe Multi Comparisons For Vitremer Specific Heat per Thickness Over Storage Time 	29
19. Fuji II LC Mean Specific Heat Results	29
 Scheffe Multi Comparisons For Fuji II LC Specific Heat per Thickness Over Storage Time 	30
21. Photac-Fil Quick Mean Specific Heat Results	30

TA	LISTS OF TABLES (Continued) BLE	PAGE
22.	Scheffe Multi Comparisons For Photac-Fil Quick Specific Heat Over Storage Time	31
23.	Vitremer Thermal Scan Mean Results	32
24.	Scheffe Multi Comparisons For Vitremer Differential Scanning Calorimetry Thermal Results	32
25.	Vitremer Thermal Scan Mean Peak Temperature Results	33
26.	Scheffe Multi Comparisons For Vitremer Differential Scanning Calorimetry Thermal Peak Results	33
27.	Scheffe Multi Comparisons For Vitremer Differential Scanning Calorimetry Thermal Peak Results	34
28.	Fuji II LC Thermal Scan Mean Results	34
29.	Scheffe Multi Comparisons For Fuji II LC Differential Scanning Calorimetry Thermal Peak Results	34
30.	Scheffe Multi Comparisons For Fuji II LC Differential Scanning Calorimetry Thermal Peak Results	35
31.	Fuji II LC Thermal Scan Mean Peak Temperature	35
32.	Scheffe Multi Comparisons For Fuji II LC Differential Scanning Calorimetry Thermal Peak Results	36
33.	Scheffe Multi Comparisons For Fuji II LC Differential Scanning Calorimetry Thermal Peak Temperature Results	36
34.	Photac-Fil Quick Mean Thermal Scan Results	37
35.	Scheffe Multi Comparisons For Photac-Fil Quick Differential Scanning Calorimetry Thermal Results	37
	Scheffe Multi Comparisons For Photac-Fil Quick Differential Scanning Calorimetry Thermal Results	38
37.	Photac-Fil Quick Thermal Scan Mean Peak Temperature Results	38
38.	Scheffe Multi Comparisons For Photac-Fil Quick Differential Scanning Calorimetry Thermal Peak Results	39
39.	Scheffe Multi Comparisons For Photac-Fil Quick Differential Scanning Calorimetry Thermal Peak Temperature Results	39
40.	Vitremer Mean Thermal Weight Change	42
41.	Scheffe Multi Comparisons For Vitremer Water Gain	43

TAI	LISTS OF TABLES (Continued) BLE	PAGE
42.	Scheffe Multi Comparisons For Vitremer Water Loss During Specific Heat Determinations	43
43.	Scheffe Multi Comparisons For Vitremer Net Water Retention During Specific Heat Determinations	44
44.	Scheffe Multi Comparisons For Vitremer DSC Scan Weight Loss	45
45.	Scheffe Multi Comparisons For Vitremer Total Weight Loss	45
46.	Scheffe Multi Comparisons For Vitremer Ratio of Total Weight Loss to Water Gain	46
47.	Scheffe Multi Comparisons For Vitremer Ratio of DSC Thermal Scan Weight Loss to Water Gain	46
48.	Fuji II LC Mean Thermal Weight Change	47
49.	Scheffe Multi Comparisons For Fuji II LC Water Gain	48
50.	Scheffe Multi Comparisons For Fuji II LC Water Loss During Specific Heat Determinations	48
51.	Scheffe Multi Comparisons For Fuji II LC Net Water Retention During Specific Heat Determinations	49
52.	Scheffe Multi Comparisons For Fuji II LC Weight Loss During DSC Thermal Scan	50
53.	Scheffe Multi Comparisons For Fuji II LC Total Weight Loss	50
54.	Scheffe Multi Comparisons For Fuji II LC Total Weight Loss To Water Gain	51
55.	Scheffe Multi Comparisons For Fuji II LC DSC Thermal Scan Weight Loss to Water Gain	51
56.	Photac-Fil Quick Mean Thermal Weight Change	52
57.	Scheffe Multi Comparisons For Photac-Fil Quick Water Gain	53
58.	Scheffe Multi Comparisons For Photac-Fil Quick Water Loss During Specific Heat Determinations	53
59.	Scheffe Multi Comparisons For Photac-Fil Quick Net Water Retention During Specific Heat Determinations	54
60.	Scheffe Multi Comparisons For Photac-Fil Quick Total Weight Loss	54
61.	Scheffe Multi Comparisons For Photac-Fil Quick DSC Thermal	55

LISTS OF TABLES (Continued)

TABLE	PAGE
62. Scheffe Multi Comparisons For Photac-Fil Quick Total Weight Loss To Water Gain	55
63. Scheffe Multi Comparisons For Photac-Fil Quick DSC Thermal Scan Weight Loss to Water Gain	56
64. Vitremer Immediate Solubility Experimental Sample Results	138
65. Vitremer Immediate Solubility Control Sample Results	139
66. Vitremer 24-Hour Solubility Experimental Sample Results	140
67. Vitremer 24-Hour Solubility Control Sample Results	141
68. Vitremer One-Week Solubility Experimental Sample Results	142
69. Vitremer One-Week Solubility Control Sample Results	143
70. Vitremer One-Month Solubility Experimental Sample Results	144
71. Vitremer One-Month Solubility Control Sample Results	145
72. Vitremer Three-Month Solubility Experimental Sample Results	146
73. Vitremer Three-Month Solubility Control Sample Results	147
74. Fuji II LC Immediate Solubility Experimental Sample Results	148
75. Fuji II LC Immediate Solubility Control Sample Results	149
76. Fuji II LC 24-Hour Solubility Experimental Sample Results	150
77. Fuji II LC 24-Hour Solubility Control Sample Results	151
78. Fuji II LC One-Week Solubility Experimental Sample Results	152
79. Fuji II LC One-Week Solubility Control Sample Results	153
80. Fuji II LC One-Month Solubility Experimental Sample Results	154
81. Fuji II LC One-Month Solubility Control Sample Results	155
82. Fuji II LC Three-Month Solubility Experimental Sample Results	156
83. Fuji II LC Three-Month Solubility Control Sample Results	157
84. Photac-Fil Quick Immediate Solubility Experimental Sample Results	158
85. Photac-Fil Quick Immediate Solubility Control Sample Results	159
86. Photac-Fil Quick 24-Hour Solubility Experimental Sample Results	160
87. Photac-Fil Quick 24-Hour Solubility Control Sample Results	161

LISTS OF TABLES (Continued)

TABLE	PAGE
88. Photac-Fil Quick One-Week Solubility Experimental Sample Results	162
89. Photac-Fil Quick One-Week Solubility Control Sample Results	163
90. Photac-Fil Quick One-Month Solubility Experimental Sample Results	164
91. Photac-Fil Quick One-Month Solubility Control Sample Results	165
92. Photac-Fil Quick Three-Month Solubility Experimental Sample Results	166
93. Photac-Fil Quick Three-Month Solubility Control Sample Results	167
94. Vitremer 24-Hour Hardness	169
95. Vitremer One-Week Hardness	170
96. Vitremer One-Month Hardness	171
97. Vitremer Three-Months Hardness	172
98. Fuji II LC 24-Hour Hardness	173
99. Fuji II LC One-Week Hardness	174
100. Fuji II LC One-Month Hardness	175
101. Fuji II LC Three-Months Hardness	176
102. Photac-Fil Quick 24-Hour Hardness	177
103. Photac-Fil Quick One-Week Hardness	178
104. Photac-Fil Quick One-Month Hardness	179
105. Photac-Fil Quick Three-Months Hardness	180
106. Vitremer Specific Heat Results	182
107. Fuji II LC Specific Heat Results	183
108. Photac-Fil Quick Specific Heat Results	184
109. Vitremer DSC Thermal Scan Results	186
110. Fuji II LC DSC Thermal Scan Results	187
111. Photac-Fil Quick DSC Thermal Scan Results	188
112. Vitremer Immediate Thermal Weight Data	190
113. Vitremer 24-hour Thermal Weight Data	191
114 Vitremer One-Week Thermal Weight Data	192

LISTS OF TABLES (Continued)

TABLE	PAGE
115. Vitremer One-Month Thermal Weight Data	193
116. Vitremer Three-Month Thermal Weight Data	194
117. Fuji II LC Immediate Thermal Weight Data	195
118. Fuji II LC 24-Hour Thermal Weight Data	196
119. Fuji II LC One-Week Thermal Weight Data	197
120. Fuji II LC One-Month Thermal Weight Data	198
121. Fuji II LC Three-Month Thermal Weight Data	199
122. Photac-Fil Quick Immediate Thermal Weight Data	200
123. Photac-Fil Quick 24-Hour Thermal Weight Data	201
124. Photac-Fil Quick One-Week Thermal Weight Data	202
125. Photac-Fil Quick One-Month Thermal Weight Data	203
126 Photac-Fil Quick Three-Month Thermal Weight Data	204

LIST OF FIGURES

FIGURE	PAGE
1. Solubility Specimen Mold Before Material Placement	12
2. Solubility Mold Material Placement	12
3. Specimen Prior to Photopolymerization	12
4. Specimen Photopolymerization	13
5. Specimen Before Storage	13
6. Vitremer Group Solubility per Storage Time	60
7. Vitremer Group Solubility per Specimen Thickness	60
8. Overlay of Vitremer Solubility Results	61
9. Fuji II LC Group Solubility per Storage Time	62
10. Fuji II LC Group Solubility per Specimen Thickness	62
11. Overlay of Fuji II LC Solubility Results	63
12. Photac-Fil Quick Solubility per Storage Time	64
13. Photac-Fil Quick Solubility per Specimen Thickness	64
14. Overlay of Photac-Fil Quick Solubility Results	65
15. Material Experimental Group Marginal Means Solubility Results as a Function of Storage Time	66
16. Vitremer Top Surface Mean KHN	72
17. Vitremer Bottom Surface Mean KHN	73
18. Fuji II LC Top Surface Mean KHN	74
19. Fuji II LC Bottom Surface Mean KHN	75
20. Photac-Fil Quick Top Surface Mean KHN	76
21. Photac-Fil Quick Bottom Surface Mean KHN	76
22. Vitremer Mean Bottom/Top KHN Ratios Over Storage Time	78
23. Fuji II LC Mean Bottom/Top KHN Ratios Over Storage Time	79
24. Photac-Fil Quick KHN Bottom/Top KHN Ratios Over Storage Time	80
25. Vitremer Mean Specific Heat Results	86
26. Fuji II LC Mean Specific Heat Results	87

LIST OF FIGURES (Continued)

FIGURE	PAGE
27. Photac-Fil Quick Mean Specific Heat Results	88
28. Vitremer Mean Thermal Scan Enthalpy Results	90
29. Vitremer 24 Hour Mean Thermal Scans	91
30. Vitremer One Week Mean Thermal Scans	91
31. Vitremer One Month Mean Thermal Scans	92
32. Vitremer Three Month Mean Thermal Scans	92
33. Vitremer Thermal Scan Mean Peak Temperature Results	93
34. Fuji II LC Mean Thermal Scan Enthalpy Results	94
35. Fuji II LC One Week Thermal Scans	95
36. Fuji II LC Three Month Thermal Scans	95
37. Fuji II LC Thermal Scan Mean Peak Temperature Results	96
38. Photac-Fil Quick Mean Thermal Scan Enthalpy Results	97
39. Photac-Fil Quick 24 Hour Mean Thermal Scans	98
40. Photac-Fil Quick One Week Mean Thermal Scans	98
41. Photac-Fil Quick One Month Mean Thermal Scans	99
42. Photac-Fil Quick Three Month Mean Thermal Scans	99
43. Photac-Fil Quick Thermal Scan Mean Peak Temperature Results	100
44. Vitremer Mean Water Loss During Specific Heat Determination	103
45. Fuji II LC Mean Water Loss During Specific Heat Determination	104
46. Photac-Fil Quick Mean Water Loss During Specific Heat Determination	105
47. Vitremer Mean Net Water Retention During Specific Heat Determination	106
48. Fuji II LC Mean Net Water Retention During Specific Heat Determination	107
49. Photac-Fil Quick Mean Net Water Retention During Specific Heat Determination	108
50. Vitremer Thermal Scan Mean Weight Loss	111
51. Fuji II LC Thermal Scan Mean Weight Loss	112
52. Photac-Fil Quick Thermal Scan Mean Weight Loss	113

LIST OF FIGURES (Continued)

FIGURE	PAGE
53. Vitremer Thermal Determinations Mean Total Weight Loss	115
54. Fuji II LC Thermal Determinations Mean Total Weight Loss	115
55. Photac-Fil Quick Thermal Determinations Mean Total Weight Loss	116
56. Vitremer Ratio Mean Thermal Scan Loss/Water Gain	118
57. Fuji II LC Ratio Mean Thermal Scan Loss/Water Gain	118
58. Photac-Fil Quick Ratio Mean Thermal Scan Loss/Water Gain	119
59. Vitremer Ratio Mean Total Weight Loss/Water Gain	120
60. Fuji II LC Ratio Mean Total Weight Loss/Water Gain	121
61. Photac-Fil Quick Ratio Mean Total Weight Loss/Water Gain	122

CHAPTER 1: INTRODUCTION AND LITERATURE REVIEW

Glass-ionomer cements were invented and developed by Wilson and Kent in the early 1970's.
Glass ionomer cements generally consist of a mixture of various polyacrylic acids and tartaric acid that react with a fluoroaluminosilcate glass. The setting reaction occurs via an acid-base reaction between the acids and the calcium and aluminum ions that are leached from the glass component.
The early glass-ionomer cements have undergone incremental changes in formulations in attempts to improve their physical properties and clinical handling characteristics. These modifications have included the use of alternative polyacids, water-activated dehydrated polyacid powders,
netal additions, and smaller mean glass particle size.**

Resin-modified glass-ionomers (RMGI) are materials that consist substantially of glass-ionomer components (water-soluble polymeric acids, ion-leachable glass, and water) combined with organic, photopolymerizable monomers and an initiation system. Although the exact composition of the RMGI materials is proprietary, the resin composition in the set material has been estimated to be approximately 4.5 to six percent. These materials were developed to hopefully improve the mechanical properties and early moisture sensitivity observed with the conventional glass-ionomer materials. The RMGI materials basically exist in two forms:

One in which part of the water is replaced by water-soluble, photopolymerizable 2-hydroxyethylmethacrylate monomer (HEMA), and one in which methacrylate groups are predominately grafted onto the polyacrylic molecules. The Fuji II LC (GC America Corp, Alsip, IL) and Photac-Fil Quick (3M ESPE, St. Paul, MN) have been identified as members of the former, while Vitremer (3M ESPE) is considered to be a member of the latter. Although HEMA is the predominate methacrylate used, Tenas some formulations may use a combination of HEMA and diurethane dimethacrylate.

Nicholson³ described the setting reaction of RMGI materials as complex compared to conventional glass ionomer materials, especially in terms of the interactions of HEMA and other added constituents. Although a required constituent, HEMA lengthens both the working and

setting times of a glass ionomer material, as well as decreases the compressive strength. ^{19,20}
HEMA also affects the configuration of the polyacrylic acid component and has been theorized to induce possible phase separations: polymerized HEMA components are water insoluble and may precipitate out of solution, establishing different domains within the material. ¹⁶
Paradoxically, some RMGI materials have been reported to not afford a long working time in the clinical environment, as the light-activated polymerization component of some products have been described as sensitive to ambient light exposure. ²¹ Some concern has been voiced regarding deeper regions of RMGI materials where light penetration may not be of sufficient concentration for photopolymerization. ³ Manufacturers have added additional chemically-activated (i.e., auto cure) methods in order to address clinical concerns. ²¹ There is some apprehension concerning the leaching of unreacted monomers, as such, RMGI materials may be the subject of future biocompatibility investigations. ²² Interestingly, a recent report details the formation of a HEMA-free RMGI material that may address future cytotoxicity concerns of HEMA in these materials. ²³

RMGI Setting Reaction

The setting reaction of RMGI materials has been studied as a function of delay or omission of visible light activation, ¹³ with the chemical setting process also investigated using infrared spectroscopy. ²⁴⁻²⁷ Li *et al.* ²⁷ investigated the setting behavior of 2-mm thick samples of RMGI materials using Fourier transform infrared (FTIR) spectroscopy and diametral tensile strength (DTS). They reported that DTS improved with time, which reflected a continuing setting reaction. FTIR analysis reported only aliphatic C=C groups were present in the RMGI liquid and resin but could not determine if the HEMA groups were free in solution or grafted on the polyacrylic acid backbone chain. In comparing the effect of delaying visible light curing on the strength of RMGI, Li *et al.* reported in the absence of light activation Vitremer obtained a higher DTS compared to Photac Fil Aplicap. ²⁷ Interestingly, delay of photoinitiation had no effect on Vitremer, whereas delay of light activation significantly lowered the DTS of Photac-Fil. ²⁷

Kakaboura *et al.*²⁸ also used FTIR spectroscopic analysis using one-millimeter-thick RMGI specimens, identifying the 1638 cm⁻¹ methacrylate C=C stretch vs. the 1712 cm⁻¹ ester C=O stretch as a means for assessing the resin reaction. For monitoring the acid/base reaction, the authors compared the 1740 cm⁻¹ ester stretch to the 1600-1500 cm⁻¹ peak of carboxylate salts formed. They reported that Fuji II LC, Photac-Fil, and Vitremer did exhibit an acid/base reaction after 20 minutes of dark storage and that the acid/base reaction was significantly slowed by the light activation process.²⁸ Of the three RMGI products investigated, Photac-Fil Aplicap (a previous formulation of Photac-Fil Quick) formed the lowest acid/base reaction products while no difference was found between Fuji II LC and Vitremer. The authors suggested that this finding was due to fact that Photac Fil Aplicap contained most of the polyacrylic acid in the powder component, and required more water infiltration for full reaction. After 20 minutes of dark storage, Fuji II LC and Photac-Fil Aplicap were found to have limited methacrylate curing efficiency, whereas Vitremer exhibited a self-curing methacrylate curing reaction that did not appear to interfere with the acid/base reaction. However, the authors could not speculate on the influence of immediate light activation on the retardation of the acid/base reaction.²⁸

de Gee *et al.*¹³ studied the effects of delay/omission of photopolymerization on the three-body wear characteristics of Fuji II LC, Photac-Fil, and Vitremer. These researchers reported that Photac-Fil Aplicap did not seem to contain an auto-polymerization catalyst for the HEMA component, and while delayed light activation favored the development of the acid-base component, total omission of light activation resulted in a relatively weak material that could leach unpolymerized HEMA.¹³ In contrast, Fuji II LC was found to contain an auto-curing HEMA mechanism, but in omission of light activation the poly-HEMA chains polymerized at a rate that retarded the acid/base reaction which led the authors to speculate that the HEMA network approached full formation within a few minutes of mixing. Vitremer, in whose structure is formed by a chemical integration of the polyalkenoate and poly-HEMA network, was found to have better physical properties and the researchers observed that the HEMA auto-cure polymerization reaction did not appear to interfere with the simultaneous acid/base reaction.¹³

Ando *et al.*²⁹ reported on the effect of delay of light curing for up to 40 minutes of a RMGI orthodontic cement (Fuji Ortho LC, GC America). These researchers found that although a trend towards reduction of bond strength was observed, there was no significant difference between tensile bond strengths of the RMGI orthodontic cement observed between 0, 5, 20, and 40 minutes delay of light exposure.²⁹

Wan et al.24 investigated Fuji II LC using FTIR spectroscopy and reported that although the photopolymerized reaction is responsible for the initial set, the detected acid-base reaction continued at a rate that leveled at 168 hours after mixing. They confirmed earlier reports³ that HEMA retards the acid-base reaction by inducing stereochemical changes to the polyacrylic polymer, while the set HEMA resin matrix established physical barriers that slowed the process of diffusion.²⁴ Watson et al.³⁰ evaluated the RMGI/dentin interface during the initial setting reaction using laser-scanning confocal microscopy and reported that water movement from the dentinal tubules into the RMGI material was observed during light polymerization and further theorized that this water movement was essential for the formation of the "absorption layer" observed with RMGI materials.30 Young31 noted the difficulty of using FTIR spectroscopy for distinguishing RMGI setting changes due to photopolymerization and those arising from the acid-base reaction. Young observed peak loss at 1322 and 1300 cm⁻¹ bands due to water (1640 cm⁻¹ stretch band) masking the 1633 cm⁻¹ C=C stretch band and reported that 90 percent of the monomer was polymerized after 20 seconds of light exposure. 31 However, the acid-base reaction of unpolymerized sections were reported to proceed slower compared to conventional glass ionomer products. 31 Young et al. 32 studied the acid/base reaction of glass ionomer cements using Raman spectroscopy and confirmed the necessity of water availability for the acid/base initial reaction to proceed.

Palmer *et al.*³³ evaluated the HEMA released as a function of curing method for two RMGI liner and two RMGI restorative materials. Both of the restorative materials polymerized by five minutes when auto-cured in the absence of light activation, but both materials released

significantly more HEMA compared to the same products that were photopolymerized. These researchers also reported that HEMA release from both RMGI materials was similar after photopolymerization, but Vitremer restorative material released significantly less HEMA than Fuji II LC, suggesting that the Fuji II LC auto-cure mechanism is less efficient compared to Vitremer restorative. ³³ Interestingly, Palmer *et al.* corroborated earlier research by Kanchanavasita *et al.* ³⁴ that Fuji II LC would benefit from a visible light photopolymerization time longer than that recommended by the manufacturer. ³³

RMGI Adhesion Studies

Resin-modified glass-ionomer products have been demonstrated to adhere to both enamel and dentin. 35-46 The adhesion mechanism was first described by Wilson *et al.* 48 through the development of ionic crosslinks at the restorative-tooth interface via the then-described "ion exchange layer." In the ion-exchange layer, Wilson *et al.* 48 postulated the formation of an ionic bond between the polyacrylic acid and the hydroxyapatite component of the mineralized tissue. This mechanism was later confirmed by Yoshida *et al.* 49 using X-ray Photoelectron Spectroscopy while the same investigators further detailed the interaction between the polyalkenoic acids and hydroxyapatite in a later report. 50 The bonding of RMGI to dentin was later found to be two-fold, being both the ionic bond between tooth hydroxyapatite and polyacrylic acids as confirmed by Yoshida *et al.* 49,50 but also by a micromechanical means by the formation of a hybrid layer, analogous to, but of smaller thickness of that seen with dentin bonding agent systems. 51-53

Sidhu and Watson⁵⁴ studied the RMGI tooth/restoration interfacial characteristics using confocal fluorescence microscopy and theorized that the seemingly structure-less, identified "absorption layer" was possibly due to water movement within the RMGI material, either by water diffusion within the setting restoration material or from the dentin and pulp of a vital tooth. The absorption layer was also observed to form over time after initial setting of the RMGI matrix and was found only in deep dentin, not being present when the RMGI material was bonded to

enamel or superficial dentin.⁵⁴ Sidhu *et al.*⁵⁵ investigated the shear bond strength of RMGI materials and reported that the absorption layer appeared to have an important role in mediating the bond between RMGI and dentin. Tay *et al.*⁵⁶ further researched the absorption layer and confirmed that it was formed only in the presence of available moisture and was a necessary feature for the development of bond strength of the RMGI to dentin. Furthermore, silver-nitrate staining techniques identified capillary pore spaces that were theorized to provide the continuous water, arising from the hydrated dentin, required for full development of the absorption layer.⁵⁶

Yiu *et al.*⁵⁷ investigated the interaction of RMGI materials with moist dentin using transmission electron microscopy (TEM) and scanning electron microscopy (SEM) ultrastructural methods and reported unique spherical bodies within the RMGI matrix only in areas that were close to the RMGI-dentin interface. Further analysis of these spherical bodies allowed the authors to theorize that the spherical body core was composed of a poly-HEMA matrix derived from water transmission from the moist dentin, while the periphery contained products from the slower ongoing acid-base reaction.⁵⁷ Pereira *et al.*⁴⁷ investigated RMGI restorative material bond strength to dentin as a function of dentin regional differences and simulated pulpal pressure, and reported that pulpal pressure had a strong influence on bond strength. These researchers also observed that although water is essential to the RMGI setting reaction, excess water may be detrimental to the early setting reaction.⁴⁷

RMGI Depth of Cure Studies

The depth of cure of RMGI materials has been reported in some studies and reviews, with various methodology and materials reported. Burke *et al.*⁵⁸ investigated the depth of cure of early RMGI lining agents using a scrape methodology and reported that the liners displayed greater hardness at increased depths when measured 12 hours after photopolymerization, suggesting that the RMGI products possessed an auto-cure acid/base reaction in addition to photopolymerization. Swift *et al.*⁵⁹ investigated the depth of cure of Fuji II LC, Photac-Fil

Aplicap, and Vitremer compared to two compomers. They photopolymerized rectangular 5 mm x 5 mm molds that were 9 mm in depth and measured the Knoop Hardness at time intervals of 10 minutes, 24 hours, and seven days after photopolymerization. For all materials tested, only the top five millimeters of the specimens had enough consistent hardness for measurement. At the end of seven days, there was no statistical difference in hardness between specimen top surfaces and five-millimeter levels. However, due to lower hardness values at increased depths, the authors recommended that it would be prudent to place and cure RMGI materials in two to three millimeter increments. One literature review²¹ reported depth of cure values of RMGI lining and restorative materials that ranged from 1.54 to 2.26 mm for lining materials and 2.68 to 2.97 for restorative materials. However, this review only referenced the methodology used for depth of cure determinations and does not identify the source of the data. Accordingly, it cannot be ascertained if this reported data has been subject to scientific peer review.

The depth of cure of RMGI materials has not received the attention that has been directed to the resin composite restorations. The main cause for this lack of interest may be due to the fact that most clinicians oversimplify RMGI materials compared to conventional glass—ionomer materials and reason that the depth of cure is not as significant, assuming that the acid/base component reaction will compensate for inadequate penetration of light.³ Research has established that the RMGI material setting reaction is indeed complex and earlier investigations regarding RMGI depth of cure may have been too simplistic in nature. For instance, early laboratory studies were not comprehensive in scope as per physical properties and no investigations to date have included *in vitro* models that investigates the *in vivo* pulpal pressure hydration effect on the physical properties of the deeper RMGI regions. Accordingly, it can be conjectured that water influx from pulpal pressure may have a dilution effect on these deeper regions that have not experienced photopolymerization, with a possible negative effect on material physical and mechanical properties.

Purpose


The purpose of this investigation was to characterize the depth of cure of resin-modified glass-ionomer restorative materials using physical property testing to include hardness, solubility, as well as thermal analysis using differential scanning calorimetry. The research hypothesis is that RMGI materials will display a depth of cure gradient that is more dependent on visible light penetration rather than what residual acid-base chemical setting reaction that may remain.

CHAPTER TWO: METHODS AND MATERIALS

The RMGI materials consisted of Fuji II LC (GC America), Vitremer (3M/ESPE), and Photac Fil Quick (3M/ESPE). All materials were prepared and mixed per manufacturers instructions in a dental materials lab fitted with special lighting that excluded the light spectrum to which photoinitiators are sensitive (F40/GO Gold, Sylvania, Danvers, MA 01923), except for specimens prepared for thermal analysis which was accomplished in a different laboratory. For precapsulated materials, mechanical mixing was accomplished with a microprocessor-controlled calibrated and programmed mixing unit (Automix, Sybron Dental Specialties, Orange CA). Materials that were not precapsulated were manually prepared per manufacturer directions, placed in a syringe delivery system (AccuDose Low Viscosity, Centrix Inc., Shelton CT), and injected into respective molds. Visible light polymerization was accomplished using a halogen-based, visible-light polymerization unit (Optilux 501, Sybron Dental Specialties, Orange CA) whose performance was assessed using a laboratory laser power meter (PowerMax 5200 with PM10 Probe, Molectron, Portland OR) before each use. When not in use, specimens were stored in dark conditions at 98 ± 2 percent humidity at 37 ± 2 °C.


For the purpose of this master's thesis, physical property and thermal analysis methods were observed, depending on the method, at intervals of immediately after fabrication, 24 hours, one week, one month, and at three months. Continued observations will continue after completion of the thesis at six, nine, and twelve months after preparation for the purpose of information gathering for submission of manuscripts.


Solubility: Solubility testing was based on ANSI/ADA Specification No. 27: Resin-Based Filling Materials (1997),⁶⁰ using Delrin molds with a diameter of 8.0 mm and depths of 2, 3, 4, and 5 mm. Ten specimens (five molds experimental and five molds control) were prepared for each material and time period evaluated. Specimens of each respective depth were placed on a mylar-covered glass slide and placed on a dental vibrator. Immediately prior to placing the material, a pre-weighed, 50 mm long piece of unwaxed dental floss was positioned so that approximately 4 mm extended into the bottom of each empty mold (Figure 1).




Specimens were filled the entire depth of the mold (Figure 2) and the top surface was then covered with a mylar strip and a glass slide with pressure placed in order to produce a flat surface (Figure 3). The experimental specimens were photopolymerized for the manufacturer recommended curing time from the top of the mold (Figure 4) and then stored (Figure 5). Control specimens were photopolymerized from both top and bottom of the molds for the manufacturer recommended curing time.

Figure 2. Material Placement







At the appointed time of evaluation, each specimen was removed from its respective mold and was weighed to the nearest 0.001 g, the weight of the floss subtracted, and this value designated as Wt1. Each specimen was placed into a bottle containing 100 mL of distilled water in a manner in which they were suspended from contact with the bottle and completely immersed in the water. The bottle was then sealed and stored for 24 hours at 37 ± 2 °C, after which the specimens were removed from the water, rinsed with distilled water, and blot dried with absorbent paper, and then weighed. The specimens were then stored in a desiccator for 24 hours at a relative humidity of 10 ± 5 percent and then weighed to the nearest 0.001 g. The specimens were weighed daily until a constant weight was reached (\pm 0.001 g). The weight of the floss was subtracted from this value with the result designated as the final weight (Wt2). The solubility (S), expressed as a percentage of weight, was calculated using the formula:

$$S = Wt1 - Wt2 X 100$$

Mean solubility was compared with analysis of variance (ANOVA) followed by Scheffé's posthoc pairwise comparison test to determine if significant differences existed as per curing method and depth of cure at the 0.05 level of significance.

Hardness: Delrin molds were used with a diameter of 5.0 mm and depths of 2, 3, 4, and 5 mm. Five molds of each depth were prepared for each time period evaluated. Specimens of each respective depth were placed on a mylar-covered glass slide and filled the entire depth of the cylinder. The top surface was then covered with a mylar strip and a glass slide with pressure placed in order to produce a flat surface. The specimens were photopolymerized from the top orientation for the time recommended by the manufacturer and then stored in dark conditions as previously specified.

At the appointed time of evaluation Knoop Hardness numbers (KHN) were determined with a hardness tester (M-400-G2, LECO, St Joseph, MI) using a 100 g load and a dwell time of 10 s. Three replicate Knoop Hardness (KH) measurements were made at the center of the top and bottom surfaces at the appointed time period for evaluation. The three values from each specimen surface were averaged to yield a top and bottom summary value for each specimen. The summary values for each of the specimen surfaces were then averaged to yield a mean and standard deviation. Values were compared with ANOVA followed by Scheffé's post-hoc pairwise comparison test to determine if significant differences existed among the depths of cure at the 0.05 level of significance.

Differential Scanning Calorimetry: RMGI materials were prepared as per manufacturer's recommendations and placed either in 40 μ I differential scanning calorimetry (DSC) aluminum crucibles (Mettler Toledo GmbH, Schwerzenbach CH) for the 2 mm depth sample, or 100 μ I aluminum DSC crucibles for the 3 and 4 mm depth samples. Five mm depth samples were not prepared for the DSC evaluation due to the lack of a suitable crucible. Each respective material was polymerized with a visible light-curing unit for the manufacturer recommended time, weighed on an analytical balance (AG245, Mettler Toledo GmbH), and then stored in dark

conditions as previously specified. An unfilled reference aluminum DSC crucible of both 40and 100 μ l was also stored with each sample group. At the appointed time of evaluation, specimens were weighed, placed into a DSC unit (DSC822°, Mettler Toledo GmbH), and were first subjected to a specific heat capacity determination thermal scan from 20 °C to 60 °C at a rate of 10 °C per minute as per the protocol of Khalil and Atkins. 61 Each scan was preceded and followed by a five-minute isothermal 20 °C period. Each specimen was subjected to three heat capacity determination scans with the weight of the specimen registered between each scan. The heat capacity of each scan was determined as the mean of the specific heat values calculated by the DSC software (Star^e, Version 8.10, Mettler Toledo GmbH) over the temperature range. The average of the three specific heat determination scans was recorded as the mean with standard deviation for each sample, as well as the weight loss using the analytical balance. After the three specific heat capacity determination scans, the specimen was then subjected to a thermal scan of 37 °C to 240 °C⁶¹ with the thermal flow characteristics of the material analyzed with thermal analysis software (Stare, Version 8.10, Mettler Toledo GmbH). Mean specific heat capacity, thermal scan findings, and weight loss was determined for each material thickness and compared with ANOVA and Scheffe' post hoc multiple comparison tests.

CHAPTER THREE: RESULTS

Solubility

The complete solubility data can be found in Tables 64-93 in Appendix A.

Time/Thickness	EXPERIMENTAL		CONTROL		t-test
	Mean	SD	Mean	SD	p =
IMMEDIATE					
2 mm	4.19	0.28	3.91	0.35	0.193
3 mm	4.70	0.47	4.27	0.40	0.160
4 mm	7.60	1.89	5.05	0.17	0.017
5 mm	5.46	0.24	4.76	0.28	0.003
24 Hour					
2 mm	5.58	0.53	5.69	0.28	0.682
3 mm	5.45	0.36	5.65	0.50	0.511
4 mm	5.79	0.34	5.04	0.46	0.020
5 mm	4.83	0.29	4.85	0.29	0.891
ONE WEEK					
2 mm	5.41	0.44	6.06	0.41	0.041
3 mm	5.47	0.61	5.52	0.41	0.874
4 mm	5.20	0.41	5.60	0.39	0.167
5 mm	5.61	0.47	6.01	0.73	0.337
ONE MONTH					
2 mm	7.05	0.78	6.52	0.31	0.195
3 mm	7.79	0.65	7.49	0.79	0.540
4 mm	6.32	0.33	7.23	0.89	0.060
5 mm	7.84	0.68	8.69	0.68	0.086
THREE MONTH					
2 mm	4.93	0.55	4.51	0.15	0.140
3 mm	5.29	0.48	5.43	0.50	0.668
4 mm	5.97	0.36	5.17	0.75	0.062
5 mm	5.75	0.52	5.90	0.84	0.743

Mean solubility results for the Vitremer samples can be seen in Table 1. Multivariate ANOVA analysis of the Vitremer solubility results as a function of specimen cure method, thickness, and storage time identified significant differences (p < 0.0001) in solubility for all factors, with an interaction

noted between the factors storage time and thickness (p < 0.0001) for both curing methods.

Within the experimental group, solubility at the one-month storage time was significantly greater (p < 0.0001) than all of the other storage times. Considering experimental group thickness, the 4 mm samples demonstrated a trend with the most solubility, exhibiting significantly more solubility than the 2 mm samples (p = 0.001) but did not exhibit significant difference than the 3 mm or 5 mm samples.

The Vitremer control group similarly exhibited significantly more solubility at the one-month storage time (p <0.0001) but there were differences noted between the other thicknesses as well. The immediate Vitremer control samples demonstrated significantly less (p \leq 0.001) solubility than the other storage times, while there were no significant differences identified between the 24-hour and one week and the one-week and one-month storage times. *Post* hoc analysis of the Vitremer control samples as a function of thickness showed that the 5 mm samples exhibited significantly more solubility (p < 0.0001) than the 2 mm samples, with no difference in solubility noted between the 2 mm, 3 mm, and 4 mm samples (p = 0.190) and the 3 mm, 4 mm, and 5 mm samples (p = 0.051).

Independent t-test analysis comparing the Vitremer experimental and control samples at equal sample thickness and storage times revealed, except for four instances, that there was no overall difference in solubility between the two curing methods. The immediate 4 mm and 5 mm as well as the 24-hour 4 mm control samples were found to exhibit significantly less solubility than the experimental samples for the same thickness and storage time. Paradoxically, the one-week 2 mm control samples were found to exhibit significantly more solubility than the experimental samples of the same thickness and storage time.

Time/Thickness	EXPERIM	MENTAL	CONT	ROL	t-test	Fuji II LC solubility
Tille/Tillckiless	Mean	SD	Mean	SD	p =	l aji ii 20 oolubiiity
IMMEDIATE					Sec. 17. 14.	samples can be
2 mm	3.89	0.31	2.61	1.03	0.029	
3 mm	3.82	1.80	1.80	0.46	<0.0001	seen in Table 2.
4 mm	4.77	1.87	2.15	0.12	0.014	N. 4. 101 1 - 1 -
5 mm	4.79	0.20	2.13	0.09	<0.0001	Multivariate
24 HOURS					A SECTION	ANOVA analysis o
2 mm	3.85	0.18	3.29	0.09	<0.0001	ANOVA arialysis o
3 mm	3.24	0.11	3.18	0.06	0.288	the results as a
4 mm	3.22	0.16	3.03	0.07	0.039	
5 mm	3.20	0.10	2.91	0.06	0.001	function of
ONE WEEK						7
2 mm	4.90	1.65	5.12	0.20	0.767	specimen cure
3 mm	4.64	0.12	4.59	0.10	0.473	
4 mm	4.79	0.20	4.69	0.22	0.485	method identified
5 mm	4.74	0.68	4.58	0.12	0.619	significant
ONE MONTH						Sigrimoant
2 mm	5.11	0.24	4.89	0.24	0.181	differences (p <
3 mm	5.17	0.13	5.01	0.11	0.076	
4 mm	5.32	0.41	4.96	0.16	0.108	0.0001) in solubilit
5 mm	5.25	0.22	4.80	0.15	0.005	
THREE MONTHS						both within and
2 mm	4.72	0.58	4.04	0.13	0.034	batana a tha tara
3 mm	5.56	0.86	4.11	0.37	0.008	between the two
4 mm	4.72	0.29	4.14	0.33	0.018	groups. In
5 mm	4.39	0.17	3.84	0.20	0.002	groups. III
= 5 values presented in	red are st	atistically	significan	t p < 0.0	5	addition, there

were differences noted between the two curing method groups when analyzed within the groups with two-way ANOVA as a function of time and thickness.

For the Fuji II LC experimental solubility group there was an identified interaction between the two factors of storage time and thickness (p = 0.044) while Scheffe *post hoc* analysis did not identify a significant difference in solubility when considering thickness (p = 0.960) for the experimental group. A significant difference (p < 0.0001) was identified within the experimental group for the factor storage time, with the one-month samples indicated with significantly more solubility (p < 0.002) than the immediate and 24-hour samples. However, the one-month samples exhibited statistically similar solubility with the one-week and three-month samples (p = 0.334).

The Fuji II LC control solubility group was identified as having significant differences within the group for both sample thickness and storage time (p < 0.001) with no interaction between the factors of thickness and storage time (p = 0.079). The Fuji II LC 2 mm samples were found to exhibit significantly more solubility than 3 mm and 5 mm (p < 0.039) but were statistically similar to the 4 mm samples (p = 0.157). *Post hoc* analysis for the Fuji II LC control samples as a function of storage time indicated that the one-month samples exhibited significantly more solubility (p < 0.0001) than the other storage times except for the one-week samples (p = 0.556).

Independent t-test analysis comparing the Fuji II LC experimental and control samples at equal sample thickness and storage times revealed that over half of the control samples exhibited significantly less solubility than the experimental samples. This was most frequently noted with the immediate, 24-hour, and three-month groups. There were no significant differences in solubility noted between paired samples within the one-week group, and the majority of the one-month group.

The results of the Photac-Fil Quick solubility can be seen in Table 3. Multivariate ANOVA analysis of the results as a function of specimen cure method identified significant differences (p < 0.0001) in solubility both within and between the two groups. Additionally, there were differences noted between the two groups when analyzed within each group using two-way ANOVA analysis as a function of time and thickness.

Similar to the Fuji II LC samples, the Photac-Fil Quick experimental solubility group exhibited an interaction between the two factors of specimen thickness and storage time (p < 0.0001) but a significant difference was not found in solubility between the different sample thicknesses (p = 0.690). A significant difference in solubility was noted when analyzed per storage time (p < 0.0001) with the one-month experimental samples demonstrating significantly more solubility

T! /Tl. ! . !	EXPERIM	MENTAL	CONT	TROL	t-test	:
Time/Thickness	Mean	SD	Mean	SD	p =	immediate
IMMEDIATE						and three
2 mm	3.89	1.10	4.00	0.51	0.842	
3 mm	4.96	1.38	3.60	0.33	0.0642	month
4 mm	5.52	0.83	3.65	0.16	0.001	
5 mm	5.64	0.70	3.29	0.17	<0.0001	samples
24 Hour						0.0001) h
2 mm	7.96	0.43	6.01	0.33	<0.0001	0.0001) b
3 mm	6.58	0.72	5.45	0.55	0.024	statistical
4 mm	6.00	0.46	5.28	0.19	0.012	Statistical
5 mm	5.75	0.37	5.25	0.85	0.257	similar to
ONE WEEK						
2 mm	6.07	1.07	5.89	0.44	0.719	24-hour a
3 mm	6.80	0.20	5.53	0.58	0.002	
4 mm	6.15	0.35	5.91	0.37	0.328	one-week
5 mm	6.31	0.20	5.71	0.41	0.019	a a manda a
ONE MONTH						samples
2 mm	7.17	0.50	6.24	0.41	0.013	0.660).
3 mm	6.30	0.34	5.69	0.14	0.006	0.000).
4 mm	6.78	0.22	5.62	0.14	<0.0001	The Phota
5 mm	6.31	0.52	5.46	0.13	<0.0001	
THREE MONTH						Fil Quick
2 mm	5.44	0.27	5.26	0.25	0.332	
3 mm	5.21	0.28	4.86	0.37	0.129	control
4 mm	5.56	0.41	5.12	0.21	0.066	a a luda ilia.
5 mm	5.46	0.13	5.09	0.46	0.204	solubility
						group wa
= 5 values presented in						'

having significant differences within the group for both sample thickness and storage time (p < 0.001) but without an identified interaction between the factors of thickness and storage time (p = 0.523). The Photac-Fil Quick 2 mm control samples were found to curiously exhibit significantly more solubility than the other thicknesses (p < 0.018) with no significant difference noted between the other thicknesses (p = 0.590). *Post hoc* analysis of the Photac-Fil Quick control samples as a function of storage time indicated that the 24-hour, one-week, and one-month samples were statistically similar (p = 0.359) and exhibited significantly more solubility than the immediate samples (p < 0.0001) and the three-month samples (p = 0.034).

Independent t-test analysis of the Photac-Fil Quick experimental and control samples at equal sample thickness and storage times revealed that over half of the control samples exhibited significantly less solubility than the experimental samples. This was most frequently noted with the 24-hour and one-month groups. There were no significant differences in solubility noted between paired samples between the samples at three months, with mixed results with the immediate and one-week groups.

Hardness

The complete data concerning hardness testing may be observed in Tables 94 – 105 in Appendix B.

The mean KHN value results for Vitremer can be observed in Table 4. The data was analyzed

	24 Hours		1 Month	
2mm Top	24.5	22.4	27.1	27.1
sd	0.4	0.3	0.4	0.5
2mm Bottom	23.1	19.9	24.4	24.9
sd	0.3	0.2	0.4	0.6
3mm Top	24.9	22.4	26.7	23.3
sd	0.3	0.3	0.3	0.5
3mm Bottom	22.7	18.0	23.9	20.0
sd	0.4	0.1	0.6	0.5
4mm Top	23.9	22.1	27.7	21.7
sd	0.4	0.2	0.4	0.4
4mm Bottom	20.5	17.3	21.4	16.6
sd	0.7	0.2	0.2	0.2
5mm Top	24.1	22.8	27.8	24.5
sd	0.6	0.2	0.5	0.4
5mm Bottom	17.1	15.4	18.5	16.8
sd	0.6	0.2	0.2	0.5

for possible differences in hardness of different layers as a function of thickness and if any possible changes over the stated storage periods. In the comparison of top surface hardness, ANOVA did not identify a significant difference (p = 0.47) in the top surface hardness values at the 24-hour storage time, but identified differences at one week (p = .0054), one-month (p =

0.0064), and three months (p = 0.01) At the one-month storage time, the only significant differences were noted between the 3 mm and 5 mm top surfaces (p = 0.045), while for the three-month aged specimens Scheffe multiple comparison test identified a difference only between the 2 mm and 4 mm top surfaces (p = 0.0107).

For the Vitremer top surface hardness compared over the three-month storage times, ANOVA identified a significant difference (p < 0.0001) within the 2 mm samples (p < 0.0001), the 3 mm samples (p = 0.0002), the 4 mm samples (p = 0.0004), and the 5 mm samples (p < 0.0001).

Table 5. Scheffe Multi Comparisons For Vitremer Top Surface **Hardness Over Storage Time Thickness** Comparisons 4 mm 5 mm 2 mm 3 mm p value 24-Hour vs. 1 Week 0.0098 0.0179 0.3286 0.0469 < 0.0001 24-Hour vs. I Month 0.2829 0.2569 0.067 0.5458 0.2663 0.8408 24-Hour vs. 3 Months 0.1663 0.0003 0.0018 < 0.0001 0.0002 1 Week vs. 1 Month 1 Week vs. 3 Months 0.0001 0.2288 0..9989 0.0086 < 0.0001 0.0208 0.0014 1 Month vs. 3 Months 0.9873 Results shaded red are considered significant at a 95% confidence

interval

The results of the post hoc multi comparison tests can be seen in Table 5.

Concerning the analysis of the Vitremer bottom KHN

over the storage times,

one-way

ANOVA
identified a
significant
difference

within the 2

mm (p <

0.0001), 3 mm (p < 0.0001),

Table 6. Scheffe Multi Comparisons For Vitremer Bottom Surface **Hardness Over Storage Time Thickness** 4 mm 5 mm 3 mm Comparisons 2 mm p value 0.258 24-Hour vs. 1 Week < 0.0001 < 0.0001 0.016 0.425 0.834 0.745 24-Hour vs. I Month 0.867 24-Hour vs. 3 Months 0.108 0.068 0.024 0.834 0.042 < 0.0001 0.003 < 0.0001 1 Week vs. 1 Month 0.997 0..713 < 0.0001 0.007 1 Week vs. 3 Months 0.003 0.004 0.282 1 Month vs. 3 Months 0.375 n = 5Results shaded red are considered significant at a 95% confidence interval

4 mm (p < 0.0001) and 5 mm (p = 0.033) sample groups. The results of the *post hoc* multi comparison tests can be seen in Table 6.

The mean bottom/top KHN ratios of the Vitremer samples by thickness and time can be seen in Table 7. For each thickness analyzed over the storage time, ANOVA did not identify a difference within the 2 mm, 4 mm, or 5 mm groups over the storage time. However, a difference (p < 0.0001) was identified within the 3 mm samples with the one week samples

	24 Hours	1 Week	1 Month	3 Months
2mm	93.0	88.7	89.8	93.4
sd	1.1	1.3	1.1	6.5
3mm	90.5	79.8	89.6	86.4
sd	1.7	1.4	2.4	3.5
4mm	84.7	77.7	77.2	80.1
sd	5.9	1.8	1.8	9.8
5mm	71.5	67.6	66.3	67.0
sd	9.9	1.2	1.3	5.6

identified as having a significantly lower mean bottom/top ratio than the other storage time periods (p < 0.0001).

When the bottom/top ratio of the entire Vitremer group was

compared with two-way ANOVA, a significant difference (p < 0.0001) was identified within the group with no interaction found between time and thickness (p = 0.258). Scheffe post hoc analysis identified that all thicknesses exhibited a significantly different mean bottom/top KHN ratio (p < 0.0001) from each other. When the KHN bottom/top ratio was analyzed as a function of time, Scheffe identified a significant difference between the 24-hour groups and one-week (p = 0.001) with no significant difference noted between any of the other storage times.

Table 8. Fuji	i II LC Me	an Hardi	ness Res	ults (KHN)	
	24 Hours	1 Week	1 Month	3 Months	
2mm Top	20.8	20.7	20.0	20.6	
sd	0.3	0.3	0.5	1.0	
2mm Bottom	18.4	18.8	18.2	18.6	
sd	0.3	0.5	0.3	0.5	
3mm Top	20.4	20.6	20.1	20.6	
sd	0.3	0.3	0.6	0.1	
3mm Bottom	15.2	18.2	16.7	19.4	
sd	0.2	0.3	0.2	0.1	
4mm Top	20.2	21.3	19.6	21.3	
sd	0.1	0.4	0.4	0.3	
4mm Bottom	11.6	15.4	11.8	13.1	
sd	0.3	0.6	0.5	0.5	
5mm Top	20.2	21.3	20.1	21.0	
sd	0.3	0.2	0.4	0.2	
5mm Bottom	9.7	12.8	10.9	10.6	
sd	0.4	0.5	0.4	0.4	
n = 5	Astronomy and the supply				

The results of the Fuji II LC hardness testing can be seen in Table 8. For the Fuji II LC samples, ANOVA did not identify a significant difference for the top surfaces at the 24-hour storage (p = 0.05), one-week (p = 0.47), one-month (p = 0.51), or the three-month (p = 0.24) storage times.

Analysis of the Fuji II LC top surfaces over the storage times did not identify a significant change in the top surface hardness with the 2 mm samples (p = 0.1320), 3 mm samples (p = 0.3659), nor the 5 mm samples (p = 0.0671). ANOVA identified a significant change within the 4 mm samples (p = 0.001) with Scheffe finding a difference between the one-week and one-month samples (p = 0.0009) and the one-month and three-month samples (p = 0.0004). For the Fuji II LC bottom surfaces over the different storage times, ANOVA did not identify a difference within the 2 mm samples (p = 0.372), but identified differences with bottom surface hardness with the 3 mm (p < 0.0001), 4 mm (p < 0.0001), and 5 mm (p < 0.0001) samples. The results of the *post hoc* analysis can be seen in Table 9.

	Thickness				
Comparisons	2 mm	3 mm	4 mm	5 mm	
	p value				
24-Hour vs. 1 Week	0.793	< 0.0001	0.001	< 0.0001	
24-Hour vs. I Month	0.919	0.001	0.629	0.254	
24-Hour vs. 3 Months	0.948	< 0.0001	0.109	0.136	
1 Week vs. 1 Month	0.425	0.001	0.010	0.002	
1 Week vs. 3 Months	0.981	0.025	0.109	0.004	
1 Month vs. 3 Months	0.648	< 0.0001	0.629	0.981	

The Fuji II LC
bottom/top KHN
ratio results can
be observed in
Table 10. When
analyzed per
thickness as a
function of
storage time,
ANOVA did not

-	24 Hours	1 Week	1 Month	3 Months
2mm	88.4	92.7	91.1	91.3
sd	1.9	9.3	2.7	4.2
3mm	74.7	88.5	82.8	94.0
sd	1.4	2.3	3.8	1.8
4mm	55.8	72.1	62.5	62.2
sd	6.1	2.0	5.8	6.2
5mm	47.6	60.1	52.9	51.4
sd	1.5	3.9	5.7	3.8

interval

identify a significant
different in change of the
KHN ratio for the 2 mm
samples (p = 0.641).
ANOVA did identify
significant changes in the
3 mm (p < 0.0001), 4 mm

(p = 0.002), and 5 mm (p = 0.001) groups. The results of the *post hoc* multi comparison test can be seen in Table 11.

	Thickness			
Comparisons	2 mm	3 mm	4 mm	5 mm
		p	value	
24-Hour vs. 1 Week	0.658	< 0.0001	0.002	0.002
24-Hour vs. I Month	0.885	0.001	0.298	0.274
24-Hour vs. 3 Months	0.857	< 0.0001	0.348	0.541
1 Week vs. 1 Month	0.972	0.019	0.082	0.077
1 Week vs. 3 Months	0.982	0.027	0.067	0.028
1 Month vs. 3 Months	0.99	< 0.0001	0.99	0.954

	24 Hours	1 Week	1 Month	3 Months
2mm Top	11.5	13.7	17.5	17.0
sd	0.4	0.3	0.5	0.5
2mm Bottom	9.4	12.0	16.3	15.4
sd	0.2	0.2	0.1	0.4
3mm Top	11.6	13.6	18.1	16.8
sd	0.4	0.2	0.6	0.2
3mm Bottom	9.2	11.3	15.6	13.5
sd	0.1	0.3	0.2	0.5
4mm Top	12.0	12.8	18.2	16.9
sd	0.2	0.3	0.2	0.1
4mm Bottom	8.4	10.5	14.5	11.2
sd	0.5	0.2	0.4	0.1
5mm Top	11.6	12.7	17.9	17.3
sd	0.3	0.2	0.6	0.4
5mm Bottom	8.7	10.2	11.2	10.7
sd	0.4	0.1	0.2	0.2

The results of the Photac-Fil Quick Knoop hardness testing can be seen in Table 12. For the analysis of the top surface KHN of the Photac-Fil Quick samples over the different storage times, no difference was noted between the top surface hardness at the 24-hour storage (p = 0.73), onemonth (p = 0.50), nor three-month (p = 0.05) storage times. However, a significant difference in top surface hardness was noted

(p = 0.008) within the one-week group, with a significant difference found between the 2 mm and 5 mm top surfaces (p = 0.0213). Regarding the changes of top surface hardness as a function of time, ANOVA identified significant differences for all of the thicknesses: 2 mm (p <

0.0001), 3 mm (p = 0.0014), 4 mm (p = 0.0004), and 5 mm (p < 0.0001). Results of the Scheffe post hoc tests can be seen in Table 13.

	Thickness			
Comparisons	2 mm	3 mm	4 mm	5 mm
		р	value	
24-Hour vs. 1 Week	< 0.0001	0.5369	0.3286	0.0469
24-Hour vs. I Month	< 0.0001	0.0317	0.067	< 0.0001
24-Hour vs. 3 Months	< 0.0001	09894	0.2663	0.8408
1 Week vs. 1 Month	< 0.0001	0.0019	0.0018	< 0.0001
1 Week vs. 3 Months	< 0.0001	0.3654	09989	0.0086
1 Month vs. 3 Months	0.2031	0.0586	0.0014	< 0.0001

Table 14. Scheffe Multi Comparisons For Photac-Fil Quick Bottom Surface Hardness Over Storage Time **Thickness** Comparisons 3 mm 4 mm 5 mm 2 mm p value 0.366 24-Hour vs. 1 Week < 0.0001 < 0.0001 0.001 24-Hour vs. I Month < 0.0001 < 0.0001 < 0.0001 0.028 0.032 24-Hour vs. 3 Months < 0.0001 < 0.0001 < 0.0001 0.491 1 Week vs. 1 Month < 0.0001 < 0.0001 < 0.0001 1 Week vs. 3 Months < 0.0001 0.629 0.524 < 0.0001 1 Month vs. 3 Months 0.337 < 0.0001 < 0.0001 0.99 Results shaded red are considered significant at a 95% confidence interval

With analysis of the Photac-Fil Quick bottom surface hardness over the storage times, ANOVA also identified a significant change in KHN with all thicknesses: 2 mm (p <0.0001), 3 mm (p < 0.0001), 4mm (p <0.0001), and 5

mm (p = 0.01).

Results of the Scheffe *post hoc* multi comparison test for each of these groups can be seen in Table 14.

The mean bottom/top KHN ratio for the Photac-Fil Quick groups over the storage time can be seen in Table 15. ANOVA analysis of each thickness over the different storage times identified a significant difference within the 2 mm (p = 0.02), 3 mm (p = 0.055), 4 mm (p < 0.0001), and 5

mm samples (p = 0.003). When the change in bottom/top KHN ratios were analyzed by two-way ANOVA as a function of specimen thickness and storage time, a significant difference was identified (p < 0.0001) with a significant interaction also identified between the effects of time

	24 Hours	1 Week	1 Month	3 Months
2mm	82.8	87.9	92.1	91.6
sd	5.9	4.7	4.3	2.9
3mm	78.8	83.6	86.8	80.9
sd	5.5	4.7	4.2	2.5
4mm	70.8	81.3	80.6	65.6
sd	2.6	2.9	4.6	4.0
5mm	75.4	77.8	61.7	62.7
sd	10.5	6.3	4.1	6.0

and thickness (p < 0.0001). Results of the Scheffe *post* hoc analysis can be seen in Table 16.

	Thickness			
Comparisons	2 mm	3 mm	4 mm	5 mm
A	Selection .	р	value	
24-Hour vs. 1 Week	0.402	0.415	0.003	0.961
24-Hour vs. I Month	0.042	0.075	0.006	0.056
24-Hour vs. 3 Months	0.058	0.901	0.201	0.083
1 Week vs. 1 Month	0.569	0.728	0.992	0.021
1 Week vs. 3 Months	0.663	0.810	< 0.0001	0.032
1 Month vs. 3 Months	0.999	0.249	< 0.0001	0.997

Differential Scanning Calorimetry

Time Period	2mm	3mm	4mm
Immediate	1.21 (0.06)	1.06 (0.03)	0.96 (0.04)
24 Hours	1.22 (0.09)	1.05 (0.05)	0.92 (0.09)
1 Week	1.28 (0.02)	1.10 (0.04)	0.98 (0.03)
1 Month	1.26 (0.02)	1.12 (0.02)	1.04 (0.05)
3 Months	1.36 (0.08)	1.23 (0.11)	1.15 (0.04)

Specific Heat

Determination: Data concerning the specific heat determinations can be found in Tables 106 – 108,

Appendix C. The mean specific heat determination results for Vitremer as per thickness and time period can be seen in Table 17. For all storage time periods, the 2-mm thickness samples exhibited higher specific heat than the 3- or 4-mm thick samples. Two-way ANOVA analysis of the Vitremer specific heat results as a function of thickness over storage time identified

	Thickness				
Comparisons	2 mm	3 mm	4 mm		
	p value				
Immediate vs. 24-Hours	0.99	0.998	0.911		
Immediate vs. I Week	0.508	0.941	0.973		
Immediate vs. 1 Month	0.870	0.673	0.328		
Immediate vs. 3 Months	0.025	0.008	0.001		
24 Hours vs. 1 Week	0.621	0.822	0.598		
24 Hours vs. 1 Month	0.935	0.483	0.078		
24 Hours vs. 3 Months	0.038	0.004	< 0.0001		
1 Week vs. 1 Month	0.967	0.977	0.656		
1 Week vs. 3 Months	0.481	0.046	0.003		
1 Month vs. 3 Months	0.180	0.147	0.108		

significant differences in both thickness and time (p < 0.0001) without a significant interaction (p = 0.690). Accordingly, Scheffe analysis found that over the storage times there were significant differences between in Vitremer specific heat values between all thicknesses (p <

0.0001) with significant differences in specific heat noted between the three-month and all of the other times (p < 0.0001). Results of the Scheffe *post hoc* testing for the Vitremer groups can be seen in Table 18. One-way ANOVA of the sample specific heat per thickness over the different storage times identified a significant increase within the 2 mm (p = 0.009), 3 mm (p = 0.001), and 4 mm (p < 0.0001) samples.

Time Period	2mm	3mm	4mm
Immediate	1.04 (0.02)	0.93 (0.03)	0.83 (0.03)
24 Hours	1.21 (0.01)	0.97 (0.02)	0.91 (0.02)
1 Week	1.17 (0.02)	1.07 (0.05)	1.02 (0.04)
1 Month	1.16 (0.08)	0.97 (0.03)	0.96 (0.03)
3 Months	1.05 (0.17)	0.99 (0.03)	0.97 (0.05)

The mean specific heat
determination results for Fuji
II LC as per thickness and
storage time can be seen in
Table 19. Two-way ANOVA

analysis of the Fuji II LC specific heat results as a function of thickness and storage time identified significant differences with both time and thickness (p < 0.0001) and identified a significant interaction between the two factors (p = 0.003). The results of the Scheffe multi comparison test concerning sample storage time can be seen in Table X. One-way ANOVA of the different thicknesses over the different storage times identified a significant change within the 2 mm (p = 0.013), 3 mm (p < 0.0001), and 4 mm samples (p < 0.0001). Results of the Scheffe *post hoc* analysis for each of the thicknesses can be seen in Table 20.

Table 20. Scheffe Multi Co Specific Heat per Thickne			
		Thickness	
Comparisons	2 mm	3 mm	4 mm
		p value	
Immediate vs. 24-Hours	0.032	0.367	0.020
Immediate vs. I Week	0.132	< 0.0001	< 0.0001
Immediate vs. 1 Month	0.184	0.562	< 0.0001
Immediate vs. 3 Months	0.999	0.137	< 0.0001
24 Hours vs. 1 Week	0.953	0.006	0.006
24 Hours vs. 1 Month	0.896	0.997	0.424
24 Hours vs. 3 Months	0.052	0.976	0.305
1 Week vs. 1 Month	0.999	0.003	0.239
1 Week vs. 3 Months	0.196	0.023	0.342
1 Month vs. 3 Months	0.265	0.884	0.999
n = 5			
Results shaded red are con confidence interval	nsidered sig	gnificant at a	95%

Time Period	2mm	3mm	4mm
Immediate		0.93 (0.05)	
24 Hours	1.19 (0.04)	1.11 (0.03)	0.89 (0.06
1 Week		0.93 (0.03)	
1 Month		1.00 (0.04)	
3 Months	1.16 (0.05)	1.02 (0.01)	0.97 (0.03

The mean specific heat
determination results for
Photac-Fil Quick can be seen
in Table 21. Two-way
ANOVA analysis of the
specific heat changes as a

function of both time and thickness identified significant specific heat changes with the Photac-Fil Quick samples for both time and thickness (p < 0.0001) with a significant interaction between the two factors (p = 0.002). The results of the Scheffe *post hoc* analysis can be seen in Table 22. One-way ANOVA did not identify a significant specific heat difference with the 2 mm samples over the storage time (p = 0.65), but identified a significant difference with the 3 mm (p < 0.0001) and 4 mm samples (p = 0.001).

	Thickness				
Comparisons	2 mm	3 mm	4 mm		
	p value				
Immediate vs. 24-Hours	0.355	< 0.0001	0.410		
Immediate vs. I Week	0.994	0.999	0.011		
Immediate vs. 1 Month	0.999	0.091	0.084		
Immediate vs. 3 Months	0.634	0.009	0.003		
24 Hours vs. 1 Week	0.182	< 0.0001	0.377		
24 Hours vs. 1 Month	0.372	0.002	0.888		
24 Hours vs. 3 Months	0.989	0.025	0.136		
1 Week vs. 1 Month	0.992	0.141	0.888		
1 Week vs. 3 Months	0.391	0.015	0.972		
1 Month vs. 3 Months	0.655	0.825	0.554		

Differential Scanning Calorimetry Thermal Scan:

The complete data detailing the differential scanning calorimetry thermal results and peak temperatures can be found in Appendix D, Tables 109 – 111.

The mean results for the Vitremer thermal scan for each thickness and storage time period can be seen in Table 23. The enthalpy results are listed as negative because it was an endothermic phenomenon.

Time Period	2mm	sd	3mm	sd	4mm	sd
Immediate	-2.97	0.75	-40.78	6.62	-47.26	8.52
24 Hours	-27.39	8.27	-47.95	4.12	-54.99	5.62
1 Week	-43.33	6.51	-64.18	6.42	-67.69	8.51
1 Month	-49.72	4.26	-88.69	10.39	-103.45	16.39
3 Months	-90.03	13.82	-103.56	20.73	-104.17	35.08

	Thickness				
Comparisons	2 mm	3 mm	4 mm		
	p value				
Immediate vs. 24-Hours	0.003	0.935	0.985		
Immediate vs. I Week	< 0.0001	0.114	0.652		
Immediate vs. 1 Month	< 0.0001	< 0.0001	0.008		
Immediate vs. 3 Months	< 0.0001	< 0.0001	0.007		
24 Hours vs. 1 Week	0.077	0.417	0.912		
24 Hours vs. 1 Month	0.007	0.002	0.025		
24 Hours vs. 3 Months	< 0.0001	< 0.0001	0.023		
1 Week vs. 1 Month	0.808	0.089	0.149		
1 Week vs. 3 Months	< 0.0001	0.002	0.136		
1 Month vs. 3 Months	< 0.0001	0.502	0.999		

Two-way ANOVA analysis of the Vitremer enthalpy results as a factor of both thickness and storage time identified significant changes concerning both thickness (p < 0.0001) and time (p < 0.0001) but did not identify any significant interactions between the two factors (p = 0.144). The Scheffe *post hoc* analysis can be seen in

Table 24. Using one-way ANOVA to analyze thermal changes of each thickness as a function of time, a significant difference was identified within the 2 mm (p < 0.0001), 3 mm (p < 0.0001), and 4 mm (p < 0.0001) samples as the respective samples aged.

The results of the mean peak temperature associated with the endothermic peak of the thermal scan for the Vitremer samples can be seen in Table 25. Two-way ANOVA analysis of the data for both thickness and storage time identified significantly different peak temperatures for both sample thickness and storage time (p < 0.0001) and indicated significant interactions between the two factors (p < 0.0001).

Time Period	2mm	sd	3mm	sd	4mm	sd
Immediate	182.94	4.80	157.83	1.99	150.72	9.17
24 Hours	108.02	2.91	130.81	8.85	133.52	2.53
1 Week	110.68	3.25	134.56	5.79	133.71	2.84
1 Month	114.63	2.65	132.05	0.96	131.49	1.21
3 Months	121.54	3.19	129.78	1.68	130.69	4.52

One-way ANOVA of the data identified significant differences for the individual thicknesses over storage time revealed significant mean peak temperature differences between the 2 mm and 3 mm (p = 0.012) as well as between the 2 mm and 4 mm specimens (p = 0.006) but not between the 3 mm and 4mm specimens (p = 0.924). Results of the Scheffe multicomparison tests can be seen in Tables 26 and 27.

	Thickness					
Comparisons	2 mm	3 mm	4 mm			
	p value					
Immediate vs. 24-Hours	< 0.0001	0.935	0.002			
Immediate vs. I Week	< 0.0001	0.114	0.002			
Immediate vs. 1 Month	< 0.0001	< 0.0001	0.001			
Immediate vs. 3 Months	< 0.0001	< 0.0001	< 0.0001			
24 Hours vs. 1 Week	0.824	0.879	0.999			
24 Hours vs. 1 Month	0.093	0.998	0.986			
24 Hours vs. 3 Months	< 0.0001	0.999	0.953			
1 Week vs. 1 Month	0.525	0.969	0.980			
1 Week vs. 3 Months	0.002	0.752	0.941			
1 Month vs. 3 Months	0.074	0.979	0.999			

		St	orage Time		
Comparisons	Immediate	24 Hour	One Week	One Month	Three Months
		米里阿里斯斯	p value		
2 mm vs. 3 mm	< 0.0001	< 0.0001	< 0.0001	< 0.0001	0.012
2 mm vs. 4 mm	< 0.0001	< 0.0001	< 0.0001	< 0.0001	0.006
3 mm vs. 4 mm	0.238	0.790	0.958	0.890	0.924

The mean enthalpy results for the Fuji II LC samples can be observed in Table 28.

Time Period	2mm	sd	3mm	sd	4mm	sd
Immediate	-15.37	12.64	-33.53	3.88	-34.48	2.48
24 Hours	-36.87	4.77	-34.65	5.90	-38.81	3.43
1 Week	-52.41	3.35	-80.03	1.46	-72.90	1.98
1 Month	-46.12	10.50	-76.58	5.54	-78.28	4.83
3 Months	-47.66	7.24	-74.94	2.71	-81.54	5.63

	Thickness					
Comparisons	2 mm	3 mm	4 mm			
	p value					
Immediate vs. 24-Hours	0.014	0.998	0.658			
Immediate vs. I Week	< 0.0001	< 0.0001	< 0.0001			
Immediate vs. 1 Month	< 0.0001	< 0.0001	< 0.0001			
Immediate vs. 3 Months	< 0.0001	< 0.0001	< 0.0001			
24 Hours vs. 1 Week	0.116	< 0.0001	< 0.0001			
24 Hours vs. 1 Month	0.570	< 0.0001	< 0.0001			
24 Hours vs. 3 Months	0.421	< 0.0001	< 0.0001			
1 Week vs. 1 Month	0.843	0.885	0.459			
1 Week vs. 3 Months	0.936	0.589	0.080			
1 Month vs. 3 Months	0.999	0.989	0.842			

Two-way ANOVA of the
Fuji II LC enthalpy results
identified a significant
difference for both
thickness and storage
time (p < 0.0001) as well
as an interaction between
the effects of time and
thickness (p < 0.0001).
One-way ANOVA of the
Fuji II LC samples as a
function of individual
thickness over time

revealed a significant difference with the 2 mm (p < 0.0001), 3 mm (p < 0.0001), and 4 mm samples (p < 0.0001). Results of the Scheffe multi comparison analysis can be seen in Table 29.

Table 30. Schef Calorimetry The			or Fuji II LC	Differential	Scanning
	SALE DE LA CONTRACTION DE LA C	S	torage Time	DE PROPERTY	7. 经公司股份
Comparisons	Immediate	24 Hour	One Week	One Month	Three Months
	建筑等		p value		
2 mm vs. 3 mm	0.011	0.799	< 0.0001	< 0.0001	< 0.0001
2 mm vs. 4 mm	0.008	0.842	< 0.0001	< 0.0001	< 0.0001
3 mm vs. 4 mm	0.982	0.471	0.003	0.941	0.238
n = 5 Results shaded	red are consid	ered signific	ant at a 95%	confidence	interval

One-way ANOVA of the Fuji II thickness per time period identified a significant difference between the different thicknesses at the immediate (p < 0.004), one-week (p < 0.0001), one-month (p < 0.0001), and three-month time groups (p < 0.0001). A significant difference was not noted in the 24-hour sample group (p = 0.470). Results of the *post hoc* tests can be seen in Table 30.

The mean results of the peak temperature associated with the endothermic peak during the thermal scan for the Fuji II LC samples can be seen in Table 31.

Time Period	2mm	sd	3mm	sd	4mm	sd
Immediate	142.43	54.93	174.70	3.13	163.61	1.22
24 Hours	126.53	17.29	146.68	9.81	141.19	2.63
1 Week	126.15	12.76	147.29	1.28	137.95	1.35
1 Month	128.58	7.66	142.89	2.23	135.63	3.98
3 Months	132.87	9.81	141.71	1.00	138.86	4.03

Two-way ANOVA of the mean Fuji II LC enthalpy peak identified significant differences within the samples (p < 0.0001) for both thickness and storage time period but did not indicate a significant interaction between the two factors (p = 0.911). One-way ANOVA of the Fuji II LC

	Thickness					
Comparisons	2 mm	3 mm	4 mm			
	p value					
Immediate vs. 24-Hours	0.926	< 0.0001	< 0.0001			
Immediate vs. I Week	0.919	< 0.0001	< 0.0001			
Immediate vs. 1 Month	0.954	< 0.0001	< 0.0001			
Immediate vs. 3 Months	0.988	< 0.0001	< 0.0001			
24 Hours vs. 1 Week	0.999	0.999	0.645			
24 Hours vs. 1 Month	0.999	0.865	0.163			
24 Hours vs. 3 Months	0.998	0.706	0.860			
1 Week vs. 1 Month	0.999	0.789	0.861			
1 Week vs. 3 Months	0997	0.611	0.995			
1 Month vs. 3 Months	0.999	0.998	0.655			

separate thicknesses as a function of storage time did not identify a significant difference in mean peak temperature with the 2 mm (p = 0.865), but identified a significant change in peak temperature within the 3 mm (p < 0.0001) and 4 mm samples (p < 0.0001) as a function of storage time. Results of

the post hoc analysis can be seen in Table 32.

One-way ANOVA analysis of the Fuji II LC thermal scan peak temperature as a function of storage time did not identify a significant difference in peak temperature with the immediate group (p = 0.301), 24-hour group (p = 0.053), or the three-month samples (p = 0.116). However, a significant change of peak enthalpy temperature was found within the one-week samples (p = 0.003) and the one-month samples (p = 0.004). Results of the Scheffe *post hoc* analysis can be seen in Table 33.

	AND	St	orage Time)	
Comparisons	Immediate	24 Hour	One Week	One Month	Three Months
		The state of the s	p value		
2 mm vs. 3 mm	0.311	0.061	0.003	0.004	0.125
2 mm vs. 4 mm	0.588	0.194	0.81	0.153	0.352
3 mm vs. 4 mm	0.861	0.772	0.184	0.139	0.776

The mean enthalp	y results of the Photac-F	il Quick samples car	be seen in Table 34.
------------------	---------------------------	----------------------	----------------------

Time Period	2mm	sd	3mm	sd	4mm	sd
Immediate	-7.76	7.66	-34.08	5.84	-43.54	2.60
24 Hours	-34.97	2.40	-70.90	3.15	-58.27	8.62
1 Week	-46.79	13.72	-67.72	4.58	-70.99	5.58
1 Month	-47.73	6.28	-74.03	2.00	-77.16	3.76
3 Months	-59.08	7.40	-92.25	2.39	-99.95	2.44

Two-way ANOVA of the mean enthalpy as a function of time and samples thickness identified a significant difference for both time and thickness (p < 0.0001) and an interaction between the two factors (p = 0.002). One-way ANOVA analysis of the individual thicknesses over the storage time revealed a significant difference with the 2 mm (p < 0.0001), 3 mm (p < 0.0001), and 4 mm samples (p < 0.0001) as they aged over the storage times. Results of the Scheffe post hoc testing can be seen in Table 35. One-way ANOVA analysis of the change in thermal

	Thickness				
Comparisons	2 mm	3 mm	4 mm		
	p value				
Immediate vs. 24-Hours	0.001	< 0.0001	0.014		
Immediate vs. I Week	< 0.0001	< 0.0001	< 0.0001		
Immediate vs. 1 Month	< 0.0001	< 0.0001	< 0.0001		
Immediate vs. 3 Months	< 0.0001	< 0.0001	< 0.0001		
24 Hours vs. 1 Week	0.319	0.849	0.041		
24 Hours vs. 1 Month	0.250	0.855	0.001		
24 Hours vs. 3 Months	0.005	< 0.0001	< 0.0001		
1 Week vs. 1 Month	0.999	0.292	0.588		
1 Week vs. 3 Months	0.283	< 0.0001	< 0.0001		
1 Month vs. 3 Months	0.358	< 0.0001	< 0.0001		

significant differences
within the different
thicknesses at the
immediate (p < 0.0001),
24-hour (p < 0.0001), oneweek (p = 0.003), onemonth (p < 001), and
three-month (p < 0.0001)
storage times. The results

of the multi comparison

enthalpy for the collective

groups as a function of

storage time identified

testing can be seen in Table 36.

confidence interval

	W. F. S. C.	St	orage Time	9	(第4) 24 图
Comparisons	Immediate	24 Hour	One Week	One Month	Three Months
	to her he		p value	BULL SHOP	
2 mm vs. 3 mm	< 0.0001	< 0.0001	0.012	< 0.0001	< 0.0001
2 mm vs. 4 mm	< 0.0001	< 0.0001	0.005	< 0.0001	< 0.0001
3 mm vs. 4	0.085	0.021	0.856	0.569	0.077

Results of the Photac-Fil Quick mean peak temperature of the DSC thermal scan can be seen in Table 37. Two-way ANOVA analysis of the mean thermal peak temperature values as a function of storage time and thickness indicated significant differences within the samples (p < 0.0001) with interactions identified between the two factors (p < 0.0001).

Time Period	2mm	sd	3mm	sd	4mm	Sc
Immediate	123.27	49.51	194.33	4.16	174.59	1.0
24 Hours	134.69	23.75	149.37	5.59	154.31	3.3
1 Week	141.22	4.42	145.34	2.08	143.29	1.7
1 Month	136.69	2.91	140.14	0.42	141.02	0.9
3 Months	129.47	4.01	137.51	0.51	137.43	1.9

One-way ANOVA of the mean peak values for Photac-Fil Quick for the individual sample thicknesses over the storage times did not identify a significant difference within the 2 mm samples (p = 0.813), but found a significant change in temperature within the 3 mm samples (p < 0.0001), and the 4 mm samples (p < 0.0001) over the storage times. Results of the Scheffe post hoc analysis can be seen in Table 38.

		Thickness	
Comparisons	2 mm	3 mm	4 mm
		p value	
Immediate vs. 24-Hours	0.968	< 0.0001	< 0.0001
Immediate vs. I Week	0.855	< 0.0001	< 0.0001
Immediate vs. 1 Month	0.944	< 0.0001	< 0.0001
Immediate vs. 3 Months	0.997	< 0.0001	< 0.0001
24 Hours vs. 1 Week	0.996	0.562	< 0.0001
24 Hours vs. 1 Month	0.999	0.015	< 0.0001
24 Hours vs. 3 Months	0.998	0.001	< 0.0001
1 Week vs. 1 Month	0.999	0.315	0.638
1 Week vs. 3 Months	0995	0.049	0.012
1 Month vs. 3 Months	0.994	0.858	0.214

One-way ANOVA of the Photac-Fil Quick mean thermal peak values for the samples as a function of storage time identified a significant difference in thermal peak values within the immediate (p = 0.006), one-month (p = 0.006), and three-month (p = 0.001) storage groups. However, a significant difference was not found at 24-hours (p =

0.121) and at one-week (p = 0.155). Results of the post hoc analysis can be seen in Table 39.

	A CANADA	St	orage Time		
Comparisons	Immediate	24 Hour	One Week	One Month	Three Months
	Share the beautiful to	er its acquising.	p value	A MILE TO THE	
2 mm vs. 3 mm	0.007	0.305	0.155	0.033	0.002
2 mm vs. 4 mm	0.047	0.139	0.589	0.009	0.002
3 mm vs. 4 mm	0.569	0.864	0.596	0.746	0.999

Differential Scanning Calorimetry Thermal Weight Change

The complete data for the weight changes during thermal analysis can be seen in Appendix E, Tables 112 – 116 for Vitremer, Tables 117 – 121 for Fuji II LC, and Tables 122 – 126 for Photac-Fil Quick.

For all tables, "H20 Gain" represents the sample weight difference noted between fabrication and the end of storage in 100% humidity. "H20 Loss" represents the amount of weight change noted during the specific heat determinations, while "Net H20" represents the difference between water gain during storage and water loss during the specific heat determinations. "Scan Loss" represents the weight change during the 37 – 240 °C thermal scan, and "Total Loss" represents the sum of weight loss during specific heat determinations and thermal scan. "Total/Gain" represents the ratio (percent) of entire weight loss to water gain, and "Scan/Gain" likewise represents the ratio of weight loss during the DSC thermal scan to water gain. The latter two analyses were undertaken to explore any possible relationships between the amount of water gained and weight loss as a function of storage time.

Some assumptions were made during the formulation of the weight determinations. First, it was assumed that any weight change during the specific heat determinations would be due to water loss, presumably of the loosely-bound nature. This assumption was reinforced that the DSC chamber is constantly purged with low-humidity nitrogen gas and that temperature span of the specific heat determination (20 º - 60 º C) would not be sufficient to cause reasonable loss of any HEMA methacrylate components (BP 120 °C62) that may be present as well as other glassionomer components. The second assumption is that the during the DSC thermal scan (37 º -240 °C) the majority of weight loss is due to water of the more firmly-bound nature. This assumption is albeit made with the knowledge that some polyacrylic acids used in resinmodified glass-ionomer restorative materials have a boiling point (~ 210-220 °C⁶³) at the higher end of the thermal scan used in this evaluation. Accordingly, any material loss associated with an endothermal peak above 100 °C and below the polyacrylic acid boiling points should be assumed to be mostly due to water. Furthermore, the second assumption was reinforced with the results of pilot studies that revealed no endothermal peaks evident during the latter end of the DSC thermal scans. A third assumption made in the DSC weight change formulations was that any RMGI polymer formation would be stable during the temperature range of the DSC

thermal scan. This last supposition was established after pilot studies revealed that RMGI materials required temperatures in excess of 400 °C to degrade during DSC thermal scans.

The mean weight change during the thermal analysis for Vitremer is presented in Table 40. Two-way ANOVA identified significant gain in water (p < 0.0001) in the Vitremer samples as a function of storage time, but a difference was not noted (p = 0.891) between the different sample thicknesses, as well as no interactions between the factors of time and thickness (p = 0.389). The immediate samples were not included in any of the water gain analyses since these samples did not have the opportunity to gain water and could abnormally skew the statistical results. The results of the *post hoc* tests concerning Vitremer water gain can be seen in Table 41.

Table 41. Scheffe Multi Comparisons For Vitremer Water Gain	
Comparisons	p value
2 mm vs. 3 mm	0.913
2 mm vs. 4 mm	0.999
3 mm vs. 4 mm	0.921
24 Hours vs. 1 Week	0.235
24 Hours vs. 1 Month	0.005
24 Hours vs. 3 Months	<0.0001
1 Week vs. 1 Month	0.393
1 Week vs. 3 Months	<0.0001
1 Month vs. 3 Months	0.002
n = 15 Results shaded red are co significant at a 95% confid interval	A The second sec

time increased. The results of the Scheffe *post* hoc analysis per sample thickness and storage time can be seen in Table 42.

Analysis of the net water retention during specific heat determination for the Vitremer samples as a function of sample thickness and storage time determined that the Vitremer samples retained significantly more during the specific heat determinations as storage time increased (p < 0.0001) and as the samples increased in thickness (p = 0.001) without an interaction in the factors of storage time and

Two-way ANOVA analysis of the Vitremer samples water loss during the specific heat determinations identified a significant difference (p < 0.0001) in water loss between the samples for both factors storage time and sample thickness with an interaction noted (p = 0.038) between the two factors. The 2 mm samples were noted to lose significantly more water than the other thicknesses and more water was lost as storage

Table 42. Scheffe Multi Comparisons For Vitremer Water Loss During Specific Heat Determinations	
Comparisons	p value
2 mm vs. 3 mm	0.006
2 mm vs. 4 mm	<0.0001
3 mm vs. 4 mm	0.014
Immediate vs. 24 Hours	0.244
Immediate vs. 1 Week	0.069
Immediate vs. 1 Month	0.047
Immediate vs. 3 Months	<0.0001
24 Hours vs. 1 Week	0.979
24 Hours vs. 1 Month	0.941
24 Hours vs. 3 Months	<0.0001
1 Week vs. 1 Month	0.999
1 Week vs. 3 Months	<0.0001
1 Month vs. 3 Months	0.001
n = 15 Results shaded red are co significant at a 95% confid interval	

specimen thickness (p = 0.456). Itemized results of the *post hoc* analysis can be seen in Table 43.

Table 43. Scheffe Multi Comparisons For Vitremer Net Water Retention During Specific Heat Determinations	
Comparisons	p value
2 mm vs. 3 mm	0.038
2 mm vs. 4 mm	0.001
3 mm vs. 4 mm	0.401
Immediate vs. 24 Hours	0.885
Immediate vs. 1 Week	0.018
Immediate vs. 1 Month	<0.0001
Immediate vs. 3 Months	<0.0001
24 Hours vs. 1 Week	0.188
24 Hours vs. 1 Month	0.001
24 Hours vs. 3 Months	<0.0001
1 Week vs. 1 Month	0.254
1 Week vs. 3 Months	<0.0001
1 Month vs. 3 Months	0.130
n = 15 Results shaded red are co significant at a 95% confid interval	

Two-way ANOVA of Vitremer weight loss during the DSC thermal scan identified significant changes in weight as a function of thickness and storage time (p < 0.0001), without an interaction between the two factors (p = 0.226). The thicker samples were noted to exhibit more weight loss and increased weight loss was also noted with longer storage times. The results of the *post hoc* analysis can be seen in Table 44.

Analysis of the total weight loss (specific heat determination loss + thermal scan loss) identified a significant increase in weight loss as a function of storage time (p < 0.0001) but did not identify a significant difference when

analyzed as per specimen thickness (p = 0.085). Furthermore, no interaction between the factors was identified (p = 0.238). The results of the multi comparison *post hoc* testing can be seen in Table 45.

Table 44. Scheffe Multi Comparisons For Vitremer DSC Scan Weight Loss		
Comparisons	p value	
2 mm vs. 3 mm	0.038	
2 mm vs. 4 mm	0.001	
3 mm vs. 4 mm	0.401	
Immediate vs. 24 Hours	0.885	
Immediate vs. 1 Week	0.018	
Immediate vs. 1 Month	<0.0001	
Immediate vs. 3 Months	<0.0001	
24 Hours vs. 1 Week	0.188	
24 Hours vs. 1 Month	0.001	
24 Hours vs. 3 Months	<0.0001	
1 Week vs. 1 Month	0.254	
1 Week vs. 3 Months	<0.0001	
1 Month vs. 3 Months	0.130	
n = 15 Results shaded red are considered significant at a 95% confidence interval		

Comparisons	p value
2 mm vs. 3 mm	0.124
2 mm vs. 4 mm	0.259
3 mm vs. 4 mm	0.925
Immediate vs. 24 Hours	0.736
Immediate vs. 1 Week	0.005
Immediate vs. 1 Month	<0.0001
Immediate vs. 3 Months	<0.0001
24 Hours vs. 1 Week	0.150
24 Hours vs. 1 Month	0.004
24 Hours vs. 3 Months	<0.0001
1 Week vs. 1 Month	0.656
1 Week vs. 3 Months	<0.0001
1 Month vs. 3 Months	<0.0001

Two-way ANOVA analysis of the ratio (expressed as a percentage) of total weight loss to water gain of the Vitremer samples identified significant differences (p < 0.0001) among both factors of specimen thickness and storage time with a significant interaction indicated between the two factors (p < 0.0001). In summary, there were significant differences between the 2 mm and 4 mm specimens and all storage time periods demonstrated significant reduction in the ratio except between the one-week and one-month storage time. Results of the *post hoc* testing can be seen in Table 46.

Analysis of the ratio of the weight loss observed during the DSC thermal scan to water gain of the Vitremer samples indicated similar results, with significant differences found among the factors of storage time and thickness (p < 0.0001) with a significant interaction between the two

factors (p < 0.0001). Scheffe *post hoc* comparison results presented similar results as that seen for the total weight loss to water gain ratio evaluation and are posted in Table 47.

Table 46. Scheffe Multi Comparisons For Vitremer Ratio of Total Weight Loss to Water Gain	
Comparisons	p value
2 mm vs. 3 mm	0.925
2 mm vs. 4 mm	0.004
3 mm vs. 4 mm	0.010
24 Hours vs. 1 Week	<0.0001
24 Hours vs. 1 Month	<0.0001
24 Hours vs. 3 Months	<0.0001
1 Week vs. 1 Month	0.236
1 Week vs. 3 Months	<0.0001
1 Month vs. 3 Months	0.033
n = 15 Results shaded red are co significant at a 95% confid interval	

	Comparisons For Vitremer Ratio of DSC Thermal Scan Weight Loss	
Comparisons	p value	
2 mm vs. 3 mm	0.218	
2 mm vs. 4 mm	<0.0001	
3 mm vs. 4 mm	<0.0001	
24 Hours vs. 1 Week	<0.0001	
24 Hours vs. 1 Month	<0.0001	
24 Hours vs. 3 Months	<0.0001	
1 Week vs. 1 Month	0.161	
1 Week vs. 3 Months	<0.0001	
1 Month vs. 3 Months	0.033	
n = 15 Results shaded red are considered significant at a 95% confidence interval		

The mean weight change during the thermal analysis for Fuji II LC is presented in Table 48.

Two-way ANOVA analysis of Fuji II LC water gain (except the immediate groups) as a function of thickness and storage time, identified significant differences (p < 0.0001) for both factors, with an interaction indicated between the two factors (p = 0.001). Interesting, the 2 mm samples demonstrated more water gain than the thicker samples and although water gain increased as storage time increased, water gain was observed to peak at the one-week storage time.

Results of the *post hoc* testing can be seen in Table 49.

Table 48. Fuji II LC Mean Thermal Weight Change (%Wt)

0.00 0.70 0.12 0.70 0.12 3.88 0.50 4.56 0.40 0.00 0.00 0.62 0.36 -0.62 0.36 4.45 0.16 5.04 0.29 0.00 0.00 0.32 0.06 -0.32 0.06 4.75 0.15 5.06 0.21 3.46 0.64 3.04 0.37 0.43 0.41 5.17 0.28 8.38 1.14 245.77 1.01 0.17 0.966 0.12 0.05 0.11 4.73 0.61 5.88 0.18 596.46 1.04 0.15 0.86 0.19 0.11 5.05 0.22 5.86 0.19 576.00 4.94 1.47 3.57 1.36 1.37 0.23 5.98 0.17 8.88 0.62 220.46 4.29 0.81 2.13 0.22 6.49 0.79 8.88 0.62 244.35 4.14 1.06 2.87 1.17			H20 Gain	SD	H20 Loss	SD	Net H20	SD	Scan Loss	OS	Total Loss	SD	Scan Loss SD Total Loss SD Total/H20 Gain	OS	Scan/H20 Gain	SD
e 3mm 0.00 0.02 0.36 -0.62 0.36 4.45 0.16 5.04 0.29 4mm 0.00 0.00 0.32 0.06 -0.32 0.06 4.75 0.15 5.06 0.27 2mm 3.mm 1.01 0.17 0.966 0.12 0.06 0.11 4.73 0.61 5.86 0.18 5.64 2mm 1.01 0.17 0.966 0.12 0.06 0.11 4.73 0.61 5.88 0.18 5.640 2mm 4.94 1.47 3.57 1.36 1.37 0.23 5.98 0.17 9.33 1.31 205.38 4 mm 3.71 1.08 1.79 0.96 1.92 0.32 6.49 0.79 8.88 0.62 220.46 2 mm 4.14 1.06 2.87 1.17 1.27 0.25 6.92 0.34 8.61 0.86 20.34 4 mm 3.22 0.08		2mm	000	000	0.70	0.12	-0.70	0.12	3.88	05.0	4.56	0.40				
4mm 0.00 0.32 0.06 -0.32 0.06 4.75 0.15 5.05 0.21 2mm 3.46 0.64 3.04 0.37 0.43 0.41 5.17 0.28 8.38 1.14 245.77 3mm 1.01 0.17 0.966 0.12 0.06 0.11 4.73 0.61 5.88 0.18 596.46 4mm 1.04 0.15 0.86 0.09 0.19 0.11 5.06 0.22 5.86 0.17 9.33 1.31 246.37 4mm 4.29 1.47 3.57 1.35 1.37 0.23 5.98 0.17 9.33 1.31 206.38 4mm 3.71 1.08 1.79 0.96 1.92 0.32 6.64 0.22 5.86 0.17 9.33 1.31 206.38 4mm 3.71 1.08 1.77 1.72 0.25 6.92 0.34 8.61 0.36 244.36 4mm	Immediate	3mm	000	000	0.62	0.36	-0.62	0.36	4.45	0.16	5.04	0.29				
2mm 3.46 0.64 3.04 0.37 0.43 0.41 5.17 0.28 8.38 1.14 245.77 3mm 1.01 0.17 0.966 0.12 0.06 0.11 4.73 0.61 5.88 0.18 596.46 4mm 1.04 0.15 0.86 0.09 0.19 0.11 4.73 0.61 5.88 0.19 576.00 2mm 4.94 1.47 3.57 1.35 1.37 0.23 5.98 0.17 9.33 1.31 206.38 2mm 4.29 0.81 2.21 0.63 2.08 0.22 6.49 0.79 8.88 0.62 220.46 2mm 4.14 1.06 2.87 1.17 1.27 0.25 6.92 0.33 1.06 244.36 4mm 3.71 1.06 2.87 1.17 1.27 0.25 6.92 0.34 861 0.86 244.36 4mm 3.22 0.08		4mm	000	0.00	0.32	90.0	-0.32	0.06	4.75	0.15	5.05	0.21				
3mm 1,01 0.17 0.966 0.12 0.05 0.11 4.73 0.61 5.88 0.18 696.46 4mm 1,04 0.15 0.86 0.09 0.19 0.11 5.06 0.22 5.85 0.19 576.00 2mm 4.94 1,47 3.57 1.36 1.37 0.23 5.98 0.17 9.33 1.31 206.38 4mm 3.71 1.08 1.79 0.96 1.92 0.32 6.64 0.22 7.91 1.06 244.35 4mm 3.71 1.08 1.79 0.96 1.92 0.32 6.64 0.22 7.91 1.06 244.35 4mm 3.71 1.06 2.87 1.17 1.27 0.25 6.92 0.34 8.61 0.36 244.35 4mm 3.22 0.06 1.04 2.18 0.07 7.11 0.18 8.51 0.30 244.84 2mm 4.08 1.03		2mm	3.46	0.64	3.04	0.37	0.43	0.41	5.17	0.28	8:38	1.14	245.77	38.34	153.78	30.35
4mm 1.04 0.15 0.86 0.09 0.19 0.11 5.05 0.22 5.85 0.19 576.00 2mm 4.94 1.47 3.57 1.35 1.37 0.23 5.98 0.17 9.33 1.31 206.38 3mm 4.29 0.81 2.21 0.63 2.08 0.22 6.49 0.79 888 0.62 220.46 2mm 4.14 1.06 2.87 1.17 1.27 0.25 6.92 0.73 8.88 0.62 220.46 3mm 3.63 0.40 1.51 0.43 2.12 0.18 7.11 0.18 8.61 0.86 223.25 4mm 3.22 0.09 1.04 2.18 0.07 7.11 0.18 8.51 0.30 244.84 3mm 4.29 1.73 3.38 0.94 0.92 0.98 8.94 0.98 11.51 2.1 21 317.83 3mm 5.08	24 Hours	3mm	10.1	0.17	996.0	0.12	90:0	0.11	4.73	0.61	5.88	0.18	596.46	88.72	499.72	80.53
2mm 4.94 1.47 3.57 1.36 1.37 0.23 5.98 0.17 9.33 1.31 206.38 3mm 4.29 0.81 2.21 0.63 2.08 0.22 6.49 0.79 8.88 0.62 200.46 2mm 4.14 1.06 2.87 1.17 1.27 0.25 5.92 0.34 8.61 0.86 223.25 3mm 3.63 0.40 1.51 0.43 2.12 0.18 7.11 0.18 8.51 0.30 244.84 4mm 3.22 0.08 1.04 0.14 2.18 0.07 7.11 0.11 8.61 0.30 244.84 2mm 4.29 1.73 3.38 0.94 0.92 0.98 8.94 0.98 11.51 2.1 302.09 3mm 5.08 1.99 3.18 1.03 1.26 9.02 1.96 11.9 2.33 302.09 4mm 5.47 2.44		4mm	1.04	0.15	0.85	0.09	0.19	0.11	5.05	0.22	5.85	0.19	976.00	84.76	493.20	79.32
3mm 4.29 0.81 2.21 0.63 2.08 0.22 6.49 0.79 8.88 0.62 220.46 4mm 3.71 1.08 1.79 0.96 1.92 0.32 6.64 0.22 7.91 1.05 244.35 2mm 4.14 1.06 2.87 1.17 1.27 0.25 5.92 0.34 8.61 0.86 223.25 3mm 3.63 0.40 1.51 0.43 2.12 0.18 7.11 0.18 8.51 0.30 244.84 4mm 3.22 0.08 1.04 0.14 2.18 0.07 7.11 0.11 8.08 0.21 259.10 2mm 4.29 1.73 3.38 0.94 0.92 0.98 8.94 0.98 11.51 2.1 317.83 3mm 5.08 1.99 3.18 1.03 1.26 9.02 1.96 11.51 2.33 302.09 4mm 5.47 2.44		2mm	4.94	1.47	3.57	1.35	1.37	0.23	96.9	0.17	9.33	1.31	205.38	29.28	130.97	36.36
4mm 3.71 1.08 1.79 0.96 1.92 0.32 6.64 0.22 7.91 1.05 244.35 2mm 4.14 1.06 2.87 1.17 1.27 0.25 5.92 0.34 8.61 0.86 223.25 4mm 3.63 0.40 1.51 0.43 2.12 0.18 7.11 0.18 8.51 0.30 244.84 4mm 3.22 0.08 1.04 0.14 2.18 0.07 7.11 0.11 8.08 0.21 259.10 2mm 4.29 1.73 3.38 0.94 0.92 0.98 8.94 0.98 11.51 2.1 317.83 3mm 5.08 1.99 3.18 1.03 1.25 9.02 1.96 11.9 2.33 302.09 4mm 5.08 1.99 3.16 1.74 3.02 1.98 10.42 12.46 1.68 29.77	One Week	3mm	4.29	0.81	221	0.63	2.08	0.22	6.49	0.79	888	0.62	220.46	31.48	167.58	37.67
2mm 4.14 1.06 2.87 1.17 1.27 0.25 5.92 0.34 8.61 0.86 223.25 nth 3mm 3.63 0.40 1.51 0.43 2.12 0.18 7.11 0.18 8.51 0.30 244.84 4mm 3.22 0.08 1.04 0.14 2.18 0.07 7.11 0.11 8.08 0.21 259.10 2mm 4.29 1.73 3.38 0.94 0.92 0.98 8.94 0.98 11.51 2.1 317.83 3mm 5.08 1.99 3.18 1.03 1.25 9.02 1.96 11.9 2.33 302.09 4mm 5.07 2.44 2.46 1.74 3.02 10.38 10.42 12.46 1.68 297.72		4mm	3.71	1.08	1.79	96.0	1.92	0.32	6.64	0.22	7.91	1.05	244.35	37.64	195.83	48.78
nth 3mm 3.63 0.40 1.51 0.43 2.12 0.18 7.11 0.18 8.51 0.30 4mm 3.22 0.08 1.04 0.14 2.18 0.07 7.11 0.11 8.08 0.21 2mm 4.29 1.73 3.38 0.94 0.92 0.98 8.94 0.98 11.51 2.1 3mm 5.08 1.99 3.18 1.03 1.25 9.02 1.96 11.9 2.38 4mm 5.47 2.44 2.46 1.74 3.02 0.98 10.38 0.42 1.54 1.68		Zmm	4.14	1.06	2.87	1.17	1.27	0.25	5.92	0.34	19:8	0.86	223.25	31.28	152.89	40.66
4mm 3.22 0.08 1.04 0.14 2.18 0.07 7.11 0.11 8.08 0.21 2mm 4.29 1.73 3.38 0.94 0.92 0.98 8.94 0.96 11.51 2.1 3mm 5.08 1.99 3.18 1.03 1.9 1.25 9.02 1.96 11.9 2.33 4mm 5.47 2.44 2.46 1.74 3.02 0.98 10.38 0.42 1.246 1.68	One Month	3mm	3.63	0.40	151	0.43	2.12	0.18	11.7	0.18	8.51	0.30	244.84	17.11	202.15	23.81
2mm 4.29 1.73 3.38 0.94 0.92 0.98 8.94 0.98 11.51 2.1 3mm 5.08 1.99 3.18 1.03 1.9 1.25 9.02 1.96 11.9 2.33 4mm 5.47 2.44 2.46 1.74 3.02 1.08 10.38 1.45 1.68		4mm	3.22	0.08	1.04	0.14	2.18	0.07	7.11	0.11	8.08	0.21	259.10	4.45	226.88	5.75
3mm 5.08 1.99 3.18 1.03 1.9 1.25 9.02 1.96 11.9 2.33 4mm 5.47 2.44 2.46 1.74 3.02 1.98 11.38 11.45 1.746 1.68	Three	Zmm	4.29	1.73	3.38	0.94	0.92	0.98	8.94	96.0	11.51	2.1	317.83	44.02	239.90	26.19
17 1 17 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	Illee Illee	3mm	90'9	1.99	3.18	1.03	1.9	1.25	9.02	1.98	11.9	233	302.09	10.65	261.45	11.27
C. 1. 2. 1. 2. 1. 2. 2. 2. 2. 2. 2. 2. 2. 2. 2. 2. 2. 2.	SIIIIOM	4mm	5.47	2.44	2.46	1.74	3.02	0.98	10.38	0.42	12.46	1.68	292.72	9.81	265.02	11.72

n = 5

Table 49. Scheffe Multi Comparisons For Fuji II I Gain	LC Water
Comparisons	p value
2 mm vs. 3 mm	0.002
2 mm vs. 4 mm	<0.0001
3 mm vs. 4 mm	0.584
24 Hours vs. 1 Week	<0.0001
24 Hours vs. 1 Month	<0.0001
24 Hours vs. 3 Months	0.115
1 Week vs. 1 Month	0.103
1 Week vs. 3 Months	<0.0001
1 Month vs. 3 Months	<0.0001
n = 15 Results shaded red are co significant at a 95% confid interval	and the second of the second o

Analysis of the net water retention during specific heat determination of the Fuji II LC samples identified significant differences for both storage time and specimen thickness (p < 0.0001) as well as an interaction between both factors (p < 0.0001). Post hoc analysis identified significant differences for all thickness combinations except between the 3 mm and 4 mm samples and all storage time periods were significantly different except between the one-week and one-month storage times. The results of the multi comparison testing can be seen in Table 51.

Two-way ANOVA analysis of the Fuji II LC water loss during specific heat determination identified significant weight loss in the samples both when analyzed per thickness and storage time (p < 0.0001) with an interaction between the two factors also identified (p = 0.012). The results of the Scheffe multi comparison testing can be seen in Table 50.

Table 50. Scheffe Multi Comparisons For Fuji II LC Water Loss During Specific Heat Determinations	
Comparisons	p value
2 mm vs. 3 mm	<0.0001
2 mm vs. 4 mm	<0.0001
3 mm vs. 4 mm	0.150
Immediate vs. 24 Hours	<0.0001
Immediate vs. 1 Week	<0.0001
Immediate vs. 1 Month	<0.0001
Immediate vs. 3 Months	0.109
24 Hours vs. 1 Week	0.003
24 Hours vs. 1 Month	0.942
24 Hours vs. 3 Months	0.308
1 Week vs. 1 Month	0.030
1 Week vs. 3 Months	<0.0001
1 Month vs. 3 Months	0.061
n = 15 Results shaded red are considered significant at a 95% confidence	

interval

Table 51. Scheffe Multi Comparisons For Fuji II LC Net Water Retention During Specific Heat Determinations	
Comparisons	p value
2 mm vs. 3 mm	<0.0001
2 mm vs. 4 mm	<0.0001
3 mm vs. 4 mm	0.129
Immediate vs. 24 Hours	<0.0001
Immediate vs. 1 Week	<0.0001
Immediate vs. 1 Month	<0.0001
Immediate vs. 3 Months	<0.0001
24 Hours vs. 1 Week	<0.0001
24 Hours vs. 1 Month	<0.0001
24 Hours vs. 3 Months	<0.0001
1 Week vs. 1 Month	0.962
1 Week vs. 3 Months	<0.0001
1 Month vs. 3 Months	<0.0001
n = 15 Results shaded red are co significant at a 95% confid interval	

Analysis of the Fuji II LC samples for weight loss during the DSC thermal scan identified significant differences for both storage time and specimen thickness (p < 0.0001) as well as an interaction between both factors (p < 0.0001). Post hoc analysis identified significant differences for all thickness combinations except between the 3 mm and 4 mm samples and all storage time periods were significantly different except between the one-week and one-month and the one-week and three-month storage times. The results of the multi comparison testing can be seen in Table 52. Analysis of the Fuji II LC samples concerning total weight loss (specific heat determination + DSC thermal scan) closely mirrored the results of the DSC thermal scan with significant differences identified (p < 0.0001) for both storage time and

specimen thickness with an interaction noted between the two factors (p < 0.0001). The results of the *post hoc* testing can be seen in Table 53.

Table 52. Scheffe Multi Comparisons For Fuji II LC Weight Loss During DSC Thermal Scan		
Comparisons	p value	
2 mm vs. 3 mm	<0.0001	
2 mm vs. 4 mm	<0.0001	
3 mm vs. 4 mm	0.183	
Immediate vs. 24 Hours	0.001	
Immediate vs. 1 Week	<0.0001	
Immediate vs. 1 Month	<0.0001	
Immediate vs. 3 Months	<0.0001	
24 Hours vs. 1 Week	<0.0001	
24 Hours vs. 1 Month	<0.0001	
24 Hours vs. 3 Months	<0.0001	
1 Week vs. 1 Month	0.148	
1 Week vs. 3 Months	0.881	
1 Month vs. 3 Months	0.012	
n = 15 Results shaded red are considered significant at a 95% confidence interval		

Table 53. Scheffe Multi Comparisons For Fuji II LC Total Weight Loss	
Comparisons	p value
2 mm vs. 3 mm	0.045
2 mm vs. 4 mm	<0.0001
3 mm vs. 4 mm	0.237
Immediate vs. 24 Hours	0.001
Immediate vs. 1 Week	<0.0001
Immediate vs. 1 Month	<0.0001
Immediate vs. 3 Months	<0.0001
24 Hours vs. 1 Week	<0.0001
24 Hours vs. 1 Month	<0.0001
24 Hours vs. 3 Months	0.164
1 Week vs. 1 Month	0.794
1 Week vs. 3 Months	<0.0001
1 Month vs. 3 Months	0.001
n = 15 Results shaded red are considered significant at a 95% confidence interval	

Two-way ANOVA analysis of the ratio (expressed as a percentage) of total weight loss to water gain of the Fuji II LC samples identified significant differences (p < 0.0001) with respect to specimen thickness and storage time and a significant interaction indicated between the two factors (p < 0.0001). There were significant differences between the all thickness pairs, except between the 3 mm and 4 mm specimens, and all storage time periods demonstrated significant reduction in the ratio except between the one-week and one-month storage time. Results of the post hoc testing can be seen in Table 54.

Analysis of the ratio of the weight loss observed during the DSC thermal scan to water gain of the Fuji II LC samples indicated similar results, with significant differences found among storage time and thickness (p < 0.0001) with a significant interaction between the two factors (p < 0.0001). Interestingly, for the 3 mm and 4 mm samples, the ratio demonstrated its highest

value at 24 hours, which then decreased at the one-week storage time. At this point, all three thicknesses of the samples demonstrated an increase in the ratio. Scheffe *post hoc* comparison results presented similar results as that seen for the total weight loss to water gain ratio evaluation, and the results are posted in Table 55.

Table 54. Scheffe Multi Comparisons For Fuji II LC Total Weight Loss to Water Gain		
Comparisons	p value	
2 mm vs. 3 mm	<0.0001	
2 mm vs. 4 mm	<0.0001	
3 mm vs. 4 mm	0.987	
24 Hours vs. 1 Week	<0.0001	
24 Hours vs. 1 Month	<0.0001	
24 Hours vs. 3 Months	<0.0001	
1 Week vs. 1 Month	0.699	
1 Week vs. 3 Months	<0.0001	
1 Month vs. 3 Months	0.005	
n = 15 Results shaded red are considered significant at a 95% confidence interval		

Table 55. Scheffe Multi Comparisons For Fuji II LC DSC Thermal Scan Weight Loss to Water Gain	
Comparisons	p value
2 mm vs. 3 mm	<0.0001
2 mm vs. 4 mm	<0.0001
3 mm vs. 4 mm	0.987
24 Hours vs. 1 Week	<0.0001
24 Hours vs. 1 Month	<0.0001
24 Hours vs. 3 Months	<0.0001
1 Week vs. 1 Month	0.335
1 Week vs. 3 Months	<0.0001
1 Month vs. 3 Months	0.004
n = 15 Results shaded red are considered significant at a 95% confidence interval	

The mean results of the Photac-Fil Quick thermal analysis weight changes can be seen in Table 56. Two-way ANOVA analysis of the Photac-Fil Quick samples for water gain (immediate groups excluded) identified a significant increase in weight as a factor of storage time (p < 0.0001) but found no significant difference in weight gain as a factor of specimen thickness (p = 0.323). There was no significant interaction between the two factors found (p = 0.106). Overall results found that the samples demonstrated significant increase in water gain with increase in storage time, with the three-month samples displaying significantly more water uptake than the other samples. Results of the *post hoc* testing can be seen in Table 57.

Table 56. Photac-Fil Quick Mean Thermal Weight Change (%Wt)

Numediate Sum 0.00	000 000	0.98	0.40	86 0,	0.40	1 22	144	CONTROL OF THE PERSON NAMED IN			STATE OF THE PARTY		
e 3mm 4mm	000 000 000	970	0.5	-	0.10	CC.4	0.41	5.26	0.30				
Zmm	000 000	200	0.04	-0.5	0.04	4.31	1.14	4.79	1.1				
Zmm	0.43	131	90.0	-0.31	90.0	5.64	0.13	5.94	0.17				
, Janes	99.0	2.02	0.22	-0.54	0.59	5.42	15.0	7.33	99.0	90'895	293.20	410.74	205.90
	0.0	1.64	0.48	0.51	0.27	7.24	0.32	8.76	0.71	441.83	112.12	363.12	106.24
4mm 1.47	0.52	1.19	0.33	0.27	0.21	6.46	0.58	7.60	0.81	554.79	105.37	470.39	99.24
Zmm 1.61	0.14	1.44	0.49	0.17	0.35	90'9	EE'0	7.43	0.27	471.06	33.74	381.99	56.34
One Week 3mm 1.66	0.16	96.0	0.12	0.7	0.09	6.52	0.15	741	0.19	457.86	41.40	398.90	41.52
4mm 1.62	0.1	0.79	0.03	0.83	0.09	6.7	0.23	7.43	0.22	466.27	21.23	416.51	19.64
2mm 1.78	0.37	1.6	0.2	0.18	0.31	6.13	98.0	7.74	0.25	459.95	111.13	366.01	98.14
One Month 3mm 1.93	0.08	96:0	0.14	86.0	0.10	6.88	0.11	7.76	0.22	408.86	11.06	389.02	10.54
4mm 1.97	0.45	0.75	0.19	1.22	0.25	7.18	0.13	7.88	0.28	419.86	76.22	379.55	78.24
Zmm	0.32	1.86	0.22	29.0	0.27	6.41	67.0	8.16	0.47	332.42	31.92	256.55	24.91
Illiet 3mm 2.78	0.18	0.99	000	1.79	110	7.61	0.12	8.52	0.17	315.32	14.74	278.82	14.64
monuns 4mm 3.64	1.07	0.91	0.12	272	0.96	8.05	0.14	8.91	0.17	334.57	158.59	299.44	140.56

n = 5

Table 57. Scheffe Multi Comparisons For Photac-Fil Quick Water Gain	
Comparisons	p value
2 mm vs. 3 mm	0.441
2 mm vs. 4 mm	0.663
3 mm vs. 4 mm	0.936
24 Hours vs. 1 Week	0.924
24 Hours vs. 1 Month	0.504
24 Hours vs. 3 Months	<0.0001
1 Week vs. 1 Month	0.189
1 Week vs. 3 Months	<0.0001
1 Month vs. 3 Months	0.001
n = 15 Results shaded red are considered significant at a 95% confidence interval	

loss appeared to peak at the 24-hour storage time and then appeared to lessen and become stable, at least for the 3 mm and 4 mm samples. Results of the Scheffe post hoc analysis are presented in Table 58.

Analysis of the net water retention of the Photac-Fil Quick samples during the specific heat determinations revealed significant differences for both sample thickness and storage time (p < 0.0001) without an interaction between the two factors (p = 0.168). Specifically the samples retained more water as sample thickness and storage

Analysis of Photac-Fil Quick water loss during the specific heat determinations indicated significant differences among sample thickness and storage time (p < 0.0001) without an interaction between the two factors (p = 0.355). The different sample thicknesses were all significantly different from each other, with the 2 mm samples displaying more water loss than the other two thicknesses (p < 0.013). Considering sample storage time, water

Table 58. Scheffe Multi Comparisons For Photac-Fil Quick Water Loss During Specific Heat Determinations		
Comparisons	p value	
2 mm vs. 3 mm	<0.0001	
2 mm vs. 4 mm	<0.0001	
3 mm vs. 4 mm	0.013	
Immediate vs. 24 Hours	<0.0001	
Immediate vs. 1 Week	<0.0001	
Immediate vs. 1 Month	<0.0001	
Immediate vs. 3 Months	<0.0001	
24 Hours vs. 1 Week	<0.0001	
24 Hours vs. 1 Month	<0.0001	
24 Hours vs. 3 Months	0.030	
1 Week vs. 1 Month	0.998	
1 Week vs. 3 Months	0.260	
1 Month vs. 3 Months	0.548	
n = 15 Results shaded red are considered significant at a 95% confidence interval		

time increased. Results of the Scheffe multi comparison testing are presented in Table 59.

Table 59. Scheffe Multi Comparisons For Photac-Fil Quick Net Water Retention During Specific Heat Determinations	
p value	
<0.0001	
<0.0001	
0.413	
0.001	
<0.0001	
<0.0001	
<0.0001	
0.001	
0.046	
<0.0001	
0.707	
<0.0001	
<0.0001	

significant at a 95% confidence

interval

Two-way ANOVA analysis of the Photac-Fil Quick samples for weight loss during the thermal DSC scan similarly to the total loss analysis identified significant differences in weight as a function of storage time (p < 0.0001) and sample thickness (p = 0.0001) with an interaction noted between the two factors (p < 0.0001). In general, total weight loss increased as storage time increased and sample thickness increased. Results of the post hoc testing can be seen in Table 61.

Two-way ANOVA analysis of the Photac-Fil Quick samples for total weight loss during specific heat determinations and thermal DSC scan identified significant differences in weight as a function of storage time (p < 0.0001) and sample thickness (p = 0.023) with an interaction noted between the two factors (p < 0.0001). In general, total weight loss increased as storage time increased and sample thickness increased. Results of the *post hoc* testing can be seen in Table 60.

Table 60. Scheffe Multi Comparisons For Photac-Fil Quick Total Weight Loss	
Comparisons	p value
2 mm vs. 3 mm	<0.0001
2 mm vs. 4 mm	<0.0001
3 mm vs. 4 mm	0.413
Immediate vs. 24 Hours	0.001
Immediate vs. 1 Week	<0.0001
Immediate vs. 1 Month	<0.0001
Immediate vs. 3 Months	<0.0001
24 Hours vs. 1 Week	0.001
24 Hours vs. 1 Month	0.046
24 Hours vs. 3 Months	<0.0001
1 Week vs. 1 Month	0.707
1 Week vs. 3 Months	<0.0001
1 Month vs. 3 Months	<0.0001
n = 15 Results shaded red are considered significant at a 95% confidence interval	

Table 61. Scheffe Multi Comparisons For Photac-Fil Quick DSC Thermal Scan Weight Loss	
Comparisons	p value
2 mm vs. 3 mm	<0.0001
2 mm vs. 4 mm	<0.0001
3 mm vs. 4 mm	0.413
Immediate vs. 24 Hours	0.001
Immediate vs. 1 Week	<0.0001
Immediate vs. 1 Month	<0.0001
Immediate vs. 3 Months	<0.0001
24 Hours vs. 1 Week	0.998
24 Hours vs. 1 Month	0.236
24 Hours vs. 3 Months	<0.0001
1 Week vs. 1 Month	0.385
1 Week vs. 3 Months	<0.0001
1 Month vs. 3 Months	0.006
n = 15	

Results shaded red are considered significant at a 95% confidence

interval

Two-way ANOVA analysis of the ratio (expressed as a percentage) of total weight loss to water gain of the Photac-Fil Quick samples identified significant difference for weight change as a function of storage time (p < 0.0001) but not for specimen thickness (p = 0.342) without an identified interaction indicated between the two factors (p < 0.903). There were significant differences noted between the 24-hour and threemonth as well as between the one-week and three-month storage time periods. Results of the post hoc testing can be seen in Table 62.

Table 62. Scheffe Multi Comparisons For Photac-Fil Quick Total Weight Loss to Water Gain		
Comparisons	p value	
2 mm vs. 3 mm	0.366	
2 mm vs. 4 mm	0.928	
3 mm vs. 4 mm	0.582	
24 Hours vs. 1 Week	0.613	
24 Hours vs. 1 Month	0.199	
24 Hours vs. 3 Months	<0.0001	
1 Week vs. 1 Month	0.868	
1 Week vs. 3 Months	0.020	
1 Month vs. 3 Months	0.129	
n = 15 Results shaded red are considered significant at a 95% confidence interval		

Analysis of the ratio of the weight loss observed during the DSC thermal scan to water gain of the Photac-Fil Quick samples indicated similar results, with significant differences found with storage time (p = 0.001) but not with sample thickness (p = 0.310). As with the analysis of total weight loss to water gain, a significant interaction between the two factors was not identified (p = 0.927). The highest ratio was noted at the 24-hour storage time, at which point the ratio of scan loss to water gain was observed to decrease with storage time. Scheffe *post hoc* comparison results are presented in Table 63.

Table 63. Scheffe Multi Comparisons Ratio For Photac-Fil Quick DSC Thermal Scan Weight Loss to Water Gain		
Comparisons	p value	
2 mm vs. 3 mm	0.992	
2 mm vs. 4 mm	0.452	
3 mm vs. 4 mm	0.382	
24 Hours vs. 1 Week	0.976	
24 Hours vs. 1 Month	0.607	
24 Hours vs. 3 Months	0.003	
1 Week vs. 1 Month	0.845	
1 Week vs. 3 Months	0.011	
1 Month vs. 3 Months	0.089	
n = 15 Results shaded red are considered significant at a 95% confidence interval		

CHAPTER FOUR: DISCUSSION In terms of setting chemistry, resin-modified glass-ionomer restorative materials are more complex than their conventional predecessors.3 Although Andrzejewska et al.64 reported that the addition of aqueous polyacids enhanced the light-induced polymerization of HEMA, the converse is not the same, as organic, photopolymerizable monomers alters the conventional acid-base setting reaction.3 Anstice and Nicholson65 reported that the addition of HEMA to a conventional restorative glass-ionomer restorative material resulted in a slower setting reaction with reduced compressive strength. The cause of this effect was assumed to be reduction of ionic reactions due to reduction of dielectric constant of HEMA.⁶⁵ conformational changes that result in a more tightly coiled thus less-reactive polyacrylic acid polymer,3 and that the aqueous HEMA solution provides less solubility for the polyacrylic acids to enter solution and react.⁶⁵ Conventional clinical wisdom states that depth of cure is not important with resin-modified glassionomer restorative materials, since clinicians seem to assume that the conventional acid-base reaction will take place at deeper areas in which light penetration is insufficient to enable the polymerization of HEMA.3 However, clinicians need to realize that these deeper unpolymerized areas represent an area in which the acid-base reaction is diluted with HEMA. Although some type of setting reaction may take place, it occurs at a much slower pace and results in a material with inferior physical properties.^{3,65}


The overall aim of this research was to identify and possibly quantify any post-visible-light irradiation reactions that may occur in resin-modified glass-ionomer restorative materials in depths of two to five millimeters. Three different areas of investigation were undertaken: Solubility, Hardness, and Thermal Analysis.


Solubility: The solubility testing method was partly based on ANSI/ADA Specification No. 27: Resin-Based Filling Materials (1997). Accordingly, the method utilized was designed to provide specimens of uniform surface area, thickness, equal exposure duration to moisture, and dehydration of the specimens until weight stability to the nearest 0.001 grams was observed. However, this evaluation differed from either the ANSI/ADA Specification No. 27 and/or clinical

conditions in three ways. The first variation is that in this investigation solubility was solely assessed as a function of percentage weight loss, whereas the ANSI/ADA specification is based on weight loss per calculated volume. This parameter was chosen in an attempt to establish a protocol that would allow expedient comparison of solubility, lessen material loss, as well as sample dehydration during extra handling required of volume measurements. Although the chosen thesis parameter, being based on weight alone, could be criticized for lacking the precision of the ANSI/ADA specification, the parameter was envisioned to allow expedient comparative solubility analysis.

The second variation is that this method involved specimens of different photopolymerization protocols with four different thicknesses per group. The "experimental" group represents a photopolymerization method that mirrors clinical placement; that is, the specimens are light activated as per manufacturer recommendations from the top surface of the specimen cylinder. This thesis "control" group involved both the specimen top and bottom surface photoactivation. Although the control group represents a clinically unrealistic situation, the control samples were designed to allow possible comparison of material solubility behavior and kinetics between the two polymerization schemes.


The third area in which this evaluation differed is that a protective varnish or resin coating was not applied to the polymerized samples. Although originally touted to be less susceptible than conventional glass-ionomer materials to the effects of early moisture contamination and later dehydration, resin-modified glass-ionomer restorative materials were found to also require a surface protecting agent. ^{10,11,66,67} Surface protectants were not applied to the specimens during this investigation as it was felt that the surface coatings could reduce the experiment sensitivity to solubility differences that could exist due to the different curing regimens.







Graphical results for the multivariate analysis of the Vitremer solubility groups as a function of curing mode, and storage time can be seen in Figure 6. The experimental groups exhibited similar solubility at all storage times except for the one month time period. The control groups exhibited increasing solubility with the increase in storage time which appeared to peak at one month and then decrease. The results for multivariate Vitremer solubility for both curing modes as a function of thickness can be seen in Figure 7.

Comparative solubility trends between the two Vitremer groups are compared with the results graphically superimposed in Figure 8. The experimental samples are represented by the dotted lines while the control samples are depicted with the solid lines. It can be readily observed that solubility over time between the two curing modes exhibits a similar trend, as independent of each thickness by curing mode pair only identified four instances of significant difference in solubility.



Graphical results for the multivariate analysis of the Fuji II LC solubility groups as a function of curing mode and storage time can be seen in Figure 9. The Fuji II LC experimental groups exhibited statistically less solubility at 24 hours but increased solubility as storage time progressed. The control samples likewise exhibited a significant increase in solubility as storage time progressed, which likewise to the Vitremer samples, appears to peak at one month storage.






Graphical results for the multivariate analysis of the Fuji II LC solubility groups as a function of curing mode and specimen thickness can be seen in Figure 10. Interestingly, the experimental groups exhibited similar solubility regardless of sample thickness while the control groups were observed to demonstrate decreased solubility as the sample increased in thickness.


An overlay of the Fuji II LC control and experimental solubility results is depicted in Figure 11.


The experimental results are annotated with dotted lines and the control results are represented

by solid lines. Independent t-test analysis of each curing mode thickness pair at the same storage time identified that the control group for the immediate, three-month, and all but the 3 mm samples of the 24-hour group exhibited significantly less solubility than the experimental samples. Similar to the Vitremer samples, the overall solubility trend of increasing solubility with increased storage time with solubility peaking at one month is evident.




Graphical results for the multivariate analysis of the Photac-Fil Quick solubility groups as a function of curing mode and storage time can be seen in Figure 12. With both curing modes there was a significant increase in solubility with increasing storage time that exhibits a significant decrease after one month of storage. When considered per specimen thickness (Figure 13), the Photac-Fil Quick experimental specimens, similar to that seen of Fuji II LC, do not differ significantly in solubility. Also similar to Fuji II LC, the Photac-Fil Quick control samples exhibit significantly less solubility with increased sample thickness, with the 2 mm specimens displaying significantly more solubility than the other thicknesses.






Photac-Fil Quick combined control and experimental solubility results are graphically depicted in Figure 14, with the dashed lines depicting the experimental samples and the solid lines annotating the control samples. Independent t-test analysis of each curing mode thickness pair at the same storage time identified that the control group samples exhibited significantly less

solubility overall than the experimental samples, except for the three-month storage time at which all exhibited similar solubility. Similar to Vitremer and Fuji II LC, the solubility of Photac-Fil Quick exhibits an increasing trend of solubility that decreases after one month of storage.



The marginal means of each material as a function of storage time are depicted in Figure 15. The experimental group is depicted as this curing mode is representative of clinical conditions. It was not the purpose of this thesis to compare the performance of each material: any such contrast is accomplished only to observe overall material trends under similar conditions. All three materials displayed a trend of increasing solubility with increase of storage time that peaked at the one-month storage interval. This solubility behavior was not totally expected, as it could be assumed that solubility would decrease as the RMGI matrix matured with increasing time. Glass-ionomer materials have been reported to exhibit continued matrix maturity with time,³ and the exact time intervals involved with RMGI maturity have not been fully revealed. Sidhu *et al.*⁶⁸ elucidated some information concerning the matrix maturity of RMGI materials,

and found that Fuji II LC, Photac-Fil Quick, and Vitremer did not exhibit matrix maturity at three months, as evidenced by all materials being susceptible to dehydration effects. At six months, Fuji II LC exhibited resistance to dehydration suggesting matrix maturity, and it required 12 months for Vitremer and Photac-Fil Quick to reach the same level. It can be inferred that the decrease in solubility noted at three months for each of these materials is associated with the matrix maturity process; perhaps this solubility behavior provides evidence that the resinmodified glass-ionomer materials are susceptible to the effects of dehydration for up to one month. However, further evaluation is required for further storage time periods before more definitive comparisons should be attempted.



For the three materials, it is interesting to also note that the trend within each material is that the solubility between the different curing methods was not overwhelmingly significantly distinct.

The experimental curing method (visible light activation at the top surface) simulates the clinical situation. The control curing method (visible light activation at the top and bottom surfaces) was designed to allow possible comparison of solubility dynamics between the clinical and a more

thoroughly-cured RMGI sample. It was expected that the additional curing of the bottom surfaces afforded to the control samples would effect almost uniformly less-soluble samples, as compared to the experimental group. Review of Tables 1, 2, and 3 will show the t-test results between individual specimens of the same thickness and storage time, and will demonstrate individual instances of non-difference between samples as per curing mode.

The lack of uniform difference between curing modes of the RMGI materials might be explained by the resin content of the RMGI materials. RMGI have been reported to be composed of approximately 4.5 – 6 percent resin in the set product. Although the presence of resin may purportedly improve some physical properties of glass ionomer materials, perhaps the amount of resin present is in insufficient quantity to affect solubility behavior. Accordingly, the solubility behavior of the RMGI materials may be more dictated by the conventional, acid-base components of the hybrid material.

Another unexpected finding under the conditions of this study was the lack of solubility differences between the different depths of samples. Additional one-way ANOVA of the Vitremer samples as a factor of thickness did not find a significant difference in solubility for the experimental samples (p = 0.141) as well as the control samples (p = 0.243). Analysis of Fuji II LC samples did not identify a significant difference in solubility within the experimental samples (p = 0.986) or the control samples (p = 0.732). Likewise, analysis of Photac-Fil Quick samples did not find a significant difference among the experimental samples (p = 0.912) nor the control specimens (p = 0.179). Although experimental conditions were followed as precisely as possible, solubility testing did result in some groups yielding variant standard deviations. However, it is thought that the wide standard deviations observed would not affect the overall results reported in this evaluation.

One explanation that may have resulted in the lack of solubility difference findings of this study could be due to the kinetics of water penetration into the samples. Water diffusion into most

dental materials follows Fick's First Law of Diffusion, which follows that for disc specimens (where edge effects can be ignored) the following equation applies:

$$M_t / M_{\infty} = 2(Dt/\pi l^2)^{1/2}$$

where

 M_t = the mass up take at time t (seconds) M_{∞} = is the equilibrium uptake I = the thickness of the specimen D = the diffusion coefficient.⁶⁹

Nicholson⁷⁰ reported the water absorption of two resin-modified glass-ionomer restorative materials (Vitremer and Baseline VLC [no longer manufactured]) and found that, based on solution used, the two materials exhibited absorption behavior that followed Fick's First Law for only an initial period of time. Nicholson identified that curing duration effected a dependent action on the diffusion coefficient, in that longer visible light polymerization times had an inverse effect on the diffusion coefficient. Nicholson reasoned that the longer curing duration polymerized more previously-soluble materials into a more insoluble state. Furthermore, an increase in the degree of polymerization would reduce the ease in which water molecules could potentially diffuse through the structure.⁷⁰ Although as water diffusion rate into the resinmodified glass-ionomer materials declined with increasing time, Nicholson reported that it still required a total of ten days to reach equilibrium with a 2.8-millimeter thick disc of Vitremer.⁷⁰

Nicholson's work was reinforced by two reports by Jevnikar *et al.*^{71,72} using micro-magnetic resonance imaging techniques. In their first report,⁷¹ the authors compared water infiltration into both conventional and resin-modified glass ionomer restorative materials for up to 192 hours and reported that in the first 24 hours water diffused 0.72 mm into the conventional material and 0.33 mm into resin-modified glass ionomer material. The conventional glass-ionomer materials required 96 hours for full water equilibrium, however; the resin-modified glass-ionomer material required 192 hours, which led the authors to conclude that resin-modified glass-ionomer restorative materials are more impervious to water infiltration than their conventional

counterparts.⁷¹ In their second report, Jevnikar *et al.*⁷² investigated the effect of surface coatings on water migration into resin-modified glass-ionomer restorative materials using similar techniques as in their previous report and further demonstrated that water diffusion into these materials are limited for the first 24 hours.

This information from Nicholson⁷⁰ and Jevnikar *et al.*^{71,72} may provide insight pertaining to why a profound expected difference was not observed in solubility between the control and experimental groups during this evaluation. Accordingly, since only a 24-hour water exposure time was used for solubility determination, the water diffusion rate into the respective materials represented only near-surface solubility behaviors. Hence, the solubility results were limited by the rate of water diffusion into the thicker materials and did not represent the solubility behavior of the entire sample. Future solubility evaluations involving resin-modified glass-ionomer restorative materials should include increased water exposure time to equilibrium conditions.

Insight into the behavior of increasing solubility over the first month could be provided by the lower-kinetic reactions that continue with glass-ionomer materials. It has been reported that with conventional glass-ionomer materials both compressive and flexural strengths decline after a few months of water storage. 73,74 Although these changes were first thought to be due to hydrolytic instability, 73,74 it was reported later by Nicholson and Abiden that these changes occurred even with anhydrous storage in oil. 75 It is thought that the observed strength degradation is due to increased material brittleness resulting from the continued slow reactions that lead to higher polymer cross-link density within the glass-ionomer material. These continued reactions in addition could create internal stresses of such magnitude that stress release is afforded by the creation of internal defects that communicate to the surface.

Although this remains conjecture, such communications could allow additional increase of surface area exposed to water and result with increased solubility. If the reactions were to cease after a certain period of time, no further stresses would be generated and solubility at that point would either plateau or decline. Decline could be possible if the proposed surface alterations could be healed, as resin-modified glass-ionomer restorative materials have been

reported to self-heal cracks upon exposure to water.^{67,76} Nevertheless, the solubility results provides strong evidence that the continuing polyalkenoate reaction can be susceptible to solution effects for up to one month after fabrication.

Microhardness

The use of hardness to observe the setting and curing behavior of resin-modified glass-ionomer materials has been investigated, ^{59,76-83} as the change in hardness may reflect a material's state of cure and possibly demonstrate the continuation of a setting reaction. ⁸¹ Xie *et al.* ⁸² evaluated the physical properties of glass-ionomer materials that included Knoop hardness after seven days storage and found that as a whole, the resin-modified glass-ionomer materials exhibited lower Knoop hardness numbers than the conventional materials.

Swift *et al.*⁵⁹ reported the depth of cure of different resin-modified glass ionomer materials using Knoop Hardness for up to seven days. These researchers evaluated the depth of cure of Fuji II LC, Photac-Fil, Vitremer, and two compomers (Geristore and Variglass) for seven days with each material evaluated in depths up to nine millimeters. Reported results stated that all materials could only be analyzed at depths five millimeters and less, as depths greater than five millimeters were too soft to obtain KHN readings. All of the resin-modified glass-ionomer materials (except Photac-Fil) exhibited an increase in hardness at all depths over the seven-day storage time. Interestingly, at the end of seven days Fuji II LC was reported to have a higher KHN at the 5 mm depth than the top surface, as with Vitremer when evaluated at one day.⁵⁹ Their results contrasted with the findings of this thesis, that the bottom hardness did not achieve the same KHN of that seen of the top surfaces. Although they reported that all of the resin-modified glass-ionomer products exhibited satisfactory depth of cure at five millimeters, it was still recommended that these products be placed in thicknesses of three millimeters or less.⁵⁹

Yap⁷⁸ compared Vickers hardness of different glass-ionomer, composite, and compomer restorative materials for up to six months. In this report, Fuji II LC and Photac-Fil demonstrated


significantly greater hardness at 24 hours than at six months, which was contrasted by Vitremer that initially demonstrated lower hardness but demonstrated increasing hardness with storage time and was comparable in hardness to Fuji II LC and Photac-Fil at the six-month evaluation period. The findings of this thesis found similar results with Vitremer whose top-surface KHN increased over a period of time, but this study contrasted in that Photac-Fil Quick top-surface KHN was demonstrated to increase to a maximum hardness at one month. Fuji II LC top-surface hardness was found to maintain a steady value over the evaluation period.

Gladys *et al.*⁷⁶ reported Vickers hardness determinations of Fuji II LC (both encapsulated and hand mixed), Vitremer, and Photac-Fil as the part of an extensive evaluation of restorative material physical property determinations. These investigators reported that the Vickers hardness of the three materials were somewhat similar, with Fuji II LC exhibiting slightly more hardness values than Vitremer, which in turn demonstrated slightly more hardness values than Photac-Fil.⁷⁶ In this thesis study, Vitremer was found to exhibit mean top-surface KHN values slightly more than Fuji II LC, which in turn demonstrated higher values than Photac-Fil Quick.

Ellakura *et al.*⁷⁹ compared the surface microhardness of resin-modified and conventional glass-ionomer materials as a function of time and water exposure. These researchers followed Vickers microhardness at 1, 7, 15, 30, 90, 180, and 365 days of Fuji II LC, Photac-Fil, and Vitremer as well as Ketac-Fil, Ketac-Molar, and Ketac-Silver. Results found that the resin-modified materials all demonstrated significantly lower microhardness than some of the conventional glass-ionomer materials (Ketac-Fil and Ketac-Molar).⁷⁹ Among the resin-modified materials, all materials demonstrated lower hardness at the 24-hour time period but then demonstrated thereafter an increase in hardness. However, Photac-Fil hardness peaked at 15 days while the hardness of Vitremer and Fuji II LC peaked at 60 days. After this time period, the surface hardness of the materials decreased in various amounts up to one year.⁷⁹ Results of this thesis corroborated with their results in that most of the RMGI materials demonstrated an

increase in top-surfaced KHN but differed in that materials demonstrated (except Fuji II LC) maximum hardness at three months.


Kanchanavasita *et al.*⁸³ evaluated the surface hardness of two restorative (Fuji II LC and Vitremer) and two liner resin-modified glass-ionomer materials for up to one year with samples stored both in distilled water and artificial saliva. In this report, both of the restorative materials stored in distilled water reached maximum hardness at 24 hours, and it was reported that the hardness of Fuji II LC and Vitremer were statistically similar. The restorative materials maintained their surface hardness for the duration of the evaluation. It is interesting to note that storage in artificial saliva had a detrimental affect on the maintenance of surface hardness for all of the materials evaluated.⁸³ As with the study by Ellakura *et al.*,⁷⁹ results of this thesis differed in that the materials evaluated reached a maximum KHN value at a later period in time.




The top surface hardness of the Vitremer samples over the storage periods are depicted

graphically in Figure 16. Within each thickness, the top surface KHN of the Vitremer samples underwent a significant decrease in KHN at one week, which was followed by a significant


increase in hardness at the one-month storage period. At three months, the top surfaces of all of the specimens were similar in hardness to that seen at 24 hours. Analysis of the Vitremer bottom surface KHN over the storage time periods are depicted in Figure 17. Similar to the top surface, the bottom surfaces of all the thicknesses demonstrated a significant decrease in hardness at one week, but all then demonstrated an increase in hardness at one month that was similar to 24 hours. All thicknesses except the 4 mm group exhibited bottom surface hardness at the three-month storage period that was statistically similar to that seen at 24 hours. All of the thickness groups demonstrated significant differences in hardness, with the 2 mm group demonstrating the highest KHN that progressively decreased to the 5 mm group, which had the lowest mean KHN.




The Fuji II LC top surface mean KHN segregated by sample thickness are depicted in Figure 18. The top surfaced displayed overall stability in hardness over the storage period, except for the 4 mm samples which demonstrated a significant decrease in hardness at one month. However, recovery of hardness was noted at three months. There was no difference in top surface KHN found between the different thicknesses over the storage time.




Similar analysis of the Fuji II LC bottom surface hardness by thickness can be observed in Figure 19. The bottom surfaces of the 2 mm samples remain stable, whereas the 3 mm samples undergo a significant increase in hardness at each observation period throughout the storage period. The 4 mm samples display a significant increase in hardness at the one-week time, but the hardness decreases significantly at one-month and then regains hardness to the one-week level at three months. The 5 mm samples also demonstrate a significant hardness increase at one week, but then decreases to the initial 24-hour level. Two-way ANOVA identified that all of the thicknesses demonstrated significant difference in hardness; as the thickness increased the magnitude of the bottom surface hardness decreased.



The Photac-Fil Quick top surface mean KHN are depicted in Figure 20. All of the top surfaces demonstrated progressive significant increases in KHN throughout the storage period, most notably between the 24-hour, one-week, and one-month storage periods for the 2 mm and 3 mm thicknesses; between one-week, one-month, and three-months storage for the 4 mm samples; and between the one-week and one-month storage for the 5 mm samples. Two-way ANOVA did not find a significant difference between the mean hardness of any of the top surfaces over the three-month storage period (p = 0.825)




The Photac-Fil Quick bottom surface mean KHN segregated by sample thickness are depicted in Figure 21.




Almost mirroring the KHN results of the top surface, all of the bottom surfaces of the various thicknesses exhibit progressive and significant increase in hardness up to one-month point and then demonstrate a hardness decrease, though not always significant. Furthermore, two-way ANOVA found significant differences between the thicknesses, with KHN decreasing as the thickness increased.

The use of surface hardness in the determination of composite resin restorative material depth of cure has largely utilized Knoop hardness, which has been reported to correlate well with methacrylate degree of conversion determined with FTIR as well as other material physical properties. Furthermore, the comparison of the top- and bottom surface Knoop hardness ratios have been used as a method for the determination of adequate visible-light-curing resin polymerization, with a bottom/top ratio of at least 0.80 suggesting an adequate cure. Prior published research has not explored the use of a similar ratio suggesting adequacy of cure with resin-modified glass-ionomer materials. Analysis of the KHN bottom/top ratios of the resin-modified glass-ionomer materials was investigated to identify any possible trends that could be associated with these materials.

Vitremer samples of all thicknesses demonstrated a diminishing trend of the bottom/top ratio (Figure 22) from 24 hours to one week after which the 2 mm and 3 mm samples experienced some increase of the ratio, nearly coinciding at one month. The 2 mm samples did not experience a significant change (p = 0.195) during the storage period whereas the 3 mm group demonstrated a significant loss (p < 0.0001) at one week but regained at one month to be significantly similar to the 2 mm samples (p = 0.829). The 4 mm and 5 mm samples, although displaying a decrease in the bottom/top ratio, did not experience a significant loss (p = 0.206, 0.500, respectively). At the end of the three-month storage period, none of the KHN bottom/top ratios were statistically similar (p < 0.0001).




Graphical depiction of the Fuji II LC bottom/top KHN ratio over time per thickness and storage time can be observed in Figure 23. The 2 mm samples displayed a stable ratio throughout the storage period, whereas the 3 mm samples displayed significant increases of the KHN ratio throughout the observation time (p < 0.0001). The 3 mm samples at one week displayed a mean KHN ratio that was statistically similar to the 2 mm samples (p = 0.357). Furthermore, the three-month KHN ratio for the 3 mm samples interestingly also increased to the level that was statistically similar (p = 0.233) to that of the 2 mm samples. The 4 mm samples demonstrated a significant increase at one week (p = 0.002) that was followed by an insignificant decline at one and three months. The 5 mm samples similarly demonstrated a significant KHN ratio increase at one week (p = 0.002) that was also followed by a diminishing of the ratio that was found to be similar to the ratio that presented at 24 hours.



At the end of the three-month storage period the 2 mm and 3 mm samples had similar KHN bottom/top ratios but were statistically greater than the 4 mm and 5 mm samples (p < 0.0001).

The Photac-Fil Quick bottom/top KHN ratios are presented in Figure 24. Interestingly, the ratio of the 2 mm and 3 mm samples are statistically similar (p = 0.304) which is maintained until one month. All the thicknesses demonstrate a trend for increase in the KHN ratio at one week, after which the 2 mm and 3 mm samples continue to increase until the one month time period. The 4 mm samples decline after one week whereas the 5 mm samples demonstrate a significant decline (p = 0.021) in the KHN ratio after one week. At one month and at three months, the 4 mm and 5 mm samples continue on a downward trend. At the end of the three-month period, the Photac-Fil Quick different thicknesses, similar to Vitremer, are not statistically similar (p < 0.0001).



The results of this study corroborate the findings of earlier studies. Similar to Yap,⁷⁸ Ellakuria *et al.*,⁷⁹ and Kanchanavasita *et al.*⁸³ the RMGI materials used in this evaluation demonstrated a variation in top surface KHN with storage time. However, where this study mainly differs from the previous studies is the time in which the KHN variation occurs. Yap reported that Fuji II LC and Photac-Fil Quick demonstrated the highest top-surface Vickers hardness at 24 hours while Vitremer had significantly higher hardness at 6 months.⁷⁸ Ellakuria *et al.* reported that Photac-Fil Quick demonstrated its maximum hardness at 15 days whereas Fuji II LC and Vitremer increased in hardness to maximum levels at two months.⁷⁹ Kanchanavasita *et al.* reported that Fuji II LC, Photac-Fil Quick, and Vitremer all exhibited maximum hardness at 24 hours.⁸³ In contrast to these reports, the findings of this study did demonstrate an increase in top KHN numbers, but the maximum hardness for both Vitremer and Photac-Fil Quick were observed to occur at one month. Also, contrary to earlier reports, Fuji II LC did not demonstrate a significant change in top surface KHN over the evaluation period, although an increasing trend was noted to occur at one week for some thicknesses. It should be noticed that all earlier studies involving

hardness of RMGI materials mainly involved analysis of development of just top-surface hardness. Only this report and the study published by Swift *et al.*⁵⁹ involved analysis of RMGI bottom surface hardness.

The overall results of the hardness testing of the resin-modified glass-ionomer restorative materials provide some evidence of a continuing setting reaction and reinforces the concept that the addition of resin does not totally inhibit the potential for the polyalkenoate reaction to continue. Young ²⁵ reported in a FTIR study of one-millimeter-thick, resin-modified glass-ionomer products that approximately 97 percent of the methacrylate reaction was complete after 150 seconds. Accordingly, this FTIR study provides emphasis that any continuing setting reactions in the upper regions of RMGI materials after this time are predominately the acid-base polyalkenoate setting reaction. The deeper regions are similarly assumed to predominately possess the polyalkenoate acid-base reaction, as all materials were observed to demonstrate an increase in bottom-surface hardness in varying degrees over the course of this study.

Both Vitremer and Photac-Fil Quick demonstrated a peak in bottom-surface KHN at one month, whereas Fuji II LC bottom-surface KHN appeared to be at its maximum (except for the 3 mm bottom surface) at one week. Thereafter, the bottom surface KHN of all the materials demonstrated a decrease; however the hardness was not observed to fall (except for the Vitremer 4 mm samples) significantly below the hardness exhibited at 24 hours. These findings do reinforce the concept that a lower-kinetic reaction does take place in the deeper regions of the resin-modified glass-ionomer materials. However, the decrease in hardness values after a certain time indicates, especially in the 4 mm and 5 mm samples, that the process reaches an endpoint.

Whereas the Vitremer and Fuji II LC top surface KHN remained stable, the top surface KHN of the Photac-Fil Quick samples demonstrated an increase of hardness with time that approximately mirrored the development of KHN in the bottom surfaces. The development of

hardness in the top surface of the 2 mm Photac specimens displayed significant correlations (α = 0.01) with the development of hardness in the bottom surfaces of the 3 mm (Pearson's 0.945, p < 0.0001), 4 mm samples (Pearson's 0.816, p < 0.0001) and 5 mm samples (Pearson's 0.686, p = 0.001). These correlations provide evidence that the setting reaction that enables the hardness development is possibly similar, as the bottom and top hardness of the Photac-Fil Quick samples appears to be developed by the same kinetics.

In spite of these signs of continued reactions, none of the materials developed a bottom surface KHN in the thicker samples (4 mm and 5 mm) that approximated the hardness observed in the 2 mm thick samples. This is made most evident in the examination of the KHN bottom/top ratios over the storage periods. Except for Vitremer, the ratios do reflect initial signs of a continued setting reaction in the 4 mm and 5 mm samples, but the trend is reversed at one week for both Fuji II LC and Photac-Fil Quick.

Using the ratio of the bottom KHN to the top KHN of the materials allows some interpretation of the overall effect of the development of KHN of the resin-modified glass-ionomer materials. As this ratio has been utilized extensively in the study of resin composites, ⁸⁴⁻⁸⁸ it has not been reported previously involving resin-modified glass-ionomer materials. Fuji II LC and Vitremer presented stable KHN ratios for the 2 mm samples over the storage time whereas Photac-Fil Quick demonstrated a progressive improvement in KHN ratio as time progressed up to the one-month storage period, after which a decrease was observed. The Photac-Fil Quick 3 mm KHN ratio mirrored the increase in KHN ratio of the 2 mm samples, although the values were significantly less. For Vitremer and Fuji II LC, the 3 mm KHN ratios overall increase with time, reaching statistical equivalence to the 2 mm groups briefly at one month for Vitremer and at one week and 3 months for Fuji II LC.

Both Photac-Fil Quick and Fuji II LC display the same trends in KHN bottom/top surface ratios of the 4 mm and 5 mm samples, demonstrating an increase in mean KHN ratios up to one week

of storage which drops afterwards. This similarity in hardness ratio behavior may suggest that these may be similar in composition and post-VLC setting reaction behavior in the deeper regions of the materials.

The results of KHN ratio with the Vitremer samples were somewhat surprising in nature. Vitremer is said to possess, in addition to a VLC and acid-base cures, a "dark cure" that is propagated by a reduction-oxidation reaction. Although some have classified this as a "tricure" curing mechanism, some authors take exception to this nomenclature, stating that the setting mechanism obtained by this method is the same obtained with photoactivation in which polymerization is initiated with a different method. Nevertheless, in deeper regions of Vitremer where light activation is insufficient to enable photoactivation, the combination of the additional "dark cure" and the acid-base polyalkenoate reaction would be expected to result in a hardness increase with longer storage time. Although an increase in hardness was observed in the bottom surfaces of the 4 mm and 5 mm Vitremer samples, the hardness that developed was not of sufficient magnitude to positively affect the KHN ratio. Accordingly, the ratio decreased from 24 hours to one week and remained relatively stable after that point.

Correlation analysis of Vitremer KHN ratios and solubility results did not find a correlation between KHN ratios and experimental solubility (p = 0.505). Correlation analysis between Fuji II LC experimental solubility results and KHN ratio results identified a correlation (Pearson's 0.318, p = 0.004) as well as with Photac-Fil Quick (Pearson's 0.245, p = 0.028). Although some correlations were reported using this set of data between the development of KHN ratios and experimental solubility, it should be remembered that the Pearson's correlation is considered to be liberal, tending to overestimate the true reliability between the data.⁸⁹ This is reinforced by that linear regression identified no strong reliability between sets of variables for both Fuji II LC and Photac-Fil Quick ($R^2 = 0.101, 0.06$, respectively).

The KHN ratio for the 4 mm and 5 mm thicknesses in all materials evaluated were observed to demonstrate a decreasing trend with an increase in storage time. This decrease was noted between 24 hours and one week for Vitremer, and the decreasing trend initiated after one week for Fuji II LC and Photac-Fil Quick. One explanation into this decreasing trend could possibly be afforded due to interruption of any potential continuing reactions due to water infiltration. Resin-modified glass-ionomer materials have been demonstrated to absorb water by both Fickian linear and non-linear diffusion behavior that required longer periods of time to reach equilibrium than conventional glass-ionomer materials. 70-72 It has been reported that RMGI materials of two to three millimeters thickness required up to eight 1 to ten days 1 for water infiltration to reach full equilibrium, and the decrease in KHN at or about this same time for the RMGI materials observed in this study could be due to water infiltration that dissolves and/or plasticizes either unreacted components or any matrix that may have formed.

In spite of evidence of a continuing setting reaction within the deeper regions of the resin-modified glass-ionomer materials, there is ample evidence to suggest that any reaction that may proceed in thicknesses greater than 3 mm is insufficient to make up for the amount of cure observed in 2 mm thicknesses. Using the KHN ratio as a guide, only the 3 mm Fuji II LC samples evidenced equivalence to the 2 mm samples at three months, and although Vitremer 3 mm KHN ratio samples were statistically similar to the 2 mm group at one month, the ratio was observed to decrease thereafter. It is only speculative at this time to assume that the KHN bottom/top ratio is a valid method in the analysis of resin-modified glass-ionomer depth of cure, and further research is needed in this area.


Thermal Analysis

The process of measuring the amount of heat involved in a chemical or physical change is known as calorimetry. 90 Using calorimetry the resin-modified glass-ionomer materials were analyzed as per specific heat capacity, thermal scan enthalpy, thermal scan temperature peak, and weight changes during calorimetry.

Specific Heat

The specific heat capacity of a material (often referred to as the "specific heat", or "heat capacity") is a measure of the amount of energy required to raise the temperature of one gram (or one mole) of the material by one degree Kelvin. 90,91 The specific heat of a given quantity of matter generally describes a property of the entire quantity of the substance. 90 The specific heat capacity of a substance can be measured from the sample mass and the difference in the heat flow rates determined during a programmed temperature sequence consisting of an isothermal period, followed by a rise in temperature, and then finally another isothermal period. Frequently this is also performed using a standard reference material for calibration. 92 Using a DSC there are usually four methods of specific heat capacity determination: Direct method, Steady State, Sapphire Method, and Alternating DSC (ADSC). The Direct Method is based on the dynamic method of specific heat and is calculated by dividing the heat flow by the heating rate and sample mass. The Direct Method can be accurate to within ± 5% and requires the shortest measurement time. 91 The Steady State is similar in computation as the Direct Method but subjects the sample to a periodic heating and cooling cycle and has an accuracy rate similar to the Direct Method. 91 The ADSC method varies the temperature sinusoidally as a function of time which is superimposed on the average heating rate. The ADSC method is slightly more accurate than the Direct and Steady State, but requires long measurement times. The most accurate specific heat capacity determination, known as the sapphire method, calibrates the thermal signals of a substance against a three-sapphire-disc calibration standard. 91 Due to software limitations with the DSC unit available for this research, the Direct Method for specific heat capacity was used during this evaluation.

The specific heat determinations for this evaluation followed a similar protocol reported by Khalil and Atkins.⁶¹ These investigators scanned a weighed sample of glass-ionomer product between 20 °C and 60 °C at a rate of 10 °C per minute and reported the observed specific heat at 37 °C. The evaluation in this thesis followed a similar protocol but recorded the specific heat as the mean of specific heat determination values at 5 °C intervals above the latent period over the thermal scan. As the specific heat capacity of a material is also dependent on the temperature recorded ^{90,92} it was felt that a mean of these values over a small thermal range would be representative. Under their protocol, Khalil and Atkins reported the specific heat capacity of Fuji II LC light cured to be 0.972 J/g°C and cured under dark conditions to be 1.267 J/g°C.⁶¹ The results for the Fuji II LC immediate samples under the present evaluation were found to be very similar to that reported by Khalil and Atkins.




The results of specific heat analysis of the Vitremer samples can be seen in Figure 25.


All Vitremer samples were observed to increase significantly in specific heat as time progressed. At the end of three-months, there was no significant difference in specific heat

between the 2 mm and 3 mm samples (p = 0.100) and between the 3 mm and 4 mm samples (p = 0.323). However, there was a significant difference found between the 2 mm and 4 mm samples (p = 0.007).

The specific heat results for Fuji II LC can be seen in Figure 26. Similar to that seen of Vitremer, the Fuji II LC samples demonstrated an overall increase in specific heat with storage time. At one week storage, there were significant differences found between the 2 mm and the other thicknesses (p < 0.005) but not between the 3 mm and the 4 mm samples (p = 0.182). This relationship continued until at three months, where there were no significant differences found between any of the thicknesses (p = 0.451).



The results of the Photac-Fil Quick specific heat can be seen in Figure 27.



As with the other two materials, Photac-Fil Quick demonstrated a generalized increase in specific heat over the storage time. At 24 hours, there was no difference noted between the 2 mm and 3 mm samples (p = 0.08) with differences noted between the 3 mm and 4 mm thicknesses (p < 0.001), while at one week specific heat values were not significantly different (p = 0.056). With passage of time the Photac-Fil Quick 2 mm samples increased in specific heat to demonstrate a significant difference with the 4 mm thickness (p = 0.019) while there remained no difference between the 2 mm and 3 mm (p = 0.265) or between the 3 mm and 4 mm (p = 0.297). At three months, the 3 mm and 4 mm samples maintained similar specific heat (p = 0.070) but the specific heat of the 2 mm samples was significantly increased compared to the 3 mm and 4 mm thicknesses (p < 0.0001).

Comparison of the specific heat results of all three materials demonstrate an increase in specific heat with increased time which suggests a change in the respective materials. At the end of the storage time all three materials did not exhibit a significant difference in the specific heats between the 2 mm and 3 mm samples as well as between the 3 mm and 4 mm samples.


Except for Fuji II LC, the samples did maintain a significant difference between the 2 mm and 4

mm samples. However, the results with Fuji II LC samples need to be interpreted with caution, as the 2 mm group demonstrated wide variation.

Thermal Scan


The differential scanning calorimetry (DSC) thermal scan method employed during this evaluation followed the same protocol as Khalil and Atkins. Accordingly, after the specific heat determination the samples were weighed and then submitted to a thermal scan of 37 °C to 240 °C at a rate of 10 °C per minute. A reverse scan from 240 °C to 37 °C was not accomplished as pilot studies with all three materials did not reveal any results. Results of the thermal scan are reported as both thermal peak temperature, i.e. the temperature at the minimum of the endothermic peak, and total enthalpy of the peak. This was accomplished by software integration, with all samples normalized for sample size. For graphical and analysis simplicity, total enthalpy is presented as a positive number, although the peaks were endothermic, meaning the value was really negative.


Graphic results for the Vitremer samples thermal scan is presented in Figure 28.





As storage time progressed all samples demonstrated an increase in thermal enthalpy that was associated with the thermal peak. One-way ANOVA of the results at each storage time revealed that there was no significant difference for enthalpy findings between the 3 mm and 4 mm samples throughout the time period. The 2 mm samples were significantly different than the 3 mm and 4 mm samples at all time periods (p < 0.004) except for at three months, where there was no significant difference noted between all three thicknesses (p = 0.665).

Representative images of the Vitremer thermal scans can be seen in Figures 29 – 32. It can be observed as the materials mature the individual peaks evolve from broad, ill-defined endothermic curves into more uniform, symmetrical curves whose peak temperatures can be readily defined. This thermal behavior is consistent with materials that are becoming more uniform. As will be discussed next, it can also be observed that the thermal peak temperatures display a trend towards migrating to a common value.










The analysis of the thermal peak temperatures for the Vitremer samples can be seen in Figure 33.




The mean thermal peak temperature of the 3 mm and 4 mm samples were statistically similar throughout the evaluation (p > 0.272). However, the immediate 2 mm samples demonstrated a significantly higher thermal peak temperature than the 3 mm and 4 mm (p < 0.0001) immediately after fabrication and significantly reduced than the other thicknesses at 24 hours (p < 0.0001). As time progressed the 2 mm samples demonstrated a rising trend that tended to approach the thermal peak temperatures of the 3 mm and 4 mm samples. However, at three months the 2 mm samples still maintained a significant difference (p < 0.012).


The mean enthalpy values of the Fuji II LC thermal scans are depicted in Figure 34. All samples demonstrated a significant increase in enthalpy as storage time progressed (p < 0.0001) with all samples demonstrating the most increase between fabrication and one week.


After one week, there was no difference in mean thermal scan enthalpy between the remaining

storage time intervals (p > 0.187). Significant differences were found between the enthalpy requirements of the 2 mm and the other samples (p = 0.002) but no difference in enthalpy was noted between the 3 mm and 4 mm samples (p = 0.975).




Representative thermal scans of representative Fuji II LC samples are depicted in Figures 35 and 36. As with the Vitremer, the Fuji II LC samples display some increase in uniformity with aging and the migration of the thermal peaks towards a common value.






The mean thermal peak temperatures of the Fuji II LC samples are depicted in Figure 37. With this material, the 3 mm depth samples displayed the highest mean thermal peak temperatures; the 2 mm samples the lowest temperature, while the 4 mm samples an intermediary temperature between the two extremes. The immediate Fuji II LC thermal peak temperatures were not significantly disparate (p = 0.301) as well as at 24 hours (p = 0.053). At one week, the 2 mm and 3 mm samples were significantly different (p = 0.003) which was maintained at one month (p = 0.004). At three months, the mean thermal peak temperatures of the different thicknesses were merging and were statistically similar (p = 0.116).





All of the Photac-Fil Quick samples (Figure 38) demonstrated a significant increase in thermal scan enthalpy requirements as storage time increased (p < 0.0001). The thermal scans of the samples demonstrated a significant increase between fabrication and 24 hours (p = 0.001) after which the loss stabilized with no significance between groups until one month (p > 0.431).


Although the samples demonstrated an increasing trend after one month, there was no difference noted between the one-month and three-month samples (p = 0.079). Throughout the storage period the 2 mm samples required significantly less enthalpy than the other thicknesses (p < 0.0001) while no difference was noted between the 3 mm and 4 mm samples (p = 0.926).




Representative thermal scans of Photac-Fil Quick can be seen in Figures 39 – 42. In contrast with the other two materials, the 3 mm and 4 mm Photac-Fil Quick samples displayed uniform endothermic peaks from the outset whose peaks became more focused with maturity. The 2 mm Photac-Fil Quick endothermic peak was nondescript at 24 hours but gained uniformity as time progressed and the peak temperature was also observed to shift towards that observed by the other two thicknesses.










The mean thermal peak temperatures of the Photac-Fil Quick samples are depicted in Figure 43. Similar to Fuji II LC, the 3 mm samples demonstrated the highest thermal peak temperature with the 2 mm samples the lowest. The 4 mm samples maintained a thermal peak temperature that was intermediary to the other thicknesses. The 2 mm samples were significantly lower than the 3 mm samples (p = 0.007) as well as the 4 mm samples (p = 0.047) but the 3 mm and 4 mm groups were statistically similar (p = 0.569). As with the other two materials, all three thicknesses of the Photac-Fil Quick samples demonstrated a merging trend of the peak temperatures with the passage of time. At 24 hours there was no difference between the samples (p = 0.121), as well as at one week (p = 0.151). However, at one month the merging trend was disrupted by a significant difference between the 2 mm and the other thicknesses (p = 0.033) which continued at three months (p = 0.002).



All three resin-modified glass-ionomer restorative materials demonstrated increasing thermal requirements for the DSC thermal scan with storage time, although all achieved such displaying distinct thermal behavior. Vitremer samples displayed an increasing rate throughout the


evaluation while Fuji II LC demonstrated a rise in thermal requirements at one week that essentially reached a plateau for the remainder of the observation time. Photac-Fil Quick reached a similar plateau at one week, which was maintained at one month but then evidenced an increasing trend for thermal requirements at one month. Collectively the data indicates different thermal behavior for all of the materials with increased storage time, which is evidence that all of the materials are exhibiting change in composition. The individual thermal behavior of these materials will demand more attention and evaluation: it remains to be determined if the early plateau reached by Fuji II LC is evidence of the formation of a mature and/or stable composition. As glass-ionomer materials are known to require considerable time to reach maturity, only longer periods of observation will be able definitively demonstrate this behavior by thermal means.

The peak temperature of the DSC thermal peak provides additional information concerning these RMGI materials. The peak temperature of a thermal measurement during DSC can provide material composition characteristics almost in the same manner as glass transition temperatures and melting points. It is interesting that all three of the resin-modified glassionomer materials demonstrated similar behavior in regards to thermal peak temperatures during the thermal scans. Each material demonstrated, especially with the 3 mm and 4 mm samples, higher mean thermal peak temperatures immediately after fabrication that slowly lessened with storage time. Moreover, as storage time increased, the temperatures also displayed a trend of merging towards the same range of temperatures, to the point that at intermediary times there was no significance between the peak temperatures of the Photac-Fil Quick samples, while at three months there was no significant difference between the mean thermal peak temperatures of the Fuji II LC samples. The merging of the mean peak temperature values provides additional evidence of material compositional change with the resin-modified glass-ionomer materials. Furthermore, as the peak temperatures approximate each other it suggests that the composition of each respective material is becoming more similar.


The RMGI materials were all evaluated for the different weight changes during all phases of thermal analysis, as investigation into these weight changes could provide information concerning the water uptake, ratio of loosely bound to firmly bound water, as well as possible polymer development in these materials.

Not surprisingly, all of the materials were observed to significantly increase in weight as storage time progressed. Any weight gain during storage in the 98 ± 2 percent humidity conditions was assumed to be due to water uptake. As with conventional glass-ionomer materials, RMGIs have been documented to absorb water, as water uptake is essential for the continuation of the polyalkenoate reaction. Furthermore, all three materials were also noted to experience significant weight loss during the specific heat determinations. Any weight loss observed during the specific heat determinations was assumed to be predominately due to water loss, presumably of the loosely bound nature. This assumption is founded on the short time duration of the specific heat determination scans (20 minutes) and the narrow temperature range involved (20 °C to 60 °C) that was decidedly below the boiling point of both water and HEMA.

An interesting finding regarding the weight loss during the specific heat determinations (presumed to be loosely-bound water) among the three materials is that each material demonstrated different behavior in the nature of the loosely-bound water loss. The Vitremer (Figure 44) 2 mm and 3 mm samples lost significantly more water at three months than at the onset (p < 0.0001) while the 4 mm samples demonstrated a similar trend towards increased water loss, but yet is was not significant. The 2 mm samples demonstrated a trend for decrease in water loss after 24 hours that continued until one month, with the 3 mm samples displaying similar behavior after one week. At one month there was no significant difference in water loss between all thicknesses, and although the 2 mm and 3 mm samples displayed a greater amount of water loss compared to the 4 mm samples at three months, no significant difference remained among the groups (p = 0.357).




The Fuji II LC samples (Figure 45) also demonstrated an increasing trend of water loss during the specific heat determinations with storage time but did so in a manner different than Vitremer.



At the end of the three months, each thickness of Fuji II LC samples did not experience a significant difference in water loss compared to the water loss with the immediate groups (p > 0.437). However, during the storage time all three thicknesses experienced a significant increase in water loss (p < 0.001) that reached a maximum at one week and then declined. Except for the immediate groups, the Fuji II LC 2 mm samples lost significantly more water than the 3 mm and 4 mm groups (p < 0.0001), while the latter two groups were not significantly different (p = 0.459).

The water loss behavior during the specific heat determinations for the Photac-Fil Quick samples is displayed in Figure 46.




The Photac-Fil Quick samples displayed an increase in water loss with time similar to Fuji II LC except that the maximum water loss peaked at the 24-hour storage time. After one week the 2 mm samples continued on an increased trend for water loss while the other thicknesses maintained relatively stable. One-way ANOVA identified that over the total storage time the 2 mm Photac-Fil Quick samples lost significantly more water during the specific heat determinations (p < 0.0001) compared to the 3 mm and 4 mm samples which were statistically similar (p = 0.219). Within the Photac-Fil Quick groups the 2 mm and 4 mm samples lost significantly more water at three months compared to the immediate samples (p < 0.002), while the 3 mm samples maintained statistical similarity (p = 0.054). However, all thicknesses demonstrated a significant increase in water loss at the 24-hour peak as that experienced by the immediate groups.

The difference between the weight gain during storage (water uptake) and weight loss during the specific heat determinations (loosely held water loss) results in the net retention of water,

and provides some insight into the development of the more firmly bound water associated with the maturation of the polyalkenoate hydrogel matrix.


All three materials demonstrated a significant change in the net water retention during the specific heat determination for the storage times evaluated. The results of all thicknesses of the



Vitremer samples are graphically depicted in Figure 47. Two-way ANOVA identified that the Vitremer samples demonstrated a significant increase in net water retention both per thickness and storage time (p < 0.0001). One-way ANOVA for net water retention between the thicknesses found no significant difference between the 2 mm and 3 mm samples (p = 0.241) but identified significantly more net water retention between the 4 mm and the 2 mm samples (p = 0.043). Although the 2 mm samples gained weight during storage time, Scheffe *post hoc* analysis did not identify any significant change in net water retention due to large standard deviations among the sample means. A significant change was identified within the 3 mm samples (p < 0.0001), with significant net retention of water increasing at all storage times except between the immediate and 24 hour samples. A significant change was also noted

within the 4 mm samples (p < 0.0001), with significant gains identified between storage times after one week. The data indicates that at approximately one week the Vitremer samples began to show evidence of some water retention within the material.


The mean net water retention results during the specific heat determination results for the Fuji II LC samples can be observed in Figure 48.



Two-way ANOVA identified significant changes in net water retention during specific heat determinations for the Fuji II LC samples for both thickness and storage time (p < 0.0001). All the thicknesses demonstrated an increasing trend for increased water retention with storage time that appeared to peak at one week or one month. Within the individual thicknesses, although the samples demonstrated a significant gain in water retention up to at least one week that was followed with a decrease at three months, the 2 mm samples still maintained a significant gain in water retention at three months as compared to the immediate group (p < 0.0001). All of the thicknesses displayed statistical similarity in net water retention immediately

after fabrication and at 24 hours, but at one week the 3 mm and 4 mm samples displayed significantly greater net water retention than the 2 mm specimens that was maintained up to three months (p < 0.001). The Fuji II LC samples demonstrated net water retention at approximately 24 hours that was earlier than that seen by the Vitremer samples.

The mean net water retention results during the specific heat determination results for the Photac-Fil Quick samples can be observed in Figure 49.



The Photac-Fil Quick samples displayed net water retention behavior during the specific heat determinations similar to that evidenced by the Vitremer samples. All samples demonstrated a significant increase in net water retention at three months (p < 0.0001) as compared to after immediate fabrication while the 3 mm and 4 mm thicknesses were found to exhibit significantly more water retention than the 2 mm samples (p < 0.006). The 3 mm and 4 mm Photac-Fil Quick samples exhibited net water retention starting at 24 hours while the 2 mm samples required one week.

Both Vitremer and Photac-Fil Quick displayed similar behavior regarding net water retention during the specific heat determinations. Fuji II LC demonstrated a different behavior marked by a trend of decreasing net water retention after a significant gain at one week. It is not known at this time if this trend would continue, but further analysis is planned to follow this material's behavior. It is this author's contention that this net water retention analysis provides some evidence of early mature hydrogel formation in these RMGI materials. Using this criterion, based on the conditions of this evaluation it appears that it requires one week for evidence of firmly-bound water in the matrix of both Vitremer and Photac-Fil Quick whereas evidence of the same with Fuji II LC is apparent at 24 hours.

In comparing the water loss and water retention behavior of the three RMGI materials, all materials displayed increasing water loss that was somewhat offset by increasing amounts of water retention as time progressed. It should also be noted that the thinner samples, especially the 2 mm thickness, appeared to demonstrate more water loss and less net water retention than the 3 mm or 4 mm samples. This behavior may be explained by the consideration that all of these samples were enclosed in aluminum calorimetry crucibles in which the RMGI materials displayed only one open surface to the calorimeter environment. With all three thicknesses having roughly the same surface area exposed, the thicker specimens would be expected to retain more water due to the time requirement for diffusion to take place.


Vitremer gave evidence of an interaction of water loss/net water retention behavior with a possible correlation between the two factors (Pearson's 0.316, p = 0.006) but this was discovered to be a weak correlation as linear regression did not identify any significant relationship ($R^2 = 0.10$). Moreover, the same behavior was seen with Fuji II LC (Pearson's 0.301, p = 0.009; $R^2 = 0.091$) but no such correlation was identified with Photac-Fil Quick (p = 0.450). It appears evident that there is no significant relationship between the water loss behavior and the amount of net water retained during the specific heat determinations. Certainly a relationship can be established between the amount of water gained during storage

and net water retention for all three materials using linear regression (Vitremer, $R^2 = 0.800$; Fuji II LC, $R^2 = 0.629$; Photac-Fil Quick, $R^2 = 0.747$) but this certainly begs the rationalization: water must be gained during the storage periods to be able for the water to be retained during the specific heat determinations.

The effect of net water content on the development of specific heat capacity must be considered, as mature glass-ionomer materials are postulated to contain up to five percent unbound water and anywhere from 18 to 28 percent bound water. Water has a reported specific heat capacity of $4.186 \text{ J/g}^{\circ}\text{C}^{90}$ and hence, more water content within a material could possibly affect its specific heat value. When analyzing any possible relationship between Vitremer water retention behavior and specific heat values, linear regression did not identify any strong relationships between net water retention ($R^2 = 0.014$) or water gain during storage ($R^2 = 0.184$). Similar results were found for Fuji II LC (net water retention, $R^2 = 0.042$; water gain, $R^2 = 0.370$) or Photac-Fil Quick (water gain, $R^2 = 0.0707$; net water retention, $R^2 = 0.011$). Interestingly, there appeared to be a stronger relationship between water loss during the specific heat determinations and specific heat values (Vitremer $R^2 = 0.540$; Fuji II LC $R^2 = 0.576$; Photac-Fil Quick $R^2 = 0.547$). The significance of these findings remains obscure at this time.


The resin-modified glass-ionomer materials were also analyzed for weight loss during the DSC thermal scan for identification of trends. In addition, the total weight loss during the specific heat determinations and the thermal scan was evaluated for similar reasons and these weight losses were compared to water gain during storage.

The Vitremer thermal scan weight loss data is depicted in Figure 50.




All thicknesses experienced significantly more weight loss during the thermal scan with increased storage time (p < 0.0001), with each thickness roughly displaying the same trend of weight loss. The 3 mm and 4 mm samples displayed parallel increasing weight loss behavior with increased storage time, but the 2 mm samples developed a decrease from one week to one month which then became almost identical to the 3 mm trend at three months. The 2 mm samples demonstrated significantly less weight loss during that thermal scan than the 4 mm samples (p = 0.006) while all other thicknesses were similar in weight loss.

The Fuji II LC thermal scan weight losses are posted in Figure 51. Considering total weight loss over the storage time, the Fuji II LC samples demonstrated a significant increase in weight loss (p < 0.0001). When evaluated between the sample thicknesses, the 4 mm groups demonstrated significantly more weight loss (p = 0.037) than the 2 mm samples, while no statistical difference was identified between the 2 mm and 3 mm as well as between the 3 mm and 4 mm samples (p > 0.153).



The results for the Photac-Fil Quick thermal scan weight loss can be seen in Figure 52.





The Photac-Fil Quick samples demonstrated a significantly greater weight loss at the end of the three month storage period as compared to the onset (p < 0.0001) with a significant difference in loss identified between the 2 mm and 3 mm samples (p = 0.017) and between the 2 mm and 4 mm samples (p = 0.001), but not between the 3 mm and 4 mm samples (p = 0.577).


All three RMGI materials displayed similar weight loss behavior during the DSC thermal process, demonstrating a significant increase in weight loss as storage time progressed. Interesting, review of the data will show that the 2 mm material samples between one week and one month either show a trend of either weight gain or at least a lessening in the rate of weight loss, but this trend is reversed after one month. As compared to the specific heat determination weight loss, it can be observed that the 2 mm thicknesses display in all materials significantly less weight loss than the 4 mm samples. This trend is the reverse from that observed during the specific heat determinations: a possible explanation can be afforded due to difference in thermal energy provided between the two methods and the development of more firmly-bound water as the RMGI polyalkenoate matrix matures. Accordingly, the enthalpy developed by the

thermal scan is much greater than that provided by specific heat determination and this increased energy could provide the impetus necessary to overcome some of the diffusion barriers that hinders the loss of loosely-bound water during the specific heat determinations. Furthermore, the 4 mm thicknesses could also be susceptible to more weight loss during the thermal scan due to loosely-bound water not removed during the specific heat determinations and the development of bound water as the polyalkenoate matrix matures. Based on some of the evidence provided by hardness testing, the deeper layers do not have the same degree of setting as the upper layers and the polyalkenoate reactions in the deeper layers have either been severely hampered or plasticized by the influx of water during storage. Hence, the polyalkenoate matrix in the deeper regions does not attain the level of maturity required to resist the forces of desiccation and is more prone to water loss.

The three resin-modified glass-ionomer materials were additionally analyzed for total weight loss experienced by the samples during the specific heat determinations combined with the DSC thermal scan, and the results for Vitremer can be seen in Figure 53. All of the Vitremer samples experienced significant weight loss between initial fabrication and three months (p < 0.0001), but there was no significant difference found between the thicknesses at the end of three months (p = 0.716). Interestingly, the results for Fuji II LC and Photac-Fil Quick were similar. All Fuji II LC samples (Figure 54) experienced significant weight loss between immediate fabrication and three months (p < 0.0001) but there was no significant difference in total weight loss between the thicknesses at three months (p = 0.254). The Photac-Fil Quick samples displayed similar weight loss behavior (Figure 55) demonstrating a significant increase in total weight loss over the total storage time (p < 0.0001) but without a significant difference identified between the different thicknesses (p = 0.548)








The total weight loss data provides some insight into the deposition and maturity of the water that may be contained within each resin-modified glass-ionomer material. The comparison of total weight loss to the thermal scan weight loss data will identify some differences, although in some instances the trends may appear to be similar. One striking difference between the two sets of data is that the thermal scan weight loss should involve only the more firmly-bound water that is incorporated in the polyalkenoate hydrogel. As this hydrogel becomes more mature, less water should theoretically be lost in the resin-modified glass-ionomer materials during the heat capacity determinations while more weight loss due to water could be expected to be lost during the thermal scan. For some cases, this has been found to be exactly the case during this three-month evaluation. However, continued evaluations for longer storage times are needed before definitive trends can be identified.


Analysis of the total weight loss of the materials will provide some identity to the nature of the total water behavior of the resin-modified glass-ionomer materials. Due to the maturation processes that occur with glass-ionomer materials in general, it would be expected that overall rate of water loss (as in weight percent loss over time) would decrease as storage time

progressed and the polyalkenoate matrix continues to mature. This behavior was observed at three months with Fuji II LC, Photac-Fil Quick demonstrated a stabilization trend of total weight loss in the evaluation midst, and Vitremer was still demonstrating an increasing weight loss at three months. It cannot be ascertained at present if the present findings will continue to display the same behavior or if the other materials will eventually provide similar total weight loss behavior as demonstrated by Fuji II LC. Further evaluations are ongoing for evaluation of weight loss for increased storage time.

To gain some insight into the effect of water gain during storage to the weight loss during the thermal determinations, the ratios of DSC thermal scan to water gain and total weight loss to water gain (expressed in percent) was analyzed. The immediate groups were excluded from this analysis. The results of the ratio between thermal scan weight loss and water gain for the Vitremer samples is depicted in Figure 56. The ratio of the weight loss during the thermal scans to water gain during storage for the Vitremer samples demonstrated a significant decrease as storage time progressed (p < 0.0001). There was no significant difference in the ratio reduction between the 2 mm and 3 mm thicknesses (p = 0.755) or between the 3 mm and 4 mm thicknesses (p = 0.189). However, a significant difference was identified between the 2 mm and 4 mm thicknesses (p = 0.042).


The results of the Fuji II LC ratio of thermal scan weight loss/water gain ratio can be seen in Figure 57.






One-way ANOVA evaluation of the Fuji II LC ratio of weight loss during the thermal scans to water gain identified a significant change (p <0.0001) within the samples when analyzed over the storage time. However, this was most punctuated by a difference noted between 24 hours and one week as seen by the 3 mm and 4 mm samples that demonstrated a precipitous and significant (p < 0.0001) decline in the ratio. However, after this observation there was no significant difference identified within the samples between one week and three months (p > 0.090). A significant difference in ratio behavior was noted between the 2 mm and 3 mm (p = 0.010) samples as well between the 2 mm and 4 mm samples (p = 0.004) but no difference was identified between the 3 mm and 4 mm samples (p = 0.941).

The Photac-Fil Quick samples (Figure 58) thermal scan weight loss to water gain ratio displays a decrease in ratio over time, similar to the case of Vitremer. A significant decrease in the overall



ratio was identified between 24 hours and three months (p = 0.002), however no significant difference was found between the thicknesses (p = 0.378).


An attempt at additional insight into the dynamics of the resin-modified glass-ionomer water behavior was done with evaluating the ratio of total thermal weight loss to water gain during storage. The analysis of Vitremer total thermal weight loss to water gain is graphically depicted in Figure 59. The behavior of this analysis is nearly the same as the thermal scan weight loss and water gain ratio.




As with the ratio of thermal scan weight loss to water gain, all samples demonstrated a significant reduction in the ratio with increased storage time (p < 0.0001). However, the total weight loss to water gain ratio differed from the thermal scan to water gain ratio analysis in that ANOVA found no significant difference in the ratio reduction behavior between the sample thicknesses (p = 0.253).

Analysis of the Fuji II LC ratio of total thermal weight loss to water gain (Figure 60) found that all samples displayed a significant difference in the ratio between 24 hours and three months storage (p < 0.0001) with the main difference noted between 24 hours and the other storage times, as there were no differences found after one week and three months (p > 0.145). As with

Vitremer there was no significant difference identified between the Fuji II LC thicknesses (p > 0.075).



As with the thermal scan loss to water gain ratio, the Photac-Fil Quick samples (Figure 61) demonstrated a significant reduction in total thermal weight loss to water gain ratio between 24 hours and three months (p < 0.0001, with no significant differences noted between the sample thicknesses (p = 0.429).



Water is an essential ingredient for all formulations of glass-ionomer materials as water is both a reactant in the early stages of the glass-ionomer setting reaction, but it also enables the formation of the silica-based hydrogel. As glass-ionomer materials mature, a certain water content becomes increasingly resistant to removal by desiccation and is referred to as being more "bound". The ratio to "bound" to "unbound" has been demonstrated to increase as glass-ionomer products mature. The exact timing of this conversion of the initial unbound water to the more removal-resistant bound water has not been ascertained, as the earliest detection of bound water has not been reported in the scientific literature.

The weight changes associated with thermal analysis of the resin-modified glass-ionomer restorative materials may provide some insight into the nature of water within these materials. If the assumptions of this thesis are correct, results of the net water retention analysis during the specific heat determinations have demonstrated evidence of net water retention as early as 24 hours for Fuji II LC and one week for both Photac-Fil Quick and Vitremer. It is unclear if this

early evidence of water retention is substantiation of the first formation of bound water in resinmodified glass-ionomer cements.

The ratios of weight loss during both the thermal scan and total weight loss to the water gain during storage was undertaken to possibly provide some erudition into relationship between the water gained during storage and the amount of weight loss during the thermal analysis. Another assumption made during the thermal analysis was that weight loss during the highertemperature thermal scans would predominately be due to the more firmly-bound water species. To assume that the weight loss would entirely be represented by water would be incorrect, as resin-modified glass-ionomer materials represent complex solutions of both polyalkenoate and methacrylate components, as well as photointiators and other modifiers.²¹ Some weight loss during the thermal scan with Fuji II LC and Photac-Fil Quick could be also attributed to unreacted methacrylate components, especially HEMA. However, due to the water-soluble nature of HEMA, it could be reasonable to assume that any unreacted HEMA remaining in the matrix might exist in solution with water during the thermal scan process. Since the methacrylate component is largely grafted upon the polyacrylic acid backbone in Vitremer, it is comfortable to assume that weight loss during the thermal scan could be more due to water loss than the other two materials. However, the reader must be reminded that these simplistic generalizations are made in order to attempt to observe generalized trends of the materials: the emphasis remains that these materials are complex entities and not mere simple mixtures of polyalkenoate components and methacrylates.

All materials demonstrated increased weight loss during the DSC thermal scan as storage time increased. This behavior is expected as all materials were demonstrated to gain weight during storage. Division of the weight loss during the DSC thermal scan by the observed water gain during storage provides some insight into the nature of the water into these materials at the time of the thermal scan. Accordingly, the RMGI materials certainly contain both unbound and bound water in their matrix: but by the time the samples were submitted to the thermal scan a

considerable component of the unbound water has been removed by the specific heat determination. This concept perhaps lends better credence to comparing total weight loss to the water gain during storage.

All three materials demonstrated, with few exceptions, an overall loss in the ratio of thermal total weight loss to water gain with storage time. If the content of RMGI materials are considered, the materials will consist on the outset a mixture of aluminofluorosilicate glasses, polyacrylic acid polymers, methacrylate monomers, water, photointiators, and a variety of other proprietary resin components. After preparation the material will contain, depending on its maturity, a mixture of methacrylate polymers, siliceous hydrogel matrix, and unreacted aluminoflurorsilicate glasses, polyacrylic acid polymers, and unreacted methacrylate components. Concerning the weight loss observed during the thermal scans, it is felt that any unreacted methacrylate components that are not grafted onto the polyacrylic acid polymer backbone may be removed during the thermal scan. However, loss of unreacted polyacrylic acid components contribution to weight loss during the thermal scan is considered minor as most polyacrylic acids used in the RMGI component demonstrate boiling points in excess of 200 °C and thermal endothermic peaks associated with this phase change was not apparent during the thermal scans. It is felt that probably that the major contributor to the weight loss during the DSC thermal scan is due to water loss, which may be accompanied by other minor components that may be concurrently in solution with the water at the time. As storage time progressed, the ratio of total weight loss during the thermal determinations to water gain decreased in all materials. This represents the transformation of the water that is gained from the unbound to the more stable bound nature in the mature hydrogel matrix.

CONCLUSIONS

Based on the conditions of this study, several encompassing summary statements may be made:

- 1. The solubility data reveals interesting material trends concerning resin-modified glass-ionomer restorative materials. All three materials demonstrated increased solubility until the one-month storage period. This reflects a continuation of the ongoing ion-exchange process of the polyalkenoate reaction and that resin-modified glass-ionomer materials may be more susceptible to the effects of solution until that time. Furthermore, solubility measurement, especially under the conditions of this study, does not provide sufficient sensitivity to provide meaningful data in the determination of resin-modified glass-ionomer depth of cure.
- 2. The analysis of the hardness data provided evidence of a continuing post-irradiation polyalkenoate reaction. Although hardness increase was noted in deeper layers of all materials, the increase in the 4 mm and 5 mm bottom surfaces was not of sufficient magnitude to reach statistical equivalence with the hardness exhibited by the 2 mm surfaces. At the end of three months, only Fuji II LC demonstrated a bottom surface hardness of the 3 mm samples that were statistically similar to the hardness exhibited by the 2 mm samples.
- 3. All materials evidenced an increase in specific heat capacity with material aging. This is suspected to be due to increased water uptake and retention, which is strong evidence of the continuation of a post-irradiation polyalkenoate reaction. At the end of the three-month storage there were no significant differences in the specific heat between the 2 mm and 3 mm thicknesses of materials. There remained significant differences between the 2 mm and 4 mm thicknesses in two of the three RMGI materials.
- 4. All materials exhibited evidence during the differential scanning calorimetry thermal scans that demonstrated a change in material compositional nature with storage time. Each material demonstrated different behavior in the thermal scan in regard to energy requirements and tendencies towards stabilization of the thermal energy requirements

- with increased storage time. Each material demonstrated mean thermal scan peak temperature behavior that also provides evidence that the material compositional change within each material is becoming more similar in nature.
- 5. Analysis of weight loss during the specific heat determinations identified evidence of the first events of net water retention with the resin-modified glass-ionomer materials. This suggests that the first evidence of bound water in Fuji II LC is at 24 hours and at one week with both Photac-Fil Quick and Vitremer.
- Analysis of thermal determinations weight loss compared to water gain revealed that
 the resin-modified glass-ionomer materials demonstrated water-retentive behavior in
 the form of unbound and bound water as that reported by conventional glass-ionomer
 materials.

The purpose of this investigation was to characterize the depth of cure of resin-modified glass-ionomer restorative materials using physical property testing with the research hypothesis that RMGI materials will display a depth of cure gradient that is more dependent on visible light penetration rather than what residual acid-base chemical setting reaction that may remain. Based on the analysis of the data that was collected under the conditions of this study the research hypothesis is accepted.

Information gained from the hardness testing and thermal analysis did give evidence of a continuing post-photopolymerization reaction, which is certainly indicative of the acid-base polyalkenoate setting reaction. However, the nature of this reaction, especially in the deeper 4 mm and 5 mm depths, did not provide sufficient change in physical properties as to be statistically equivalent to the 2 mm depths. There were instances in which material properties at the 3 mm depth was equivalent to that seen of 2 mm depths, however this was not consistently seen within all materials during all storage times. It is the conclusion from the conditions of this study that for best material properties resin-modified glass-ionomer restorative products should be photopolymerized in thicknesses of no more than two millimeters.

References

- 1. Wilson AD, Kent BE. The glass-ionomer cement, a new translucent dental filling material. J Appl Chem 1971;21:313.
- 2. Wilson AD. Developments in glass-ionomer cements. Int J Prosthodont 1989;2:438-46.
- 3. Nicholson JW. Chemistry of glass-ionomer cements: a review. Biomaterials 1998;19:485-94.
- 4. Prosser HJ, Powis DR, Brant PJ, Wilson AD. Characterization of glass-ionomer cements. The physical properties of current materials. J Dent 1984;12:231-40.
- 5. McLean JW, Gasser O. Glass cermet cements. Quint Int 1985;16:333-43.
- 6. Simmons JJ. The miracle mixture: glass-ionomer and alloy powder. Texas Dent J 1983;100:6-12.
- 7. Williams JA, Billington RW, Pearson GJ. The comparative strengths of commercial glassionomer cements with and without metal additions. Br Dent J 1992;172:279-82.
- 8. Guggenberger R, May R, Stefan KP. New trends in glass-ionomer chemistry. Biomaterials 1998;19:479-483.
- 9. McLean JW, Nicholson JW, Wilson AD. Proposed nomenclature for glass-ionomer cements and related materials. Quint Int 1994;25:587-9.
- 10. Sidhu SK, Watson TF. Resin-modified glass ionomer materials. A status report for the American Journal of Dentistry. Am J Dent 1995;8:59-67.
- 11. Davidson CL, Mjör IA. Advances in glass-ionomer cements. Chicago: Quintessence Publishing Co., Inc.; 1999.
- 12. Wilson AD. Resin-modified glass-ionomers. Int J Prosthodont 1990;3:425-429.
- 13. deGee, AJ, Leloup G, Werner A, Vreven J, Davidson CL. Structural integrity of resin-modified glass ionomers as affected by the delay or omission of light activation. J Dent Res 1998;77:1658-63.

- 14. Hammesfahr PD. Developments in resionomer systems. In: Glass ionomers: The next generation. Proceedings of the 2nd International Symposium on Glass Ionomers, June 1994, Philadelphia, PA. Hunt P, editor. Philadelphia PA: International Symposia in Dentistry, PC,pp.47-55.
- 15. Eliades G, Palaghias G. In vitro characterization of visible light-cured glass ionomers liners. Dent Mater 1993:9:198-203.
- 16. Culbertson BM. Glass-ionomer dental restoratives. Prog Polym Sci 2001;26:577-604.
- 17. Fuji II LC MSDS, available at http://www.gcamerica.com; accessed 05 Feb 2005.
- 18. Vitremer and Photac-Fil Quick MSDS, available at http://www3.3m.com/search/ww/en001/msdssearchform.do, accessed 05 Feb 2005.
- 19. Anstice HM, Nicholson JW. Studies in the setting of polyelectrolyte materials. 2: The effect of organic compounds on a glass polyalkenoate cement. J Mater Sci Mater Med 1994;5:299-302.
- 20. Nicholson JW, Anstice HM. The physical chemistry of light-cured glass-ionomers. J Mater Sci Mater Med 1994;S:I 19-22.
- 21. McCabe JF. Resin-modified glass ionomers. Biomaterials 1998;19:521-527.
- 22. Geurtsen W. Substances released from dental resin composites and glass ionomer cements. Eur J Oral Sci 1998;106:687-95.
- 23. Xie D, Chong IL, Wu W, Mays J. Synthesis and evaluation of HEMA-free glass-ionomer cements for dental applications. Dental Materials 2004; 20: 470-8.
- 24. Wan ACA, Yap AUJ, Hastings GW. Acid-base complex reactions in resin-modified and conventional glass ionomer cements. J Biomed Mater Res (Appl Biomater) 1999;48:700-4.
- 25. Young AM. FTIR investigation of polymerization and polyacid neutralization kinetics in resinmodified glass-ionomer dental cements. Biomaterials 2002;23:3289-95.

- 26. Young AM, Rafeeka SA, Howlett JA. FTIR investigation of monomer polymerization and polyacid neutralization kinetics and mechanisms in various aesthetic dental restorative materials. Biomaterials 2004;25:823-33.
- 27. Li J, von Beetzen M, Sundström F. Strength and setting behavior of resin-modified glass ionomers cements. Acta Odontol Scand 1995;53:311-7.
- 28. Kakaboura A, Eliades G, Palaghias G. An FTIR study on the setting mechanism of resin-modified glass ionomers restoratives. Dent Mater 1996;12:173-8.
- 29. Ando F, Komori A, Kojima I. Tensile bond strength of light-cured resin-reinforced glass ionomers cement with delayed light exposure. Odontology 2001;89:45-8.
- 30. Watson TF, Pagliari D, Sidhu SK, Naasan MA. Confocal microscopic observation of structural changes in glass-ionomer cements and tooth surfaces. Biomaterials 1998;19:581-8.
- 31. Young AM. FTIR investigation of polymerization and polyacid neutralization kinetics in glass-ionomer cements. Biomaterials 2002;23:3289-95.
- 32. Young AM, Sherpa A, Pearson G, Schottlander B, Waters DN. Use of Raman spectroscopy in the characterization of the acid-base reaction in glass-ionomer cements. Biomaterials 2001;21:1971-9.
- 33. Palmer G, Anstice HM, Pearson GJ. The effect of curing regime on the release of hydroxyethyl methacrylate (HEMA) from resin-modified glass-ionomer cements. J Dent 1999;27:303-11.
- 34. Kanchanavista W, Pearson GJ, Anstice HM. Factors contributing to the temperature rise during polymerization of resin-modified glass-ionomer cements. Biomaterials 1996;17:2305-2312.
- 35. Mitra SB. Adhesion to dentin and physical properties of a light cured glass-ionomer liner/base. J Dent Res 1991;70:72-4.
- 36. Hotta M, Aono M. Adaptation to the cavity floor of the light cured glass ionomer cement base under a composite restoration. J Oral Rehabil 1993:21:679 85.

- 37. Swift EJ, Pawlus MA, Vargas MA. Shear bond strengths of resin-modified glass-ionomer restorative materials. Oper Dent 1995; 20: 138-43.
- 38. Hinoura K, Miyazaki M, Onose H. Dentin bond strength of light cured glass ionomer cements. J Dent Res 1991:70:1542-44.
- 39. Friedl KH, Powers JM, Hiller KA. Influence of different factors on bond strength of hybrid ionomers. Oper Dent 1995:20:74-80.
- 40. Lin A, McIntyre NS, Davidson RD. Studies on the adhesion of glass-ionomer cements to dentin. J Dent Res 1992;71:1836-41.
- 41. Bourke AM, Walls AWG, McCabe JF. Bond strength to dentine of two light-activated glass polyalkenoate (lonomer) cements. J Dent 1994:22:242-6.
- 42. Chadwick RG, Woolford MJ. A comparison of the shear bond strengths to a resin composite of two conventional and two resin-modified glass polyalkenoate (ionomer) cements. J Dent 1993:21:111-6.
- 43. Sidhu SK, Sherriff F, Watson TM. Failure of resin-modified glass-ionomers subjected to shear loading. Journal of Dentistry 27 (1999) 373–381.
- 44. Kerby RE, Knobloch L. The relative shear bond strength of visible light-curing and chemically curing glass-ionomer cement to composite resin. Quintessence International 1992;23:641–644.
- 45. Compton AM, Meyers CE, Hondrum SO, Lorton LE. Comparison of the shear bond strength of a light-cured glass ionomer and a chemically cured glass ionomer for use as an orthodontic bonding agent. American Journal of Orthodontics and Dentofacial Orthopedics 1992;101:138–144.
- 46. Davis EL, Yu X, Joynt RB, Wieczkowski Jr. G, Giordano L. Shear strength and microleakage of light-cured glass-ionomers. American Journal of Dentistry 1993;6:127–129.
- 47. Pereira PNR, Sano H, Ogata M, Zheng L, Nakajima M, Tagami J, Pashley DH. Effect of region and dentin perfusion on bond strengths of resin-modified glass ionomer cements. Journal of Dentistry 28 (2000) 347–354.

- 48. Wilson AD, Prosser HJ, Powis DM. Mechanism of adhesion of polyelectrolyte cements to hydroxyapatite. J Dent Res 1983;62:590-2.
- 49. Yoshida Y, Van Meerbeek B, Nakayama Y, Snauwaert J, Hellemans L, Lambrechts P, Vanherle G, Wakasa K. Evidence of chemical bonding at biomaterial-hard tissue interfaces. J Dent Res 2000;79:709-14.
- 50. Yoshida Y, Van Meerbeek B, Nakayama Y, Yoshioka M, Snauwaert J, Abe Y, Labrechts P, Vanherle G, Okazaki M. Adhesion to and decalcification of hydroxyapatite by carboxylic acids. J Dent Res 2001;80:1565-9.
- 51. Tay FR, Smales RJ, Ngo H, Pashley DH. Effect of different conditioning protocols on adhesion of a GIC to dentin. J Adhes Dent 2001;3:153-67.
- 52. Yoshida Y, Van Meerbeek B, Snauwaert J,Hellemans L,Lambrechts P,Vanherle G, Wakasa K Pashley DH. A novel approach to AFM characterization of adhesive tooth–biomaterial interfaces. J Biomed Mater Res, 47, 85–90, 1999.
- 53. Van Meerbeek B, De Munck J, Yoshida Y, Inoue S, Vargas M, Vijay P, Van Landuyt K, Lambrechts P, Van Herle G. Adhesion to enamel and dentin: current status and future challenges. Oper Dent 2003;28:215-35.
- 54. Sidhu SK, Watson TF. Interfacial characteristics of resin-modified glass ionomer materials: A study on fluid permeability using confocal fluorescence microscopy. J Dent Res 1998;77:1749-59.
- 55. Sidhu SK, Sherriff M, Watson TF. Failure of resin-modified glass-ionomers subject to shear loading. J Dent 1999;27:373-81.
- 56. Tay FR, Sidhu SK, Watson TF, Pashley DH. Water-dependent interfacial transition zone in resin-modified glass-ionomer cement/dentin interfaces. J Dent Res 2004p83:644-49.
- 57. Yiu CKY, King NM, Pashley DH, Carvalho RM, Carriho MRO. Interaction of resin-modified glass-ionomer cements with moist dentin. J Dent 2004;32:521-30.

- 58. Burke FM, Hamlin PD, Lynch EJ. Depth of cure of light-cured glass-ionomer cements. Quintessence Int 1990;21:977-81.
- 59. Swift EJ, Pawlus MA, Vargas MA, Fortin D. Depth of cure of resin-modified glass ionomers. Dent Mater 1995;11:196-200.
- 60. ANSI/ADA Specification No. 27 Resin-Based Filling Materials: 1997.
- 61. Khalil SKH, Atkins EDT. Investigation of glass-ionomer cements using differential scanning calorimetry. J Mater Science Mater in Med 1998;9:529-533.
- 62. 2-Hydroxyethyl methacrylate MSDS, Sigma-Aldrich Corporation, St. Louis, Mo.
- 63. Polyacrylic Acid MSDS, American Polymer Standards Corporation, Mentor OH.
- 64. Andrzejewska A, Andrzejewska M, Socha E, Zych-Tomkowiak D. Effect of polyacid aqueous solutions on photocuring of polymerizable components of resin-modified glass ionomer cements. Dent Mater 2003;19;501-509.
- 65. Anstice HM, Nicholson JW. Studies in the setting of polyelectrolyte materials Part II: The effect of organic compounds on a glass poly (alkenoate) cement. J Mater Sci Mat Med 1994;5:299-302.
- 66. Chuang S, Jin Y, DDS, Tsai P, Wong T. Effect of various surface protections on the margin microleakage of resin-modified glass ionomer cements. J Prosthet Dent 2001;86:309-14.
- 67. Sidhu SK, Pilecki P, Sherriff M, Watson TF. Crack closure on rehydration of glass-ionomer materials. Eur J Oral Sci 2004;112:465-469.
- 68. Sidhul SK, Sherriff M, Watson TF. The Effects of Maturity and Dehydration Shrinkage on Resin-modified Glass-ionomer Restorations J Dent Res 1997; 76:1495-1501.
- 69. Deb S, Braden M, Bonfield W. Water absorption characteristics of modified hydroxyapatite bone cements. Biomater 1995;16:1095-1100.
- 70. Nicholson JW. The physics of water sorption by resin-modified glass-ionomer dental cements. J Mater Sci Mater Med 1997;8:691-695.

- 71. Jevnikar P, Jarh O, Sepe A, Pintar MM, Fundak N. Micro magnetic resonance imaging of water uptake by glass ionomer cements. Dent Mater 1997;13:20-23.
- 72. Jevnikar P, Sersa I, Sepe A, Jarh O, Funduk N. Effect of surface coating on water migration into resin-modified glass ionomer cements: a magnetic resonance micro-imaging study. Magn Reson Med 2000;44:686-691.
- 73. Williams JA, Billington RW. Changes in the compressive strength of glass-ionomer restorative materials with respect to time periods of 24 h to 4 months. J Oral Rehab 1991;18:163-168.
- 74. Pearson GJ, Atkinson AS. Long-term flexural strength of glass-ionomer cements. Biomaterials 1991;12:655-660.
- 75. Nicholson JW, Abiden F. Changes in compressive strength on ageing in glass polyalkenoate ('glass ionomer') cements prepared from acrylic/maleic acid copolymers. Biomaterials 1997;18:59-62.
- 76. Gladys S, Van Meerbeek B, Braem M, Lambrechts P, Vanherle G. Comparative physico-mechanical characterization of new hybrid restorative materials with conventional glass-ionomer and resin composite restorative materials. J Dent Res 1997;76:883-894.
- 77. Burke FM, Hamlin PD, Lynch EJ. Depth of cure of light-cured glass-ionomer cements. Quintessence Int 1990;21:977-981.
- 78. Yap AUJ. Post-irradiation hardness of resin-modified glass ionomer cements and a polyacid-modified composite resin. J Mater Sci Mat Med 1997;8:413-416.
- 79. Ellakuria J, Triana R, Minguez N, Soler I, Ibaseta G, Maza J, Gardia-Godoy F. Effect of one-year water storage on the surface microhardness of resin-modified versus conventional glass-ionomer cements. Dent Mater 2003;19:286-290.
- 80. Li Jianguo, von Beetzen M, Sundstrom F. Strength and setting behavior of resin-modified glass ionomer cements. Acta Odontol Scand 1995;53:311-317.

- 81. Bourke AM, Walls AW, McCabe JF. Light-activated glass polyalkenoate (ionomer) cements: the setting reaction. J Dent 1992;20:115-120.
- 82. Xie D, Brantley WA, Culbertson BM, Wang G. Mechanical properties and microstructures of glass-ionomer cements. Dent Mater 2000;16:129-138.
- 83. Khanchanavasita W, Anstice HM, Pearson GJ. Long-term surface micro-hardness of resin-modified glass ionomers. J Dent 1998;26:707-712.
- 84. DeWald JP, Ferracane JL. A comparison of four modes of evaluating depth of cure of light-activated composites. J Dent Res 1987;66:727-730.
- 85. Rueggeberg FA, Craig RG. Correlation of parameters used to estimate monomer conversion in a light-cured composite. J Dent Res 1988;67:932-937.
- 86. Watts DC, Amer O, Comb EC. Characteristics of visible light-activated composite systems. Br Dent J 1984;156:209-215.
- 87. Yearn JA. Factors affecting cure of visible light activated composites. Int Dent J 1985;35:218-225.
- 88. Pilo R, Cardash HS. Post-irradiation polymerization of different anterior and posterior light-activated resin composites. Dent Mater 1992;8:299-304.
- 89. Brunette DM. Critical thinking understanding and evaluating dental research. Chicago: Quintessence, 1996.
- 90. Holtzchaw HF Jr., Robinson WR. General Chemistry 8th ed. Lexington: D.C. Heath an^d Company, pp. 72-92.
- 91. Mettler Toledo USER COM #7. Mettler-Toledo GMBH, June 1998.
- 92. Hatakeyama T, Zhenhai L. Handbook of Thermal Analysis. Chichester: John Wiley & Sons, Inc. 1998;pp.155-6.
- 93. Wilson AD, Nicholson JW. Acid-base cements their biomedical and industrial applications. New York: Cambridge University Press, 1993: 48-51.

APPENDIX A

SOLUBILITY DATA

Table 64. Vitremer Immediate Solubility Experimental Sample Results (grams)

	•	A CHILL	Shair of Direct	des Costs				Dre H20	noct HOO						
		O I Idus Aildingucai i	alytical i lu	orano o											
am Groun								3-Jan-06 1	4-Jan-06	15-Jan-06	17-Jan-06	13-Jan-06 14-Jan-06 15-Jan-06 17-Jan-06 18-Jan-06 19-Jan-06	19-Jan-06	20-Jan-06 21-Jan-06	21-Jan-06
Sample	-	2	3	4	5	-		0.2155 0.2238		0.2082	0.207	0.2064	0.2062	0.2061	0.2064
Start weight	0.2155	0.2324	0.2169	0.2162	0.2268	0.2268 Floss total 2		0.2324 0.2425		0.2259	0.2245	0.2241	0.2236	0.2235	0.2239
Floss Weight	0 0111	0 0103	0.0106	0.01	0.0104	0.0524 3		0.2169 0.261		0.2099	0.2085	0.2079	0.2078	0.2076	0.2078
W1 (start floss)	0 2044	0 2221	0 2063	0.2062	0.2164			0.2162 0.2249			0.2086				0.208
Final constant weight	0 2064	0 2239	0 2078	0.208	0.2175			0.2268 0.2355			0.2182	0.2176	_		0.2175
Floss Weight	0.0111	0.0103	0.0106	0.01	0.0104	W	Wt-floss	1.0554 1			1.0144				1.0112
W2 (final-floss)	0.1953	0.2136	0.1972	0.198	0.2071	Me	Mean	0.21108	0.22706	0.2042	0.20288	0.20234	0.2021	0.20184	0.20224
			THE RESERVE OF THE PERSON NAMED IN COLUMN TWO IS NOT THE PERSON NAMED IN THE PERSON NAMED IN THE PERSON NAMED IN THE PERSON NAMED IN THE PERSON NA			ps		0.007618	0.015229	0.007726	0.007653	0.00771	0.007554	0.007747	0.007653
SOLUBILITY	4.452055	3.827105	4.411052	3.976722	4.297597			Mean solubility	ility	4.192906	0.2767				
nm Group															
Sample	-	2			2	-		0.288 0.2987				0.2747	0.2746		
Start weight	0.288	0.2723	0.2763		0.2804	F		0.2723 0.2832				0.2604	0.2602		
Floss Weight	0.0101	0.0089	0.0102	0.0088	0.0095	0.0475		0.2763 0.2869		0.267	0.2651	0.2643	0.2641		
W1 (start-floss)	0.2779	0.2634	0.2661	0.2609	0.2709			0.2697 0.2823		0.2587	0.2569	0.2559	0.2556		-
Final constant weight	0.2746	0.2602	0.2641	0.2556	0.2693			0.2804 0.2905		0.2721	0.2703	0.26	0.2693		
Floss Weight	0.0101	0.0089	0.0102	0.0088	0.0095		Wt - floss	1.3392	1.3941	1.2908			1.2763	-0.0475	
W2 (final-floss)	0.2645	0.2513	0.2539	0.2468	0.2598		Mean		0.27882	0.25816		183		986	
						ps		0.007209	0.006651	0.007302		0.007351	0.007459	#DIV/0i	
% SOLUBILITY	4.821878	4.593774	4.584743	5.404369	4.097453			Mean solubility	illity	4.700443	0.473978				
nm Groun															
Sample	-	2	3	4	5	-		0.3254 0.3408		0.3106	0.3074	0.3062	0.3062	0.3057	0.3059
Start weight	0.3254	0.3788	0.3765	0.3623	0.4152	0.4152 Floss total 2		0.3788 0.3912		0.3618	0.358	0.3567	0.3564	0.3563	0.3565
Floss Weight	0.0109	0.01		0.0087	0	0.0394 3		0.3765 0.39		0.3553	0.3518	0.3499	0.3493	0.3493	0.3498
W1 (start-floss)	0.3145	0.3688		0.3536	0.4152	4		0.3623 3752		0.3406	0.3368	0.3353		0.3344	0.3348
Final constant weight	0.3059	0.3565		0.3348	0.3707			0.4152 0.4122		0.3795	0.3	0.3	0.3709	0.3703	0.3707
Floss Weight	0.0109	0.01	0.00	0.0087	0		Wt-floss	1.8188	3753.495	1.7084				1.6766	1.6783
W2 (final-floss)	0.295	0.3465	0.34	0.3261	0.3707		Mean	0.36376	750.699	0.34168					0.33566
						ps		0.032395	1677.774	0.025871	0.02512	0.024811	0.024594	0.024632	0.024688
% SOLUBILITY	6.200318	6.046638	7.281156	7.777149	10.71773			Mean solubility	oillity	7.604598	1.885815				
mm Group															
Sample	-	2	3	4	5	-		0.4481 0.4624		0.4297	0.4265	0.4251	0.4249	0.4247	0.4251
Start weight	0.4481	0.4145	0.4833	0.5052	0.4921	0.4921 Floss total 2		0.4145 0.4305		0.3979	0.395	0.3938	0.3928	0.4433	0.3931
Floss Weight	0.0107	0.0097		0.0089	0.0099	0.0486		0.4833 0.4966		0.464	0.4603	0.4589	0.4579	0.4847	0.4581
W1 (start-floss)	0.4374	0.4048		Ğ				0.5052 0.5181		0.4834	0.4792	0.4778	0.4768	0.4664	0.477
Final constant weight		0.3931						0.4921 0.5017	1	0.471	0.4	0.4	0.46	0.4	0.4643
Floss Weight	0.0107	0.0097	00.0		0.0099		Wt-floss	2.2946	2.3607	2.1974					2.169
W2 (final-floss)	0.4144	0.3834	0.4487	0.4681	0.4544		Mean	0.45892	0.035145	0.43948	0.43584	5 0 034276	0.4336	0.025297	0.034084
						200		0.00002	0.000.4	0.001011		4			

Table 65. Vitremer Immediate Solubility Results Control Sample Results (grams)

12-Dec-05 13-Dec-05 14-Dec-05 15-Dec-05 17-Dec-05 14-Dec-05 15-Dec-05 17-Dec-05 12-Dec-05 12-D	START: 12 Dec 2005		WHairs Anah	vtical Plus	olone				OCT	noet H20					
1			ומתם היים		Calle				DIe UZO	חסמו וידים					
1	dnos www								12-Dec-05	13-Dec-05	14-Dec-05	15-Dec-05	16-Dec-05	17-Dec-05	18-Dec-05
Colores Colo	Sample	-	2	3	4	5		-	0.2333	0.2395	0.2273	0.2259		0.2252	
Columbia Columbia		0 2333	0.2177	0.2151	0.2352	0.2373	Floss total	2	0.2177	0.225	0.2109	0.2098		0.209	
1		0 0125	0 0125	0.0098	0.0113	0.0125	0.0586	3	0.2151	0.2216	0.2081	0.2072		0.2062	
Colorest Colorest	iel	0 2208	0 2052	0.2053	0 2239	0 2248		4	0.2352	0.242	0.2293	0.228		0.2271	
1,112 1,112 1,045 1,04		0 2252	0 200	0 2062	0 2271	0 229		2	0 2373	0 2434	0 2312	0 2299		0 229	
Table Tabl	Thee Weight	0.0125	0.203	0 0098	0 0113	0.0125		Wt.floss	1 08	1 1129	1 0482	1 0422		1 0379	
1	1035 Weight	0 2127	0 1965	0.0000	0.2158	0.2165		Mean	Hill		0 2096		0808	0 20758	
The color of the	vz (iiiiai-ii0ss)	17170	0.130	55.	7.5130	20170		ps	0 01047053	0 01020931	0.0109589		0 0	0 01081018	
1	COLUDII ITA	267	NC N	131	3 62	3 60		8	Mean colith	Ility	3 4				
1	SOLUBILITY	3.01	+7.4	4.54	3.05	0.00			anios impour	6	2:0				
17 2 3 4 5 5 5 6 6 6 6 6 6 6	mem Craston				ľ										
Color Colo	Cample		2	~	4	5		-	0 276	0 2848	0 2673	0 2658		0 2647	0.2647
Control Cont	balliple tot woight	A7C 0	0 269	0 3235	0 3047	0 2863	Flace total	2	0 269	0 2783	0 2594	0.2579			0 2567
tweight 0.2666 0.2682 0.3148 0.2963 0.2753 4 0.3047 0.3144 0.2943 0.2923 ss) 0.2666 0.2687 0.3168 0.2936 0.2738 5 0.2663 0.2767 0.2725 0.2765 ss) 0.2647 0.2459 0.3014 0.2849 0.3014 0.2849 0.2721 0.2766 0.2725 0.2765 ss) 0.2543 0.2459 0.3014 0.2643 0.2721 0.2769 0.2725 0.2765 r 4.25 4.76 4.76 3.75 4.54 4.54 4.444 7.361 1.3625 r 4.25 4.76 4.76 4.76 4.54 4.54 4.444	lose Wainht	0 0104	0.203	0 0097	0 0087	0 011	0 0506	3 1	0 3235	0 3321	0.3139	0 3121			0.3108
0.2647 0.2567 0.3108 0.2366 0.2738 5 0.2863 0.2959 0.2767 0.2755 0.0104 0.0108 0.0097 0.0087 0.011 Wt. floss 1.4089 1.4549 1.361 1.3525 0.0104 0.0108 0.0097 0.0087 0.014 0.014 0.02018 0.20098 0.2702 0.2706 1 4.25 4.06 3.75 4.54 4.54 4.27 0.4049 0.563 0.381 0.381 0.3746 0.3923 0.3457 0.3629 0.3779 0.4059 0.07 0.0109 0.0093 0.0111 0.009 0.0603 0.3447 0.3483 floss total 0.3954 0.3469 0.3574 0.3469 0.3674 0.3674 0.3469 0.3674 0.3674 0.3499 0.3674 0.3489 0.3779 0.4059 0.3674 0.3489 0.3779 0.4059 0.3679 0.3679 0.3679 0.3679 0.3679 0.3679 0.3679 0.3674	If fetart floeel	0 2656	0 2582	0 3138	0 296	0 2753	Water Control of the	4	0 3047	0 3144	0 2943	0 2923			0.2906
0.0104 0.1068 0.097 0.087 0.011 Wr. floss 14549 14549 1361 1.362 0.2543 0.2459 0.3011 0.2849 0.2628 Mean 0.28178 0.2722 0.2722 0.2705 4.25 4.76 4.05 3.75 4.54 5 Mean 0.0219786 0.229082 0.21278 0.2722 0.2725 1 2 3 4 5 4.54 5 4.278 0.217648 1.3649 0.2722 0.2725 0.2726 <t< td=""><td>inal constant weight</td><td>0.2647</td><td>0 2567</td><td>0 3108</td><td>92620</td><td>0 2738</td><td></td><td>5</td><td>0 2863</td><td>0 2959</td><td>0 2767</td><td>0.275</td><td></td><td></td><td>0.2738</td></t<>	inal constant weight	0.2647	0 2567	0 3108	92620	0 2738		5	0 2863	0 2959	0 2767	0.275			0.2738
1	loce Weight	0 0104	0.0108	0 0097	0 0087	0 011		Wt . floss	1 4089					3459	1.346
1	D /final floes	0 2543	0.2459	0 3011	0.2849	0 2628		Mean	0.28178		•			0 26918	0 2692
The color of the	formula 7							ps	0.02219786	0.02208201	0.0219278		0.02164943	0.0	0.02164548
1	SOLUBILITY	4.25	4.76	4.05	3.75	4.54			Mean solub	illity	4.2				
1		ı													
Court Cour	nm Group		2		4	5			0.363	0.3746	0.3493	0.3466		0.3451	0.3451
12 12 13 14 15 15 15 15 15 15 15	tart weight	0 363	0.381	0.3954	0.3847	0.3483	Floss total	2	0.381	0.3923	0.3674	0.3649		0.3631	0.3632
tweight 0.3451 0.3701 0.3861 0.3735 4 0.3847 0.396 0.3698 0.3666 ss) tweight 0.3451 0.3632 0.3763 0.3365 0.3305 5 0.3483 0.3592 0.3353 0.3323 0.3323 ss) 0.01 0.0109 0.0093 0.0111 0.009 Wt-floss 1.8221 1.8777 1.7524 1.738 r/ 0.357 0.3542 0.3216 Mean 0.01872023 0.01859503 0.0187022 0.01807933 r/ 5.07 4.81 4.95 5.19 5.22 Mean solubility 5.05 0.147 r/ 5.07 4.81 5.22 Mean solubility 5.05 0.147 r/ 5.07 4.83 0.5015 6.22 0.4664 0.481 0.4456 0.4456 r/ 0.0125 0.0101 0.0102 0.011 0.0517 0.4664 0.481 0.4456 0.4456 r/ 0.4734	oss Weight	0.01	0.0109	0.0093	0.0111	0.00	0.0503	3	0.3954	0.4059	0.3809	0.3779	0.3773	0.3763	0.3763
0.3451 0.3632 0.3763 0.3365 0.3306 5 0.3483 0.3592 0.3353 0.3233 0.01 0.0109 0.0093 0.0111 0.009 Wt-floss 1 8221 1 8777 1.7524 1.738 0.3351 0.3523 0.3216 Mean 0.36442 0.37554 0.35048 0.3476 1 2 3.3 4.8 5.2 Mean 0.01872023 0.01859503 0.0180702 0.01807935 1 2 3.4 5.2 Mean solubility 5.05 0.1476 0.1478 0.0125 0.0109 0.0101 0.0102 0.011 0.0647 0.4664 0.481 0.4489 0.4456 0.0125 0.0109 0.0101 0.0102 0.011 0.0647 0.0507 0.5182 0.4856 0.4456 0.0126 0.0109 0.0101 0.0102 0.011 0.0647 0.507 0.5182 0.4891 0.4456 0.0126 0.0109 0.0101 0	// (start-floss)	0.353	0.3701	0.3861	0.3736	0.3393		4	0.3847	0.396	0.3698	0.3666	0.3663	0.3652	0.3653
0.01 0.0109 0.0093 0.0111 0.009 Wt-floss 1 8221 1 8777 1 7524 1 738 5.07 4.81 4.95 5.19 5.21 Mean 0.01872023 0.01859503 0.01810202 0.01807935 5.07 4.81 4.95 5.19 5.22 Amean 0.01872023 0.01859503 0.01810202 0.01807935 1 2 3 4 5 1 0.4839 0.4864 0.4819 0.4489 0.4456 0.0125 0.0109 0.0101 0.0102 0.011 0.0547 3 0.507 0.4839 0.4456 0.4456 0.0125 0.0109 0.0101 0.0102 0.011 0.0547 3 0.507 0.5182 0.4456 0.4456 0.0126 0.0109 0.0101 0.0102 0.011 0.0547 3 0.507 0.5182 0.4856 0.4856 0.0126 0.0109 0.0101 0.0102 0.011 0.0547 0.0501	nal constant weight	0.3451	0.3632	0.3763	0.3653	0.3306		5	0.3483	0.3592	0.3353	0.3323	0.3316	0.3306	0.3306
Columbia Columbia	loss Weight	0.01	0.0109	0.0093	0.0111	0.00		Wt-floss	1.8221		1.752		1,7354	1.73	1,7302
5.07 4.81 4.95 5.19 5.22 addition 5.05 0.01872023 0.01859503 0.01810202 0.01807935 1 2 3 4 5 4 5 1 0.4839 0.4961 0.4658 0.4455 0.4455 0.0125 0.0109 0.0101 0.0102 0.011 0.0547 3 0.507 0.4839 0.4456 0.4456 0.4564 0.4555 0.4963 0.5015 Floss total 2 0.4664 0.4489 0.4456 0.4456 0.4714 0.4555 0.4963 0.5017 0.5077 0.5182 0.4456 0.4456 0.4564 0.4014 0.4893 0.5077 0.5182 0.4456 0.4456 0.4564 0.4964 0.4964 0.4964 0.4964 0.4456 0.4456 0.4564 0.4964 0.4964 0.4964 0.4489 0.4456 0.4456 0.4664 0.4483 0.5074 0.4893 0.5142 0.4486	72 (final-floss)	0.3351	0.3523	0.367	0.3542	0.3216		Mean	0.36442		0.3504		3 0.34708	0.346	0.34604
5.07 4.81 4.95 5.19 5.22 Mean solubility 5.05 0.17 1 2 3 4 5 1 0.4839 0.4961 0.4658 0.4625 0.4155 0.0125 0.0109 0.0101 0.0102 0.011 0.0547 3 0.507 0.4839 0.4861 0.4856 0.4456 0.4714 0.4555 0.4964 0.481 0.481 0.4489 0.4456 0.4504 0.4714 0.4855 0.4964 0.481 0.481 0.4871 0.4564 0.4964 0.481 0.4864 0.481 0.4489 0.4456 0.4714 0.4865 0.4967 0.6477 0.577 0.5182 0.4811 0.4881 0.4564 0.4964 0.4964 0.4964 0.4967 0.4881 0.4881 0.4881 0.4564 0.4964 0.4883 0.5014 0.41912 0.4881 0.4881 0.012 0.010 0.010 0.010 0.010 <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td>ps</td> <td>0.01872023</td> <td>0.01859503</td> <td>0.0181020</td> <td></td> <td>0.01809677</td> <td>0.0181045</td> <td>0.01812691</td>								ps	0.01872023	0.01859503	0.0181020		0.01809677	0.0181045	0.01812691
1 2 3 4 5 1 0.4839 0.4961 0.4658 0.4625 0.4839 0.4864 0.507 0.4893 0.5015 Floss total 0.4664 0.481 0.4899 0.4456 0.0125 0.0109 0.0101 0.0102 0.011 0.0547 3 0.507 0.5182 0.4903 0.4871 0.4774 0.4555 0.4969 0.4779 0.4905 4 0.4893 0.5014 0.4719 0.4685 0.4564 0.4736 0.4796 5 0.507 0.5182 0.4903 0.4871 0.4564 0.4736 0.4796 5 0.507 0.5182 0.4719 0.4685 0.6125 0.0109 0.0101 0.0102 0.011 0.4796 5 0.5015 0.5186 0.4851 0.4813 0.0125 0.0109 0.0101 0.011 0.011 0.011 0.011 0.011 0.011 0.011 0.011 0.011 0.011 0.011 0.011 </td <td>SOLUBILITY</td> <td>2.07</td> <td>4.81</td> <td>4.95</td> <td>5.19</td> <td>5.22</td> <td></td> <td></td> <td>Mean solut</td> <td>illity</td> <td>5.0</td> <td></td> <td></td> <td></td> <td></td>	SOLUBILITY	2.07	4.81	4.95	5.19	5.22			Mean solut	illity	5.0				
1 2 3 4 5 1 0.4839 0.4961 0.4568 0.4625 0.4839 0.4664 0.507 0.4893 0.5015 Floss total 2 0.4664 0.481 0.4489 0.4456 0.0125 0.0109 0.0101 0.0102 0.011 0.0647 3 0.507 0.5182 0.4089 0.4871 0.4714 0.4555 0.4969 0.4791 0.4905 4 0.4893 0.5014 0.4719 0.4685 0.4794 0.4906 0.4796 5 0.5015 0.5136 0.4871 0.4885 0.0125 0.0109 0.0101 0.0102 0.011 Whtfloss 2.3934 2.4556 2.3073 2.2903 0.4479 0.4562 0.4686 Mean 0.4786 0.46112 0.46142 0.46142 0.46142 0.46142	mm Group														
0.4839 0.4664 0.507 0.4839 0.5015 Floss total 2 0.4664 0.507 0.4869 0.4459 0.4459 0.4459 0.4871 0.0125 0.0109 0.0101 0.0102 0.011 0.057 0.507 0.5182 0.4903 0.4871 0.4604 0.4433 0.4905 4 0.4893 0.5014 0.4719 0.4685 0.0125 0.0109 0.0101 0.0102 0.011 When 0.4786 0.4586 0.45806 0.45806 0.4479 0.4324 0.4562 0.4686 Mean 0.4786 0.45806 0.45806	Sample	-	2	3	4	5		- 0	0.4839	0.4961	0.4658	0.4625	0.461/	0.4504	0.4604
0.0125 0.0109 0.0101 0.0102 0.011 0.0547 3 0.507 0.5182 0.4903 0.4871 0.4714 0.4555 0.4969 0.4791 0.4905 4 0.4893 0.514 0.4719 0.4885 0.4604 0.4433 0.4847 0.4664 0.4796 5 0.5015 0.5136 0.4861 0.4813 0.0125 0.0109 0.0101 0.0102 0.011 Wirfloss 2.3934 2.4556 2.3073 2.2903 0.4479 0.4324 0.4562 0.4686 Mean 0.47868 0.46112 0.46102 0.46102	tart weight	0.4839	0.4664	0.507	0.4893	0.5015	rioss total	7	0.4664	0.461	0.4469	0.4450	0.4440	0.4434	0.4433
0.4714 0.4555 0.4969 0.4791 0.4995 4 0.4883 0.5014 0.4719 0.4065 0.4604 0.4433 0.4847 0.4664 0.4796 5 0.5015 0.5136 0.4851 0.4813 0.0125 0.0109 0.0101 0.0102 0.011 Wtrfloss 2.3934 2.4556 2.3073 2.2903 0.4479 0.4324 0.4766 0.4686 Mean 0.47868 0.46112 0.46102 0.46102	loss Weight	0.0125	0.0109	0.0101	0.0102	0.011		e .	0.507	0.5182	0.4903	0.4871	0.4861	0.4848	0.4847
0.4604 0.4433 0.4847 0.4664 0.4796 5 0.5015 0.5136 0.4851 0.4813 0.0125 0.0109 0.0101 0.0102 0.011	V1 (start-floss)	0.4714	0.4555	0.4969	0.4791	0.4905		4	0.4893	0.5014	0.4/19	0.4685	0.46/6	0.4664	0.4664
85) 0.0125 0.0109 0.0101 0.0102 0.011 Wr.floss 2.3934 2.4556 2.3073 2.2903 0.4179 0.4324 0.4746 0.4562 0.4686 Wean 0.47868 0.49112 0.46146 0.45806	inal constant weight	0.4604	0.4433	0.4847	0.4664	0.4796		2	0.5015	0.513	0.485	0.481	0.48	0.479	0.4796
ss) 0.4479 0.4324 0.4746 0.4562 0.4686 Mean 0.47868 0.49112 0.46146 0.45806	loss Weight	0.0125	0.0109	0.0101	0.0102	0.011		Wt-floss	2.3934						2.2797
	V2 (final-floss)	0.4479	0.4324	0.4746	0.4562	0.4686		Mean	0.47868	0.4477467	0.4614	6 0.45800	0.45724	0.45604	0.45594
0.01595101 0.01477457 0.01641465 0.01654746					00.7		THE STREET	DS.	0.0159310	0.01477457	0.0164146	3 U.U163474			

Table 66. Vitremer 24-Hour Solubility Results Experimental Sample Results (grams)

		O'Haus Ana	O'Haus Analytical Plus Scale	Scale			pre H20	post H20					
um Groun							25-Jan-06	26-Jan-06	27-Jan-06	28-Jan-06	29-Jan-06	30-Jan-06	31-Jan-06
Sample	-	2	3	4	2	-	4	0.2116	0.1985	0.1979			
Start weight	0.2094	0.2271	0.2214	0.2245	0.2416	0.2416 Floss total 2	0.2271 0.2296	3.2296	0.2171	0.2168	0.2167		
Floss Weight	9600.0	0.0092	0.0107	0.0094	0.0121	0.0513	0.2214 0.2239	0.2239	0.2104	0.2099	0.2098		
W1 (start-floss)	0.1998	0.2179	0.2107	0.2151	0.2295	4	0.2245 0.2267	7,2267	0.2131	0.2125	0.2125		
Final constant weight	0.1978	0.2167	0.2098	0.2125	0.2273	5	0.2416 0.2434	0.2434	0.2279	0.2273	0.2273		
Floss Weight	9600.0	0.0092	0.0107	0.0094	0.0121	Wt-floss	1.073	1.0842	1.016		1.0131		
W2 (final-floss)	0.1882	0.2075	0.1991	0.2031	0.2152	Mean	0.2146	0.21684	0.2032	0.20268	0.20262		
						ps	0.0115795	0.0114299	0.0106658	0.0106856	0.0107185		
SOLUBILITY	5.81	4.77	5.51	5.58	6.23		Mean solubility	llity	5.58	0.53			
nm Group													
Sample	-	2	3	4	2	-	0.3127 0.3185	0.3185		0.2963	0.2961		
Start weight	0.3127	0.316	0.2836	0.2691	0.3099	0.3099 Floss total 2	0.316 0.322	0.322		0.3007	0.3004		
Floss Weight	0.011	0.0102	0	0.0099	0.0098	0.0409 3	0.2836 0.2909	0.2909	0.2681	0.2672	0.267		
W1 (start-floss)	0.3017	0.3058	0.2836	0.2592	0.3001	4	0.2691 0.2759	0.2759	0.2551	0.2544	0.2542		
Final constant weight	0.2961	0.3004	0.267	0.2542	0.2947	9	0.3099 0.3156	0.3156	0.2957	0.2948	0.2947		
Floss Weight	0.011	0.0102	0	0.0099	0.0098	Wt - floss		1.482		1.3725			
W2 (final-floss)	0.2851	0.2902	0.267	0.2443	0.2849	Mean	0.29008	0.2964		0.2745			
						ps	0.0207601	0.0201764	0.0207181	0.0205943	0.0205873		
SOLUBILITY	5.50	5.10	5.85	5.75	90.9		Mean solubility	ility	5.45	0.36			-
om Groun													
Sample		2	3	4	5		0 4022 0 4105	0.4105	0.3818	0.3804	0.3796	0.3791	0.3791 0.3795
Start woight	0.4022	0 3973	0 39	0 3961	0 3823	0 3823 Floss total 2	0 3973 0 4035	0 4035	0.3782	0 3771	0.3765	0.3758	0.3758 0.3762
Flore Weight	0 0107	0 012	0 0101	6600 0	0 0109	0.0536 3	0.39	0.39 0.3975	0.3692	0.3677	0.3672	0.3662	0.3662 0.3668
W1 (start.floss)	0 3915	0 3853	0.3799	0.3862	0.3714	4	0.3961 0.4038	0.4038	0.3747	0.3734	0.3728	0.3722	0.3722 0.3724
Final constant weight	0 3795	0.3762		0 3724	0.3622	. 2	0.3823 0.3901	0.3901	0.3639	0.3629	0.3624	0.3617	0.3617 0.3622
Floss Weight	0 0107	0.012		0.0099	0.0109	Wt-floss	1.9143	1.9518		1.8079			1.8035
A7 (final-floss)	0 3688	0.3642		0.3625	0.3513	Mean	0.38286	0.39036	0.36284	0.36158	0.36098	0.36028	0.3607
							0.0076575	0.0076709	0.0071248	0.0070601	0.0069534	0.0070643	0.0069901
SOLUBILITY	5.80	5.48	6.11	6.14	5.41		Mean solubility	ility	5.79	0.34			
Sample		2	3	P	5	-	0 5015 0 5083	0.5083	0.4813	0.4792	0.4784	0.4774	0.4774 0.4781
Start weight	0 5015		0	0.5022	0.5121	0.5121 Floss total 2	0.5052 0.5098	0.5098	0.4854	0.4839	0.4831	0.4825	0.4825 0.4829
Floss Weight	0 011			0.0112	0.0101	0.0546 3	0.4808 0.4871	0.4871	0.4611	0.4601	0.459	0.4586	0.4586 0.4592
W1 (start-floss)	0.4905	0.4941		0.491	0.502		0.5022 0.5084	0.5084	0.4801	0.4783	0.4776	0.477	0.477 0.4773
Final constant weight	0.4781	0.4829	0.4592	0.4773	0.4861	5	0.5121 0.5182	0.5182	0.4892	0.48	0.48	0.4859	0.486
Floss Weight	0.011	0.0111	0.0112	0.0112	0.0101	Wt-floss	2.4472	2.4772	2.3425	2.3345		2.3268	2.329
W2 (final-floss)	0.4671	0.4718	0.448	0.4661	0.476	Mean	0.48944	0.49544					
						ps	0.0117104	0.0115205	0.0108516	0.0105918	0.0106732	0.0105545	0.0104327
The second secon													

Table 67. Vitremer 24-Hour Solubility Results Control Sample Results (grams)

		Julians Alla	Undus Analytical Files Scale	Scale			DIE 1120	post 1120					
mm Group							21-Jan-06	21-Jan-06 22-Jan-06	23-Jan-06	24-Jan-06	25-Jan-06 26-Jan-06	26-Jan-06	27-Jan-06
Sample	-	2	3	4	2	-	0.2085	0.2085 0.2107	0.1988	0.1971	0.1971		
Start weight	0.2085	0.2213	0.2317	0.2236	0.2436	0.2436 Floss total 2	0.2213	0.2213 0.2237	0.2118	0.2094	0.2091		
Floss Weight	0.0107	0.012	0.0115	0.012	0.0117	0.0579 3	0.2317	0.2317 0.2348	0.2224	0.2204	0.2201		
W1 (start-floss)	0.1978	0.2093	0.2202	0.2116	0.2319	4	0.2236	0.2236 0.2269	0.2129	0.2107	0.2109		
Final constant weight	0.1971	0.2091	0.2201	0.2109	0.2306	5	0.2436	0.2436 0.2463	0.2334	0.2307	0.2306		
Floss Weight	0.0107	0.012	0.0115	0.012	0.0117	Wt-floss	1.0708	1.0845		1.0104	1.0099		
W2 (final-floss)	0.1864	0.1971	0.2086	0.1989	0.2189	Mean	0.21416	0.2169	0.20428	0.20208	0.20198		
						ps	0.013	0.013219	0.01291	0.012619	0.01256		
SOLUBILITY	5.76	5.83	5.27	00.9	5.61	Charles of the same	Mean solubility	pillity	5.69	0.28			
Create Contract													
Sample		2	2	4	5	-	0 297	0 297 0 3044	0 2844	0 2819	0 2806		
Start weight	0 297	0 2968	0.2656	0.2686	0.2855	0 2855 Floss total 2	0.2968	0.2968 0.3032	0.2854		0.282		
Floss Weight	0.0124	0.0108	0.0115	0.0104	0.0119	0.057 3	0.2656	0.2656 0.2713	0.2553	0.2528	0.2523		
W1 (start-floss)	0.2846	0.286	0.2541	0.2582	0.2736		0.2686 0.276	0.276	0.2553	0.2522	0.252		
Final constant weight	0.2806	0.282	0.2523	0.252	0.2701	5	0.2855	0.2855 0.2919	0.2734	0.2705	0.2701		
Floss Weight	0.0124	0.0108	0.0115	0.0104	0.0119	Wt - floss	s 1.3565	1.3898	1.2968	1.2825	1.28		
W2 (final-floss)	0.2682	0.2712	0.2408	0.2416	0.2582	Mean	0.2713	0.27796	0.25936	0.2565	0.256		
						ps	0.01502	0.015239	0.014878	0.014823	0.014662		
% SOLUBILITY	5.76	5.17	5.23	6.43	5.63		Mean solubility	pillity	5.65	0.50			
mm Group													
Sample	-	2	3	4	2	-	0.3631	0.3631 0.3798	0.3483	0.3447	0.3437	0.3433	0.3433 0.3432
Start weight	0.3631	0.4162	0.412	0.4313	0.4301	0.4301 Floss total 2	0.4162	0.4162 0.4209	0.4007	0.3965	0.3953	0.3948	0.3948 0.3946
Floss Weight	0.012	0.0112	0	0.0126	0.0118	0.0476 3	0.412	0.412 0.4177	0.3984	0.3943	0.3933	0.3928	0.3928 0.3926
W1 (start-floss)	0.3511	0.405	0.412	0.4187	0.4183	4	0.4313	0.4313 0.4362	0.4159	0.4123	0.4113	0.4107	0.4107 0.4105
Final constant weight	0.3432	0.3946	0.3926	0.4105	0.4111		0.4301	0.4	0.4	0.4	0.4	0.4112	0.4
Floss Weight	0.012	0.0112	0	0.0126	0.0118	Wt-floss	2.0051		1.9324			1.9052	
W2 (final-floss)	0.3312	0.3834	0.3926	0.3979	0.3993	Mean	0.40102					0.38104	
						ps	0.027831	0.022885	0.05	0.027891	0.027865	0.027785	0.027762
% SOLUBILITY	2.67	5.33	4.71	4.97	4.54		Mean solubility	pillity	5.04	0.46			
mm Groun													
Sample	-	2	3	4	5	-	0.5098	0.5098 0.5156	0.4948	0.489	0.4879	0.4873	0.4873 0.4871
Start weight	0.5098	0.5004	0.4859	0.4903	0.5003	0.5003 Floss total 2	0.5004	0.5004 0.5056	0.4829	0.478	0.4773	0.4776 0.477	0.477
Floss Weight	0.0101	0.0105	0.0098	0.0107	0.011	0.05213	0.4859	0.4859 0.4918	0.4705	0.4651	0.4643	0.4633 0.463	0.463
W1 (start-floss)	0.4997	0.4899	0.4761	0.4796	0.4893	4	0.4903	0.4903 0.4966	0.4743	0.4691	0.4685	0.4677	0.4677 0.4673
Final constant weight	0.4871	0.477	0.463	0.4673	0.4742	9	0.5003	0.5003 0.5064	0.4824	0.4764	0.4755	0.4744	0.4744 0.4742
Floss Weight	0.0101	0.0105	0.0098	0.0107	0.011	Wt-floss	2.4346		2.3528	2.3255	2.3214	2.3182	2.3165
W2 (final-floss)	0.477	0.4665	0.4532	0.4566	0.4632	Mean	0.48692		1000				
						ps	0.009405	0.009264	0.00	0.00	0.009059	0.009276	0.009303
W. COLINEILITY	A 5.4	4 70	4 24	A 80	E 33		Moan collinity	thillips.	484	DC U			

Table 68. Vitremer One-Week Solubility Results Experimental Sample Results (grams)

1 0.2085 0.0107 s) 0.1978 t weight 0.1971 ss) 0.1864 77 5.76 10.0124 is) 0.2846 t weight 0.2806	Haus Anal	U'Haus Analytical Pius Scale	ocale				DZH ad	DOST HZU					
	2												
	2						21-Jan-06 22-Jan-06	22-Jan-06	23-Jan-06	24-Jan-06	25-Jan-06 26-Jan-06		27-Jan-06
		3	4	2	-		0.2085	0.2085 0.2107	0.1988	0.1971	0.1971		
	0.2213	0.2317	0.2236	0.2436	0.2436 Floss total 2		0.2213 0.2237	0.2237	0.2118	0.2094	0.2091		
	0.012	0.0115	0.012	0.0117	0.0579 3		0.2317	0.2317 0.2348	0.2224	0.2204	0.2201		
	0.2093	0.2202	0.2116	0.2319	4		0.2236	0.2236 0.2269	0.2129	0.2107	0.2109		
	0.2091	0.2201	0.2109	0.2306	5		0.2436	0.2436 0.2463	0.2334	0.2307	0.2306		
	0.012	0.0115	0.012	0.0117	S	Wt-floss	1.0708	1.0845	1.0214	1.0104	1.0099		
	0.1971	0.2086	0.1989	0.2189	2	Mean	0.21416	0.2169	0.20428	0.20208			
					ps	P	0.013	0.013219	0.01291	0.012619	0.01256		
	5.83	5.27	00.9	5.61	September 1		Mean solubility	billity	5.69	0.28			
	2	3	4	2	•		0.297	0.297 0.3044	0.2844	0.2819	0.2806		
	0.2968	0.2656	0.2686	0.2855	0.2855 Floss total 2		0.2968	0.2968 0.3032	0.2854	0.2821	0.282		
	0.0108	0.0115	0.0104	0.0119	0.057 3		0.2656	0.2656 0.2713	0.2553	0.2528	0.2523		
	0.286	0.2541	0.2582	0.2736	4		0.2686 0.276	0.276	0.2553	0.2522	0.252		
	0.282	0.2523	0.252	0.2701	5		0.2855	0.2855 0.2919	0.2734	0.2705	0.2701		
Floss Weight 0.0124	0.0108	0.0115	0.0104	0.0119	>	Wt - floss	1.3565	1.3898		1.2825	1.28		
W2 (final-floss) 0.2682	0.2712	0.2408	0.2416	0.2582	2	Mean	0.2713	0.27796	0.25936	0.2565	0.256		
					ps	P	0.01502	0.015239	0.014878	0.014823	0.014662		
% SOLUBILITY 5.76	5.17	5.23	6.43	5.63	A STATE OF THE PARTY OF THE PAR		Mean solubility	pillity	5.65	0.50			
Amm Grain													
Sample 1	2	3	4	5	-		0.3631	0.3631 0.3798	0.3483	0.3447	0.3437	0.3433 0.3432	0.3432
Start weight 0.3631	0.4162	0.412	0.4313	0.4301	0.4301 Floss total 2		0.4162	0.4162 0.4209	0.4007	0.3965	0.3953	0.3948 0.3946	0.3946
#	0.0112	0	0.0126	0.0118	0.0476 3		0.412	0.412 0.4177	0.3984	0.3943	0.3933	0.3928 0.3926	0.3926
W1 (start-floss) 0.3511	0.405	0.412	0.4187	0.4183	4		0.4313	0.4313 0.4362	0.4159	0.4123	0.4113	0.4107 0.4105	0.4105
Final constant weight 0.3432	0.3946	0.3926	0.4105	0.4111	5		0.4301	0.4301 0.4352	0.4167	0.4129	0.4118	0.4112 0.4111	0.4111
Floss Weight 0.012	0.0112	0	0.0126	0.0118		Wt-floss	2.0051					1.9052	1.9044
W2 (final-floss) 0.3312	0.3834	0.3926	0.3979	0.3993	=	Mean	0.40102					0.38104	0.38088
					_{(h}	ps	0.027831	0.022885	0.05	0.02	0.027865	0.027785	0.027762
% SOLUBILITY 5.67	5.33	4.71	4.97	4.54	San	Sales a sign	Mean solubility	pillity	5.04	0.46			
Smm Groun													
Sample 1	2	3	4	5	-		0.5098	0.5098 0.5156	0.4948	0.489	0.4879	0.4873 0.4871	0.4871
Start weight 0.5098	0.5004	0.4859	0.4903	0.5003	0.5003 Floss total 2		0.5004	0.5004 0.5056	0.4829	0.478	0.4773	0.4776 0.477	0.477
1	0.0105	0.0098	0.0107	0.011	0.0521 3		0.4859	0.4859 0.4918	0.4705	0.4651	0.4643	0.4633 0.463	0.463
W1 (start-floss) 0.4997	0.4899	0.4761	0.4796	0.4893	4		0.4903	0.4903 0.4966	0.4743	0.4691	0.4685	0.4677 0.4673	0.4673
Final constant weight 0.4871	0.477	0.463	0.4673	0.4742	9		0.5003	0.5003 0.5064	0.4824	0.4764	0.4755	0.4744 0.4742	0.4742
Floss Weight 0.0101	0.0105	0.0098	0.0107	0.011		Wt-floss	2.4346	2.4639		2.3255		2.3182	2.3165
W2 (final-floss) 0.477	0.4665	0.4532	0.4566	0.4632		Mean	0.48692	0.49278	0.47056			0.46364	0.4633
					co.	ps	0.009405	0.009264	0.009369	0.00	0.009059	0.009276	0.009303
% SOLUBILITY 4.54	4.78	4.81	4.80	5.33			Mean solubility	billity	4.85	0.29			

Table 69. Vitremer Restorative One-Week Solubility Control Sample Results (grams)

11 2 3 4 5 5 5 5 5 5 5 5 5	VITREMER ONE WEEK SOLUBILITY (CONTROL OR	SOLUBIL	III (CONI	NOT OUR	(nnr)					000						
Colores Colo	Weight in grams							olamo,	pre H20	post HZU						
1	mm Group							Sample	5-Jan-06	6-Jan-06	7-Jan-06	8-Jan-06	9-Jan-06	10-Jan-06	11-Jan-06	11-Jan-06 12-Jan-06
Color Colo	Sample	-	2	3	4	5		-			0.182	0.1811	0.1802	0.1799	0.1794	0.1794
0.01726 0.01039 0.01019 0.01119 0.01119 0.01119 0.01119 0.01119 0.01119 0.01170 0.01120 0.0120 0.0120 0.0120 0.0120 0.	tart weight	0.1913	0.2282	0.212	0.2005	0.2337	Floss total	2			0.2173	0.2162	0.2161	0.2157	0.2152	0.2152
0.1756 0.2163 0.2063 0.1866 0.2223 0.1844 0.1849<	loss Weight	0.0128	0.0099	0.0118	0.0109	0.0117	0.0571	m			0.2023	0.2013	0.2011	0.2007	0.2003	0.2005
1,1,2 2,1,2 2,0,00.0. 0,1886 0,121 0,1386 0,1236 0,1239 0,121 0,121 0,139 0,131 0,139 0,131 0,139 0,131 0,139 0,131 0,139 0,131 0,139 0,131 0,139 0,131 0,139 0,131 0,139 0,131 0,139 0,131 0,139 0,131 0,139 0,131 0,139 0,131 0,139 0,131 0,139 0,131 0,139 0,13	/1 (start-floss)	0.1785	0.2183	0.2002	0.1896	0.222		4						0.1888	0.1886	0.1886
Column C	inal constant weight	0.1794	0.2152	0.2005	0.1886	0.2211		c,						0.2213	0.2209	0.2211
Fig. 16. Fig. 16.	loss Weight	0.0128	0.0099	0.0118	0.0109	0.0117		Wt-floss	1.0086	1.022	0.9578	0.9526	0.9514	0.949	0.947	0.947
1	(final-floss)	0.1666	0.2053	0.1887	0.1777	0.2094		Mean	0.20		0.19					
Column C								ps	0.02	0.02	0.02			0.02	0.02	0.02
1	SOLUBILITY	29.9	5.96	5.74	6.28	5.68			Mean solu	ps/killid	90'9	0.41				
1	11111		*													
0.2916 0.309 0.3097 0.3202 0.2916 Floss total 2 0.399 0.3147 0.2946 0.293 0.2939 0.3147 0.2946 0.2939 0.2784 0.2929 0.2993 0.2789 0.0789 0.0789 <td>mm Group</td> <td></td> <td>2</td> <td></td> <td>4</td> <td>2</td> <td></td> <td>-</td> <td></td> <td>0.295</td> <td>0 2798</td> <td>0 2786</td> <td>0 2784</td> <td>0 2775</td> <td>0 2772</td> <td>0.277</td>	mm Group		2		4	2		-		0.295	0 2798	0 2786	0 2784	0 2775	0 2772	0.277
0.2415 0.2983 0.2787 0.3147 0.2973 0.2965 2959 0.2415 0.2983 0.2887 0.2887 0.2789 0.2922 0.2923 0.2771 0.2973 0.2764 0.2771 0.2415 0.2983 0.2887 0.2887 0.2789 0.0275 0.441 1.4628 0.2771 0.2771 0.2773 0.2773 0.2774 0.2771	ample tot woight	n 2916	0 309	1 3097	CCPC D	n 2916	Floss total	2		0.3147	0 2946	0.293	0 2928	0 2919	0.2914	0.2913
0.2815 0.2933 0.29467 0.28678 0.27898 0.27789 4 0.2925 0.29633 0.27894 0.27898 0.27898 0.27898 0.27898 0.27899 0.47899	nee Weight	0.2310	0.000	0.000	0.0095	0.0127	0.063	ı m			0.2973	0.296	2959	0.295	0.2944	0.2945
Colored Colo	// (start-floss)	0.2815	0.2993	0.2987	0.2827	0.2789		4			0.2793	0.2781	0.2778	0.2771	0.2767	0.2768
Color Colo	nal constant weight	0.277	0.2913	0.2945	0.2768	0.2749		5			0.2771	0.2763	0.2754	0.2754	0.2749	0.2749
Mean 0.269 0.2816 0.2815 0.2673 0.2622 Mean 0.01 0.01 0.01 1323 18 Mean 0.01	oss Weight	0.0101	0.0097	0.011	0.0095	0.0127		Wt-floss			1.3751			1.3639	1.3616	1.3615
1,000,000,000,000,000,000,000,000,000,0	2 (final-floss)	0.2669	0.2816	0.2835	0.2673	0.2622		Mean	0.29		0.28					
6.19 5.91 5.02 5.45 5.99 Mean solubility/sd 5.52 0.41 1 2 3 4 5 1 0.3826 0.3877 0.3641 0.3622 0.3617 0.3826 0.4341 0.4018 0.3825 0.4038 0.0148 0.0497 3 0.4018 0.4174 0.4144 0.4144 0.4144 0.4144 0.4144 0.4144 0.4144 0.4144 0.4144 0.4144 0.4144 0.4144 0.4144 0.4144 0.4144 0.4148 0.4147 0.4144								sq	0.01	0.01	0.01			0.01	0.01	0.01
1 2 3 4 5 1 0.3826 0.3877 0.3641 0.3622 0.3617 0.3826 0.4341 0.4018 0.0108	SOLUBILITY	5.19	5.91	5.09	5.45	5.99			Mean solu	ps/killidi	5.52					
1 2 3 4 5 1 0.3826 0.3877 0.3641 0.3622 0.3617 0.3826 0.4341 0.4018 0.3825 0.4036 Floss total 2 0.4341 0.43769 0.4163 0.4147 0.4144 0.0082 0.009 0.0198 0.0118 0.0109 0.0497 3 0.4018 0.4078 3 0.4178 0.3827 3 0.4018 0.4078 0.3823 0.3817 0.3814 0.0082 0.0099 0.0118 0.0109 Wrtloss 1.9561 1.98019 1.864 1.864 0.3824 0.0082 0.0099 0.0118 0.0109 Wrtloss 1.9561 1.98019 1.864 1.864 1.864 0.3824 0.324 0.0082 0.0099 0.0109 Wrtloss 0.029 0.020 0.02 0.02 0.02 0.02 0.02 0.02 0.02 0.02 0.02 0.02 0.02 0.02 0.02 0.02 0.02	nm Group															
0.3826 0.4341 0.4018 0.3825 0.4038 Floss total 2 0.4376 0.43769 0.4163 0.4147 0.4144 0.0082 0.0098 0.0118 0.0109 0.0497 3 0.4018 0.4078 0.3833 0.3817 0.3811 0.3744 0.4251 0.392 0.310 0.0497 3 0.4018 0.4078 0.3827 3 0.4078 0.3829 0.4018 0.4078 0.3829 0.4018 0.4078 0.3829 0.4018 0.4078 0.3829 0.4018 0.4019 0.3829 0.4019 0.4078 0.3829 0.4019 0.4078 0.3829 0.4019 0.4078 0.3829 0.4028 0.4019 0.3829 0.4028 0.4019 0.3829 0.4028 0.4019 0.3829 0.4028 0.4028 0.3829 0.4028 0.4029 0.4029 0.4029 0.4029 0.4029 0.4029 0.4029 0.4029 0.4029 0.4029 0.4029 0.4029 0.4029 0.4029 0.4029 <td>ample</td> <td></td> <td>7</td> <td>3</td> <td>4</td> <td>5</td> <td></td> <td>-</td> <td>0.3826</td> <td>0.3877</td> <td>0.3641</td> <td>0.3622</td> <td>0.3617</td> <td>0.3606</td> <td>0.3601</td> <td>0.36</td>	ample		7	3	4	5		-	0.3826	0.3877	0.3641	0.3622	0.3617	0.3606	0.3601	0.36
0.0082 0.0099 0.0118 0.0109 0.0497 3 0.4018 0.4078 0.3833 0.3817 0.3811 0.3744 0.4251 0.392 0.3707 0.3929 4 0.3825 0.3866 0.3642 0.3628 0.3628 0.3744 0.4251 0.3793 0.361 0.3827 4 0.3826 0.3466 0.3628 0.3628 0.3628 0.3628 0.3628 0.3628 0.3628 0.3628 0.3628 0.3628 0.3628 0.3628 0.3628 0.3628 0.3649 0.3649 0.3788 0.3640 0.3788 0.3640 0.3788 0.3788 0.3788 0.3788 0.407 0.02 0	art weight	0.3826	0.4341	0.4018	0.3825	0.4038	Floss total	2	0.4341	0.43769	0.4163	0.4147	0.4144	0.4131	0.4127	0.4127
6.04 6.3744 0.4251 0.3929 4 0.3825 0.3865 0.3642 0.3628 0.3624 0.36 0.4127 0.3733 0.361 0.3827 5 0.4038 0.4101 0.3868 0.3642 0.3649 0.3846 0.3846 0.3841 0.3841 0.3841 0.3841 0.3846 0.3846 0.3841 0.3841 0.3841 0.3841 0.3841 0.3841 0.3841 0.3842 0.3842 0.3842 0.3841 0.3841 0.3841 0.3841 0.3841 0.3841 0.3441 0.3742 0.3742 0.3742 0.3742 0.3742 0.3742 0.3742 0.3742 0.3742 0.02	oss Weight	0.0082	0.009	0.0098	0.0118	0.0109		m	0.4018		0.3833	0.3817	0.3811	0.38	0.3796	0.3793
0.36 0.4127 0.3793 0.361 0.3827 5 0.4038 0.4101 0.3868 0.3845 0.3841 0.3841 0.3841 0.3841 0.3841 0.3841 0.3841 0.3841 0.3841 0.3842 0.3841 0.0109 0.0109 0.0109 0.0109 0.0109 0.0109 0.0109 0.0109 0.0109 0.0109 0.0109 0.020 0.02 <td>// (start-floss)</td> <td>0.3744</td> <td>0.4251</td> <td>0.392</td> <td>0.3707</td> <td>0.3929</td> <td></td> <td>4</td> <td>0.3825</td> <td></td> <td>0.3642</td> <td>0.3628</td> <td>0.3624</td> <td>0.3613</td> <td>0.361</td> <td>0.361</td>	// (start-floss)	0.3744	0.4251	0.392	0.3707	0.3929		4	0.3825		0.3642	0.3628	0.3624	0.3613	0.361	0.361
Course C	nal constant weight	0.36	0.4127	0.3793	0.361	0.3827		2	0.40	0	0.3858	0.38	98	0.3	0.3	88.0
6.04 5.03 5.74 5.80 5.37 Mean solubility/sd of 2002 0.03 0.04 <th< td=""><td>loss Weight</td><td>0.0082</td><td>0.009</td><td>0.0098</td><td>0.0118</td><td>0.0109</td><td></td><td>Wt-floss</td><td>-</td><td><u>.</u></td><td></td><td>2000</td><td></td><td>_</td><td></td><td>and a second</td></th<>	loss Weight	0.0082	0.009	0.0098	0.0118	0.0109		Wt-floss	-	<u>.</u>		2000		_		and a second
6.04 5.03 5.74 5.80 5.37 San Diagram 0.02	72 (final-floss)	0.3518	0.4037	0.3695	0.3492	0.3/18		Mean	9.39					0.00	0.00	0.37
1 2 3 4 5 1 0.4378 0.4381 0.4109 0.4092 0.4086 0.4378 0.461 0.4646 0.4611 0.4209 Floss total 2 0.461 0.4637 0.4372 0.4364 0.4362 0.0103 0.0101 0.0088 0.0094 0.0093 0.0499 3 0.4496 0.4529 0.4265 0.4265 0.4261 0.4077 0.4116 4 0.4811 0.465 0.4265 0.4266 0.4261 0.4069 0.4242 0.4242 0.3664 0.0394 0.0394 0.0394 0.4811 0.485 0.4269 0.4568 0.0103 0.0101 0.0088 0.0094 0.0093 Writioss 2.2025 2.2176 2.0907 2.0791 0.3966 0.4245 0.4856 0.3894 Mean 0.44 0.42 0.42 0.42 0.3966 0.4256 0.4256 0.360 2.090 2.090 0.02 0.3967 </td <td>SOLIBILITY</td> <td>F 04</td> <td>5113</td> <td>574</td> <td>580</td> <td>537</td> <td></td> <td>S</td> <td>Mean solu</td> <td>bility/sd</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td>	SOLIBILITY	F 04	5113	574	580	537		S	Mean solu	bility/sd						
1 2 3 4 5 1 0.4378 0.4381 0.4109 0.4092 0.4096 0.4096 0.4096 0.4096 0.4096 0.4096 0.4096 0.4096 0.4096 0.4096 0.4096 0.4096 0.4096 0.4096 0.4097 0.4096 0.4496 0.4372 0.4354 0.4354 0.4352 0.01013 0.01014 0.0098 0.0094 0.0093 0.0496 3 0.4496 0.4595 0.4295 0.4286 0.4261 0.4069 0.4242 0.4549 0.3987 4 0.4811 0.485 0.4296 0.4591 0.4568 0.4006 0.4003 0.0103 0.0101 0.0088 0.0093 Wthloss 2.2055 2.2176 2.0907 2.0907 2.0791 0.3966 0.4253 0.4156 0.4856 0.3894 Mean 0.4 0.42 0.42 0.403 0.02 0.02 0.02 0.02 0.02 0.02 0.02																
0.4378 0.461 0.4496 0.4811 0.4209 Floss total 2 0.4637 0.4537 0.4362 0.4362 0.4362 0.4362 0.4362 0.4362 0.4362 0.4362 0.4362 0.4362 0.4362 0.4362 0.4362 0.4362 0.4362 0.4362 0.4263 0.4263 0.4263 0.4263 0.4263 0.4263 0.4263 0.4268 0.4268 0.4268 0.4268 0.4368 </td <td>mm Group</td> <td></td> <td>,</td> <td>~</td> <td>,</td> <td>ıc</td> <td></td> <td>-</td> <td>0.4378</td> <td>0.4381</td> <td>n 4109</td> <td>0.4092</td> <td>O 4086</td> <td>0.4077</td> <td>0.4069</td> <td>0.4069</td>	mm Group		,	~	,	ıc		-	0.4378	0.4381	n 4109	0.4092	O 4086	0.4077	0.4069	0.4069
0.0103 0.0101 0.0088 0.0094 0.0093 0.0479 3 0.4496 0.4529 0.4285 0.4285 0.4265 0.4261 0.4268 0.4269 0.4269 0.4568 0.4568 0.4568 0.4568 0.4568 0.4568 0.4568 0.4568 0.4568 0.4568 0.4568 0.4568 0.4008 0.4003 0.4003 0.428 0.4242 0.4549 0.3987 5 0.4209 0.4258 0.4024 0.4006 0.4003 0.0103 0.0101 0.0088 0.0094 0.0093 Wr.Hoss 2.2025 2.2176 2.0907 2.0791 0.3966 0.4233 0.4154 0.4455 0.3894 Mean 0.44 0.44 0.42 0.42 0.42 0.3966 0.4233 0.4154 0.4455 0.3894 Mean 0.04 0.04 0.02 0.02 0.02 0.02 0.02 0.02 0.02 0.02 0.02 0.02 0.02 0.02 0.02 0.02 0.02 </td <td>tart weight</td> <td>0.4378</td> <td>0.461</td> <td>0.4496</td> <td>0.4811</td> <td>0.4209</td> <td>Floss total</td> <td></td> <td>0.461</td> <td>0.4637</td> <td>0.4372</td> <td>0.4354</td> <td>0.4352</td> <td>0.4338</td> <td>0.4332</td> <td>0.4334</td>	tart weight	0.4378	0.461	0.4496	0.4811	0.4209	Floss total		0.461	0.4637	0.4372	0.4354	0.4352	0.4338	0.4332	0.4334
0.4275 0.4509 0.4408 0.4717 0.4116 4 0.4811 0.485 0.4569 0.4569 0.4568 0.4069 0.4334 0.4242 0.4549 0.3987 5 0.4209 0.4288 0.4006 0.4006 0.4003 0.0103 0.0101 0.0088 0.0093 Wr.floss 2.2025 2.2176 2.0902 2.0807 2.0791 0.3966 0.4233 0.4154 0.4455 0.3894 Mean 0.44 0.44 0.42 0.42 0.42 sq 0.020 0.02 0.02 0.02 0.02 0.02 0.02 0.02 0.02	loss Weight	0.0103	0.0101	0.0088	0.0094	0.0093			0.4496	0.4529	0.4285	0.4265	0.4261	0.4249	0.4242	0.4242
0.4069 0.4334 0.4242 0.3987 5 0.4209 0.4268 0.4006 0.4006 0.4003 0.0103 0.0101 0.0088 0.0094 0.0093 Wr.floss 2.2025 2.2176 2.0902 2.0807 2.0791 0.3966 0.4233 0.4154 0.4455 0.3894 Mean 0.44 0.44 0.42 0.42 0.42 sd 0.02 0.02 0.02 0.02 0.02 0.02 0.02 0.02 0.02	V1 (start-floss)	0.4275	0.4509	0.4408	0.4717	0.4116		4	0.4811	0.485	0.4591	0.4569	0.4568	0.4555	0.4548	0.4549
0.0103 0.0101 0.0088 0.0094 0.0093 Wt-floss 2.2025 2.2176 2.0902 2.0807 2. 0.3966 0.4233 0.4154 0.4455 0.3894 Mean 0.44 0.44 0.42 0.42 0.42 sd 0.02 0.02 0.02	inal constant weight	0.4069	0.4334	0.4242	0.4549	0.3987		5	0.4209	0.4258	0.4024	0.4006	0.4003	0.3992	0.3985	0.3987
ss) 0.3966 0.4233 0.4154 0.4455 0.3894 Mean 0.44 0.44 0.42 0.42 0.42 sd 0.02 0.02 0.02	loss Weight	0.0103	0.0101	0.0088	0.0094	0.0093		Wt-floss				2.	2	2	2	2
0.02 0.02 0.02	VZ (final-floss)	0.3966	0.4233	0.4154	0.4455	0.3894		Mean	0.44							
Manuschillianiad C 04	, mail 100	1	5	22.2	22.2	00.3		sq	0.02	0.02				0.02	0.02	70.U

Table 70. Vitremer One-Month Solubility Experimental Sample Results (grams)

VIIKEMER ONE MONIN SOLUBILITY (EAFENIMENTAL GROOT	3050	-												
		O'Haus An	O'Haus Analytical Plus Scale	Scale				pre H20	post HZO					
nm Group								6-Jan-06	7-Jan-06	8-Jan-06	9-Jan-06	10-Jan-06	11-Jan-06	12-Jan-06
Sample	-	2	3	4	2	7		0.2196	0.221	0.2065	0.2062	0.2054	0.2049	0.2051
Start weight	0.2196	0.2451	0.2364	0.2899	0.2261	0.2261 Floss total 2		0.2451	0.2453	0.23	0.2293	0.2284	0.2279	0.2281
Floss Weight	0.0134	0.013	0.0105	0.0112	0.0095	0.0576 3		0.2364	0.2363	0.2237	0.2231	0.2222	0.2219	0.2221
W1 (start-floss)	0.2062	0.2321	0.2259	0.2787	0.2166	4	4	0.2899	0.2898	0.2692	0.2685	0.2672	0.2669	0.267
Final constant weight	0.2051			0.267	0.2123	4)	5	0.2261	0.2269	0.2139	0.2134	0.2124	0.2121	0.2123
Floss Weight	0.0134			0.0112	0.0095	_	At-floss	1.1595	1.1617	1.0857	1.0829	1.078	1.0761	1.077
W2 (final-floss)	0 1917		0.2116	0.2558	0.2028		Mean	0.2319	0.23234	0.21714	0.21658	0.2156	0.21522	0.2154
						0)	ps	0.0277483	0.0272944	0.0243865	0.0242606	0.0277483 0.0272944 0.0243865 0.0242606 0.02408676	0	0.0240878
% SOI UBILITY	7.03	7.32	6.33	8.22	6.37			Mean solubility	pility	7.05	0.78			
2000000		ı	ı											
mm Group														
Sample	-	2	3	4	5		_	0.2947	0.295	0.272	0.2715	0.2702	0.2694	0.2703
Start weight	0.2947	0.3365	0.2917	0.3057	0.3051	0.3051 Floss total 2	2	0.3365	0.3374	0.3156	0.3149	0.3131	0.3126	0.3129
Floss Weight	0.0113	0.0101	0.0114	0.0105	0.0103	0.0536	3	0.2917	0.2922	0.2745	0.2736	0.2724	0.2718	0.2719
W1 (start-floss)	0.2834		0.2803	0.2952	0.2948	4	4	0.3057	0.306	0.2838	0.2829	0.2817	0.2812	0.2814
Final constant weight	0.2703		0.2719	0.2814	0.2821		5	0.3051	0.3055	0.2844	0.2836	0.2822	0.2819	0.2821
Floss Weight	0.0113		0.0114	0.0105	0.0103		Wt - floss	1.4801	1,4825	1.3767	1.3729	1.366	1.3633	1.365
W2 (final-floss)	0 259		0.2605	0.2709	0.2718		Mean	0.29602	0.2965	0.27534	0.27458	0.2732	0.27266	0.273
							ps	0.0177524	0.0177524 0.0179597 0.0174054	0.0174054		0.0174079 0.01717781		0.0172503 0.01717038
% SOLUBILITY	8.61	7.23	7.06	8.23	7.80			Mean solubility	billity	7.79	0.65			
nm Group														
Sample	-	2	3	4	2		-	0.4048	0.4072	0.3848	0.3833	0.3817	0.381	0.3811
Start weight	0.4048	,	0.4092	0.4106	0.3358	0.3358 Floss total 2	2	0.4392	0.44	0.4175	0.4162	0.4141	0.4133	0.4135
Floss Weight	0.0113		0.0125	0.0115	0.0105	0.0567	3	0.4092	0.4102	0.3859		0.3826	0.382	0.3822
W1 (start-floss)	0.3935		0.3967	0.3991	0.3253		4	0.4106	0.4113	0.3886	0.3871	0.3852	0.3847	0.3848
Final constant weight	0.3811		0.3822	0.3848	0.3153		5	0.3358	0.3361	0.3182	0.3171	0.3156	0.31	0.3153
Floss Weight	0.0113		0.0125	0.0115	0.0105		Wt-floss	1.9429	1.9481	1.8383	1.8316	1.8225	1.8195	1.8202
W2 (final-floss)	0.3698	0.4026	0.3697	0.3733	0.3048		Mean	0.38858	0.38962	0.36766	0.36632	0.3645	0.3639	0.36404
							ps	0.038335	0.0386076	0.0365797	0.0364758	0.038335 0.0386076 0.0365797 0.0364758 0.03626407	0.0361308	0.03616182
% SOLUBILITY	6.02	00.9	6.81	6.46	6.30			Mean solubility	billity	6.32	0.33			
nm Group														
Sample	-	2	3	4	5		-	0.4609	0.466	0.4342	0.4329	0.4308	0.4299	0.4303
Start weight	0.4609	0.4934	0.4585	0.4726	0.4475	0.4475 Floss total 2	2	0.4934	0.4958	0.4618	0.46	0.4577	0.4571	0.4571
Floss Weight	0.0116		0.0117	0.0106		0.0447	3	0.4585	0.4595	0.4264	0.4243	0.4223	0.4219	0.4218
W1 (start-floss)	0.4493	0.4826		0.462	0.4475		4	0.4726	0.4728	0.4378	0.4361	0.4338	0.4331	0.4332
Final constant weight	0.4303			0.4332	0.4111		5	0.4475	0.4	0.4	0.4	0.41	0.41	0.4111
Floss Weight	0.0116		0.0117	0.0106			Wt-floss	2.2882						2.1088
W2 (final-floss)	0.4187	0.4463	0.4101	0.4226	0.4111		Mean	0.45764	0.45932	0.42628	0.42448	0.42232	SEL	0.42176
							ps	0.0174438	0.0180585	0.0170679	0.0171906	0.0174438 0.0180585 0.0170679 0.0171906 0.01709248		0.0171229 0.01707879
% COLUBILITY	E 04	7 63	0 04	0 53										

Table 71. Vitremer One-Month Solubility Control Sample Results (grams)

START: 12 Dec 05		O'Haus Ans	O'Haus Analytical Plus Scale	Scale			pre H2O post H2O	0					
m Group							12-Jan-06 13-Jan-06	36 14-Jan-06	36 15-Jan-06	17-Jan-06 18-Jan-06	18-Jan-06	19-Jan-06 20-Jan-06	20-Jan-06
Sample	-	2	3	+	5	-	0.2239 0.2243	0.2123	0.2118	0.2112			
Start weight	0.2239	0.2277	0.2288	0.225		Floss	0.2277 0.2279	0.2142	0.2138	0.2131			
Floss Weight	0.014	0.0108	0.0131	0.0115	0.0106	0.06 3	0.2288 0.2292	0.2149	0.2145	0.214			
W1 (start-floss)	0.2099	0.2169	0.2157	0.2135	0.2227	4	0.225 0.2246	0.2118	0.2116	0.211			
Final constant weight	0.2112	0.2131	0.214	0.211	0.219	5	0.2333 0.233	0.2198	0.2194	0.219			
Floss Weight	0.014	0.0108	0.0131	0.0115	0.0106	Wt-floss	1.0787 1.079	1.013	1.0111	1.0083			
W2 (final-floss)	0.1972	0.2023	0.2009	0.1995	0.2084	Mean	0.21574 0.2158	58 0.2026	26 0.20222	2 0.20166			
						ps	0.003684 0.003588	88 0.003179	79 0.003155	5 0.003243			
SOLUBILITY	6.05	6.73	98.9	95.9	6.42		Mean solubility	9	6.52 0.31				
nm Group													
Sample	-	2	3	+	2	-	0.2798 0.2814	0.2615	0.2609			0.2595	
Start weight	0.2798	0.2916	0.3165	0.3069	0.3114	읎	0.2916 0.29	0.2678	0.2673	0.2663	0.2659	0.2659	
Floss Weight	0.0099	0	0.01	0.0117	0.0105	0.0421 3	0.3165 0.3167	0.2964	0.2959	0.2949	0.2945	0.2942	
W1 (start-floss)	0.2699	0.2916	0.3065	0.2952	0.3009	4	0.3069 0.3074	0.289	0.2885	0.2875	0.2869	0.2869	
Final constant weight	0.2595	0.2659	0.2942	0.2869	0.2901	5	0.3114 0.3121	0.2921	0.2916	0.2907	0.2903	0.2901	
Floss Weight	0.0099	0	0.01	0.0117	0.0105	Wt - floss	s 1.4641 1.4655	355 1.3647	1.3621	1.3574	1.355	1.3545	
W2 (final-floss)	0.2496	0.2659	0.2842	0.2752	0.2796	Mean	0.29282 0.2931	331 0.27294	94 0.27242	2 0.27148	0.271	0.2709	
							0.015172 0.01512	512 0.015638		7 0.015624	0.015663	0.015556	
% SOLUBILITY	7.52	8.81	7.28	6.78	7.08		Mean solubility	7	7.49 0.79	-			
mm Group				-			2000 0 2000 0	0000	0 2000				000
Sample	-			4	0	THE OWNER OF THE OWNER OWNER OF THE OWNER OWNE	1.39U5 U.392/	0.3639	0.3691				0.355
Start weight	0.3905		0.412	0.4036	0.3899	Floss total	0.4142 0.4153	0.3876	0.3868				0.383/
Floss Weight	0.0091			0.0106	0.0101	0.0528	0.412 0.4133	0.3904	0.3895		1		0.3866
W1 (start-floss)	0.3814			0.393	0.3798		0.4036 0.4051	0.3783	0.3775				0.3749
Final constant weight	0.366	0.3837	0.3866	0.3749	0.3576	5	0.3899 0.3933	0.3611	0.3603	0.3593	0.3586		0.3576
Floss Weight	0.0091	0.0116	0.0114	0.0106	0.0101	Wt-floss	1.9574 1.9669	569 1.8345	345 1.8304		1.821	1.81607	1.816
N2 (final-floss)	0.3569	0.3721	0.3752	0.3643	0.3475	Mean	0.39148 0.39336	338 0.3669	90996.0 698	3 0.36488	0.3642	0.363214	0.363
						ps	0.011511 0.010695	395 0.012204	204 0.012178	8 0.012099	0.012065	0.011933	0.012107
% SOLUBILITY	6.42	7.58	6.34	7.30	8.50	The second second	Mean solubility	7	7.23 0.89	6			
ım Groun													
Sample	-	2	3	4	2	-	0.466 0.4662	0.4314	0.4302	0.4288	0.428	0.4272	0.4271
Start weight	0.466	0.4662	0.4547	0.4669	0.4561	0.4561 Floss total 2	0.4662 0.4669	0.4344	0.4323	0.4307	0.4296	0.4292	0.4289
Floss Weight	0.0126	0	0.0102	0.0102	0.01	0.043 3	0.4547 0.4531	0.4198	0.4179	0.4164	0.4154	0.4144	0.4148
W1 (start-floss)	0.4534	0.4662	0.4445	0.4567	0.4461	4	0.4669 0.4667	0.434	0.4325	0.4312	0.4302	0.4307	0.4296
Final constant weight	0.4271	0.4289	0.4148	0.4296	0.4128	5	0.4561 0.4543	0.4178	0.4162	0.4148	0.4139	0.4134	0.4128
Floss Weight	0.0126	0	0.0102	0.0102	0.01	Wt-floss	2.2669 2.2642		344 2.0861	1 2.0789	2.0741	2.0719	2.0702
V2 (final-floss)	0.4145	0.4289	0.4046	0.4194	0.4028	Mean	men.		200		100		0.4140
						ps	0.006036 0.007083	0.0	8038 0.008079	9 0.008085	0.008064	0.008389	0.008152
% COLINEI TO	0 50												

Table 72. Vitremer Three-Month Solubility Experimental Sample Results (grams)

Sample S	VITREMER 5 MORITI SOLUBILITY (EAPERIMENTAL GROOT		A CHILL	Ation Divis	Cool														
1			IN SIDELLO	alytical mius	arga					nost H20									
C C C C C C C C C C	mm Group								2	8-Oct-05	9-0ct-05	10-0ct-05	11-0ct-05	12-0ct-05	13-0ct-05	14-Oct-05	15-Oct-05	16-Oct-05	17-Oct-05
Courty C	Sample	-	2	3	4	5		-	0.2072	0.2071	0.1981	0.198	0.1977	0.1959	0.1955				
Colorio Colo	start weight	0 2072	0.215	0.2202	0.2385	0.2352	Floss total	2	0.215	0.2149	0.2073	0.2065	0.2062	0.2057	0.2052				
1. 1. 1. 1. 1. 1. 1. 1.	loss Weight	0.0081	0.0084	0.007	0.0098	0.0084	0.0417		0.2202	0.2217	0.213	0.2118	0.2114	0.2103	0.21				
Vision V	V1 (start-floss)	0.1991	0.2066	0.2132	0.2287	0.2268			0.2305	0.2385	0.231	0.23	0.2301	0.2289	0.2284				
	inal constant weight	0.1955	0.2052	0.21	0.2284	0.2242		5	0.2352	0.2351	0.2265	0.2256	0.2256	0.2245	0.2242				
	loss Weight	0.0081	0.0084	0.007	0.0098	0.0084		Wt-Floss	1.0664	1.0756	1.0342	1.0302	1.0293	1.0236	1.0216				
1	V2 (final-floss)	0.1874	0.1968	0.203	0.2186	0.2158		Mean	0.21328	0.21512	0.20684	0.20604	0.20586	0.20472	0.20432				
The color of the								ps	0.011369	0.013283	0.013574	0.013294	0.01348	0.013583	0.013584				
1	SOLUBILITY .	5.88	4.74	4.78	4.42	4.85			Mean solu	billity	4.93	0.55							
1	Course Course				I		ı												
c) 114 c) 1206 c) 278 c) 284	Samole	-	2	3	4	2		-	0.314	0.3166	0.3042	0.3022	0.3019	0.3004	0.2997	0.3008	0.2991	0.2984	0.2983
1	tart weight	0 314	0 3205	0 2784	0.3154	0.2874	Floss total		0.3205	0.3219	0.3107	0.3089	0.3088	0.3072	0.3067	0.3079	0.3061	0.3052	0.3052
9j 0.3665 0.3365 0.3365 0.3365 0.3365 0.3365 0.3365 0.2374 0.2274	loss Weight	0.0075	0.0074	0.0068	0.0073	0.0071	0.0361		0.2784	0.2789	0.2698	0.2683	0.2681	0.2668	0.2669	0.268	0.2664	0.2658	0.265
National Carrell Carrell C	V1 (start-floss)	0.3065	0.3131	0.2716	0.3081	0.2803			0.3154	0.3137	0.3028	0.3008	0.3005	0.2991	0.2983	0.2995	0.2977	0.297	0.2968
Court Cour	inal constant weight	0.2983	0.3052	0.265	0.2968	0.2721		5	0.2874	0.2871	0.2765	0.2749	0.2747	0.2741	0.2735	0.2747	0.2731	0.2723	0.2721
1	loss Weight	0.0075	0.0074	0.0068	0.0073	0.0071		Wt-Floss	1.4796	1.4821	1.4279	1.419	1.4179	1.4115	1.409	1.4148	1.4063	1.4026	1.4013
1	V2 (final-floss)	0.2908	0.2978	0.2582	0.2895	0.265		Mean	0.29592	0.29642	0.28558	0.2838	0.28358	0.2823	0.2818	0.28296	0.28126	0.28052	0.28026
512 3 4.93 6.04 5.46 Mean solubility 5.29 0.48 3569 0.3756 0.3756 0.3757 0.3427 0.3427 0.3426 0.3427 0.3427 0.3427 0.3427 0.3427 0.3427 0.3427 0.3457 0.3575 0.3411 0.3389 0.3375 0.3411 0.3587 0.3547 <th< td=""><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td>ps</td><td>0.018904</td><td>0.019289</td><td>0.018338</td><td>0.018141</td><td>0.018137</td><td>0.017866</td><td>0.017628</td><td>0.017644</td><td>0.017552</td><td>0.017493</td><td>0.017747</td></th<>								ps	0.018904	0.019289	0.018338	0.018141	0.018137	0.017866	0.017628	0.017644	0.017552	0.017493	0.017747
1 2 3 4 5 1 0.3579 0.3457 0.3427 0.3423 0.3406 0.3397 0.3411 0.3389 0.3379 0.3779 0.3779 0.3770 0.3779 0.3471 0.3779 0.3779 0.3779 0.3779 0.3774 0.3779 0.3779 0.3779 0.3779 0.3770 0.0770 0.0770 0.0770 0.07710 0.07710 0.07710 0.07710 0.07710 0.07710 0.07710 0.07710 0.07710 0.07710 0.07710 0.	SOLUBILITY	5.12	4.89	4.93	6.04	5.46			Mean solu	billity	5.29	0.48							
1	mm Grain				ı		ı		Carried States and Assessment States and Ass										
0.3569 0.3756 0.3756 0.3756 0.3757 0.3692 0.3577 0.3587<	Sample	-	2	3	4	5		-	0.3579	0.3591	0.345	0.3427	0.3423	0.3406	0.3397	0.3411	0.3389	0.3379	0.3375
0.0078 0.0067 0.0071 0.0071 0.0072 0.0369 0.3794 0.3817 0.365 0.3627 0.3622 0.3604 0.3595 0.3594 0.3703 0.3685 0.3703 0.3703 0.3703 0.3685 0.3262 0.3703 0.0703 0.07170 0.017707 0.017707 0.017707 0.017707 0.017707 0.017707 0.017707 0.017707 0.017707 0.017707 0.017707 0.017707<	tart weight	0.3569	0.3756	0.3794	0.345	0.3891	Floss tota		0.3756	0.3765	0.3621	0.3597	0.3592	0.3575	0.3567	0.3579	0.3557	0.3547	0.3545
0.3491 0.3689 0.3723 0.3816 0.345 0.345 0.3472 0.3282 0.3252 0.3267 0.3267 0.3733 0.3708 0.3703 0.3685 0.3677 0.3697 0.3703 0.3686 0.3703 0.3687 0.3703 0.3703 0.3687 0.3703 0.3708 0.3703 0.3687 0.3708 0.3703 0.3687 0.3676 0.3703 0.3703 0.3867 0.3703 0.3708 0.3703 0.3883 0.3703 0.3687 0.3703 0.3703 0.3867 0.3703 0.3708 0.3703 0.3883 0.3703 0.3687 0.3703 0.3703 0.3867 0.3703 0.3708 0.3703 0.3867 0.3703 <td>loss Weight</td> <td>0.0078</td> <td>0.0067</td> <td>0.0071</td> <td>0.0071</td> <td>0.0076</td> <td>0.0363</td> <td></td> <td>0.3794</td> <td>0.3817</td> <td>0.365</td> <td>0.3627</td> <td>0.3622</td> <td>0.3604</td> <td></td> <td>0.3608</td> <td>0.3584</td> <td>0.3574</td> <td>0.3573</td>	loss Weight	0.0078	0.0067	0.0071	0.0071	0.0076	0.0363		0.3794	0.3817	0.365	0.3627	0.3622	0.3604		0.3608	0.3584	0.3574	0.3573
Court Cour	V1 (start-floss)	0.3491	0.3689	0.3723	0.3379	0.3815		4	0.345		0.3305	0.3282	0.3279	0.3259		0.3267	0.3242	0.3236	0.3232
Court Cour	inal constant weight	0.3375	0.3545	0.3573	0.3232	0.3655			0.3891	0.39	0.3733	0.3708	0.3703	0.3685		0.3691	0.3667	0.3658	0.3655
Secondary Colored	loss Weight	0.0078	0.0067	0.0071	0.0071	0.0076		Wt-floss	1.8107	1.8188	1.7396	1.7278	1.7256		98	1,7193	1.7076	1.7031	1.7017
The color of the	N2 (final-floss)	0.3297	0.3478	0.3502	0.3161	0.3579		Mean	0.047707	0.017216	0.047246	0.047462	0.017066	600	0 017097	0 047027	0.047064	0.046027	0.04034
1 2 3 4 5 1 0.4492 0.4463 0.4311 0.4286 0.4279 0.4259 0.4249 0.4263 0.4277 0.4163 0.4277 0.4277 0.4473 0.4277 0.4278 0.4277 0.4278 0.4277 0.4473 0.4477 0.4473 0.4474 0.4474 0.4474 0.4479 0.4479 0.4479 0.4479 0.4479 0.4479 0.4479 0.4479 0.4479 0.4479 0.4479 0.4479 0.4479 0.4479 0.4479 0.4479 0.4479 0.4479 0.4479	SOI UBILITY	99 9	5.72		6.45	6.19		200	Mean solu	billity	5.97	0.36							
1 2 3 4 5 1 0.4492 0.4453 0.4311 0.4286 0.4259 0.4259 0.4249 0.4269 0.4451 0.4286 0.4259 0.4249 0.4269 0.4451 0.4461 0.4483 0.4481							ı												
1 2 3 4 5 1 0.4492 0.4413 0.4218 0.4279 0.4259 0.4249 0.4263 0.4217 0.4227 0.4182 0.4219 0.4229 0.4249 0.4268 0.4219 0.4219 0.4186 0.4219 0.4419	imm Group																		
1 0.4492 0.4456 0.4456 0.4457 0.4583 0.4454 0.4454 0.4452 0.4450 0.4457 0.4452 0.4458 0.4417 0.4452 0.4458 0.4417 0.4452 0.4458 0.4417 0.4452 0.4458 0.4417 0.4452 0.4458 0.4417 0.4452 0.4458 0.4417 0.4452 0.4458 0.4417 0.4452 0.4458 0.4417 0.4452 0.4458 0.4417 0.4452 0.4458 0.4417 0.4452 0.4458 0.4417 0.4452 0.4459 0.4554 0.4458 0.4417 0.4452 0.4459 0.4554 0.4458 0.4417 0.4452 0.4459 0.4554 0.4473	Sample	-	2	3	4	5			0.4492		0.4311	0.4286	0.4279	0.4259		0.4263	0.4237	0.4227	0.4224
colory 0.0077 0.0072 0.007266 0.43376 0.43128 0.43674 0.42674 0.42674 0.42674 0.42674 0.42674 0.42674 0.42674 0.42674 0.42674 0.42674 0.42674 0.42674 0.42674 0.42674 0.42674 0.42674 <td>Start weight</td> <td>0.4492</td> <td>0.4436</td> <td>0.46/5</td> <td>0.4563</td> <td>0.475</td> <td>rioss tota</td> <td></td> <td>0.4430</td> <td></td> <td>0.4242</td> <td>0.4210</td> <td>0.4212</td> <td>0.4194</td> <td></td> <td>0.42</td> <td>0.417</td> <td>0.4100</td> <td>0.410</td>	Start weight	0.4492	0.4436	0.46/5	0.4563	0.475	rioss tota		0.4430		0.4242	0.4210	0.4212	0.4194		0.42	0.417	0.4100	0.410
0.4413 0.4589 0.4581 0.4583 0.4586 0.4442 0.4417 0.4419 0.4389 0.4381 0.4386 0.4386 0.4450 0.4586 0.4586 0.4587 0.4586 0.4586 0.4587 0.4586 0.4587 0.4586 0.4587 0.4586 0.4587 0.4587 0.4587 0.4587 0.4587 0.4587 0.4587 0.4474 0.4474 0.4587 0.4587 0.4587 0.4587 0.4474 0.4474 0.4474 0.4587 0.4587 0.4587 0.4474 0.4474 0.4587 0.4587 0.4587 0.4587 0.4587 0.4474 0.4474 0.4487 0.4587 0.4587 0.4587 0.4587 0.4587 0.4587 0.4587 0.4587 0.4587 0.4587 0.4584 0.4584 0.4587 0.01328 0.01328 0.01348 0.01348 0.01328 0.01328 0.01328 0.01328 0.01328 0.01328 0.01328 0.01328 0.01328 0.01328 0.01328 0.01328 0.01328 0.01328	Floss Weight	0.0079	0.0077		0.0072	0.0077	0.0382		0.4675		0.4515	0.4491	0.4483	0.447		0.4468	0.4441	0.4432	0.442
0.4224 0.415 0.4429 0.4354 0.4473 5 0.475 0.456 0.4554 0.4556 0.4557 0.4454 0.456 0.	W1 (start-floss)	0.4413	0.4359		0.4511	0.4673		4	0.4583		0.4442	0.4417	0.4409	0.4389		0.4393	0.4366	0.4358	0.4354
0.0079 0.0077 0.0077 0.0077 Wtfloss 2.2554 2.2601 2.1688 2.1564 2.1527 2.1437 2.1394 2.1454 2.1316 2.1272 ssj 0.0175 0.0077 0.00	Final constant weight	0.4224	0.416		0.4354	0.4473			0.475		0.456	0.4534	0.4526	0.4507		0.4512	0.4483	0.4474	0.447
0.4145 0.4083 0.4352 0.4282 0.4396 Mean 0.45108 0.45202 0.43376 0.43128 0.43054 0.42874 0.42788 0.42908 0.42632 0.42544 0.42874 0.43788 0.43998 0.43632 0.42544 0.42874 0.43788 0.43998 0.4328	Floss Weight	0.0079	0.0077	0.0077	0.0072	0.0077		Wt-floss	2.2554		2.1688	2.1564	2.1527	2.1437	2.1394	2.1454	2.1316	2.1272	2.1258
sd 0.01286 0.013816 0.01345 0.01339 0.01348 0.013509 0.013286 0.013286 0.013288	N2 (final-floss)	0.4145	0.4083	0.4352	0.4282	0.4396		Mean	0.45108	0.45202	0.43376	0.43128	0.43054	0.42874	0.42788	0.42908	0.42632	0.42544	0.42516
								ps	0.01286	0.013916	0.013462	0.013422	0.013339	0.013438	0.013509	0.013286	0.013275	0.013268	0.01331

Table 73. Vitremer Three-Month Solubility Control Sample Results (grams)

			2211									
								14-0ct	15-Oct	16-Oct	17-Oct	18-Oct
mm Group							Sample	pre H20	post H20			
Sample	-	2	3	4	2		-	0.2158	0.2159	0.2035	0.2032	
Start weight	0.2158	0.1883	0.2117	0.2105	0.2265 F	0.2265 Floss total	2	0.1883	0.1887	0.1795	0.1792	
Floss Weight	0.0073	0.0074	0.0071	6900.0	0.008	0.0367	3	0.2117	0.2117	0.1993	0.199	
W1 (start-floss)	0.2085	0.1809	0.2046	0.2036	0.2185		4	0.2105	0.2111	0.2005	0.2003	
Final constant weight	0.2132	0.1872	0.2089	0.2079	0.2248		5	0.2265	0.2268	0.2163	0.2161	
Floss Weight	0.0073	0.0074	0.0071	0.0069	0.008		Wt-floss	1.0161	1.0175	0.9624	0.9611	
N2 (final-floss)	0.2059	0.1798	0.2018	0.201	0.2168		Mean	0.20322	0.2035	0.19248	0.19222	
							ps	0.013949	0.013886	0.013217	0.01325	
% SOLUBILITY	4.59	4.51	4.68	4.51	4.28		Mean solubility	billity	4.51	0.15		
mm Group												
Sample	-	2	3	4	5		-	0.2871	0.2884		0.272	
Start weight	0.2871	0.3305	0.3202	0.3048	0.3079 F	0.3079 Floss total	2	0.3305	0.3319	0.3155	0.315	
Floss Weight	9900.0	9.0000	0.0064	0.0062	6900.0	0.0337	3	0.3202	0.3215	0.3038	0.3034	
W1 (start-floss)	0.2805	0.3229	0.3138	0.2986	0.301		4	0.3048	0.3064	0.2891	0.2887	
Final constant weight	0.272	0.315	0.3034	0.2887	0.2892		5	0.3079	0.3084	0.2897	0.2892	
Floss Weight	9900.0	0.0076	0.0064	0.0062	6900.0		Wt-floss	1.5168	1.5229	1.4369	1.4346	
N2 (final-floss)	0.2654	0.3074	0.297	0.2825	0.2823		Mean	0.30336	0.30458	0.28738	0.28692	
							ps	0.016432	0.016471	0.016302	0.016309	
% SOLUBILITY	5.38	4.80	5.35	5.39	6.21		Mean solubility	pillity	5.43	0.50		
Imm Group												
Sample	-	2	3	4	5		-	0.386	0.3885	0.3693	0.3685	0.3684
Start weight	0.386	0.3907	0.3852	0.3911	0.3788 F	0.3788 Floss total	2	0.3907	0.3928	0.3685	0.3675	0.367
Floss Weight	0.0067	9900.0	0.0064	0.0061	0.0067	0.0325	3	0.3852	0.3875	0.3693	0.3687	0.3686
W1 (start-floss)	0.3793	0.3841	0.3788	0.385	0.3721		4	0.3911	0.3932	0.3732	0.3723	0.3722
Final constant weight	0.3684	0.367	0.3686	0.3722	0.3575		5	0.3788	0.3799	0.3586	0.3577	0.3575
Floss Weight	0.0067	9900.0	0.0064	0.0061	0.0067		Wt-floss	1.8993		1.8064	1.8022	1.8012
W2 (final-floss)	0.3617	0.3604	0.3622	0.3661	0.3508		Mean	0.37986	0.38188	0.36128	0.36044	0.36024
							ps	0.004998	0.005373	0.005449	0.005476	0.005511
% SOLUBILITY	4.64	6.17	4.38	4.91	5.72		Mean solubility	pillity	5.17	0.75		
Comple		2	2	,	4		•	DACA O	0.4283	0 3953	0 3938	A595 A
Start weight	0 4244	0 4837	0 4942	0.4816	0.4733 F	0.4733 Floss total	2	0.4837		0.4586	0.4576	0.4573
Floss Weight	0.0076	0.0071	0.0077	0.0075	0.0072	0.0371	3	0.4942	0.4967	0.4682	0.4669	0.4667
W1 (start-floss)	0.4168	0.4766	0.4865	0.4741	0.4661		4	0.4816	0.4831	0.4561	0.455	0.455
Final constant weight	0.3936	0.4573	0.4667	0.455	0.4485		5	0.4733		0.4497	0.4486	0.4485
Floss Weight	9.0000	0.0071	0.0077	0.0075	0.0072		Wt-floss	2.3201	2.3337	2.1908		2.184
W2 (final-floss)	0.386	0.4502	0.459	0.4475	0.4413		Mean	0.46402	0.46674	B	0.43696	0.4368
							ps	0.02733	0.026681	0.028883	0.02903	0.02904
% SOLUBILITY	7.39	5.54	5.65	5.61	5.32		Mean solubility	pillity	5.90	0.84		

Table 74. Fuji II LC Immediate Solubility Experimental Sample Results (grams)

FUJI II LU IMMEDIATE SULUDILI II (EAPERIMENIAL GROUP)	20000			-																	
		O'Haus Analytical Plus Scale	lytical Plus	Scale																	
Zmm Group								4-0ct-05				8-0ct-05	9-0ct-05 1	0-0ct-05 1	1-Oct-05 1	2-0ct-05 1	3-0ct-05 1	4-0ct-05 1	9-0ct-05 10-0ct-05 11-0ct-05 12-0ct-05 13-0ct-05 14-0ct-05 15-0ct-05 16-0ct-05 17-0ct-05	-0ct-05 17	7-Oct-05
Sample	-	2	3	4	2		-	0.2466	0.2376	0.2372	0.237		-								
Start weight	0.2466	0.2597	0.2572	0.2493	0.2621 F	1 Floss total	2	0.2597	0.2502	0.2494	0.249							-			
Floss Weight	0.0103	0.0101	0.00	0.0089	0.00	0.0473	3	0.2572	0.2484	0.2489	0.2486			-					and the same and the same	-	
W1 (start-floss)	0.2363	0.2496	0.2482	0.2404	0.2531		4	0.2493	0.2404	0.24	0.2402										
Final constant weight	0.237	0.249	0.2486	0.2402	0.2524		2	0.2621	0.253	0.2526	0.2524										
Floss Weight	0.0103	0.0101	0.009	0.0089	0.00		Wt-floss	1.2276	1.1823	1.1808	1.1799										
W2 (final-floss)	0.2267	0.2389	0.2396	0.2313	0.2434		Mean	0 24552			0.23598										
							DS .	0.006/15 0.006601	1099		0.006515				-	-					
% SOLUBILITY	4.06	4.29	3.46	3.79	3.83		The second second	Mean solubility	oillity	3.89	0.31										
Jame Grass																					
Cample		2	2	,	4		-	0 336	0 3246	D 3234	0 3231	-									
Sample	0 326	0 2533	3535 0	0 3623	N 3556	Fines Intal	0	0.353	0.2414	0 3407	0 3403										
Start weight	0.330	0.0000	0.0000	0.0000	0.000	0.0547	2 6	0.3636	0 3545	0 3508	0 3606										
My form Boar	0.0100	0.0100	0.0105	0.0035	0.0103	0.0103	. 4	0.3623	0.306	0.3300	0.3386						İ				
WI (Statt-Huss)	0.3534	0.046.0	0.000	0 2400	0.2437			O SEEE	0 2437	0 242	0 2427										
Final constant weight	0.3237	0.3403	0.3500	0.3400	0.0400		20	4 7404	4 5504	4 666	4 5525					ı	İ				
Floss Weight	0.0106	80L0.0	CULU.U	0.0035	0.0103		WI-HOSS	1617.1	1.0001	0000	00000					-					
W2 (final-floss)	0.3125	0.3295	0.3401	0.339	0.3324		Mean	0.34382			0.3307										
							sq	0.011045	6384		0.010865										
% SOLUBILITY	3.96	3.80	3.68	3.91	3.74			Mean solubility	pility	3.82	0.12								+	1	
4mm Group												,000				1					
Sample	-	2	3	4	2		-	0.4378	0.4238	0.423	0.4229	0.4221		1							
Start weight	0.4378	0.4584	0.4538	0.4609	0.4677	0.4677 Floss total		0.4584	0.4439	0.443	0.4426	0.4221									
Floss Weight	0.0113	0.0104	0.0116	0.011	0.0104	0.0547		0.4538	0.4389	0.438	0.4376	0.437									
W1 (start-floss)	0.4265	0.448	0.4422	0.4499	0.4573		4	0.4609	0.4439	0.443	0.4427	0.442									
Final constant weight		0.4221	0.437	0.442	0.449		2	0.4677	0.4513	0.4499	0.4494	0.449									
Floss Weight	0.0113	0.0104	0.0116	0.011	0.0104		Wt-floss	2 2239	2.1471	2.1422	2.1405	2.1175									
W2 (final-floss)	0.4108	0.4117	0.4254	0.431	0.4386		Mean	0.444/8	638	0.42844	318	0.4235				-	Ì				
							bs	0.01120/	1970	0.010088		0.012044	the last contract of the	and the second second	Action was a second		The second second second second				
% SOLUBILITY	3.68	8.10	3.80	4.20	4.09	No. of Street, or other Persons and the street, or other persons are street, or other persons and the street, or other persons and the street, or other persons and the street, or other persons are street, or other persons and the street, or other persons are street, or other persons and the street, or other persons are street, or other persons and the street, or other persons are street, or other persons and the street, or other persons are street, or other persons and the street, or other persons are street, or other persons and the street, or other persons are street, or other persons and the street, or other persons are street, or other persons and the street, or other persons are STATE OF THE PARTY	Mean solubility	pility	4.71	1.87											
Smm Groun																					
Sample		,	2	Y	4			0.51	0.4958	0 4949	0 4943	0 4936	0 4921	0 491	0 4904	0 4892	0 4486	0 4893	0.4878	0 4872	0.487
Start weight	0.51	0 5404	0 5437	0 5481	0 5453	0 5453 Floss total	2	0 5404	0.5245	0.5236	0.523	0.5221	0.5206	0.5195	0.519	0.5175	0.5169	0.5178	0.5161	0.5156	0.5153
Floor Moinh	0.0107	0.0104	0.0104	0 0 100	0.0106	0.0623	3	0 5437	0 528	0 527	D 5264	0 5255	0 5242	0 523	0 5225	0 521	0 5205	0 5214	0 5197	0 5191	0 519
Willess meight	0 4993	0.63	0 5333	0 5379	0 5347			0 5481	0 5307	0 5303	0 5291	0 5282	0 527	0 5258	0 5253	0 524	0 5233	0.5241	0 5224	0 522	0 522
MI (statt-lines)		0 5453	0.000	0.530	0 540	-	· u	O EAE2	0.5974	0 6268	0 5254	0 527A	0.5227	0 5216	0 5243	0 5198	0 5192	0 5202	0 5183	0 518	0 518
Final constant weight		0.0103	610.0	770.0	0.010		200	0.0400	0.0271	0.2530	0.02531	0.0214	2 5242	2002	2 5253	2 5400	2010.0	2020	2 543	2 5005	2 500
Floss Weight	/0L0.0	0.0104	0.0104	0.0102	0.0100		VVI-IIOSS	70007	00007	2.3433	0000	C++C-7	C+CC.7	0076.7	70707	7510.7	70/4/7	CU2C.2	710.7	2.3030	2.303
W2 (final-floss)	0.4763	0.5049	9805 0	0.5118	0.5074		Mean	0.52/04	200	an.	887	0.5089		90I -	506		58D.	800		589	BLOCK O
							ps	0.015623 0.014382	4382	0.014415		0.014591	0.01429	0.014261	0.014324	0.014226	0.032004	0.014302	0.014194	0.014263	0.014323
% COLLIBILITY	461	474	4 63	4.85	5.11			Mean solubility	hility	4 79	0.00										

Table 75. Fuji II LC Immediate Solubility Control Sample Results (grams)

	i.O	O'Haus Analytical Plus Scale	Haus Analytical Plus	Scale								
m Group								11-Nov-05	post H2O 12-Nov-05 13-Nov-05 14-Nov-05 15-Nov-05	13-Nov-05	14-Nov-05	15-Nov-05
Sample	-	2	3	4	5		-	0.2233	0.2298	0.2194	0.2187	0.2187
Start weight	0.2233	0.2194	0.2368	0.2365	0.2321	0.2321 Floss total	2	0.2194	0.2263	0.2156	0.2151	0.2152
Floss Weight	0.0103	0.0095	0.0098	9600.0	0.0101	0.0493	3	0.2368	0.2446	0.2325	0.2318	0.2317
W1 (start-floss)	0.213	0.2099	0.227	0.2269	0.2274		4	0.2365	0.244	0.232	0.2314	0.2315
Final constant weight	0.2187	0.2152	0.2317	0.2315	0.2274		5	0.2321	0.2394	0.2277	0.2272	0.2274
Floss Weight	0.0103	0.0095	0.0098	9600.0	0.0101		Wt-floss	1.0988	1.1348	1.0779	1.0749	1.0752
W2 (final-floss)	0.2084	0.2057	0.2219	0.2219	0.2173		Mean	0.21976	0.22696	0.21558	0.21498	0.21504
							ps	0.007897	0.008347	0.007604	0.007576	0.007559
SOLUBILITY	2.16	2.00	2.25	2.20	4.44		ALLES STORY	Mean solubility	hillity	2.61	1.03	
ım Group							•	2000		0000	,000	0000
Sample	-	2	3	4	2000		- 0	0.3027	0.3121	0.2976	0.2964	0.2968
Start weight	0.3027	0.3425	0.3381	0.3326	0.3531	0.3531 Floss total	7	0.3425	0.3558	0.3397		0.3391
Floss Weight	0.0099	0.0089	0.0112	0.0106	0.0107	0.0513	3	0.3381	0.3462	0.3316		0.331
W1 (start-floss)		0.3336	0.3269	0.322	0.3424		4	0.3326	0.3424	0.3278		0.3269
Final constant weight		0.3391	0.331	0.3269	0.3462		5	0.3531	0.3627	0.3468		0.3462
Floss Weight	0.0099	0.0089	0.0112	0.0106	0.0107		Wt-floss	1.6177	1.6679	1.5922		1.5887
W2 (final-floss)	0.2869	0.3302	0.3198	0.3163	0.3355		Mean	0.32354	0.33358	0.31844	0.31766	0.31774
							ps	0.018943	0.019454	0.018875	0.019106	0.018961
% SOLUBILITY	2.02	1.02	2.17	1.77	2.02	STATE OF STA	The second	Mean solubility	ubility	1.80	0.46	
im Groun												
Sample		2	3	4	2		-	0 4359	0 4464	0.4284	0 4275	0 4276
Start weight	0.4359	0.4333	0.4652	0.4508	0.4247	0.4247 Floss total	2	0.4333		0.425		0.424
Floss Weight	0.0083	0.0107	0.0108	0.0082	0	0.038	3	0.4652		0.4562	0	0.455
W1 (start-floss)	0.4276	0.4226	0.4544	0.4426	0.4247		4	0.4508		0.4421	0.441	0.441
Final constant weight		0.424	0.455	0.441	0.4155		5	0.4247		0.4165	0	0.4155
Floss Weight	0.0083	0.0107	0.0108	0.0082	0		Wt-floss	2.1719	2.224	2.1302	2.1246	2.1251
W2 (final-floss)	0.4193	0.4133	0.4442	0.4328	0.4155		Mean	0.43438	0.4448	0.42604	0.42492	0.42502
							ps	0.01603	0.016916	0.015624	0.015577	0.015522
SOLUBILITY	1.94	2.20	2.24	2.21	2.17			Mean solubility	ubility	2.15	0.12	
ım Group												
Sample	-	2	3	4	2		-	0.533	0.5457	0.5241	0.5224	0.5223
Start weight	0.533	0.5363	0.541	0.538	0.5444	Floss total	2	0.5363	0.5472	0.5269	0.5255	0.5255
Floss Weight	0.0121	0.0102	0.0097	0.0091	0.0103	0.0514	3	0.541	0.5558	0.5317	0.5299	0.5298
W1 (start-floss)	0.5209	0.5261	0.5313	0.5289	0.5341		4	0.538	0.5485	0.5282	0.5266	0.5265
Final constant weight	t 0.5223	0.5255	0.5298	0.5265	0.5323		9	0.5444	0.554	0.5339	0.5323	0.5323
Floss Weight	0.0121	0.0102	0.0097	0.0091	0.0103		Wt-floss	2.6413	2.6998	2.5934	2.5853	2.585
W2 (final-floss)	0.5102	0.5153	0.5201	0.5174	0.522		Mean	0.52826	0.53996	0.51868	0.51706	0.517
							ps	0.00437	0.004414	0.003884	0.003857	0.003877
VTI HOH ICA	200	200	***	-								

Table 76. Fuji II LC 24-Hour Solubility Experimental Sample Results (grams)

ruji II CC 24 Hour s	Columnia	solubility Laper Illiental	TO TO	1	-		-		-	-	-			
,		O'Haus Analytical Plus Scale	lytical Plus	Scale	-									-
							_		post H20					
nm Group								29-Nov	30-Nov	1-Dec	1-Dec 2-Dec-05	3-Dec-05 4-Dec-05	4-Dec-05	5-Dec-05
Sample		2	3	4	2		-	0.2206	0.2244	0.2146	0.2141	0.2136		
Start weight	0 2206	0 2378	0 2382	0.224	0.2459	0.2459 Floss total	2	0.2378	0.2419	0.2312	0.2308	0.2303		
Flore Weight	0.0119	0 0102	0 0101	0.0103	0 0104	0 0529	3	0 2382	0 2418	0 2316	0 231	0 2305		
Mrt Jetart floor	0.2087	0 2276	0 2281	0 2137	0 2355		4	0 224	0 2275	0 2181	0 2176	0 2173		
i (stait-iinss)	0.200	0 2202	0 2205	0 2472	0 2204		· u	0 2450	0.00	0 2304	0 2306	0 2304		
Final constant weight	0.2130	0.2303	0.6200	0.0173	0.2301		0	4 4435	4 4207	4 0047	4 0700	4 0750		
Floss Weight	0.0119	0.0102	0.0101	0.0103	0.0104		Wt-floss	1.1136	1.132/	/180.T	1.0792	1.0769		
W2 (final-floss)	0.2017	0.2201	0.2204	0.207	0.2277		Mean	0.22272	0.22654	0.21634	0.21584	0.21538		
							ps	0.010616	0.010782	0.010224	0.010224	0.010181		
SOLUBILITY	3.35	3.30	3.38	3.14	3.31			Mean solubility	oillity	3.29	0.09			
nm Group														
Sample	-	2	3	4	5		-	0.322	0.3274	0.3135	0.3126	0.3122		
Start weight	0.322	0.3534	0.3342	0.3499	0.3357	0.3357 Floss total	2	0.3534	0.3593	0.3435	0.3429	0.3424		
Floss Weight	0.0103	0.0109	0.0111	0.0117	0.011	0.055	3	0.3342	0.3391	0.3256	0.3247	0.3241		
W1 (start-floss)	0.3117	0.3425	0.3231	0.3382	0.3247		4	0.3499	0.3555	0.3408	0.3399	0.3393		
Final constant weight	0.3122	0.3424	0.3241	0.3393	0.3251		5	0.3357	0.3418	0.3265	0.3256	0.3251		
Floss Weight	0.0103		0.0111	0.0117	0.011		Wt-floss	1.6402	1.6681	1.5949	1.5907	1.5881		
N2 (final-floss)	0.3019	0.3315	0.313	0.3276	0.3141		Mean	0.32804	0.33362	0.31898	0.31814	0.31762		
							ps	0.012739	0.012931	0.012275	0.012359	0.012312		
% SOLUBILITY	3.14	3.21	3.13	3.13	3.26			Mean solubility	billity	3.18	90.0			
										-	-			
mm Group		,			4		-	0.4522	0.4642	0.4478	0 4462	0 4454	0 AAE1	0.445
Sample	0 4502	0 4649	0 4507	0.4672	0 4357	0 4357 Flace total		0.4518	0.4681	0.4470	0.4480		0.4478	0.4478
Start weight	0.4303	0	0.430	0.400	0 0007	O OCAA	2 0	0.4507	0.4679	0 4306	0 4303		0.4373	0.4373
FIOSS Weight	0.0113		0.011	0.0102	0.003		0 4	0.4500	0.4577	0.4350	0.4303	0.4370	0.437	0.436
WI (start-rioss)	0.445		0.4373	0.4436	0.420		t u	0.4367	0.446	0.4758	0.4243	0.4736	0.4233	0.4433
Final Constant weight	0.443	0	0.4313	0.0100	0.4532		Wt floed	2 2003	2 2395	2 1549	2 1478	2 1441	2 1425	2 1424
AP Itinal Floral	0.4337	•	0 AZE	0.4333	0.4135		Mean	0.44186	0.4479	0 43098	0.42956	0.42882	0.4285	0 42848
(scort-iniii) 7/	1004.0	2004-0	0.7500	2001	200		ps	0.010337	0.010372	0.009779	0.009827	0.009854	0.009781	0.009818
% SOLUBILITY	2.98	3.11	3.05	3.06	2.93			Mean solubility	billity	3.03	0.07			
mm Groun														
Sample	-	2	3	4	5		-	0.538	0.5461	0.526	0.5242	0.5232	0.5227	0.5227
Start weight	0.538	0.5498	0.5189	0.5247	0.5474	0.5474 Floss total	2	0.5498	0.5562	0.5377	0.5358	0.5348	0.5345	0.5345
Floss Weight	0.0101	0.0109	0.0089	0.0105	0.0118	0.0522	3	0.5189	0.5246	0.5067	0.5051	0.5042	0.5038	0.5038
W1 (start-floss)	0.5279	0.5389	0.51	0.5142	0.5356		4	0.5247	0.53	0.5125	0.5105	0.5096	0.5094	0.5094
Final constant weight	0.5227		0.5038	0.5094	0.532		9	0.5474	0.5536	0.5351	0.5334	0.5323	0.5321	0.532
Floss Weight	0.0101		0.0089	0.0105	0.0118		Wt-floss	2.6266	2.6583	2.5658	2.5568	2.5519	2.5503	2.5502
W2 (final-floss)	0.5126	0.5236	0.4949	0.4989	0.5202		Mean	0.52532	0.53166	0.51316	0.51136	0.51038	0.51006	0.51004
							ps	0.01364	0.01414	0.013653	0.013628	0.013556	0.0135785	0.0135572
VIII II II II II II II	000	100	000	000	-						-			

Table 77. Fuji II LC 24-Hour Solubility Control Sample Results (grams)

FUJI II 24-HK SOLUBILI IT (CONTROL GROUP)	NOO LIT	NOL GRO	1 11					Ī	-				
		O'Haus Analytical Plus Scale	ytical Plus	Scale				pre H20	post H20				
nm Group								10-Jan-06	11-Jan-06	12-Jan-06	13-Jan-06	14-Jan-06	15-Jan-06
Sample	-	2	3	4	2	-		0.2482 0.2504	0.2504	0.24	0.2394	0.2389	0.239
Start weight	0.2482	0.233	0.2409	0.2507	0.2533	0.2533 Floss total 2		0.233 0.2355	0.2355	0.2249	0.2244	0.2239	0.2239
Floss Weight	0.0114	0.0106	0.0126	0.0116	0.0103	0.0565 3		0.2409 0.2444	0.2444	0.2329	0.2321	0.2319	0.2319
W1 (start-floss)	0.2368	0.2224	0.2283	0.2391	0.243	4		0.2507 0.2542	0.2542	0.2429	0.2422	0.2418	0.2418
Final constant weight	0.239	0.2239	0.2319	0.2418	0.2445	5		0.2533 0.256	0.256	0.2455	0.2446	0.2442	0.2445
Floss Weight	0.0114	0.0106	0.0126	0.0116	0.0103	Wt	Wt-floss	1.1696	1.184	1.1297	1.1262	1.1242	1.1246
W2 (final-floss)	0.2276	0.2133	0.2193	0.2302	0.2342	Mean	=	0.23392	0.2368	0.22594	0.22524	0.22484	0.22492
						ps		0.00824906	0.0083271	0.00824906 0.0083271 0.00835153 0.00825033 0.00825245 0.00833469	0.00825033	0.00825245	0.0083346
SOLUBILITY	3.89	4.09	3.94	3.72	3.62			Mean solubility	llity	3.85	0.18		
um Groun													
Sample		2	3	4	5	-		0.3593 0.3644	0.3644	0.3492	0.348	0.3474	0.3474
Start weight	0 3593	0.3617	0.3584	0.3644	0.3654	0.3654 Floss total 2		0.3617 0.3684	0.3684	0.3521	0.3508	0.35	0.3502
Floss Weight	0.0097	0.0119	0.0106	0.0102	0.0087	0.0511 3		0.3584 0.3644	0.3644	0.3488	0.3482	0.3474	0.3474
W1 (start-floss)	0.3496	0.3498	0.3478	0.3542	0.3567	4		0.3644 0.3696	0.3696	0.3547	0.3539	0.3531	0.3532
Final constant weight	0.3474	0.3502	0.3474	0.3532	0.3541	9		0.3654 0.3706	0.3706	0.3561	0.3549	0.3542	0.3541
Floss Weight	0.0097	0.0119	0.0106	0.0102	0.0087		Wt - floss	1.7581	1.7863	1.7098	1.7047	1.701	1.7012
W2 (final-floss)	0.3377	0.3383	0.3368	0.343	0.3454	Mean	ue.	0.35162	0.35726	0.34196	0.34094	0.3402	
						ps		0.00306317	0.0029175	0.0029175 0.00324145 0.00317695	0.00317695	0.00315785	0.00314452
% SOLUBILITY	3.40	3.29	3.16	3.16	3.17			Mean solubility	ility	3.24	11.0		
um Group													
Sample	-	2	3	4	2	-		0.4681 0.4733	0.4733	0.4553	0.4539	0.4527	0.4528
Start weight	0.4681	0.4536	0.4775	0.4561	0.4526	Floss total		0.4536 0.4607	0.4607	0.4428	0.4412	0.4402	0.4403
Floss Weight	0.011	0.0109	0.0103	0.011	0.0095	0.0527 3		0.4775 0.4833	0.4833	0.4656	0.464	0.4628	0.4628
W1 (start-floss)	0.4571	0.4427	0.4672	0.4451	0.4431	4		0.4561 0.4607	0.4607	0.4434	0.442	0.4409	0.4409
Final constant weight	0.4528	0.4403	0.4628	0.4409	0.4384			0.4526 0.458	0.458	0.44	0.43	0.43	
Floss Weight	0.011	0.0109	0.0103	0.011	0.0095		Wt-floss	2.2552	2.2833				
W2 (final-floss)	0.4418	0.4294	0.4525	0.4299	0.4289	Mean	an	0.45104	9/8	0.43906	0.43758		
						ps		0.01084329		0.0107884 0.01059316 0.01055021	0.01055021	0.01046217	0.01048012
SOLUBILITY	3.35	3.00	3.15	3.41	3.20			Mean solubility	ility	3.22	0.16		
nm Group													
Sample	-	2	3	4	2	-		0.5434 0.5501	0.5501	0.5289	0.5271	0.5258	0.5259
Start weight	0.5434	0.5321	0.5479	0.5304	0.5438	0.5438 Floss total 2		0.5321 0.5402	0.5402	0.5193	0.5175	0.5164	0.5162
Floss Weight	0.0109		0.0102	0.0089	0.0098	0.0508 3		0.5479 0.5533	0.5533	0.5337	0.532	0.5306	0.5306
W1 (start-floss)	0.5325	0.5211	0.5377	0.5215	0.534	4		0.5304 0.5357	0.5357	0.5162	0.5143	0.5133	0.5132
Final constant weight	0.5259	0.5162	0.5306	0.5132	0.527	5		0.5438 0.5498	0.5498	0.5299	0.5284	0.527	0.527
Floss Weight	0.0109	0.011	0.0102	0.0089	0.0098		Wt-floss	2.6468	2.6783				2.5621
W2 (final-floss)	0.515	0.5052	0.5204	0.5043	0.5172	Mean	an	0.52936	0.53566	0.51544	0.5137	0.51246	
						ps		0.00777541	0.0074851	0.0074851 0.00746726 0.00756987 0.00739135	0.00756987	0.00739135	0.0074761
VTI IIGII ITV	-		-										

Table 78. Fuji II LC One-Week Solubility Experimental Sample Results (grams)

Б	٥	"Haus Anal	O'Haus Analytical Plus Scale	Scale										
mm Group								pre H20 28-Nov-05	29-Nov-05	30-Nov-05	1-Dec-05	2-Dec-05	3-Dec-05	4-Dec-05
Sample	-	2	3	4	5									
Start weight	0.2134	0.224	0.2288	0.2287	0.2296	0.2296 Floss total	-	0.2134		0.2042	0.2036			
Floss Weight	0.0114	0.0112	0.0105	0.0105	0.0089	0.0525	2	0.224	0.2301	0.2191	0.2186	0.2184		
W1 (start-floss)	0.202	0.2128	0.2183	0.2182	0.2207		3	0.2288	0.2288	0.2184	0.2178	0.2179		
Final constant weight	0.2035	0.2184	0.2179	0.2185	0.2135		4	0.2287		0.219	0.2184			
Floss Weight	0.0114	0.0112	0.0105	0.0105	0.0089		5	0.2296	0.224	0.214	0.2135	0.2135		
W2 (final-floss)	0.1921	0.2072	0.2074	0.208	0.2046		Wt-floss	1.072	1.0734	1.0222	1.0194	1.0193		
							Mean	0.2144	0.21468	0.20444	0.20388	0.20386		
							ps	0.00679706	0.00667623	0.0063638	0.0063767	0.0064162		
% COLUBII ITY	4 9n	263	99 A	467	7 29			Mean solubility	lity	4.90	1.65			
	200	50.7	200	10.1	27.1			2000	6	2	3			
mm Groun														
Sample		2	3	4	2									
Start weight	0 3639	0 3242	0 352	0 3483	0.3581	0.3581 Floss total	-	0.3639	0.3648	0 3498	0.3485	0 3485	0 3482	0 3481
Floce Wainht	0.0102	0 0119	0.0121	0 0118	0 0102	0 0562	2	0 3242		0.3108	0.3102		0 3098	0 3097
Wil fetart florel	0 3537	0 3123	0 3366	0 3365	0 3479			0 352		0 3375		0	0 3363	0 3361
Final constant wortht	0 3481	0 3097	0 3361	0.3327	0 3414		4	0 3483		0 3341			0.3329	0 3327
Floe Weight	0.0102	0.0119	0.0121	0.0118	0 0102		. 2	0 3581	0.3599	0 3429	0.342		0 3414	0 3414
My Minal floor	0.3370	0 2078	N 324	0 3200	0 3312		Wt.flose	1 6903		1 6189	1 6141		16124	1 6118
VZ (IIIIIaI-11035)	0.000	0.63.0	V.354	0.0500	71000		Magan	0 33806	STATE OF	0 22278	0 33282	0 32278	0 32248	92662
							Po	0.01523893	0.01523893 0.01544225	0.014783	0 0145579	0	0.0145615	0.0145719
K. COLUBILITY	4.47	A 6A	4 68	4 64	4 80		2	Mean solubility	lity	4 64	0.12			
3OCODICIL I	17.1	-	20.7		20.						•			
mm Group														
Sample	-	2	3	4	5									
Start weight	0.4851	0.4358	0.4401	0.4337	0.4458	0.4458 Floss total	-	0.4851	0.4857	0.4646	0.4633	0.4629	0.4626	0.4624
Floss Weight	0.0104	0.0115	0.0119	0.0094	0.0104	0.0536	2	0.4358	0.4363	0.4167	0.4154	0.415	0.4148	0.4145
W1 (start-floss)	0.4747	0.4243	0.4282	0.4243	0.4354		3	0.4401	0.4424	0.4229			0.4209	0.4209
Final constant weight	0.4624	0.4145	0.4209	0.4134	0.4245		4	0.4337	0.4367	0.4154	0.4143		0.4136	0.4134
Floss Weight	0.0104	0.0115	0.0119	0.0094	0.0104		9	0.4458		0.4268		0.4251	0.4248	0.4245
W2 (final-floss)	0.452	0.403	0.409	0.404	0.4141		Wt-floss	2.1869		2.0928		2.0848	2.0831	2.0821
							Mean	0.43738	0.43936	0.41856	0.41736	0.41696	0.41662	0.41642
							ps	0.02119516	0.02059908	0.02028095	0.0202281	0.0202236	0.0202314	0.0202345
6 SOLUBILITY	4.78	5.02	4.48	4.78	4.89			Mean solubility	llity	4.79	0.20			
mm Group			•											
Sample	1	7	200	4	30 5174		•	0 5450	0 5505	02020	0 5353	0 5345	0 5343	A C 2 A
Start Weight	0.5400	0.0446	0.000	0.040	0.047	0.047 FIUSS LUCA	- 6	0.5400	0.5303	0.3260				0.324
Mid fetert floor	0.00	0.6326	0.6413	0.5324	0 5361	0.0000	3 6			0 5312				0 5281
Final constant woight	0 5241	0 5219	0.5281	0.5147	0 5169		4			0 5174		0 5153	0 5148	0 5147
IIIai colistalit weigitt	0.0241	0.0415	0.020	2010	440			0.5474		0 5407		0 5474		0 6160
Floss Weight	1110.0	0.0115	0.0105	0.0106	0.011		184 6122	0.5471		0.5197		9 5524		2010.0
WZ (final-floss)	0.5124	0.5104	0.3170	0.3041	ecoc.o		WI-HOSS	0 53548		0.3040	0 6113	0.51068	0 51016	0.5304
							Po Po		0 0	0.005531	0 0	0	0 0054623	0 0054137

Table 79. Fuji II LC One-Week Solubility Control Sample Results (grams)

O'Haus Analytical Plus		"Haus Ana	O'Haus Analytical Plus Scale	Scale											
								pre H20	post H20						
m Group									6-Dec-05 7-Dec-05	7-Dec-05	8-Dec-05	9-Dec-05	10-Dec-05	8-Dec-05 9-Dec-05 10-Dec-05 11-Dec-05 12-Dec-05	12-Dec-05
Sample	-	2	3	4	5		-	0.2244	0.2244	0.2137	0.2136	0.2137			
Start weight	0.2244	0.2284	0.2293	0.2305	0.2275	0.2275 Floss total	2	0.2284	0.228	0.2181	0.2179	0.2178			
Floss Weight	0.0119	0.0122	0.01	0.0108	0.0103	0.0552	3	0.2293	0.2288	0.2177	0.2176	0.2178			
W1 (start-floss)	0.2125	0.2162	0.2193	0.2197	0.2172		4	0.2305	0.2297	0.2186	0.2185	0.2186			
Final constant weight	0.2137	0.2178	0.2178	0.2186	0.2166		5	0.2275	0.2272	0.2168	0.2166	0.2166			
Floss Weight	0.0119	0.0122	0.01	0.0108	0.0103		Wt-floss	1.0849	1.0829	1.0297	1.029	1.0293			
W2 (final-floss)	0.2018	0.2056	0.2078	0.2078	0.2063		Mean	0.21698	0.21658	0.20594	0.2058	0.20586			
							ps	0.002308	0.002025	0.001949	0.001937	0.001926			
% SOLUBILITY	5.04	4.90	5.24	5.42	5.02		Mean	solubility		5.12	0.20				
m Group															
Sample	-	2	3	4	5		-	0.3728	0.3735	0.3569	0.3563	0.3563			
Start weight	0.3728	0.3267	0.3449	0.3427	0.3776	0.3776 Floss total	2	0.3267	0.3265	0.3124	0.3118	0.312			
Floss Weight	0.0102	0.0106	0.0084	0.01	0.0106	0.0498	3	0.3449	0.3453	0.3295	0.329	0.329			
W1 (start-floss)	0.3626	0.3161	0.3365	0.3327	0.367		4	0.3427	0.343	0.3279	0.3273	0.3275			
Final constant weight	0.3563	0.312	0.329	0.3275	0.3612		9	0.3776	0.3782	0.3619	0.3612	0.3612			
Floss Weight	0.0102	0.0106	0.0084	0.01	0.0106		Wt-floss	1,7149	1.7167	1.6388	1.6358	1.6362			
W2 (final-floss)	0.3461	0.3014	0.3206	0.3175	0.3506		Mean	0.34298	0.34334	0.32776	0.32716	0.32724			
							ps	0.021567	0.021888	0.020962	0.020924	0.02084			
SOLUBILITY	4.55	4.65	4.73	4.57	4.47		Mean	Mean solubility		4.59	0.10				
m Group		•					•	0 400	70010	0 4740	0 4700	0 4000	1001	0 4600	
Sample	-	7	5	4	0		- (0.4922	0.4937	21.14.0	0.4702	0.4030	0.4004	0.4002	
Start weight	0.4922	0.4615	0.4658	0.4533	0.4433	0.4433 Floss total	7	0.4615	0.4642	0.4438	0.4429	0.4424	0.4414		
Floss Weight	0.0097	0.0116	0.0097	0.0098	רטרטיט	60000	2	0.4656	0.4678	0.4483	0.4473		0.4456		
W1 (start-floss)	0.4825	0.4499	0.4561	0.4435	0.4332	-	4	0.4533	0.4545	0.4349	0.4338		0.4323		
Final constant weight	0.4682	0.4411	0.4456	0.4321	0.4227	-	9	0.4433		0.4252	0.4242	0.4238	0.4229	0.4227	
Floss Weight	0.0097	0.0116	0.0097	0.0098	0.0101		Wt-floss	2.2652	- 1	2.1725	2.1675	2.1651	2.1599	2.1588	
W2 (final-floss)	0.4585	0.4295	0.4359	0.4223	0.4126		Mean	0.45304	0.45486	0.4345	0.4335	0.43302	0.43198	0.43176	
							ps	0.018335	0.018319	0.017265	0.017278	0.017197	0.017107	0.017108	
SOLUBILITY	4.97	4.53	4.43	4.78	4.76		Mean	solubility		4.69	0.22				
-															
Sample		2		4	5		-	0.5629	0.5655	0 5419	0 5407	0 5399	0 5387	0 5384	0 5382
Start weight	0.5629	0.5409	0.5467	0.5598	0.5597	0.5597 Floss total	2	0.5409	0.5431	0.5206	0.5191		0.5174		0.5167
Floss Weight	0.0138	0.0129	0.0114	0.0124	0.012	0.0625	3	0.5467	0.5491	0.5259	0.5246	0.5239	0.5228	0.5223	0.5222
W1 (start-floss)	0.5491	0.528	0.5353	0.5474	0.5477		4	0.5598	0.5616	0.5372	0.5359	0.5353	0.534	0.5336	0.5336
Final constant weight	0.5382	0.5167	0.5222	0.5336	0.5352		9	0.5597	0.5624	0.5388	0.5375	0.5369	0.5357	0.5353	0.5352
Floss Weight	0.0138	0.0129	0.0114	0.0124	0.012		Wt-floss	2.7075	2.7192	2.6019	2.5953	2.592	2.5861	2.5839	2.5834
W2 (final-floss)	0.5244	0.5038	0.5108	0.5212	0.5232		Mean	0.5415	0.54384	0.52038	0.51906	74D0	0.51722		0.51668
							ps	0.009621	969600.0	0.009146	0.009238	0.009208	0.009143	0.009254	0.009241
WITH HOUSE W															

Table 80. Fuji II LC One-Month Solubility Experimental Sample Results (grams)

TUJI II ONE MOINTII SOLUDICII I (EXI ENIMEINIAE GINO)	(1000				0011				
	3	O'Haus Analytical Plus Scale	ytical Plus	Scale			Dre HZU	DOST HZU				
nm Group							3-Dec-05 4	4-Dec-05	5-Dec-05	6-Dec-05	7-Dec-05	8-Dec-05
Sample	-	2	3	4	5	-	0.2274 0.2291	2291	0.2182	2172	0.2166	0.2168
Start weight	0.2274	0.2416	0.2387	0.3087	0.2401	0.2401 Floss total 2	0.2416 0.2421	2421	0.2308	0.2294	0.2292	0.2292
Floss Weight	0.0105	0.0113	0.0114	0.0116	0.0119	0.0567 3	0.2387 0.2383	.2383	0.2282	0.227	0.2267	0.2273
W1 (start-floss)	0.2169	0.2303	0.2273	0.2971	0.2282	4	0.3087 0.3093		0.2949	0.2931	0.2284	0.2928
Final constant weight	0.2168	0.2292	0.2273	0.2928	0.2289	2	0.2401 0.2407	2407	0.2303	0.2291	2928	0.2289
Floss Weight	0.0105	0.0113	0.0114	0.0116	0.0119	Wt-floss	1.1998 1	1.2028	1.1457	2172.922	0.2289	1.1383
W2 (final-floss)	0.2063	0.2179	0.2159	0.2812	0.217	Mean	0.23996	0.24056	0.22914	434.5844	0.04578	0.22766
						ps	0.032574	0.032485	0.030847	971.2385	1309.341	0.030505
% SOLUBILITY	4.89	5.38	5.02	5.35	4.91		Mean solubility	illity	5.11	0.24		
	ı	ı			ı							
nm Group							0 2000	0,000	2000	07000	2000	1000
Sample	-	7	5	4	0		0.3211 0.3213	3213	0.3067	0.3043	0.3045	0.3045
Start weight	0.3211	0.3551	0.3527	0.3444	0.3365	FIOS	0.3551 0.3556	3556	0.3402	0.338	0.3377	0.3377
Floss Weight	0.0083	0.0085	0.0089	0.0088	0.0085	0.043	0.3527 0.3519		0.3375	0.3353	0.3349	0.3348
W1 (start-floss)	0.3128	0.3466	0.3438	0.3356	0.328	4	0.3444 0.3457		0.3302	0.3281	0.3275	0.3275
Final constant weight	0.3045	0.3377	0.3348	0.3275	0.3192	5	0.3365 0.3377	3377	0.3214	0.3197	0.3193	0.3192
Floss Weight	0.0083	0.0085	0.0089	0.0088	0.0085	Wt - floss	1.6668	1.6692	1.593	1.583	1.5809	1.5807
N2 (final-floss)	0.2962	0.3292	0.3259	0.3187	0.3107	Mean	0.33336	0.33384	0.3186	0.3166	0.31618	0.31614
						ps	0.013774	0.013626	0.013583	0.013385	0.013398	0.01339
SOLUBILITY	5.31	5.05	5.21	5.04	5.27		Mean solubility	illity	5.17	0.13		
			ı	ı	ı							
mm Group			6		u		0.4446.0.444		0 3059	0 3030	A 202A	0 3031
Sample	- 0	7	0 4420	2000	0 4052		0 4246 0 4240		0.0000	0.3332	0.002	0.4440
Start weight	0.4146	0.4316	0.4132	0.4320	0.0000	FIOSS TOTAL	0.4310 0.4319	4313	0.4 140	0.412	0.4112	0.4112
Floss Weight	0.0068	0.0000	0.0031	0.0077	0.0009	0.045	0.4132 0	1014	0.3304	0.3333	0.3333	0.000
W1 (start-floss)	0.4058	0.423	0.4041	0.4249	0.4764		0.4326 0.4332	4332	0.4124	0.41	0.4093	0.4093
Final constant weight	0.3921	0.4112	0.3932	0.4093	0.4012		0.4253 0.425	425	0.4042	0.4018	0.4013	0.4012
Floss Weight	0.0088	0.0086	0.0091	0.0077	0.0089		2.0742	2.0767	1.9803		1.9644	1.9639
W2 (final-floss)	0.3833	0.4026	0.3841	0.4016	0.3923	Mean	0.41484	0.41534	0.39606		0.39288	0.39278
						ps	0.009178	0.008901	0.008741	0.0	0.008732	0.008835
% SOLUBILITY	5.54	4.82	4.95	5.48	5.79		Mean solubility	illity	5.32	0.41		
,		I	I		ı							
mm Group							0.000		0 1010	0 7.40	0.5450	0 540
Sample	-	7	3	4	c		0.5452 0.5444	2444	0.5212	0.518	0.5169	0.5168
Start weight	0.5452	0.5473	0.5351	0.5581	0.5456	E	0.5473 0.5492	.5492	0.5225		0.5178	0.51//
Floss Weight	0.0082	0.0087	0.0106	0.0094	0.0089	0.0458 3	0.5351 0.5354		0.5141		9609.0	0.5094
W1 (start-floss)	0.537	0.5386	0.5245	0.5487	0.5367	4	0.5581 0.5572		0.5335	0.5302	0.5292	0.5292
Final constant weight	0.5168	0.5177	0.5094	0.5292	0.5171	5	0.5456 0.5465	.5465	0.5212	0.518	0.5172	0.5171
Floss Weight	0.0082	0.0087	0.0106	0.0094	0.0089	Wt-floss	2.6855	2.6869	2.5667	2.55	2.5449	2.5444
W2 (final-floss)	0.5086	0.509	0.4988	0.5198	0.5082	Mean	0.5371	0.53738	0.51334	0.51	0.50898	0.50888
						ps	0.008176	0.007896	0.006981	0.007004	0.007032	0.007102
VI IIII III	E 20	6 50	4 90	5.27	531		Moan collubility	illity	525	000		

Table 81. Fuji II LC One-Month Solubility Control Sample Results (grams)

FUJI II ONE MONTH SOLUDILITY (CONTROL GROOF)				-					2						
START: 06 Dec 2006		O'Haus Analytical Plus Scale	ytical Plus	Scale				pre HZU post HZU	Q						
mm Group								6-Dec-05 7-Dec-05	9	9-Dec-05	10-Dec-05	11-Dec-05	12-Dec-05	13-Dec-05	14-Dec-05
Sample	-	2	3	4	2		_	0.2118 0.2122		0.2031	0.2027	0.2023			
Start weight	0.2118	0.2363	0.2414	0.24	0.2431	0.2431 Floss total 2	2	0.2363 0.2372		0.2259	0.2253	0.2253			
Floss Weight	0.0099	0.0088	0.0105	0.0114	0.0121	0.0527 3	3	0.2414 0.2415		0.2315	0.2309	0.2307			
W1 (start-floss)	0.2019	0.2275	0.2309	0.2286	0.231		4	0.24 0.2405		0.2289	0.2284	0.2282			
Final constant weight	0.2023	0.2253	0.2307	0.2282	0.2313		5	0.2431 0.2435	0.2325	0.2321	0.2317	0.2313			
Floss Weight	0.0099	0.0088	0.0105	0.0114	0.0121		Wt-floss	1.1199 1.1222	1.071	1.0688	1.0663	1.0651			
W2 (final-floss)	0.1924	0.2165	0.2202	0.2168	0.2192		Mean	0.22398 0.2	0.22444 0.2142	12 0.21376	5 0.21326	5 0.21302			
							ps	0.0129455 0.0129367		0.012109 0.0121021 0.01205653	0.01205653	3 0.0121185			
% SOLUBILITY	4.71	4.84	4.63	5.16	5.11		September 1	Mean solubility	4.89	89 0.24					
	I	I					ı								
imm Group		(0 2500 0 254		0 3353	2000		0000	10000	2000
Sample	-	7	3	4	0 2000			0.3509 0.3512	0.3365	0.3353	0.3346	0.3341	0.3339	0.3334	0.3335
Start weight	0.3509	0.3278	0.3394	0.3336	0.3420	0.3426 FIOSS TOTAL 2	7	0.3270 0.3201		0.3130	0.3020	0.3122	0.3124	0.3119	0.000
Floss Weight	0.0092	0.0085	0.0094	CULU.U	0.0034	0.047.3	2	0.3394 0.3394		0.3248	0.3241			0.3229	0.3229
W1 (start-floss)	0.3417	0.3193	0.33	0.3231	0.3332		4	0.3336 0.3337		0.3186	0.3179			0.3166	0.317
Final constant weight	0.3335	0.3119	0.3229	0.317	0.3264		2	2	0.329	0.32	0.321	0.32		0.3264	0.3264
Floss Weight	0.0092	0.0085	0.0094	0.0105	0.0094		Wt - floss						1.5669	1.5644	1.5647
W2 (final-floss)	0.3243	0.3034	0.3135	0.3065	0.317		Mean	0.32946 0.3	0.32966 0.3156	2201		2 0.31344		0.31288	
							ps	0.0087925 0.0088189	0.008	26 0.0083909		0.00849688 0.00843759	0.00831998	0.00834847	0.00834943
% SOLUBILITY	5.09	4.98	2.00	5.14	4.86	Salam Malah	A STATE OF	Mean solubility	5.	5.01 0.11					
2000							ı								
Sample		2		4	4		-	0 4492 0 4484	0 4304	0 4289	0 4277	0.427	0 4269	0 4263	0 4264
Start weight	0 4492	0.4519	0 4606	0 439	0 4621	0 4621 Floss total 2	2	0.4519 0.4522		0.4319	0.431	0.4303	0.4299	0.4294	0.4295
Flose Weight	0 0089	0 0094	9600 0	0000	0 0084	0.0453 3	3	0.4606 0.4613		0.4409	0.4398	0.4391	0.4389	0.4384	0.4385
W1 (start-floss)	0 4403	0.4425	0.451	0.43	0.4537		4	0.439 0.4404		0.4203	0.4193	0.4186	0.4185	0.418	0.418
Final constant weight	0 4264	0.4295	0.4385	0.418	0.4405		5	0.4621 0.4628		0.4433	0.4421	0.4414	0.4412	0.4406	0.4405
Floss Weight	0.0089	0.0094	9600.0	0.00	0.0084		Wt-floss	2.2175 2	2.2198 2.1277	77 2.12	2 2.1146	5 2.1111	2.1101	2.1074	2.1076
W2 (final-floss)	0.4175	0.4201	0.4289	0.409	0.4321		Mean	0.4435 0.4	0.44396 0.42554	54 0.424	4 0.42292	2 0.42222	0.42202	0.42148	0.42152
							ps	0.0093709 0.0092937	2937 0.0095454	54 0.0093246	6 0.0092686		0.0092686 0.00922887	0.00921293	0.00919277
% SOLUBILITY	5.18	90.9	4.90	4.88	4.76			Mean solubility	4	4.96 0.16	9				
Simm Group	ı														
Sample	-	2	3	*	5		-	0.546 0.5466	0.5255	0.5235	0.5221	0.5212	0.5209	0.5201	0.52
Start weight	0.546	0.5701	0.5644	0.4795	0.5151	0.5151 Floss total 2	2	0.5701 0.5758	0.5492	0.5469	0.5456	0.5446	0.5444	0.5437	0.5436
Floss Weight	0.0099	0.0097	0.0098	0.0087	0.0086	0.0467 3	3	0.5644 0.5668	0.5437	0.5414	0.5401	0.5391	0.5389	0.5381	0.5382
W1 (start-floss)	0.5361	0.5604	0.5546	0.4708	0.5065		4	0.4795 0.4813	0.4601	0.4586	0.4573	0.4566	0.4565	0.4559	0.4588
Final constant weight	0.52	0.5436	0.5382	0.4558	0.4915		5	0.5151 0.5159	0.496	0.4943	0.4932	0.4925	0.4924	0.4914	0.4915
Floss Weight	0.0099	0.0097	0.0098	0.0087	0.0086		Wt-floss				8 2.5116			2.5025	2.5054
W2 (final-floss)	0.5101	0.5339	0.5284	0.4471	0.4829		Mean	0.52568 0	0.52794 0.50556	96 0.5036		2 0.50146	0.50128	0.5005	
							ps	0.0377281 0.0	0.038836 0.0370086	0.036		0.0366249 0.03648294 0.03642619	0.03642619	0.03640038	0.03530938
% COLUBILITY	A 9.5	173	4 70	503	33 Y			Manage and Laboratory		400	-				

Table 82. Fuji II LC Three-Month Solubility Experimental Sample Results (grams)

Fuji II THREE MONTH SOLUBIL		III (EAPERIMENTAL GROUP)	MENIAL	L LOON								-
		O'Haus Analytical Plus Scale	ytical Plus	Scale			pre H20	post H20				
m Group			Ī				3-Nov-05	4-Nov-05	7-Nov-05	8-Nov-05	9-Nov-05	10-Nov-05
Sample	-	2	3	4	5	-	0.213	0.2137 0.2145	0.2051	0.2502		
Start weight	0.2137	0.221	0.2265	0.2196	0.2204	0.2204 Floss total 2	0.22	0.221 0.2217	0.2123	0.2123		
Floss Weight	0.0071	0 008	0.0072	0.0079	0.0079	0.03813	0.226	0.2265 0.2277	0.2173	0.2173		
W1 (start-floss)	0 2066	0.213	0.2193	0.2117	0.2125		0.219	0.2196 0.2218	0.2113	0.2114		
Final constant weight	0 2052	0 2123	0 2173	0.2114	0.212	5	0.220	0.2204 0.2221	0.2121	0.212		
Flore Weight	0 0071	0 008	0 0072	0 0079	0 0079	Wt-floss	1.06	1 0697	1.02	1.0651		
W7 (final floce)	0 1981	0 2043	0 2101	0 2035	0 2041	Mean		2 0 2 1 3 9 4	0.204			
(contraction)						ps	0.0045544	4 0.0046881	0.0046881 0.0043465	0		
% SOLUBILITY	411	4.08	4.20	3.87	3.95		Mean solubility	ubility	4.04			
m Group												
Sample	-	2	3	4	5	-	0.323	0.3234 0.3256	0.3108	0.3108	0.3109	
Start weight	0.3234	0.2963	0.3074	0.3058	0.3111	0.3111 Floss total 2	0.296	0.2963 0.2995	0.2864	0.2861	0.286	
Floss Weight	0.0081	0.0079	0.0077	0.0074	0.0078	0.0389 3	0.307	0.3074 0.3092	0.2939	0.2937	0.2934	
W1 (start-floss)	0.3153	0 2884	0.2997	0.2984	0.3033	4	0.305	0.3058 0.3086	0.2934	0.2936	0.293	
Final constant weight	0.3109	0 286	0.2934	0.293	0.2985	5	0.311	0.3111 0.313	0.299	0.2988	0.2985	
Floss Weight	0.0081	0.0079	0.0074	0.0074	0.0078	Wt - floss	oss 1.5051	1 1.517	1,4446	1,4441	1.4429	
W2 (final-floss)	0.3028	0.2781	0.286	0.2856	0.2907	Mean	0.30102	2 0.3034	0 28892	0.28882	0.28858	
						ps	0.0098191	1 0.0094637	0.0090681	0.0091397	0.0092651	
% SOLUBILITY	3.96	3.57	4.57	4.29	4.15		Mean solubility	ubility	4.11	0.37		
	ŀ											
Cample		6	3	4	2	-	0 414	0 4148 0 4191	0 399	0 3988	0 3982	0.398
Start weight	0 4148	0 4221	0.4262	0.383	0.4325	Floss total 2	0.422	0.4221 0.4275	0.4086	0.4083	0.4076	0.4072
Floce Weight	0 0083	0 0075	0 0071	9900 0	0 0086	0 0086 0 0381 3	0.426	0.4262 0.4293	0.4097	0.4097	0.4089	0.4083
W1 letart-floed	0.4065	0 4146	0 4191	0.3764	0 4239	4	0.38	0.383 0.3852	0.3678	0.3677	0.3671	0.367
Final constant weight	0.398	0.4072	0.4083	0.367	0.4136		0.432	0.4325 0.4345	0.4144	0.4145	0.4138	0.4136
Floss Weight	0.0083		0.0071	9900.0	0.0086	Wt-floss		5 2.0575	1.9614	1.9609	1.9575	1.956
W2 (final-floss)	0.3897		0.4012	0.3604	0.405	Mean	0.4081	11 0.4115	0.39228	0.39218		0.3912
						ps	0.0193878	8 0.0197545	0.0187963	0.0188266	0.0187696	0.01865294
% SOLUBILITY	4.13	3.59	4.27	4.25	4.46		Mean solubility	ubillity	4.14	0.33		
m Group				ľ								
Sample	-	2	3	4	5	-	0.519	0.5194 0.5253	0.501	0.5008	0.5003	0.4999
Start weight	0.5194	0.4841	0.5421	0.4894	0.5092	0.5092 Floss total 2	0.484	0.4841 0.487	0.4678	0.4674	0.4669	0.4666
Floss Weight	0.0084	0.0073	0.0078	0.008	0.0083	0.0398 3	0.542	0.5421 0.545	0.5208	0.5207	0.5201	0.5197
W1 (start-floss)	0.511	0.4768	0.5343	0.4814	0.5009	4	0.489	0.4894 0.493	0.4733	0.4727	0.4723	0.4715
Final constant weight	0.4998	0.4665	0.5197	0.4714	0.4904	9	0.50	0.5092 0.5133	0.4917	0.4915	0.49	0.4904
Floss Weight	0.0084	0.0073	0.0078	0.008	0.0083	Wt-floss	ss 2.5044	4 2.5238	2.4148			2.4083
W2 (final-floss)	0.4914	0.4592	0.5119	0.4634	0.4821	Mean	0.50088	18 0.50476	and a		3048	Name of Street
				-		ps	0.0234939	9 0.0237223	0.0214478	0.021	0.0215625	0.02161798
W COLLIBILITY	200	000	4 40		-		****					

Table 83. Fuji II LC Three-Month Solubility Control Sample Results (grams)

START: 08 Aug 2005		O'Haus Analytical Plus Scale	ytical Plus	Scale			pre H20	post H20				
m Group							8-Nov-05	9-Nov-05	10-Nov-05	10-Nov-05 11-Nov-05 12-Nov-05	12-Nov-05	13-Nov-05
Sample	-	2	3	4	5	-	0.2224	0.2224 0.2222	0.212	0.2069	0.2117	0.2119
Start weight	0.2224	0.23	0.233	0.2096	0.2395	0.2395 Floss total 2	0.23	0.23 0.2312	0.2199	0.2193	0.2189	0.2188
Floss Weight	0.0074	0.0075	0.0075	0.0094	0.0087	0.0405 3	0.233	0.233 0.2316	0.222	0.2212	0.2211	0.2212
W1 (start-floss)	0.215	0.2225	0.2255	0.2002	0.2308	4	0.2096	0.2096 0.2107	0.2006	0.2003	0.2001	0.2001
Final constant weight	0.2119	0.2188	0.2212	0.2001	0.2287	9	0.2395	0.2395 0.2402	0.2296	0.2289	0.2286	0.2287
Floss Weight	0.0074	0.0075	0.0075	0.0074	0.0087	Wt-floss	1.094	1.0954	1.0436	1.0361	1.0399	1.0402
W2 (final-floss)	0 2045	0 2113	0 2137	0.1927	0.22	Mean	0.2188	0.21908	0.20872	0.20722	0.20798	0.20804
						ps	0.011458	0.011198	0.011022	0.011525	0.010776	0.01079
SOLUBILITY	4.88	5.03	5.23	3.75	4.68		Mean solubility	billity	4.72	0.58		
ım Group												
Sample	-	2	3	4	2	-	0.2973	0.2973 0.2979	0.2838	0.2832	0.2828	0.2826
Start weight	0.2973	0.2945	0.2936	0.3028	0.2946	0.2946 Floss total 2	0.2945 0.294	0.294	0.2793	0.2787	0.2784	0.2782
Floss Weight	0.0078	0.008	0.0077	9900.0	0.0067	0.0368 3	0.2936	0.2936 0.2943	0.2816	0.2813	0.2809	0.2806
W1 (start-floss)	0.2895	0.2865	0.2859	0.2962	0.2879	4	0.3028	0.3028 0.2995	0.2838	0.283	0.2826	0.2825
Final constant weight	0.2826	0.2782	0.2806	0.2825	0.2784	5	0.2946	0.2946 0.2951	0.2794	0.2789	0.2785	0.2784
Floss Weight	0.0078	0.008	0.0077	9900.0	0.0067	Wt - floss	s 1.446	1.444	1.3711	1.3683	1.3664	1.3655
W2 (final-floss)	0.2748	0.2702	0.2729	0.2759	0.2717	Mean	0.2892	0.2888	0.27422	0.27366	0.27328	0.2731
						ps	0.003753	0.00242	0.002225	0.002158	0.002131	0.002128
% SOLUBILITY	90.5	69.6	4.55	6.85	5.63		Mean solubility	pillity	5.56	0.86		
ım Group												
Sample	-	2	3	4	5	-	0.454	0.454 0.4553	0.4351	0.4337	0.4331	0.4331
Start weight	0.454	0.3436	0.4194	0.427	0.4257	0.4257 Floss total 2	0.3436	0.3436 0.3432	0.328	0.3272	0.3267	0.3264
Floss Weight	9.0000	0.0077	0.0072	0.0077	0.0082	0.0384 3	0.4194	0.4194 0.4208	0.4012	0.3999	0.3995	0.3993
W1 (start-floss)	0.4464	0.3359	0.4122	0.4193	0.4175	4	0.427	0.427 0.4274	0.4093	0.4083	0.4076	0.407
Final constant weight	0.4331	0.3264	0.3993	0.4077	0.4076	5	0.4257	0.4257 0.4259	0.4093	0.4081	0.4077	0.4076
Floss Weight	9/00.0	0.0077	0.0072	0.0077	0.0082	Wt-floss	2.0313	2.0342		1.9388		1.935
W2 (final-floss)	0.4255	0.3187	0.3921	0.4	0.3994	Mean	0.40626					0.387
						ps	0.041507	0.042084	0.04041	0.040206	0.040177	0.040244
% SOLUBILITY	4.68	5.12	4.88	4.60	4.34		Mean solubility	pillity	4.72	0.29		
nm Group												
Sample	-	2	3	4	5	-	0.5361	0.5361 0.538	0.5164	0.5148	0.5142	0.514
Start weight	0.5361	0.5183	0.522	0.5312	0.5405	0.5405 Floss total 2	0.5183	0.5183 0.5196	0.4972	0.4955	0.4948	0.4947
Floss Weight	0.00	0.0081	0.0077	0.0071	0.007	0.0389 3	0.522	0.522 0.5258	0.5013	0.4996	0.4889	0.4989
W1 (start-floss)	0.5271	0.5102	0.5143	0.5241	0.5335	4	0.5312	0.5312 0.5329	0.511	0.5089	0.5084	0.5085
Final constant weight	0.514	0.4947	0.4989	0.5085	0.5175	5	0.5405	0.5405 0.5414	0.5202	0.5182	0.5176	0.5175
Floss Weight	0.00	0.0081	0.0077	0.0071	0.007	Wt-floss	2.6092	2.6188			2.485	
W2 (final-floss)	0.505	0.4866	0.4912	0.5014	0.5105	Mean	0.52184	0.52376			0.497	0.49892
						ps	0.009342	0.008893	0.009779	0.009697	0.012429	0.009721
VTI IIGII ITV	** *	-	** *					-				

Table 84. Photac-Fil Quick Immediate Solubility Experimental Sample Results (grams)

Filorac Illillieurate Sociolici I pro								001	0011					
		O'Haus Analytical Plus Scale	ytıcal Plus	Scale				DZH aud	DOST HZO					
m Group								19-Jan-06	17-Dec-05	18-Dec-05	19-Jan-06 17-Dec-05 18-Dec-05 19-Dec-05 20-Dec-05 21-Dec-05 22-Dec-05	20-Dec-05	21-Dec-05	22-Dec-06
Sample	-	2	3	4	5		-	0.2564	0.2564 0.2508	0.2318	0.2307	0.2305		
Start weight	0.2564	0.2526	0.2643	0.2601	0.2705	0.2705 Floss total	2	0.2526 0.246	0.246	0.2275	0.265	0.2262		
Floss Weight	0.0117	0.0116	0.0118	0.0101	0.0113	0.0565	3	0.2643 0.2584	0.2584	0.2425	0.2413	0.2412		
W1 (start-floss)	0.2447	0.241	0.2525	0.25	0.2592		4	0.2601 0.2662	0.2662	0.2478	0.2466	0.2463		
Final constant weight	0.2459	0.2409	0.253	0.255	0.2607		5	0.2705 0.2662	0.2662	0.2447	0.2438	0.2434		
Floss Weight	0.0117	0.0116	0.0118	0.0101	0.0113		Wt-floss 1.2474		1.2311	1.1378	1.1709	1.1311		
N2 (final-floss)	0.2342	0.2293	0.2412	0.2449	0.2494		Mean	0.24948	0.24622	0.22756	0.23418	0.22622		
							ps	0.006954	0.009074	0.008749	0.01246	0.008698		
% SOLUBILITY	4.29	4.85	4.48	2.04	3.78			Mean solubility	billity	3.89	1.10			
ım Group														
Sample	-	2	3	4	2		-	0.3465	0.3465 0.3344	0.3125	0.311	0.3106	0.31	
Start weight	0.3465	0.3573	0.3622	0.3482	0.3494	0.3494 Floss total	2	0.3573	0.3573 0.3495	0.3266	0.3251	0.3245	0.3241	
Floss Weight	0.0123	0.0122	0.0126	0.0117	0.011	0.0598	3	0.3622	0.3622 0.3552	0.3332	0.3313	0.3307	0.3302	
W1 (start-floss)	0.3342	0.3451	0.3496	0.3365	0.3384		4	0.3482	0.3482 0.3572	0.3337	0.332	0.3316	0.331	
Final constant weight	0.3285	0.3394	0.348	0.3372	0.3261		5	0.3494	0.3494 0.3698	0.3459	0.3441	0.3435	0.3431	
Floss Weight	0.0123	0.0122	0.0126	0.0117	0.011		Wt-floss	1.7038	1.7063	1.5921	1.5837	1.5811	1.5786	
W2 (final-floss)	0.3162	0.3272	0.3354	0.3255	0.3151		Mean	0.34076	0.34126	0.31842	0.31674	0.31622	0.31572	
							ps	0.006727	0.012871	0.012184	0.012053	0.011993	0.012047	
% SOLUBILITY	5.39	5.19	4.06	3.27	6.89		RAME AND ADDRESS OF THE PERSON NAMED IN COLUMN TWO IN COLU	Mean solubility	billity	4.96	1.38			
din Group		2		•	4		•	0.4637	0 4637 0 4387	0.4083	0.4062	0.4054	0.4049	0 4045
Salliple	0 4627	0.4602	0 A626	0.4741	0 AED3	0 AERS Fines total		0.4603	0.4603 0.4521	0.4238	0.4221	0.4211	0.4208	0.4205
Start weight	0.400	0.400	0.4020	0.114	0.400	0.067	3 6	0.4626	0.4626 0.4638	0.4226	0.4303	0.4292	0.4287	0.4285
rioss weight	0.0135	0.0123	0.0133	0.014	0.0123	ğ	7	0.4023 0.403	0.4030	0.432	0.450	0.4204	0.440	0.4404
W1 (start-floss)	0.4502	0.44/4	0.4400	0.4530	0.440		4	0.4741	0.4741 0.455	0.450	0.4421	0.4204	0.4150	0.4134
rinal constant weight	0.435	0.4376	0.44	0.4510	0.434		With floors	2 2630	2 2222				2 0637	
Floss weignt	0.0135	0.0123	0.0139	0.0144	0.0123		VVI-IIOSS	2.2333		1		No.		1
W2 (final-floss)	0.4195	0.4249	0.4261	0.43/4	0.4217		Mean	0.005735	0.01599	0.015575	0.015499	0.015251		0.041046
% SOLUBILITY	6.82	5.03	5.02	4.85	5.87		8	Mean solubility	billity					
nm Group					4		-	0 5581	0 5581 0 5536	0.516	0 5135	0.5125	0 5116	0.5116
Start weight	0 5581	0 565	0 556	0 5543	0 5536	0 5536 Floss total	2	0.565	0.565 0.5544	0.5161	0.5139	0.5127	0.5121	0.512
Floss Weight	0 0113		0.0118	0.0107	0.0112	0.056	3	0.556	0.556 0.5682	0.5339	0.5313	0.5301	0.5293	0.529
W1 (start-floss)	0.5468		0.5442	0.5436	0.5424		4	0.5543	0.5543 0.5772	0.5421	0.5399	0.5386	0.5379	0.5379
Final constant weight	0.526	0.5377	0.5214	0.5206	0.5272		5		0.5536 0.5616	0.5249	0.5226	0.521	0.5207	0.5206
Floss Weight	0.0113	0.011	0.0118	0.0107	0.0112		Wt-floss	2.731	2.759	2.577		2.5589		
W2 (final-floss)	0.5147	0.5267	9605.0	0.5099	0.516		Mean	0.5462	0.5518	0.5154	0.51304			0.51102
							ps	0.00459	0.009912	0.011391	0.011402	0.011337	0.01133	0.01131

Table 85. Photac-Fil Quick Immediate Solubility Control Sample Results (grams)

		O'Haus Ana	Analytical Plus Scale	Scale				Och ese	0001 1000					
om Grossa								16-Dec-05	17-Dec-05	18-Dec-05	19-Dec-05	16-Dec-05 17-Dec-05 18-Dec-05 19-Dec-05 20-Dec-05 21-Dec-05	21-Dec-05	22-Dec-05
amalo		2	2	V	4		-	0 2403	0 2508	0 2318	0 2307	0 2305		
Sample	0 0 0 0	7367	00700	0 2546	0 2524	0 2524 Elana tatal		0 2367	0.246	0 2276		0 2262		
Start weight	0.2403	0.2307	0.2499	0.0000	0.000	PIOSS LOCAL	7 (0.2307	0 0 0 0	0.0405		0.2202		
Floss Weight	1110.0	1710.0	1110.0	0.0036	COLO.O	0.0550	2	0.2433	0.2564	0.2423	0.2413			
W1 (start-floss)	0.2286	0.2246	0.2382	0.245	0.2426		4	0.2546	0.2662	0.2478	0.2466			
Final constant weight	0.2305	0.2262	0.2412	0.2463	0.2434		5	0.2531	0.2662	0.2447	0.2438	0.2434		
Floss Weight	0.0117	0.0121	0.0117	9600.0	0.0105		Wt-floss	1.179	1.232	1.1387	1.1333	1.132		
W2 (final-floss)	0.2188	0.2141	0.2295	0.2367	0.2329	(A)(B)	Mean	0.2358	0.2464	0.22774	0.22666	0.2264		
							ps	0.007974	0.009074	0.008749	0.008715	0.008698		
SOLUBILITY	4.29	4.67	3.65	3.39	4.00			Mean solubility	oility	4.00	0.51			
am Group														
Sample	-	2	3	4	5		-	0.3214	0.3344	0.3125	0.311	0.3106	0.31	
Start weight	0.3214	0.336	0.3436	0.3425	0.3539	0.3539 Floss total	2	0.336	0.3495	0.3266	0.3251	0.3245	0.3241	
Floss Weight	0.0112	0.0122	0.011	0.0109	0.01	0.0553	3	0.3436	0.3552	0.3332	0.3313	0.3307	0.3302	
W1 (start-floss)	0.3102	0.3238	0.3326	0.3316	0.3439		4	0.3425	0.3572	0.3337	0.332	0.3316	0.331	
Final constant weight	0.31	0.3241	0.3302	0.331	0.3431		5	0.3539	0.3698	0.3459	0.3441	0.3435	0.3431	
Floss Weight	0.0112	0.0122	0.011	0.0109	0.01		Wt-floss	1.6421	1.7108	1.5966	1.5882		1.5831	
W2 (final-floss)	0.2988	0.3119	0.3192	0.3201	0.3331		Mean	0.32842	0.34216	0.31932	0.31764	0.31712	0.31662	
							ps	0.011969	0.012871	0.012184	0.012053	0.011993	0.012047	
% SOLUBILITY	3.68	3.68	4.03	3.47	3.14			Mean solubility	billity	3.60	0.33			
				ı										
Sample		2	3	4	5		-	0.4201	0.4387	0.4083	0.4062	0.4054	0.4049	0.4045
Start weight	0 4201	0 4352	0 4451	0.4344	0.4618	0.4618 Floss total	2	0.4352	0.4521	0.4258	0.4221	0.4211	0.4208	0.4205
Floss Weight	0 0103	9600 0	0.0119	0.0112	0.0112	0.0542	3	0.4451	0.4638	0.4325	0.4303	0.4292	0.4287	0.4285
W1 (start-floss)	0 4098	0 4256	0 4332	0 4232	0.4506		4	0.4344	0.453	0.423			0.4196	0.4194
Final constant weight	0 4045	0.4205	0.4285	0.4194	0.4465		5	0.4618	0.4816	0.4508	0.4486	0.4472	0.4467	0.4465
Floss Weight	0.0103	9600.0	0.0119	0.0112	0.0112		Wt-floss	2.1424	2.235	2.0862	2.074	2.0691	2.0665	2.0652
W2 (final-floss)	0.3942	0.4109	0.4166	0.4082	0.4353		Mean	0.42848	0.447	0.41724	0.4148	0.41382	0.4133	0.41304
							ps	0.015406	0.01599	0.015476	0.015499	0.015251	0.015262	0.0153304
% SOLUBILITY	3.81	3.45	3.83	3.54	3.40			Mean solubility	billity	3.61	0.20			
am Grain														
Sample	-	2	3	4	2		-	0.5293	0.5536	0.516	0.5135	0.5125	0.5116	0.5116
Start weight	0.5293	0.5301	0.5461	0.5546	0.5381	0.5381 Floss total	2	0.5301	0.5544	0.5161	0.5139		0.5121	0.512
Floss Weight	0.0115	0.0114	0.0097	0.0086	0.0101	0.0513	3	0.5461	0.5682	0.5339	0.5313	0.5301	0.5293	0.529
W1 (start-floss)	0.5178	0.5187	0.5364	0.546	0.528		4	0.5546	0.5772	0.5421	0.5399		0.5379	0.5379
Final constant weight	0.5116	0.512	0.529	0.5379	0.5206		5	0.5381	0.5616	0.5249	0.5226	0.5211	0.5207	0.5206
Floss Weight	0.0115	0.0114	0.0097	0.0086	0.0101		Wt-floss	2.6469	2.7637	2.5817	2.5699		2.5603	2.5598
W2 (final-floss)	0.5001	0.5006	0.5193	0.5293	0.5105		Mean	0.52938	0.55274	0.51634	0.51398		0.51206	0.51196
							sq	0.010792	0.009912	0.011391	0.011402	0.011333	0.01133	0.0113103
% COLUMNITY	2 40	2 40								-	-			

Table 86. Photac-Fil Quick 24-Hour Solubility Experimental Sample Results (grams)

									001						
		O'Haus Analytical Plus Scale	lytical Plus	Scale				DZH evd	DOST HZO						
mm Group							4	4-Feb-06 5	5-Feb-06	90-qa-1-9	7-Feb-06	8-Feb-06	9-Feb-06	10-Feb-06	11-Feb-06
Sample	-	2	3	4	2	-	_	0.2257 0	0.2267	0.2088	0.2074	0.2072	0.2072		
Start weight	0.2257	0.2298	0.2356	0.2414	0.2554	0.2554 Floss total 2	_	0.2298 0	0.2313	0.2136	0.2123	0.212	0.212		
Floss Weight	0.0097	0.0111	0.0092	0.0091	0.0114	0.0505 3	_	0.2356 0	0.2367	0.2206	0.2192	0.2189	0.2189		
W1 (start-floss)	0.216	0.2187	0.2264	0.2323	0.244	4	_	0.2414 0	0.2426	0.225	0.2234	0.2232	0.2232		
Final constant weight	0.2072	0.212	0.2189	0.2232	0.2361	5	_	0.2554 0	0.257	0.2379	0.2366	0.2363	0.2361		
Floss Weight	0.0097	0.0111	0.0092	0.0091	0.0114	3	Wt-floss '	1.1374	1.1438	1.0554	1.0484	1.0471	1.0469		
W2 (final-floss)	0.1975	0.2009	0.2097	0.2141	0.2247	M	Mean	0.22748	0.22876	0.21108	0.20968	0.20942	0.20938		
						ps		0.0115945	0.0117551	0.0112384	0.0112424	0.0112231	0.0115945 0.0117551 0.0112384 0.0112424 0.0112231 0.0111484		
% SOLUBILITY	8.56	8.14	7.38	7.83	7.91		STATE OF	Mean solubility	ility	7.96	0.43				
Smm Group															
Sample	-	2	3	4	2	-		0.367 0.371		0.3454	0.3432	0.3426	0.3423		
Start weight	0.367	0.3693	0.3472	0.3045	0.3615	0.3615 Floss total 2		0.3693 0.3724		0.3512	0.349	0.3485	0.3483		
Floss Weight	0.0108	0.0104	0.0118	0.0102	0	0.0432		0.3472 0.35		0.3273	0.3253	0.3246	0.3242		
W1 (start-floss)	0.3562	0.3589	0.3354	0.2943	0.3615			0.3045 0.3065		0.2851	0.2834	0.2829	0.2826		
Final constant weight	0.3423	0.3483	0.3242	0.2826	0.3405	5		0.3615 0.3641		0.3435	0.3413	0.3	0.3		
Floss Weight	0.0108	0.0104	0.0118	0.0102	0	3	Wt - floss	1.7063	1.7208	1.6093	1.599				
W2 (final-floss)	0.3315	0.3379	0.3124	0.2724	0.3405	Z	Mean	0.34126	0.34416	0.32186	0.3198				
						ps		0.0267945 0.0273616 0.0268835 0.0266665	0.0273616	0.0268835	0.0266665	0.0266709	0.0266812		
% SOLUBILITY	6.93	5.85	6.86	7.44	5.81		10 St. 15 C	Mean solubility	ility	6.58	0.72				
				ı											
4mm Group					4	•		2320 0 4540		0.4274	0 4356	0 4349	0 4345	0.4330	0 4337
Sample	- 0		0 000			- 0		0.4545	I	0.4314	0.4330	0.4540	0.4545	0.4303	0.4300
Start weight	0.4595		0.4612		0.4444	ĕ		0.45346 0.4623		0.434	0.4352	0.4310	0.4313	0.4307	0.4300
Floss Weight	0.0096	0.0034	2010.0		1010.0	0.0469		0.4625 0.4625		0.4373	0.4357	0.435	0.4340	0.4545	0.434
W1 (start-floss)	0.4499	0.4472	0.451		0.4343	4		0.4402 0.4318		0.4008	0.3992	0.3987	0.3984	0.3978	0.3975
Final constant weight	0.4337		0.434		0.4198			0.4501 0.4522		0.4233	0.4	0.42	0.47	0.41	0.4198
Floss Weight	9600.0		0.0102		0.0101		Wt-floss	2.1934	2.2264	2.0839	1			0000000	2.0669
W2 (final-floss)	0.4241	0.4214	0.4238	0.3879	0.4097		Mean	0.43868 0.44528 0.41678	0.44528	0.0455444	0.41500	0.04144		0.0153002 0.01530	0.41330
W. COLLIBILITY	573	577	E 03	6 78	5.66	Pa .	1000	Mean soluhility	IIItv	6.00	0.0154635			2001	
N SOCOULILL	9.50		20.0		200	I									
5mm Group															
Sample	-	2	3	4	2	-		0.5246 0.5331	5331	0.5008	0.4988	0.498	0.4971	0.4965	0.4969
Start weight	0.5246	0.5206	0.55	0.5201	0.5284	0.5284 Floss total 2		0.5206 0.526		0.4937	0.4917	0.4911	0.4901	0.4895	0.4899
Floss Weight	0.0116	0.0097	0.0108	0.0085	0.0105	0.0511 3		0.55	0.55 0.5577	0.5254	0.5234	0.5227	0.5218	0.521	0.5213
W1 (start-floss)	0.513	0.5109	0.5392	0.5116	0.5179	4		0.5201 0.5285	.5285	0.4939	0.4919	0.4912	0.4903	0.4895	0.4898
Final constant weight	0.4969	0.4899	0.5213	0.4898	0.4967	9		0.5284 0.5372		0.501	0.4988	0.4981	0.4372	0.4963	0.4967
Floss Weight	0.0116	0.0097	0.0108	0.0085	0.0105		Wt-floss	2.5926	2.6314	2.4637	2.4535	2.45	2.3854	2.4417	2.4435
W2 (final-floss)	0.4853	0.4802	0.5105	0.4813	0.4862		Mean	0.51852	0.52628	0.49274	0.4907	0.49	0.47708		
						sq		0.0123506	0.0126069	0.0130374	0.0130452	0.0123506 0.0126069 0.0130374 0.0130452 0.0130318 0.0308753	0.0308753	0.0130103	0.0129847
WITH HOLLING W															

Table 87. Photac-Fil Quick 24-Hour Solubility Control Sample Results (grams)

THOUSE STEEDINGS SOCIETY (SOUTHERS STORE)		A	0	-				OCT CA	OCT TOOK							
	0	UHaus Analytical Plus Scale	Aical Plus	Scale				DZU ald	OZU 180d							
Imm Group								1-Feb-06	2-Feb-06	3-Feb-06	4-Feb-06	5-Feb-06	90-qa-J-9	7-Feb-06	8-Feb-06	9-Feb-06
Sample	-	2	3	4	2	-		0.238	0.2403		0.2238	0.2234				
Start weight	0.238	0.232	0.2372	0.2433	0.2626	0.2626 Floss total 2		0.232	0.2346	0.2204	0.2193	0.219				
Floss Weight	0.0102	0.0094	0.0106	0.0095	0.0103	0.05 3		0.2372	0.2389	0.2246	0.2233	0.223				
W1 (start-floss)	0.2278	0.2226	0.2266	0.2338	0.2523	4		0.2433	0.2446	0.231	0.2296	0.2294				
Final constant weight	0.2234	0.219	0.223	0.2294	0.2485	5		0.2626	0.2644	0.2504	0.2489	0.2485				
Floss Weight	0.0102	0.0094	0.0106	0.0095	0.0103	3	Wt-floss	1.1631	1.1728	1.1014	1.0949	1.0933				
W2 (final floss)	0.2132	0.2096	0.2124	0.2199	0.2382	M	Mean	0.23262	0.23456	0.22028	0.21898	0.21866				
						ps		0.011866	0.0116513	0.011866 0.0116513 0.0118647 0.0117259 0.0116961	0.0117259	0.0116961				
% SOLUBILITY	6.41	5.84	6.27	5.95	5.59			Mean solubility	billity	6.01	0.33					
3mm Group	,		,			•		0 2004	2020	0 2000	2010	0070	20100			
Sample	-	7	3	4	C			0.3681	0.3681 0.3737	0.3523	0.3497	0.3489	0.3487			
Start weight	0.3681	0.3683	0.3673	0.3438	0.3677	0.3677 Floss total 2		0.3683	0.3683 0.3/43		0.3516	0.3508	0.3507			
Floss Weight	9600.0	0.0097	0.0097	0.0113	0.0107	0.0513		0.3673	0.3673 0.3733		0.3503	0.3495	0.349			
W1 (start-floss)	0.3585	0.3586	0.3576	0.3325	0.357	4		0.3438	0.3438 0.3486		0.3264	0.3259	0.3256			
Final constant weight	0.3487	0.3507	0.349	0.3256	0.345	5		0.3677	0.3677 0.3724	0.3489	0.3461	0.34	0.345			
Floss Weight	9600.0	0.0097	0.0097	0.0113	0.0107	S	Wt - floss	1.7642	1.7913		1.6731					
W2 (final-floss)	0.3391	0.341	0.3393	0.3143	0.3343	2	Mean	0.35284	0.35826	0.33704	0.33462	0.3339	0.3336			
						ps		0.010762	0.0111235	0.010762 0.0111235 0.0106348	0.010497	0.010497 0.0103685	0.0103844			
% SOLUBILITY	5.41	4.91	5.12	5.47	6.36			Mean solubility	billity	5.45	0.55				-	and the same of th
		ı	ı													
4mm Group				The same of the sa												
Sample	-	2	3	4	2	-		0.4549	0.4549 0.4621	0.4378	0.4346	0.4332	0.4325	0.4309	0.4312	0.4311
Start weight	0.4549	0.4346	0.4625	0.4402	0.4501	0.4501 Floss total 2		0.4346	0.4346 0.4407	0.4184	0.4156	0.4149	0.4144	0.4138	0.4134	0.4132
Floss Weight	0.0095	0.0105	0.0104	0.0093	0.0091	0.0488 3		0.4625	0.4625 0.4708	0.4452	0.4422	0.4414	0.4408	0.4398	0.4396	0.4394
W1 (start-floss)	0.4454	0.4241	0.4521	0.4309	0.441	4		0.4402 0.449	0.449	0.4222	0.4194	0.4185	0.418	0.4174	0.4169	0.4167
Final constant weight	0.4311	0.4132	0.4394	0.4167	0.4261	5		0.4501	0.4	0.43	0.4288	0.4	0.4	0.4	0.4261	0.4261
Floss Weight	0.0095	0.0105	0.0104	0.0093	0.0091	S	Wt-floss	2.1935	- 1		2.0918					2.0777
W2 (final-floss)	0.4216	0.4027	0.429	0.4074	0.417	2	Mean	0.4387	0.44604	0.4213	0.41836	0.41742	0.41684	0.4159	0.41568	0.41554
						ps	_	0.011198	0.0116114	0.0110063	0.0108799	0.0107832	0.011198 0.0116114 0.0110063 0.0108799 0.0107832 0.0107208 0.0104503 0.0106303	0.0104503	0.0106303	0.0106473
% SOLUBILITY	5.34	5.05	5.11	5.45	5.44		-	Mean solubility	billity	5.28	0.19					
Same Grain																
Sample		2	~	4	2	-	-	0.5364	0 5364 0 5453	0 5195	0 5157	0 5147	0 5139	0 5127	0.5125	0.5122
Start weight	0 5364	0 5222	0 5309	0 5402	0.5579	0 5579 Floss total 2		0 5222	0.5222 0.5301	0.5014	0.4982	0.4971	0.4964	0.4953	0.4948	0.4948
Floss Weight	0 0093	0.0102	0 0087	0.0103	0.0092	0.0477 3		0.5309	0.5309 0.5388	0.5127	0.5091	0.5081	0.5074	0.5065	0.5059	0.5058
W1 (start-floss)	0.5271	0 512	0.5222	0.5299	0.5487	4		0.5402	0.5402 0.5464	0.5214	0.5181	0.517	0.5164	0.5154	0.5149	0.5148
Final constant weight	0 5122	0 4948	0 5058	0.5148	0.5212	9		0.5579	0.5579 0.5621	0.5295	0.5253	0.5242	0.523	0.5219	0.5214	0.5212
Floss Weight	0.0093	0.0102	0.0087	0.0103	0.0092		Wt-floss	2.6399	2.675	2.5368	2.5187	2.5134	1 2.5094	2.5041	2.5018	2.5011
W2 (final-floss)	0 5029	0.4846	0.4971	0.5045	0.512		Mean	0.52798	0.535	0.50736	0.50374	0.50268	3 0.50188	0.50082	0.50036	0.50022
						ps	-	0.013251		0.011767 0.0105316		0.0102238	0.010228 0.0102238 0.0100847 0.0100716 0.0100998 0.0100214	0.0100716	0.0100998	0.0100214

Table 88. Photac-Fil Quick One-Week Solubility Experimental Sample Results (grams)

O'Haus Analytical		O'Haus Analytical	Wical Plus	Plus Scale										
								pre H20	post H20					
nm Group										10-Dec-05	11-Dec-05	12-Dec-05	13-Dec-05	14-Dec-05
Sample	-	2	3	4	5		-	0.2246	0.2248	0.2116	0.2104	0.2102		
Start weight	0.2246	0.2365	0.2523	0.2331	0.2541	0.2541 Floss total	2	0.2365	0.237	0.2231	0.2223	0.222		
Floss Weight	0.0099	0.0113	0.0106	0.0106	0.0119	0.0543	3	0.2523	0.2566	0.2436	0.2426	0.2422		
W1 (start-floss)	0.2147	0.2252	0.2417	0.2225	0.2422		4	0.2331	0.2347	0.2196	0.2188	0.2184		
Final constant weight	0.2102	0.222	0.2422	0.2184	0.2385		2	0.2541	0.255	0.2395	0.2386	0.2385		
Floss Weight	0.0099	0.0113	0.0106	0.0106	0.0119		Wt-floss	1.1463	1.1538	1.0831	1.0784	1.077		
W2 (final-floss)	0.2003	0.2107	0.2316	0.2078	0.2266		Mean	0.22926			0.21568	0.2154		
							ps	0.0127186	0.0137438	0.01358113	0.01361756	0.01361756		
SOLUBILITY	6.71	6.44	4.18	6.61	6.44			Mean solubility	oility	6.07	1.07			
Cample		2	~	7	2		-	0.3519	0 3524	0 3319	0.3303	0 3299	0.3292	0.3294
Start weight	0 3519	0 3174	0 3438	0 3207	0 328	0.328 Floss total	2	0.3174		0.2989	0.2977	0.2972	0.2965	0.2967
Flose Weight	0 0111	0 0001	0.0101	0.011	0.0101	0.0514	3	0.3438	0.344	0.3242	0.3228		0.3215	0.3215
W1 (start-floss)	0.3408	0.3083	0.3337	0.3097	0.3179		4	0.3207	0.3218	0.3014	0.2998	0.2993	0.2986	0.2988
Final constant weight	0.3294	0.2967	0.3215	0.2988	0.3059		5	0.328	0.329	0.3087	0.307	0.3066	0.3059	0.3059
Floss Weight	0.0111	0.0091	0.0101	0.011	0.0101		Wt-floss	1.6104	1.6143	1.5137	1.5062	1.5038	1.5003	1,5009
W2 (final-floss)	0.3183	0.2876	0.3114	0.2878	0.2958		Mean	0.32208	0.32286	0.30274	0.30124	0.30076	0.30006	0.30018
							ps	0.0149286	0.0145684	0.01444393 0.01439364	0.01439364	0.01439906	0.01439906	0.01437613
SOLUBILITY	09.9	6.71	6.68	7.07	6.95			Mean solubility	billity	08.9	0.20			
mm Group														
Sample	-	2	3	4	2		-	0.4372	0.4378	0.4143	0.4122	0.4115	0.4105	0.4107
Start weight	0.4372	0.4091	0.4685	0.4119	0.4684	0.4684 Floss total	2	0.4091	0.4108	0.386	0.3841	0.3835	0.3826	0.3828
Floss Weight	0.0118	0.0124	0.0118	0.0117	0.0126	0.0603	3	0.4685	0.4709	0.4472	0.4445		0.4427	0.4427
W1 (start-floss)	0.4254	0.3967	0.4567	0.4002	0.4558		4	0.4119		0.3915	0.3891	0.3886	0.3876	0.3876
Final constant weight	0.4107	0.3828	0.4427	0.3876	0.4403		5	0.4684		0.4444	0.4419		0.4402	0.4403
Floss Weight	0.0118	0.0124	0.0118	0.0117	0.0126		Wt-floss	2.1348		2.0231	2.0115		2.0033	2.0038
W2 (final-floss)	0.3989	0.3704	0.4309	0.3759	0.4277		Mean	0.42696	800				0.40066	0.40076
							ps	0.0290098	0.0290053	0.02863926	0.02839574	0.02832467	0.02827149	0.02823875
% SOLUBILITY	6.23	6.63	59.65	6.07	6.16	0		Mean solubility	pility	6.15	0.35			
mm Grain														
Sample	-	2	3	4	5		-	0.5571	0.5588	0.5293	0.5264	0.5253	0.5239	0.5238
Start weight	0.5571	0.5784	0.5353	0.5461	0.5541	0.5541 Floss total	2		0.5802	0.5477	0.5445	0.5433	0.542	0.542
Floss Weight	0.0102	0.0122	0.0102	0.0111	0.0114	0.0551	3	0.5353	0.5377	0.5071	0.5045	0.5035	0.5021	0.5021
W1 (start-floss)	0.5469	0.5662	0.5251	0.535	0.5427		4	0.5461		0.5165			0.5108	0.5109
Final constant weight	0.5238	0.542	0.5021	0.5109	0.5208		5	0.5541	0.5551	0.5262	0.5231	0.5221	0.5209	0.5208
Floss Weight	0.0102	0.0122	0.0102	0.0111	0.0114		Wt-floss	2.7159	2.7235	2.5717	2.5567	2.5513	2.5446	2.5445
W2 (final-floss)	0.5136	0.5298	0.4919	0.4998	0.5094		Mean	0.54318	0.5447	0.51434	0.51134		0.50892	0.5089
							ps		0.0159172	0.0159443 0.0159172 0.01522524 0.01506745	0.01506745	0.0150024	0.01504204	0.01501889
W COLUBII ITY	000		000											

Table 89. Photac-Fil Quick One-Week Solubility Control Sample Results (grams)

FINAL ONE WEEK SULUDICITY (COMBO) GROOT	,				-	The second secon	001	0011							
		O'Haus Analytical Plus Scale	ytical Plus	Scale			pre H20	post HZO							
mm Group							4-Jan-06	5-Jan-06	6-Jan-06	7-Jan-06	8-Jan-06	9-Jan-06	10-Jan-06	11-Jan-06	12-Jan-06
Sample	-	2	3	4	5	-	0.2236	0.2255	0.2138	0.2126	0.2124	0.2123			
Start weight	0.2236	0.2463	0.2553	0.2339	0.2295	0.2295 Floss total 2	0.2463	0.248	0.2345	0.2331	0.2327	0.2326			
Floss Weight	0.0126	0.0111	0.0103	0.0099	9600.0	0.0535 3	0.2553	0.2558	0.2438	0.2415	0.2415	0.2413			
W1 (start-floss)	0.211	0.2352	0.245	0.224	0.2199	4	0.2339	0.2343	0.2227	0.2209	0.2205	0.2205			
Final constant weight	0.2123	0.2326	0.2413	0.2205	0.2151	5	0.2295	0.2291	0.2175	0.2157	0.2152	0.2151			
Floss Weight	0.0126	0.0111	0.0103	0.0099	9600.0	Wt-floss	1,1351	~	1.0788	1.0703	•	-			
W2 (final-floss)	0.1997	0.2215	0.231	0.2106	0.2055	Mean	0.22702	0.22784	0.21576	0.21406	0.21376	0.21366			
						ps	0.0128861	0.0128861 0.0128893		0.0121925	0.0122997	0.01245 0.0121925 0.0122997 0.0122572			
% SOLUBILITY	5.36	5.85	5.71	5.98	6.55		Mean solubility	billity	5.89	0.44					
mm Group															
Sample	-	2	3	4	5	-	0.3546	0.3572		0.33/4	0.3367	0.3366			
Start weight	0.3546	0.3305	0.321	0.33/1	0.3304	은	0.3305	0.3338			0.3145	0.3143		and the same decision in the same	
Floss Weight	0.0111	0.0093	0.0103	0.0092	9600.0	0.0495	0.321	0.3214			0.3023	0.3021			
W1 (start-floss)	0.3435	0.3212	0.3107	0.3279	0.3208		0.3371	0.3386		0.321	0.3206	0.3205			
Final constant weight	0.3366	0.3143	0.3021	0.3205	0.3104		0.3	0.3	0.3138	0.3112	0.3106	0.3			
Floss Weight	0.0111	0.0093	0.0103	0.0092	9600.0					1.5382	1.5352				
W2 (final-floss)	0.3255	0.305	0.2918	0.3113	0.3008	Mean	0.32482	0.32468	0.31024	0.30764	0.30704	0.30688			
						ps	0.0125059	0.0146878	0.0125059 0.0146878 0.0130335 0.0129073 0.0128873 0.0129328	0.0129073	0.0128873	0.0129328			
% SOLUBILITY	5.24	5.04	80.9	90.9	6.23		Mean solubility	billity	5.53	0.58					
4mm Group								0.55		****	*00*	0000		02070	0 4020
Sample	-	7	3	4	C	T TOTAL CONTROL OF THE PARTY OF	0.4544	0.45/8		0.431	0.4301	0.4298		0.4270	0.4279
Start weight	0.4544	0.4463	0.4434	0.4371	0.4633	Ë	0.4463	0.4481			0.4226	0.4224			0.4204
Floss Weight	0.01	0.0101	0.0102	0.0103	0.0106	0.0512 3	0.4434	0.4434			0.4198				0.4175
W1 (start-floss)	0.4444	0.4362	0.4332	0.4268	0.4527	4	0.4371	0.4372		0.4128	0.4121	0			0.41
Final constant weight	0.4279	0.4204	0.4175	0.41	0.4392		0.4633	0.4	0.4	0.4421	4.0	0.40		0.4392	0.4392
Floss Weight	0.01	0.0101	0.0102	0.0103	0.0106	Wt-floss	2.1933			2.0784			2.0659	2.0641	2.0638
W2 (final-floss)	0.4179	0.4103	0.4073	0.3997	0.4286	Mean	0.43866	0.44034	0.4197	0.41568	0.41494	0.4081	0.41318	0.41282	0.41276
						ps	0.0101/18	0.0116096	0.0111901	0.0111412	0.0110/64	0.0101718 0.0116096 0.0111901 0.0111412 0.0110764 0.0086283 0.01109739	0.01109/39	0.011031/	0.011031/ 0.011095/2
% SOLUBILITY	96.9	5.94	5.98	6.35	5.32		Mean solubility	pillity	5.91	0.37					
Smm Group											1.				
Sample	-	2	3	4	5	-	0.5461	0.5499	0.5266	0.5219	0.5205	0.5199		0.5179	0.5178
Start weight	0.5461	0.5513	0.5194	0.5498	0.5602	0.5602 Floss total 2	0.5513	0.5547	0.5295	0.525	0.5236		0.5217	0.5211	0.5211
Floss Weight	0.0128	0.0113	9600.0	0.0097	0.0127	0.0561 3	0.5194	0.522	0.4969	0.4922	0.4898	0.4896		0.4875	0.4875
W1 (start-floss)	0.5333	0.54	0.5098	0.5401	0.5475	4	0.5498	0.5539	0.5288	0.5246	0.5232	0.5229		0.5207	0.5206
Final constant weight	0.5178	0.5211	0.4875	0.5206	0.5273	9	0.5602	0.5619	0.5358	0.5313	0.5302	0.5298	0.5282	0.5274	0.5273
Floss Weight	0.0128	0.0113	9600.0	0.0097	0.0127	Wt-floss	2.6707						2.5219	2.5185	2.5182
W2 (final-floss)	0.505	0.5098	0.4779	0.5109	0.5146	Mean	0.53414	0.53726	0.5123	0.50778	0.50624	0.50588	0.50438	0.5037	0.50364
						ps	0.015408	0.0154221	0.0152688	0.0153729	0.0158688	0.0154085 0.0154221 0.0152688 0.0153729 0.0158688 0.0157913 0.01571798	0.01571798	0.0157169	0.0157169 0.01568321
VEN POLITICAL INC.	-	-	-	-											

Table 90. Photac-Fil Quick One-Month Solubility Experimental Sample Results (grams)

O'Haus Analytical Plus	0	O'Haus Analytical Plus		Scale				pre H20	post H20						
Durm Genun								7-Dec-05	8-Dec-05	9-Dec-05	10-Dec-05	11-Dec-05	12-Dec-05	13-Dec-05	14-Dec-05
Sample		2	3	4	5	-		S	0.2273	0.2129	0.2118				
Start weight	0.2266	0.2244	0.216	0.2258	0.2265	0.2265 Floss total 2		0.2244 0.254	0.254	0.2098	0.2089	0.2088			
Floss Weight	0.0101	0.0098	0.0104	0.0107	0.0102	0.0512 3		0.216	0.216 0.2167	0.2032	0.2021	0.202			
W1 (start-floss)	0.2165	0.2146	0.2056	0.2151	0.2163			0.2258 0.2264	0.2264	0.2122	0.2112	0.211			
Final constant weight	0 2117	0.2088	0.202	0.211	0.2092	9		0.2265 0.2247	0.2247	0.2103	0.2093	0.2092			
Floss Weight	0.0101	0.0098	0.0104	0.0107	0.0102	>	Wt-floss	1.0681	1.0979	0.9972	0.9921	0.9915			
Wy (final floce)	0 2016	0 199	0 1916	0 2003	0 199	-	Mean	0.21362	0.21958	0.19944	0.19842	0.1983			
feed and the same		3				65	ps	0.0044808	0.0141516	0.0038441	0.0141516 0.0038441 0.00386691 0.00385071	0.00385071			
% SOLUBILITY	6.88	7.27	6.81	6.88	8.00			Mean solubility	oillity	71.17	0.50				
3mm Group															
Sample	-	2	3	4	2			0.3437	0.3437 0.3445	0.3262	0.3247			0.3232	0.3232
Start weight	0.3437	0.3389	0.3444	0.336	0.3581	0.3581 Floss total 2		0.3389 0.3387	0.3387	0.3219	0.3204			0.319	0.319
Floss Weight	0.0091	0.0098	0.0083	0.0078	0.0081	0.04313		0.3444 0.3447	0.3447	0.3273	0.3257	0.3146	0.3248	0.3242	0.3224
W1 (start-floss)	0.3346	0.3291	0.3361	0.3282	0.35	4		0.336	0.336 0.3353	0.3165	0.3151	0.3379	0.3144	0.3137	0.3138
Final constant weight	0.3232	0.319	0.3224	0.3138	0.3371	5		0.3581 0.358	0.358	0.3403	0.3386	0.3396	0.3375	0.3369	0.3371
Floss Weight	0.0091	0.0098	0.0083	0.0078	0.0081	-	Wt . floss	1.678	1.6781	1.5891	1.5814	1.5932	1.5768	1.5739	1.5724
W/2 /final flocel	0 3141	0 3092	0.3141	0 306	0.329	-	Mean	0.3356	0.33562	0.31782	0.31628	0.31864	0.31536	0.31478	0.31448
						0,	ps	0.0084975	0.0086636	0.0088362		0.01105319	0.00872726 0.01105319 0.00859343 0.00860785	0.00860785	0.00865737
% SOLUBILITY	6.13	6.05	6.55	97.9	00.9			Mean solubility	oility	6.30	0.34				
4mm Group		1								0000					, , ,
Sample	-	2	3	4	2			0.4645	0.4645 0.4632	0.4379	0.4356	0.4346	0.4342	0.4334	0.4334
Start weight	0.4645	0.438	0.4342	0.428	0.426	E S		0.438	0.438 0.4373	0.413	0.4108	0.4099	0.4095	0.4086	0.4087
Floss Weight	9600.0	0.0098	0.0103	0.0108	0.00	0.0495 3	_	0.4342	0.4342 0.4331	0.4115	0.4091	0.4083	0.4079	0.4069	0.407
W1 (start-floss)	0.4549	0.4282	0.4239	0.4172	0.417	4		0.428	0.428 0.4263	0.4029	0.4007	0.3999	0.3994	0.3988	0.3988
Final constant weight	0.4334	0.4087	0.407	0.3988	0.3976			0.426	0.42	0.40	0.399	0.3985	0.3982	0.3975	0.3976
Floss Weight	9600.0	0.0098	0.0103	0.0108	0.009		Wt-floss	2.1412					1.9997	1.9957	1.996
W2 (final-floss)	0.4238	0.3989	0.3967	0.388	0.3886		Mean	0.42824	0.42676	0.40344	0.40118	0.40034	0.39994		
						o,	ps	0.0154967	0.0154967 0.0158162	0.014648		0.01461256 0.01450889	0.01450665	0.01445867	0.01443433
% SOLUBILITY	6.84	6.84	6.42	7.00	6.81			Mean solubility	oillity	6.78	0.22				
	ı		ı												
Sample		2	3	4	2			0.5384	0.5384 0.5362	0.5067	0.5036	0.5024	0.5022	0.5011	0.5015
Start weight	0 5384	0.5457	0.5343	0.5504	0.5619	0.5619 Floss total 2	01	0.5457	0.5457 0.5429	0.5169	0.514	0.513	0.5125	0.5116	0.5118
Floss Weight	0 0107	0.0105	0.0087	0.0085	0.0087	0.04713	_	0.5343	0.5343 0.5329	0.5107	0.5075	0.5065	0.5059	0.5049	0.5052
W1 (start-floss)	0.5277	0.5352	0.5256	0.5419	0.5532		4	0.5504	0.5504 0.5493	0.5207	0.5175	0.5164	0.5161	0.5152	0.5159
Final constant weight	0.5015	0.5118	0.5052	0.5159	0.5269	4,		0.5619 0.56	0.56	0.5328	0.5293	0.5282	0.5278	0.5268	0.5269
Floss Weight	0.0107	0.0105	0.0087	0.0085	0.0087		Wt-floss	2.6836	2.6742	2.5407	2.5248	2.5194	2.5174		
W2 (final-floss)	0.4908	0.5013	0.4965	0.5074	0.5182		Mean	0.53672	0.53484	0.50814	0.50496	0.50388	0.50348	0.5025	
							ps	0.0108029	0.0108029 0.0108283		0.010094 0.00995023 0.00995942	0.00995942	0.00994862	0.00998284	0.00991731
W COLLIDII ITV	000	6 33	777	C 37	CC 3			A. L. L. 1154.		6 24	0.83				

Table 91. Photac-Fil Quick One-Month Solubility Control Sample Results (grams)

THOUSE ONE MOUNT SOCIOUETT CONTINUE GIVES				-				0						
		O'Haus Anal	Analytical Plus Scale	Scale				pre HZU	DOST HZU					
								30 acl 5	A lan Of	Sh nel A	Si lan Oc	7. Inn 06	S lan Of	90 ncl 0
nm Group	1			1		1	-	Dani-No	4.0010	Double of	o orte	DOLLO O	O Date	STAILT
Sample	-	2	3	4	2			0.2276	0.2276	0.2156	0.2153	0.2146	0.2145	
Start weight	0.2276	0.2287	0.2245	0.2103	0.2578	0.2578 Floss total 2		0.2287	0.2287	0.2174	0.2169	0.2164	0.2162	
Floss Weight	0.0107	0.0117	0.0111	0.0109	0.0121	0.0565 3		0.2245	0.2243	0.2128	0.2123	0.2117	0.2115	
W1 (start-floss)	0.2169	0.217	0.2134	0.1994	0.2457	4		0.2103	0.2093	0.1984	0.1979	0.1975	0.1972	
Final constant weight	0.2145	0.2162	0.2115	0.1972	0.2412	5		0.2578	0.2574	0.2425	0.2418	0.2411	0.2412	
Floss Weight	0.0107	0.0117	0.0111	0.0109	0.0121		Wt-floss	1.0924	1.0908	1.0302	1.0277	1.0248	1.0241	
W2 (final floce)	0.2038	0 2045	0 2004	0 1863	0 2291	2	Mean	0.21848	0 21816	0.20604	0 20554	0.20496	0.20482	
(continue)	200					Ö	ps	0.0173089	12	_	0.0158486		_	
SOLUBILITY	6 04	97.9	60 9	6.57	6.76			Mean solubility	bility	6.24	0.41			
am Group														
Sample	-	2	3	4	2	-		0.3446	0.3446 0.3443	0.3279	0.3266	0.3255	0.3252	
Start weight	0 3446	0.3621	0 3467	0.3503	0.3595	0.3595 Floss total 2		0.3621	0.3621 0.3614	0.3448	0.3435	0.3426	0.3424	
Floss Weight	0 0109	0.0115	0.0101	0.0099	0.0106	0.053 3		0.3467 0.346	0.346	0.3305	0.3291	0.3282	0.3282	
W1 (start-floss)	0 3337	0 3506	0 3366	0 3404	0.3489	NAME OF TAXABLE PARTY AND TAXABLE PARTY.		0.3503 0.35	0.35	0.3327	0.3316	0.3306	0.3305	
Final constant weight	0 3252	0.3424	0.3282	0.3305	0.3396			0.3595	0.3595 0.3601	0.3423	0.3407	0.3396	0.3396	
Floss Weight	0.0109	0.0115	0.0101	0.0099	0.0106		Wt . floss	1.7102	1.7088	1.6252	1.6185	1.6135	1.6129	
W2 (final-floss)	0.3143	0.3309	0.3181	0.3206	0.329		Mean	0.34204	0.34176	0.32504	0.3237		0.32258	
							ps	0.0077774	0.0079469	0.0077774 0.0079469 0.0074705	0.007403	0.0074216	0.007448	
% SOLUBILITY	5.81	5.62	5.50	5.82	5.70			Mean solubility	billity	5.69	0.14			
im Group														
Sample	-	2	3	4	2	-		0.467	0.467 0.4686	0.4455	0.4431	0.4415	0.4412	0.4412
Start weight	0.467	0.4688	0.4723	0.4549	0.4722	Flo		0.4688	0.4688 0.4688	0.4466	0.4445	0.443	0.4425	0.4425
Floss Weight	0.0093	9600.0	0.0105	0.0107	0.0111	0.0512		0.4723	0.4723 0.4722	0.4496	0.4472	0.446	0.4459	0.4457
W1 (start-floss)	0.4577	0.4592	0.4618	0.4442	0.4611	4		0.4549	0.4549 0.4557	0.4337	0.432	0.4305	0.4303	0.4302
Final constant weight	0.4412	0.4425	0.4457	0.4302	0.4472			0.4722	4	0.4	0.4	0.4	0.44	4
Floss Weight	0.0093	9600.0	0.0105	0.0107	0.0111		Wt-floss	2.284						
N2 (final-floss)	0.4319	0.4329	0.4352	0.4195	0.4361		Mean	0.4568	0.45736	0.43512	ener-	0.43146	0.4312	0.43112
							ps	0.0071549	0.0069141	0.0071549 0.0069141 0.0069277	0.006	0.0066951	0.0066951 0.0066991	0.006686
% SOLUBILITY	5.64	5.73	92.9	99.9	5.42		S. S. S. S. S. S. S. S. S. S. S. S. S. S	Mean solubility	billity	5.62	0.14			
Ann Grains			ı	ľ										
Sample		2	3	4	5	-		0.5675	0.5675 0.5671	0.5436	0.5405	0.5388	0.5382	0.5381
Start weight	0.5675	0.5616	0.5478	0.5466	0.5714	0.5714 Floss total 2		0.5616	0.5616 0.5611	0.5362	0.5334	0.5318	0.5315	0.5313
Floss Weight	0.0103	9600.0	0.0087	0.0087	0.0099	0.0472 3		0.5478 0.547	0.547	0.5225	0.5199	0.5182	0.5179	0.5175
W1 (start-floss)	0.5572	0.552	0.5391	0.5379	0.5615	4		0.5466	0.5466 0.5465	0.5216	0.5191	0.5175	0.5169	0.5169
Final constant weight	0.5381	0.5313	0.5175	0.5169	0.541	9		0.5714	0.5714 0.571	0.5461	0.5436	0.5416	0.5412	0.541
Floss Weight	0.0103	9600.0	0.0087	0.0087	0.0099		Wt-floss	2.7477	2.7455		2.6093	3 2.6007	2.5985	2.5976
W2 (final-floss)	0.5278	0.5217	0.5088	0.5082	0.5311		Mean	0.54954	0.5491	0.52456	0.52186	5 0.52014	0.5197	0.51952
						S	. ps	0.0113134	0.011327	0.0113134 0.0113271 0.0115046	0.0113923	0.0112899	0.0112899 0.0112833 0.0112999	0.0112999
ATT HOLLING W	-												The same of the sa	

Table 92. Photac-Fil Quick Three-Month Solubility Experimental Sample Results (grams)

	0	O'Haus Analytical Plus Scale	O'Haus Analytical Plus Scale	Scale			pre H20	post H20					
nm Group							20-Oct-05	9	22-Oct-05	24-Oct-05	25-Oct-05	26-Oct-05	27-Oct-05
Sample	-	7	3	4	c		0.2315 0.2317		1617.0	0.2182	0.2163		
Start weight	0.2315	0.2251	0.223	0.2346	0.2328	0.2328 Floss total 2	0.2251 0.257		0.2147	0.2139	0.2138		
Floss Weight	0.0078	0.0082	0.0082	0.0093	0.0082	0.0417 3	0.223 0.2236		0.2125	0.2116	0.2115		
W1 (start-floss)	0.2237	0.2169	0.2148	0.2253	0.2246	4	0.2346 0.2349		0.2233	0.2226	0.2225		
Final constant weight	0.2183	0.2138	0.2115	0.2225	0.2208	5	0.2328 0.2335		0.2219	0.2208	0.2208		
Floss Weight	0.0078	0.0082	0.0082	0.0093	0.0082	Wt-floss	1.1053	1.139	1.0498	1.0454	1.0452		
W2 (final-flose)	0 2105	0 2056	0 2033	0.2132	0.2126	Mean	0.22106	0.2278	0.20996	5 0.20908	0.20904		
		No. of the last of				ps	0.00506113		0.012456 0.00461519	9 0.00461324	0.0046408		
SOLUBILITY	6 9	521	5.35	5.37	5.34	TO SHARE STREET, SALES	Mean solubility	lity	5.44	1 0.27			
100000													
m Group													
Sample		2	3	4	5	-	0.353 0.3522		0.3362	0.3347	0.3344		
Start weight	0.353	0.3303	0.3285	0.3323	0.3447 F	Floss total 2	0.3303 0.3299		0.315	0.3139	0.3138		
Floss Weight	0.0083	0.0079	0.0075	0.0073	0.0081	0.0391 3	0.3285 0.3289		0.3119	0.3107	0.3105		
W1 (start-floss)	0.3447	0.3224	0.321	0.325	0.3366	4	0.3323 0.332		0.3175	0.3164	0.3161		
Final constant weight	0.3344	0.3138	0.3105	0.3161	0.328	2	0.3447 0.3466		0.3295	0.3283	0.328		
Floss Weight	0 0083	0 0079	0.0075	0.0073	0.0081	Wt - floss	1.6497	1.6505	1.571	1.5649	1,5637		
W2 (final-floss)	0.3261	0.3059	0.303	0.3088	0.3199	Mean	0	0.3301	0.3142	Ĭ	Ĭ		
				THE OUT OF THE OUT OUT OF THE OUT OF THE OUT OUT OF THE OUT OUT OUT OUT OF THE OUT OUT OUT OUT OUT OUT OUT OUT OUT OUT	A STREET, STATE OF THE STREET, STATE OF THE STREET, STATE OF THE STREET, STATE OF THE STREET, STATE OF THE STREET, STATE OF THE STREET, STATE OF THE STREET, STATE OF THE STREET, STATE OF THE STREET, STATE OF THE STREET, STATE OF THE STATE	ps	0.01062582	0.010724	0.01035794	0.0102279	0.01016971		
% SOLUBILITY	5.40	5.12	5.61	4.98	4.96		Mean solubility	lity	5.21	1 0.28			
ım Group													
Sample	-	2	3	4	2	-	0.4536 0.4532	2	0.4293	0.4277	0.4272		0.427
Start weight	0.4536	0.4501	0.4402	0.4382	0.4396	0.4396 Floss total 2	0.4501 0.449	0.449	0.4268	0.4252	0.4248		0.4245
Floss Weight	0.0075	0.0074	6900.0	0.0074	0.0074	0.0366 3	0.4402 0.4413		0.4212	0.4195	0.4192		0.4188
W1 (start-floss)	0.4461	0.4427	0.4333	0.4308	0.4322	4	0.4382 0.438		0.4159	0.4142	0.4136		0.4134
Final constant weight	0.427	0.4245	0.4188	0.4134	0.4164	2	0.4396 0.4394	0.4	0.4188	0.417	0.416	0.4163	0.4164
Floss Weight	0.0075	0.0074	6900.0	0.0074	0.0074	Wt-floss	2.1851	- 1		-			2.0635
W2 (final-floss)	0.4195	0.4171	0.4119	0.406	0.409	Mean	0.43702	0.43686		90.0	MY A	DUES.	0.4127
						ps s	0.00700414	0.006594	0.00555923	0.0056	0.00563294	0.00565553	0.00563933
% SOLUBILITY	96.5	5.78	4.94	5.76	5.37	Section of the second	Mean solubility	llity	5.56	6 0.41			
Cample		6	3	V	4	-	0 5608 0 5601	0 5601	0 5339	0 5318	0 5314	0 5311	0 5311
Start weight	0 5608	0.5569	0.5606	0.5547	0.5542	0.5542 Floss total 2	0.5569 0.5551		0.5299	0.5277	0.5273		0.527
Floss Weight	69000	0.0073	0.0075	0.0077	0.0073	0.0367 3	0.5606 0.562		0.5344	0.5319	0.5315	0.5312	0.531
W1 (start-floss)	0.5539	0.5496	0.5531	0.547	0.5469	4	0.5547 0.5546		0.5268	0.5247	0.5241	0.5238	0.5236
Final constant weight	0.5311	0.527	0.531	0.5236	0.5244	5	0.5542 0.554	0.554	0.5277	0.5253	0.5249	0.5246	0.5244
Floss Weight	6900.0	0.0073	0.0075	0.0077	0.0073	Wt-floss	2.7505		2.616		2.6025	2.601	2.6004
W2 (final-floss)	0.5242	0.5197	0.5235	0.5159	0.5171	Mean	0.5501	0.54982	0.5232	2 0.52094	0.5205	0.5202	0.52008
						ps	0.00314531		0.003635 0.00348755	5 0.00344703	0.00349971	0.00349971 0.00354429	0.0035442

Table 93. Photac-Fil Quick Three Month Solubility Control Sample Results (grams)

	0	Hans Analy	O'Haus Analytical Plus Scale	Scale			pre H20	post H20				
mm Group							2-Nov-05	3-Nov-05	4-Nov-05	7-Nov-05	8-Nov-05	9-Nov-05
Sample	-	2	3	4	5	-	0.2297 0.229	0.229	0.2185	0.2178	0.2178	
Start weight	0.2297	0.2266	0.2326	0.2339	0.2357	0.2357 Floss total 2	0.2266 0.2274	0.2274	0.2166	0.2157	0.2157	
Floss Weight	600.0	0.0086	0.0084	0.0077	0.0083	0.042 3	0.2326 0.2333		0.222	0.2212	0.2213	
W1 (start-floss)	0.2207	0.218	0.2242	0.2262	0.2274	4	0.2339 0.2352	0.2352	0.2221	0.2212	0.2212	
Final constant weight	0.2178	0.2157	0.2213	0.2212	0.2237	9	0.2357 0.2359		0.2245	0.2237	0.2237	
Floss Weight	0.009	0.0086	0.0084	0.0077	0.0083	Wt-floss	1.1165	1.1188	1.0617	1.0576	1.0577	
W2 (final-floss)	0.2088	0.2071	0.2129	0.2135	0.2154	Mean	0.2233	0.22376		0.21152	MARK	
						ps	0.0035937	0.0037806	0.0031517	0.0031571	0.0031675	
% SOLUBILITY	5.39	9.00	5.04	5.61	5.28	or services and the	Mean solubility	ility	5.26	0.25		
	ı											
Smm Group		٠	•		4		0 2642 0 265	1 266	0 2525	0.2542	0.0543	
Sample	0 2642	0 207	0 3248	0 2020	O 3506	0 3506 Floor total 2	0 307 0 3000	1 3292	0.2323	0.3116	0.2313	
Start weight	0.0042	0.000	0.000	0.077	0.000	0.039 01812	0 3318 0 3322		0.3159	0.3145	0 3145	
M/1 (ctart floes)	0.0033	0.0001	0 3236	0 3155	0 344	4	0 3232 0 324		0 3108	0 3094	0.3093	
Final constant weight	0.2513	0.3114	0.3145	0.3093	0.3348	2	0.3506 0.3527	7	0.3362	0.3349	0.3348	
Floss Weight	0.0093	0.0081	0.0082	0.0077	0.0066	Wt - floss		1.5632	1.489	1.4818	1,4814	
W2 (final-floss)	0.242	0.3033	0.3063	0.3016	0.3282	Mean	0.31138	0.31264	0.2978	0.29636	0.29628	
						ps	0.0325823	0.0325823 0.0329464	0.0314299	0.0313396	0.0312997	
% SOLUBILITY	90.9	4.89	5.35	4.41	4.59		Mean solubility	illity	4.86	0.37		
Imm Groun		I										
Sample	-	2	3	4	5	-	0.4623 0.4623	0.4623	0.4418	0.4396	0.4393	0.439
Start weight	0.4623	0.4256	0.4255	0.435	0.438	Floss total 2	0.4256 0.426	0.426	0.4065	0.4048	0.4046	0.4045
Floss Weight	0.0087	9.0000	0.0071	0.0077	0.0075	0.0075 0.0386 3	0.4255 0.4259	0.4259	0.4072	0.4051	0.4049	0.4046
W1 (start-floss)	0.4536	0.418	0.4184	0.4273	0.4305	4	0.435	0.435 0.4344	0.4141	0.4121	0.4118	0.4116
Final constant weight	0.439	0.4045	0.4046	0.4116	0.4167	5	0.438	0.438 0.4389	0.4191	0.4171	0.4169	0.4167
Floss Weight	0.0087	9.0000	0.0071	0.0077	0.0075	Wt-floss	2.1478	2.1489	2.0501		2.0389	2.0378
W2 (final-floss)	0.4303	0.3969	0.3975	0.4039	0.4092		0.42956		MOR_	0.40802		0.40756
						DS .	0.010070	0.0143451	0.014		0.0142303	0.014210
% SOLUBILITY	5.14	5.05	2.00	5.48	4.95	I	Mean solubility	illity	21.6	0.21		
Smm Group												
Sample	-	2	3	4	2		0.5636 0.567	0.567	0.5413	0.5379		0.5376
Start weight	0.5636	0.5327	0.4998	0.5031	0.5183	E	0.5327 0.5336	0.5336	0.5093	9905.0		0.5062
Floss Weight	0.0065	0.0077	0.0078	0.0079	0.0078	0.0377	0.4998 0.5042	0.5042	0.4765	0.474		0.4736
W1 (start-floss)	0.5571	0.525	0.492	0.4952	0.5105		0.5031 0.5025	0.5025	0.4828	0.4804		0.48
Final constant weight	0.5376	0.5062	0.4736	0.48	0.489		0.5183 0.516	0.516	0.49	0.48	0.4892	0.489
Floss Weight	0.0065	0.0077	0.0078	0.0079	0.0078		2.5798	2.5856				2.4487
W2 (final-floss)	0.5311	0.4985	0.4658	0.4721	0.4812		0.51596	0.51712	0.49284	800	100	0.48974
						ps	0.0259689	0.0259689 0.0267267	0.0259978	0.0256342	0.0256741	0.0256572
% SOI HRII ITY	A K7	505	6 33	J V CC	E 7.4		Moon on hilling	ility.	5 00	JV V		

APPENDIX B HARDNESS TESTING DATA

Table 94. Vitremer 24-Hour Hardness (KHN)

2mm		建物、建筑						
2,11111	1	2	3	4	5		VIII. (1)	
	29.8	24.8	22.9	24.7	25.0		No.	
top	29.1	24.5	24.3	25.2	25.9	avg	24.5	
	27.2	23.3	25.5	22.4	25.6	st dev	0.4	_
avg	28.7	24.2	24.2	24.1	25.5			
st dev	1.3	0.8	1.3	1.5	0.5			% bottom/top 94.28766
								34.20100
	1	2	3	4	5			
	25.5	21.3	22.1	24.1	23.1			1
oottom	24.7	23.0	22.8	23.9	24.9	avg	23.1	
	25.5	20.8	23.5	22.4	25.4	st dev	0.3	
avg	25.2	21.7	22.8	23.5	24.5			
st dev	0.5	1.2	0.7	0.9	1.2			
3mm	d elle	STATE OF THE PARTY		NAME OF			P 15 19 29	
	1	2	3	4	5			
	25.2	26.4	26.4	24.8	24.6			
top	28.0	24.6	24.0	24.6	23.6	avg	24.9	
	26.4	26.0	24.4	26.2	22.9	st dev	0.3	
avg	26.5	25.7	24.9	25.2	23.7			
st dev	1.4	0.9	1.3	0.9	0.9			% bottom/top
								91.12228
	1	2	3	4	5			
	23.7	23.8	22.0	24.3	22.0			
ottom	22.9	21.4	24.0	21.9	20.0	avg	22.7	
	23.6	23.6	23.3	23.0	22.7	st dev	0.4	
avg	23.4	22.9	23.1	23.1	21.6			
st dev	0.4	1.3	1.0	1.2	1.4			
								-
4mm			MUT I					
	1	2	3	4	5			
	26.9	23.7	25.4	25.8	22.9		A CONTRACTOR OF THE PARTY OF TH	
top	25.3	23.0	23.8	23.9	22.2	avg	23.9	
	26.0	24.8	25.5	22.8	22.4	st dev	0.4	
avg	26.1	23.8	24.9	24.2	22.5			0/ 1-1/
st dev	0.8	0.9	1.0	1.5	0.4			% bottom/top
								85.884
	1	2	3	4	5			
	21.8	24.2	21.6	19.1	17.9	1		
oottom	20.3	21.8	19.4	22.7	18.2	avg	20.5	
Cuom	20.7	20.8	20.3	22.1	17.7	st dev	0.7	
avg	20.9	22.3	20.4	21.3	17.9			
st dev	0.8	1.7	1.1	1.9	0.3			
5mm								
	1	2	3	4	5			
	25.8	22.9	26.2	23.2	23.9			
top	22.8	23.0	26.0	23.4	22.1	avg	24.1	
	25.9	24.1	25.2	24.2	24.7	st dev	0.6	
avg	24.8	23.3	25.8	23.6	23.6			0/ 1 - 1/ 5
st dev	1.8	0.7	0.5	0.5	1.3			% bottom/to 70.88958
								10.00330
	1	2	3	4	5			
	20.7	14.8	20.6	13.6	17.7			
bottom	17.5	13.3	21.4	13.9	19.4	avg	17.1	
	17.6	16.4	19.8 20.6	14.0	19.9 19.0	st dev	0.6	
	18.6	14.8						

Table 95. Vitremer One-Week Hardness (KHN)

	Week Ha	nuness						
mm	94.6					I SAME TO A		
	1	2	3	4	5			
	22.9	22.6	23.0	22.4	23.1			
ор	22.7	22.5	21.0	22.7	21.9	avg	22.4	
-	22.3	23.3	22.4	22.0	22.4	st dev	0.3	
avg	22.6	22.8	22.1	22.4	22.5			
dev	0.3	0.4	1.0	0.4	0.6			% bottom/top
dev	0.0	0.4	1.0	0.4	0.0			88.8
	1	2	3	4	5			
-	1 20.0	19.4	19.9	19.1	19.4	1		-
+	20.0			19.1	19.4		19.9	
ttom	19.9	20.0	20.5	20.0		avg	0.2	
	20.1	20.7	20.0		20.7	st dev	0.2	
dev	20.0 0.1	20.0 0.7	20.1	19.6 0.5	19.9			
dev	0.1	0.1	0.0	0.0				
nm		VEIDEN SEUT						
A CHARLES	1	2	3	4	5			
	22.9	21.4	22.8	21.6	22.4			
ор	22.4	22.2	21.6	22.1	23.4	avg	22.4	
-P	22.7	22.5	23.8	23.1	22.3	st dev	0.3	
avg	22.7	22.0	22.7	22.3	22.7	0. 001	0.0	
dev	0.3	0.6	1.1	0.8	0.6			% bottom/top
asv	0.5	0.0	I.I	0.0	0.0			80.1
							-	00.1
	1	2	3	4	5			
-	18.7	18.4	18.7	18.6	18.5			
ttom		16.9	19.1	17.7	17.8	avg	18.0	
rom	17.3 17.6	16.9	17.8	17.7	17.0	st dev	0.1	
				18.1	17.2	St dev	0.1	
vg	17.9	17.4 0.9	18.5	0.5	0.7			
dev	0.7	0.9	0.7	0.5	0.7			
nm	1	2	3	4	5			ALL SECTION AND ADDRESS.
-		22.3	22.5	22.0	22.2			
<u> </u>	22.0		21.8	22.0	22.7	avg	22.1	
ор	22.9	22.2		22.1	21.5	st dev	0.2	
	21.8	21.0	22.1			St dev	0.2	
avg	22.2	21.8	22.1	22.2	22.1			% hattam#a
dev	0.6	0.7	0.4	0.2	0.6	-		% bottom/top 78.4
-								76.4
		2	3	4	5			
			J					
	17 N		18.3	17.8	18.0			
ttom	17.0	17.2	18.3	17.8 16.9	18.0	avo	17.3	
tom	17.0 16.5	17.2 16.9	17.6	16.9	17.3	avg	17.3	
	17.0 16.5 16.5	17.2 16.9 16.8	17.6 17.0	16.9 16.7	17.3 17.1	avg st dev	17.3 0.2	
/g	17.0 16.5 16.5 16.7	17.2 16.9 16.8 17.0	17.6 17.0 17.6	16.9 16.7 17.1	17.3 17.1 17.5			
/g	17.0 16.5 16.5	17.2 16.9 16.8	17.6 17.0	16.9 16.7	17.3 17.1			
vg dev	17.0 16.5 16.5 16.7	17.2 16.9 16.8 17.0	17.6 17.0 17.6	16.9 16.7 17.1	17.3 17.1 17.5			
vg dev	17.0 16.5 16.5 16.7 0.3	17.2 16.9 16.8 17.0 0.2	17.6 17.0 17.6 0.7	16.9 16.7 17.1 0.6	17.3 17.1 17.5 0.5			
⁄g dev	17.0 16.5 16.5 16.7 0.3	17.2 16.9 16.8 17.0 0.2	17.6 17.0 17.6 0.7	16.9 16.7 17.1 0.6	17.3 17.1 17.5 0.5			
dev	17.0 16.5 16.5 16.7 0.3	17.2 16.9 16.8 17.0 0.2	17.6 17.0 17.6 0.7	16.9 16.7 17.1 0.6	17.3 17.1 17.5 0.5	st dev	0.2	
vg dev	17.0 16.5 16.5 16.7 0.3	17.2 16.9 16.8 17.0 0.2 2 22.3 23.5	17.6 17.0 17.6 0.7 3 22.7 23.7	16.9 16.7 17.1 0.6	17.3 17.1 17.5 0.5 5 23.1 22.1	st dev	0.2	
vg dev nm op	17.0 16.5 16.5 16.7 0.3 1 22.5 22.9 22.4	17.2 16.9 16.8 17.0 0.2 2 22.3 23.5 22.1	17.6 17.0 17.6 0.7 3 22.7 23.7 22.4	16.9 16.7 17.1 0.6 4 22.8 22.2 23.5	17.3 17.1 17.5 0.5 5 23.1 22.1 22.6	st dev	0.2	
vg dev	17.0 16.5 16.5 16.7 0.3 1 22.5 22.9 22.4 22.6	17.2 16.9 16.8 17.0 0.2 2 22.3 23.5 22.1 22.6	17.6 17.0 17.6 0.7 3 22.7 23.7 22.4 22.9	16.9 16.7 17.1 0.6 4 22.8 22.2 23.5 22.8	17.3 17.1 17.5 0.5 5 23.1 22.1 22.6 22.6	st dev	22.8 0.2	% hattam/ta
nm op	17.0 16.5 16.5 16.7 0.3 1 22.5 22.9 22.4	17.2 16.9 16.8 17.0 0.2 2 22.3 23.5 22.1	17.6 17.0 17.6 0.7 3 22.7 23.7 22.4	16.9 16.7 17.1 0.6 4 22.8 22.2 23.5	17.3 17.1 17.5 0.5 5 23.1 22.1 22.6	st dev	22.8 0.2	% bottom/to 67.6
nm op	17.0 16.5 16.5 16.7 0.3 1 22.5 22.9 22.4 22.6 0.3	17.2 16.9 16.8 17.0 0.2 2 22.3 23.5 22.1 22.6 0.8	17.6 17.0 17.6 0.7 3 22.7 23.7 22.4 22.9 0.7	16.9 16.7 17.1 0.6 4 22.8 22.2 23.5 22.8 0.7	17.3 17.1 17.5 0.5 5 23.1 22.1 22.6 0.5	st dev	22.8 0.2	
avg dev	17.0 16.5 16.5 16.5 16.7 0.3 1 22.5 22.9 22.4 22.6 0.3	17.2 16.9 16.8 17.0 0.2 2 22.3 23.5 22.1 22.6 0.8	3 22.7 23.7 22.4 22.9 0.7	16.9 16.7 17.1 0.6 4 22.8 22.2 23.5 22.8 0.7	17.3 17.1 17.5 0.5 5 23.1 22.1 22.6 22.6 0.5	st dev	22.8 0.2	
mm top	17.0 16.5 16.5 16.5 16.7 0.3 1 22.5 22.9 22.4 22.6 0.3	17.2 16.9 16.8 17.0 0.2 2 22.3 23.5 22.1 22.6 0.8	3 22.7 23.7 22.4 22.9 0.7	16.9 16.7 17.1 0.6 4 22.8 22.2 23.5 22.8 0.7	17.3 17.1 17.5 0.5 5 23.1 22.1 22.6 22.6 0.5	st dev	22.8 0.2	
avg dev mm	17.0 16.5 16.5 16.5 16.7 0.3 1 22.5 22.9 22.4 22.6 0.3	17.2 16.9 16.8 17.0 0.2 2 22.3 23.5 22.1 22.6 0.8	3 22.7 23.7 22.4 22.9 0.7	16.9 16.7 17.1 0.6 4 22.8 22.2 23.5 22.8 0.7	17.3 17.1 17.5 0.5 5 23.1 22.1 22.6 22.6 0.5	avg st dev	22.8 0.2	
nm op dev	17.0 16.5 16.5 16.5 16.7 0.3 1 22.5 22.9 22.4 22.6 0.3	17.2 16.9 16.8 17.0 0.2 2 22.3 23.5 22.1 22.6 0.8	3 22.7 23.7 22.4 22.9 0.7	16.9 16.7 17.1 0.6 4 22.8 22.2 23.5 22.8 0.7	17.3 17.1 17.5 0.5 5 23.1 22.1 22.6 22.6 0.5	st dev	22.8 0.2	

Table 96. Vitremer One-Month Hardness (KHN)

mm								
	1	2	3	4	5			
	25.5	28.0	27.3	26.7	27.2			
top	25.8	25.8	26.2	27.1	26.3	avg	27.1	
	25.4	28.3	26.1	28.2	28.2	st dev	0.4	
avg	25.6	27.4	26.5	27.3	27.2			
t dev	0.2	1.4	0.7	0.8	1.0	-		% bottom/to
	V.2		0.1	0.0				89.92
	1	2	3	4	5			
	23.8	25.3	25.4	24.6	25.5			
	22.1	24.2	23.1	25.3	23.3		24.4	
ottom						avg	Name and Address of the Owner, where	
	22.7	24.4	23.0	25.1	23.4	st dev	0.4	
avg	22.9	24.6	23.8	25.0	24.1			
t dev	0.9	0.6	1.4	0.4	1.2			
mm	1	2	3	4	5			
	28.0	27.9	25.1	27.6	27.0			-
ton +	27.1	26.9	27.4	26.1	27.3	avg	26.7	
top							THE RESERVE OF THE PARTY OF THE	
	26.3	26.6	26.8	27.0	24.7	st dev	0.3	
avg	27.1	27.1	26.4	26.9	26.3			V b - 11 - 1
t dev	0.9	0.7	1.2	0.8	1.4			% bottom/to
								89.45
	1	2	3	4	5			
	24.4	24.3	22.8	23.8	25.9			
ttom	24.3	25.1	23.1	22.6	23.7	avg	23.9	
LOIN	24.5	23.5	22.5	25.5	23.8	st dev	0.6	
	24.5	24.3	22.8	24.0	24.5	St dev	0.0	
avg			0.3	1.5	1.2			
dev	0.1	0.8	0.3	1.5	1.2			
							100000000000000000000000000000000000000	
mm	1	2	3	4	5			
	27.8	27.7	28.4	27.8	26.4			
100			28.8	28.2	26.1	2000	27.7	
ор	27.1	28.1				avg et dev	0.4	
	28.2	26.8	28.1	27.6	28.5	st dev	0.4	
avg	27.7	27.5	28.4	27.9	27.0			V h-44 - #
t dev	0.6	0.7	0.4	0.3	1.3			% bottom/to
								77.26
	1	2	3	4	5			
	21.3	21.2	21.4	22.3	22.8			
ttom	22.3	21.3	22.7	21.1	21.2	avg	21.4	
	20.4	19.5	21.3	21.3	20.8	st dev	0.2	
avg	21.3	20.7	21.8	21.6	21.6			
dev	1.0	1.0	0.8	0.6	1.1			
mm								
	1	2	3	4	5	-		
	28.5	28.8	27.6	27.7	27.9			
top	27.8	28.5	26.1	28.1	27.3	avg	27.8	
	27.4	26.2	28.3	28.1	28.8	st dev	0.5	
avg	27.9	27.8	27.3	28.0	28.0			
dev	0.6	1.4	1.1	0.2	0.8			% bottom/to
						-		66.50
	1	2	3	4	5			
	17.4	17.9	18.3	19.2	17.5			
ottom	19.2	19.2	18.1	18.6	18.8	avg	18.5	
	18.2	18.9	19.0	18.4	17.8	st dev	0.2	
		18.7	18.5	18.7	18.0			
avg	18.3							

Table 97. Vitremer Three-Month Hardness (KHN)

-	-			and the state of t	Marie de marie	NAME OF TAXABLE PARTY.	STORY WATER	
mm	MARK AND						19 Feb. 52	
	1	2	3	4	5			
	27.5	26.9	29.5	25.6	27.4		07.4	
top	28.0	26.8	28.7	26.9	26.0	avg	27.1	-
12002	25.6	27.1	27.0	27.2	25.7	st dev	0.5	-
avg	27.0	26.9	28.4	26.6	26.4			~
t dev	1.3	0.2	1.3	0.9	0.9			% bottom/to
								91.81034
	1	. 2	3	4	5			-
	28.0	24.0	24.6	25.4	26.9			
ottom	25.9	26.8	22.6	25.4	23.3	avg	24.9	B
	26.6	25.3	23.0	25.2	25.7	st dev	0.6	
avg	26.8	25.4	23.4	25.3	25.3			
t dev	1.1	1.4	1.1	0.1	1.8			
mm						1		
	1	2	3	4	5			
	26.8	22.4	22.4	24.7	23.4			
ор	26.9	25.0	22.8	22.8	23.3	avg	23.3	
	28.5	26.2	21.3	23.4	21.7	st dev	0.5	
avg	27.4	24.5	22.2	23.6	22.8	5. 55.	0.0	
t dev	1.0	1.9	0.8	1.0	1.0			% bottom/to
uev	1.0	1.3	0.0	1.0	1.0			85.89835
								00.0000
	4	2	3	4	5			+
	22.7		19.2	19.3	19.1			+
		18.5					20.0	
ottom	24.5	21.0	20.2	20.3	19.7	avg		
	25.0	21.4	21.6	20.5	19.2	st dev	0.5	_
avg	24.1	20.3	20.3	20.0	19.3			-
t dev	1.2	1.6	1.2	0.6	0.3			
							THE RESERVE OF THE PERSON NAMED IN	
mm								
	1	2	3	4	5			
	23.6	15.9	25.1	20.3	26.1			
top	22.6	16.1	24.3	20.7	25.7	avg	21.7	
	22.3	17.8	23.3	21.6	23.1	st dev	0.4	
avg	22.8	16.6	24.2	20.9	25.0			
st dev	0.7	1.0	0.9	0.7	1.6			% bottom/to
								76.57692
					-			
	1	2	3	4				
	20.9	13.5	3 18.6	17.6	16.0			
ottom	20.9	13.5	18.6	17.6	16.0	avo	16.6	
ottom	20.9 19.9	13.5 14.0	18.6 17.5	17.6 17.3	16.0 16.4	avg st dev	16.6	
	20.9 19.9 20.7	13.5 14.0 15.3	18.6 17.5 18.1	17.6 17.3 18.0	16.0 16.4 16.8	avg st dev	16.6	
avg	20.9 19.9 20.7 20.5	13.5 14.0 15.3 14.3	18.6 17.5 18.1 18.1	17.6 17.3 18.0 17.6	16.0 16.4 16.8 16.4		-	
avg	20.9 19.9 20.7	13.5 14.0 15.3	18.6 17.5 18.1	17.6 17.3 18.0	16.0 16.4 16.8		-	
avg	20.9 19.9 20.7 20.5	13.5 14.0 15.3 14.3	18.6 17.5 18.1 18.1	17.6 17.3 18.0 17.6	16.0 16.4 16.8 16.4	st dev	-	
avg t dev	20.9 19.9 20.7 20.5	13.5 14.0 15.3 14.3	18.6 17.5 18.1 18.1	17.6 17.3 18.0 17.6	16.0 16.4 16.8 16.4		-	
avg st dev	20.9 19.9 20.7 20.5 0.5	13.5 14.0 15.3 14.3 0.9	18.6 17.5 18.1 18.1 0.6	17.6 17.3 18.0 17.6 0.4	16.0 16.4 16.8 16.4 0.4	st dev	-	
avg st dev	20.9 19.9 20.7 20.5 0.5	13.5 14.0 15.3 14.3 0.9	18.6 17.5 18.1 18.1 0.6	17.6 17.3 18.0 17.6 0.4	16.0 16.4 16.8 16.4 0.4	st dev	-	
avg t dev	20.9 19.9 20.7 20.5 0.5	13.5 14.0 15.3 14.3 0.9	18.6 17.5 18.1 18.1 0.6	17.6 17.3 18.0 17.6 0.4	16.0 16.4 16.8 16.4 0.4	st dev	0.2	
avg t dev	20.9 19.9 20.7 20.5 0.5	13.5 14.0 15.3 14.3 0.9	18.6 17.5 18.1 18.1 0.6	17.6 17.3 18.0 17.6 0.4	16.0 16.4 16.8 16.4 0.4	st dev	24.5	
avg t dev	20.9 19.9 20.7 20.5 0.5	13.5 14.0 15.3 14.3 0.9	18.6 17.5 18.1 18.1 0.6	17.6 17.3 18.0 17.6 0.4	16.0 16.4 16.8 16.4 0.4	st dev	0.2	
avg t dev	20.9 19.9 20.7 20.5 0.5 1 26.0 24.7 25.4	13.5 14.0 15.3 14.3 0.9 2 25.0 24.9 25.1	18.6 17.5 18.1 18.1 0.6 3 25.2 25.0 27.1	17.6 17.3 18.0 17.6 0.4 4 22.8 24.1 24.0	16.0 16.4 16.8 16.4 0.4 5 24.0 23.1 23.3	st dev	24.5	
avg st dev	20.9 19.9 20.7 20.5 0.5 1 26.0 24.7 25.4 25.4	13.5 14.0 15.3 14.3 0.9 2 25.0 24.9 25.1 25.0	18.6 17.5 18.1 18.1 0.6 3 25.2 25.0 27.1 25.8	17.6 17.3 18.0 17.6 0.4 4 22.8 24.1 24.0 23.6	16.0 16.4 16.8 16.4 0.4 5 24.0 23.1 23.3 23.5	st dev	24.5	% bottom/to
avg st dev 5mm top avg	20.9 19.9 20.7 20.5 0.5 1 26.0 24.7 25.4	13.5 14.0 15.3 14.3 0.9 2 25.0 24.9 25.1	18.6 17.5 18.1 18.1 0.6 3 25.2 25.0 27.1	17.6 17.3 18.0 17.6 0.4 4 22.8 24.1 24.0	16.0 16.4 16.8 16.4 0.4 5 24.0 23.1 23.3	st dev	24.5	
avg st dev	20.9 19.9 20.7 20.5 0.5 1 26.0 24.7 25.4 25.4	13.5 14.0 15.3 14.3 0.9 2 25.0 24.9 25.1 25.0	18.6 17.5 18.1 18.1 0.6 3 25.2 25.0 27.1 25.8	17.6 17.3 18.0 17.6 0.4 4 22.8 24.1 24.0 23.6	16.0 16.4 16.8 16.4 0.4 5 24.0 23.1 23.3 23.5	st dev	24.5	% bottom/to 68.46049
avg st dev 5mm top avg	20.9 19.9 20.7 20.5 0.5 1 26.0 24.7 25.4 25.4	13.5 14.0 15.3 14.3 0.9 2 25.0 24.9 25.1 25.0 0.1	18.6 17.5 18.1 18.1 0.6 3 25.2 25.0 27.1 25.8 1.2	17.6 17.3 18.0 17.6 0.4 4 22.8 24.1 24.0 23.6 0.7	16.0 16.4 16.8 16.4 0.4 5 24.0 23.1 23.3 23.5 0.5	st dev	24.5	
avg st dev 5mm top avg	20.9 19.9 20.7 20.5 0.5 1 26.0 24.7 25.4 0.7	13.5 14.0 15.3 14.3 0.9 2 25.0 24.9 25.1 25.0 0.1	18.6 17.5 18.1 18.1 0.6 3 25.2 25.0 27.1 25.8 1.2	17.6 17.3 18.0 17.6 0.4 4 22.8 24.1 24.0 23.6 0.7	16.0 16.4 16.8 16.4 0.4 5 24.0 23.1 23.3 23.5 0.5	st dev	24.5	% bottom/to 68.46049
5mm top avg st dev	20.9 19.9 20.7 20.5 0.5 1 26.0 24.7 25.4 25.4 0.7	13.5 14.0 15.3 14.3 0.9 2 25.0 24.9 25.1 25.0 0.1	18.6 17.5 18.1 18.1 0.6 3 25.2 25.0 27.1 25.8 1.2	17.6 17.3 18.0 17.6 0.4 22.8 24.1 24.0 23.6 0.7	16.0 16.4 16.8 16.4 0.4 5 24.0 23.1 23.3 23.5 0.5	avg st dev	24.5 0.4	
avg st dev 5mm top avg st dev	20.9 19.9 20.7 20.5 0.5 1 26.0 24.7 25.4 25.4 0.7	13.5 14.0 15.3 14.3 0.9 2 25.0 24.9 25.1 25.0 0.1	18.6 17.5 18.1 18.1 0.6 3 25.2 25.0 27.1 25.8 1.2	17.6 17.3 18.0 17.6 0.4 22.8 24.1 24.0 23.6 0.7	16.0 16.4 16.8 16.4 0.4 5 24.0 23.1 23.3 23.5 0.5	avg st dev	0.2 24.5 0.4	
avg st dev 5mm top avg st dev	20.9 19.9 20.7 20.5 0.5 1 26.0 24.7 25.4 25.4 0.7	13.5 14.0 15.3 14.3 0.9 2 25.0 24.9 25.1 25.0 0.1	18.6 17.5 18.1 18.1 0.6 3 25.2 25.0 27.1 25.8 1.2 3 16.0 16.9 15.7	17.6 17.3 18.0 17.6 0.4 22.8 24.1 24.0 23.6 0.7	16.0 16.4 16.8 16.4 0.4 5 24.0 23.1 23.3 23.5 0.5 5 16.2 15.8 16.7	avg st dev	24.5 0.4	
avg st dev 5mm top avg st dev	20.9 19.9 20.7 20.5 0.5 1 26.0 24.7 25.4 25.4 0.7	13.5 14.0 15.3 14.3 0.9 2 25.0 24.9 25.1 25.0 0.1	18.6 17.5 18.1 18.1 0.6 3 25.2 25.0 27.1 25.8 1.2	17.6 17.3 18.0 17.6 0.4 22.8 24.1 24.0 23.6 0.7	16.0 16.4 16.8 16.4 0.4 5 24.0 23.1 23.3 23.5 0.5	avg st dev	0.2 24.5 0.4	

Table 98. Fuji II LC 24-Hour Hardness (KHN)

2mm			S					
mm								
	1	2	3	4	5			
	21.3	20.4	21.0	20.1	21.0			
top	22.2	21.4	20.5	21.8	21.3	avg	20.8	
	20.3	21.0	19.7	21.1	20.7	st dev	0.3	
avg	21.3	20.9	20.4	21.0	21.0			
t dev	1.0	0.5	0.7	0.9	0.3			% bottom/top
								88.2
					-			
	1	2	3	4	5	-		
	19.0	18.0	17.9	17.7	18.6		70.4	
ottom	17.9	18.7	18.9	17.3	18.5	avg	18.4	
	19.8	19.2	18.4	18.6	18.8	st dev	0.3	
avg	18.9	18.6	18.4	17.9	18.6			-
t dev	1.0	0.6	0.5	0.7	0.2			
mm								
	1	2	3	4	5			
	19.5	19.9	20.8	19.5	19.6			
top	19.8	20.0	21.2	21.1	21.3	avg	20.4	
	20.1	20.5	19.9	20.3	20.5	st dev	0.3	
avg	19.8	20.1	20.6	20.3	20.5			
t dev	0.3	0.3	0.7	0.8	0.9			% bottom/top
								74.7
	1	2	3	4	5			-
	14.9	14.6	14.8	15.8	14.6			1
ttom	14.4	15.5	15.4	15.2	15.9	avg	15.2	
LLUITI	15.1	15.3	15.4	15.7	14.3	st dev	0.2	
oua	14.8	15.1	15.3	15.6	14.9	31 064	0.2	
avg dev	0.4	0.5	0.4	0.3	0.9			+
uev	0.4	0.5	0.4	0.5	0.5			1
mm								
	1	2	3	4	5			-
	19.8	19.6	20.0	19.4	21.1		20.0	
	20.4	20.5	19.3	20.5	21.2	avg	20.2	
top	20.1	20.5						
	21.2	20.8	19.8	19.7	20.2	st dev	0.1	
avg	21.2 20.4	20.8 20.3	19.8 19.7	19.9	20.8	st dev		
avg	21.2	20.8	19.8			st dev		
avg	21.2 20.4	20.8 20.3	19.8 19.7	19.9	20.8	st dev		% bottom/to 57.7
avg	21.2 20.4 0.7	20.8 20.3 0.6	19.8 19.7 0.4	19.9 0.6	20.8 0.6	st dev		
avg	21.2 20.4 0.7	20.8 20.3 0.6	19.8 19.7 0.4	19.9 0.6	20.8 0.6	st dev		
avg t dev	21.2 20.4 0.7 1 9.5	20.8 20.3 0.6	19.8 19.7 0.4 3 9.9	19.9 0.6 4 11.9	20.8 0.6 5 12.5			
avg t dev	21.2 20.4 0.7 1 9.5 9.4	20.8 20.3 0.6 2 13.6 12.4	19.8 19.7 0.4 3 9.9 11.8	19.9 0.6 4 11.9 11.4	20.8 0.6 5 12.5 10.3	avg	11.6	
avg dev	21.2 20.4 0.7 1 9.5 9.4 10.6	20.8 20.3 0.6 2 13.6 12.4 13.5	19.8 19.7 0.4 3 9.9 11.8 11.3	19.9 0.6 4 11.9 11.4 10.8	5 12.5 10.3 10.2			
avg t dev	21.2 20.4 0.7 1 9.5 9.4 10.6 9.8	20.8 20.3 0.6 2 13.6 12.4 13.5	19.8 19.7 0.4 3 9.9 11.8 11.3	19.9 0.6 4 11.9 11.4 10.8 11.4	20.8 0.6 5 12.5 10.3 10.2 11.0	avg	11.6	
avg dev	21.2 20.4 0.7 1 9.5 9.4 10.6	20.8 20.3 0.6 2 13.6 12.4 13.5	19.8 19.7 0.4 3 9.9 11.8 11.3	19.9 0.6 4 11.9 11.4 10.8	5 12.5 10.3 10.2	avg	11.6	% bottom/toj 57.7
avg t dev ottom avg t dev	21.2 20.4 0.7 1 9.5 9.4 10.6 9.8	20.8 20.3 0.6 2 13.6 12.4 13.5	19.8 19.7 0.4 3 9.9 11.8 11.3	19.9 0.6 4 11.9 11.4 10.8 11.4	20.8 0.6 5 12.5 10.3 10.2 11.0	avg	11.6	
avg t dev ottom avg t dev	21.2 20.4 0.7 1 9.5 9.4 10.6 9.8 0.7	20.8 20.3 0.6 2 13.6 12.4 13.5 13.2 0.7	3 9.9 11.8 11.3 11.0	19.9 0.6 4 11.9 11.4 10.8 11.4 0.6	20.8 0.6 5 12.5 10.3 10.2 11.0 1.3	avg	11.6	
avg t dev ottom avg t dev	21.2 20.4 0.7 1 9.5 9.4 10.6 9.8 0.7	20.8 20.3 0.6 2 13.6 12.4 13.5 13.2 0.7	19.8 19.7 0.4 3 9.9 11.8 11.3 11.0 1.0	19.9 0.6 4 11.9 11.4 10.8 11.4 0.6	20.8 0.6 5 12.5 10.3 10.2 11.0 1.3	avg	11.6	
avg t dev ottom avg t dev	21.2 20.4 0.7 1 9.5 9.4 10.6 9.8 0.7	20.8 20.3 0.6 2 13.6 12.4 13.5 13.2 0.7	19.8 19.7 0.4 3 9.9 11.8 11.3 11.0 1.0	19.9 0.6 4 11.9 11.4 10.8 11.4 0.6	20.8 0.6 5 12.5 10.3 10.2 11.0 1.3	avg st dev	11.6	
avg t dev ottom avg t dev	21.2 20.4 0.7 1 9.5 9.4 10.6 9.8 0.7	20.8 20.3 0.6 2 13.6 12.4 13.5 13.2 0.7	19.8 19.7 0.4 3 9.9 11.8 11.3 11.0 1.0	19.9 0.6 4 11.9 11.4 10.8 11.4 0.6 4 21.7 20.0	20.8 0.6 5 12.5 10.3 10.2 11.0 1.3 5 20.7 20.5	avg st dev	11.6 0.3	
avg t dev ottom avg t dev	21.2 20.4 0.7 1 9.5 9.4 10.6 9.8 0.7 21.8 21.0	20.8 20.3 0.6 2 13.6 12.4 13.5 13.2 0.7 2 20.2 19.6 19.8	19.8 19.7 0.4 3 9.9 11.8 11.3 11.0 1.0	19.9 0.6 4 11.9 11.4 10.8 11.4 0.6 4 21.7 20.0 19.7	20.8 0.6 5 12.5 10.3 10.2 11.0 1.3 5 20.7 20.5 19.9	avg st dev	11.6	
avg t dev ottom avg t dev top	21.2 20.4 0.7 1 9.5 9.4 10.6 9.8 0.7 21.8 21.0 21.2	20.8 20.3 0.6 2 13.6 12.4 13.5 13.2 0.7 2 20.2 19.6 19.8	19.8 19.7 0.4 3 9.9 11.8 11.3 11.0 1.0	19.9 0.6 4 11.9 11.4 10.8 11.4 0.6 4 21.7 20.0 19.7 20.5	20.8 0.6 5 12.5 10.3 10.2 11.0 1.3 5 20.7 20.5 19.9 20.4	avg st dev	11.6 0.3 20.2 0.3	57.7
avg t dev ottom avg t dev top avg	21.2 20.4 0.7 1 9.5 9.4 10.6 9.8 0.7 21.8 21.0	20.8 20.3 0.6 2 13.6 12.4 13.5 13.2 0.7 2 20.2 19.6 19.8	19.8 19.7 0.4 3 9.9 11.8 11.3 11.0 1.0	19.9 0.6 4 11.9 11.4 10.8 11.4 0.6 4 21.7 20.0 19.7	20.8 0.6 5 12.5 10.3 10.2 11.0 1.3 5 20.7 20.5 19.9	avg st dev	11.6 0.3 20.2 0.3	57.7 % bottom/to
avg t dev ottom avg t dev top avg	21.2 20.4 0.7 1 9.5 9.4 10.6 9.8 0.7 21.8 21.0 21.2	20.8 20.3 0.6 2 13.6 12.4 13.5 13.2 0.7 2 20.2 19.6 19.8	19.8 19.7 0.4 3 9.9 11.8 11.3 11.0 1.0	19.9 0.6 4 11.9 11.4 10.8 11.4 0.6 4 21.7 20.0 19.7 20.5	20.8 0.6 5 12.5 10.3 10.2 11.0 1.3 5 20.7 20.5 19.9 20.4	avg st dev	11.6 0.3 20.2 0.3	57.7
ottom avg st dev 5mm top avg	21.2 20.4 0.7 1 9.5 9.4 10.6 9.8 0.7 21.8 21.0 21.2 0.6	20.8 20.3 0.6 2 13.6 12.4 13.5 13.2 0.7 2 20.2 19.6 19.8 19.9 0.3	19.8 19.7 0.4 3 9.9 11.8 11.3 11.0 1.0 3 20.3 20.6 19.8 20.2 0.4	19.9 0.6 4 11.9 11.4 10.8 11.4 0.6 4 21.7 20.0 19.7 20.5 1.1	20.8 0.6 5 12.5 10.3 10.2 11.0 1.3 5 20.7 20.5 19.9 20.4	avg st dev	11.6 0.3 20.2 0.3	57.7 % bottom/to
avg ottom avg st dev	21.2 20.4 0.7 1 9.5 9.4 10.6 9.8 0.7 21.8 21.0 21.2 0.6	20.8 20.3 0.6 2 13.6 12.4 13.5 13.2 0.7 2 20.2 19.6 19.8 19.9 0.3	19.8 19.7 0.4 3 9.9 11.8 11.3 11.0 1.0 3 20.3 20.6 19.8 20.2 0.4	19.9 0.6 4 11.9 11.4 10.8 11.4 0.6 4 21.7 20.0 19.7 20.5 1.1	20.8 0.6 5 12.5 10.3 10.2 11.0 1.3 5 20.7 20.5 19.9 20.4 0.4	avg st dev	11.6 0.3 20.2 0.3	57.7 % bottom/to
avg ottom avg t dev top avg st dev	21.2 20.4 0.7 1 9.5 9.4 10.6 9.8 0.7 21.8 21.0 21.2 0.6	20.8 20.3 0.6 2 13.6 12.4 13.5 13.2 0.7 2 20.2 19.6 19.8 19.9 0.3	19.8 19.7 0.4 3 9.9 11.8 11.3 11.0 1.0 3 20.3 20.6 19.8 20.2 0.4	19.9 0.6 4 11.9 11.4 10.8 11.4 0.6 4 21.7 20.0 19.7 20.5 1.1	20.8 0.6 5 12.5 10.3 10.2 11.0 1.3 5 20.7 20.5 19.9 20.4 0.4 5 9.7	avg st dev	11.6 0.3 20.2 0.3	57.7 % bottom/to
avg t dev ottom avg t dev top avg	21.2 20.4 0.7 1 9.5 9.4 10.6 9.8 0.7 21.8 21.0 21.2 0.6	20.8 20.3 0.6 2 13.6 12.4 13.5 13.2 0.7 2 20.2 19.6 19.8 19.9 0.3 2 8.7 10.8	19.8 19.7 0.4 3 9.9 11.8 11.3 11.0 1.0 3 20.3 20.6 19.8 20.2 0.4	19.9 0.6 4 11.9 11.4 10.8 11.4 0.6 4 21.7 20.0 19.7 20.5 1.1 4 9.2 9.6	20.8 0.6 5 12.5 10.3 10.2 11.0 1.3 5 20.7 20.5 19.9 20.4 0.4 5 9.7 9.6	avg st dev	11.6 0.3 20.2 0.3	57.7 % bottom/to
avg t dev ottom avg t dev top avg t dev	21.2 20.4 0.7 1 9.5 9.4 10.6 9.8 0.7 21.8 21.0 21.2 0.6	20.8 20.3 0.6 2 13.6 12.4 13.5 13.2 0.7 2 20.2 19.6 19.8 19.9 0.3	19.8 19.7 0.4 3 9.9 11.8 11.3 11.0 1.0 3 20.3 20.6 19.8 20.2 0.4	19.9 0.6 4 11.9 11.4 10.8 11.4 0.6 4 21.7 20.0 19.7 20.5 1.1	20.8 0.6 5 12.5 10.3 10.2 11.0 1.3 5 20.7 20.5 19.9 20.4 0.4 5 9.7	avg st dev	11.6 0.3 20.2 0.3	57.7 % bottom/to

Table 99. Fuji II LC One-Week Hardness (KHN)

			District Co.				11 10 10 10	200-170-1-12
mm	4	2	3	4	5			
-	1 20.0	2						
	20.9 19.3	21.5 22.0	21.6 19.8	20.7	20.0 19.5	21.00	20.7	
top						avg	0.3	
	20.7	22.2	19.3	20.3	19.1	st dev	0.3	
avg	20.3	21.9	20.2	21.0	19.5			
t dev	0.9	0.4	1.2	0.9	0.5			% bottom/to
								90.77
								-
	1	2	3	4	5			
	20.8	20.7	19.1	17.0	21.7		40.0	
ottom	20.7	19.3	18.9	16.4	18.3	avg	18.8	
	19.2	19.2	19.8	15.7	19.0	st dev	0.5	
avg	20.2	19.7	19.3	16.4	19.7			
t dev	0.9	0.8	0.5	0.7	1.8			
	1000			STATE OF THE PARTY				September 1
mm		THE REAL PROPERTY.		Maria Maria	THE PERSON NAMED IN	100000		CONTRACTOR OF THE PARTY.
	1	2	3	4	5			
	20.5	20.1	19.9	21.3	20.6		20.0	
top	21.9	20.2	20.6	21.2	21.1	avg	20.6	
	22.0	20.1	20.9	20.2	20.7	st dev	0.3	
avg	21.5	20.1	20.5	20.9	20.8			
t dev	0.8	0.1	0.5	0.6	0.3			% bottom/to
								88.50
	1	2	3	4	5			
	18.0	17.6	17.9	19.4	18.0			
ottom	19.8	18.0	17.7	19.6	17.1	avg	18.2	
	19.3	18.6	18.1	18.3	18.2	st dev	0.3	
avg	19.0	18.1	17.9	19.1	17.8			
t dev	0.9	0.5	0.2	0.7	0.6			
					CAROLINA CONTRACTOR		ALCOHOLD IN SECTION	COLORODO PARA MARIA
lmm								
	1	2	3	4	5			
	19.8	20.1	20.0	20.9	22.6			
top	21.0	21.2	21.6	21.0	22.5	avg	21.3	
	21.9	21.0	21.4	20.7	22.0	st dev	0.4	
avg	20.9	20.8	21.0	20.9	22.4			
t dev	1.1	0.6	0.9	0.2	0.3			% bottom/to
								72.59
	1	2	3	4	5			
	14.1	13.8	14.5	15.7	15.8			
ottom	14.4	16.1	15.9	15.6	13.9	avg	15.4	
	15.5	16.6	15.0	14.8	17.4	st dev	0.6	
	13.3							
avo			15.1	15.4	15./			
avg	14.7	15.5	0.7	0.5	15.7			
			0.7	0.5	15.7			
avg st dev	14.7	15.5						
t dev	14.7	15.5						
	14.7	15.5						
t dev	14.7	15.5	3	0.5	1.8			
t dev	14.7 0.7 1 20.9	15.5 1.5	3 21.6	0.5 4 20.3	1.8 5 22.2	avo	21.3	
t dev	14.7 0.7 1 20.9 21.4	2 2 22.1 22.2	3 21.6 21.4	4 20.3 20.0	5 22.2 21.8	avg st dev	21.3	
5mm top	14.7 0.7 1 20.9 21.4 20.6	2 22.1 22.2 22.6	3 21.6 21.4 20.3	4 20.3 20.0 19.5	5 22.2 21.8 21.8	avg st dev	21.3 0.2	
5mm top	14.7 0.7 1 20.9 21.4 20.6 21.0	2 22.1 22.2 22.6 22.3	3 21.6 21.4 20.3 21.1	4 20.3 20.0 19.5	5 22.2 21.8 21.8 21.9		0.2	% hattam/ta
5mm top	14.7 0.7 1 20.9 21.4 20.6	2 22.1 22.2 22.6	3 21.6 21.4 20.3	4 20.3 20.0 19.5	5 22.2 21.8 21.8		0.2	% bottom/to
5mm top	14.7 0.7 1 20.9 21.4 20.6 21.0	2 22.1 22.2 22.6 22.3	3 21.6 21.4 20.3 21.1	4 20.3 20.0 19.5	5 22.2 21.8 21.8 21.9		0.2	% bottom/to 60.24
5mm top	14.7 0.7 1 20.9 21.4 20.6 21.0 0.4	2 22.1 22.2 22.6 22.3 0.3	3 21.6 21.4 20.3 21.1 0.7	4 20.3 20.0 19.5 19.9 0.4	5 22.2 21.8 21.8 21.9		0.2	
5mm top	14.7 0.7 1 20.9 21.4 20.6 21.0 0.4	2 22.1 22.2 22.6 22.3 0.3	3 21.6 21.4 20.3 21.1 0.7	4 20.3 20.0 19.5 19.9 0.4	5 22.2 21.8 21.8 21.9 0.2		0.2	
6mm top avg	14.7 0.7 1 20.9 21.4 20.6 21.0 0.4	2 22.1 22.2 22.6 22.3 0.3	3 21.6 21.4 20.3 21.1 0.7	4 20.3 20.0 19.5 19.9 0.4	5 22.2 21.8 21.8 21.9 0.2	st dev	0.2	
5mm top	14.7 0.7 1 20.9 21.4 20.6 21.0 0.4	2 22.1 22.2 22.6 22.3 0.3	3 21.6 21.4 20.3 21.1 0.7	4 20.3 20.0 19.5 19.9 0.4 4 13.0 11.9	5 22.2 21.8 21.8 21.9 0.2 5 13.1	st dev	12.8	
6mm top avg	14.7 0.7 1 20.9 21.4 20.6 21.0 0.4	2 22.1 22.2 22.6 22.3 0.3	3 21.6 21.4 20.3 21.1 0.7	4 20.3 20.0 19.5 19.9 0.4	5 22.2 21.8 21.8 21.9 0.2	st dev	0.2	

Table 100. Fuji II LC One-Month Hardness (KHN)

BERTHAND TO STATE OF THE PARTY			No. of Concession, Name of Street, or other Designation, Name of Street, Name		and the second second	AND PROPERTY OF STREET	MARKATANA MARANA	DECEMBER OF THE PARTY OF THE PA
mm					SELECTION OF SELECTION			
	1	2	3	4	5			
	18.7	18.9	19.0	18.8	21.2		ATTOCKED TO A SECTION	
top	19.3	21.3	18.8	20.3	21.1	avg	20.0	
	19.0	20.7	19.3	21.2	19.9	st dev	0.5	
avg	19.0	20.3	19.0	20.1	20.7			
t dev	0.3	1.2	0.3	1.2	0.7			% bottom/top
								90.94
	1	2	3	4	5			
	17.2	17.3	19.0	18.2	18.5			
ottom	17.5	17.5	17.9	18.3	18.8	avg	18.2	
	17.5	18.5	17.3	18.7	18.7	st dev	0.3	
avg	17.4	17.8	18.1	18.4	18.7			
t dev	0.2	0.6	0.9	0.3	0.2			
mm								
	1	2	3	4	5			
	20.3	19.2	21.8	22.4	20.4			
top	20.6	19.4	20.5	19.3	19.8	avg	20.1	
	21.0	18.5	19.2	20.9	20.3	st dev	0.6	
avg	20.6	19.0	20.5	20.9	20.2		1	
dev	0.4	0.5	1.3	1.6	0.3			% bottom/top
	0.4	0.0	1.0					82.91
	1	2	3	4	5			
	17.0	17.1	16.0	17.8	15.8			
ottom	16.3	16.5	15.7	17.9	16.5	avg	16.7	
	17.5	16.7	16.2	17.3	16.9	st dev	0.2	
	16.9	16.8	16.0	17.7	16.4	J. 007	V.E	
avg	0.6	0.3	0.3	0.3	0.6			-
t dev	0.0	0.3	0.3	0.5	0.0			
mm			10 BK 15 B				100	
and the same	1	2	3	4	5			
	20.5	20.8	18.8	20.1	19.2			
ton	20.9	19.8	18.6	20.9	19.0	avg	19.6	
top			18.0	19.9	21.2	st dev	0.4	
	18.3	18.7		20.3	19.8	St dev	0.4	
avg	19.9	19.8	18.5		19.8			% bottom/to
t dev	1.4	1.1	0.4	0.5	1.2			% bottom/to
						-		00.49
		^	2		5			
	1	2	3	4				
	14.4	12.2	12.3	11.2	11.1		44.0	
ottom	12.5	12.3	11.8	12.7	11.0	avg	11.8	
	15.1	13.2	11.9	9.9	12.5	st dev	0.5	
vg	14.0	12.6	12.0	11.3	11.5			
dev	1.3	0.6	0.3	1.4	0.8			-
								-
	TO SERVICE THE PARTY OF THE PAR							THE RESERVE AND ADDRESS OF THE PERSON NAMED IN COLUMN 1
mm				ALL SHEET AND				
	1	2	3	4	5			-
	20.8	20.2	20.0	19.9	20.2		00	
top	19.2	18.0	20.3	20.9	20.5	avg	20.1	E
	20.5	20.0	19.9	21.8	19.9	st dev	0.4	
avg	20.2	19.4	20.1	20.9	20.2			
t dev	0.9	1.2	0.2	1.0	0.3			% bottom/to
								54.23
-								
	1	2	3	4	5			
				0.4	11.5			
		10.0	11.9	9.4	11.3			
ottom	9.1			9.4	11.6	avg	10.9	
ottom	9.1 9.7	11.0	9.9	9.7	11.6		10.9	
ottom	9.1					avg st dev		

Table 101. Fuji II LC Three-Month Hardness (KHN)

	3 Month H	uruncoo						
2mm								
	1	2	3	4	5			
	20.3	20.7	22.3	21.7	20.9		20.0	
top	20.4	19.6	20.7	18.0	18.3	avg	20.6	
	21.1	20.1	21.5	23.7	20.2	st dev	1.0	1/
avg	20.6	20.1	21.5 0.8	21.1	19.8			% bottom/top
st dev	0.4	0.6	0.0	2.9	1.3			90.23
	1	2	3	4	5			
	19.8	18.9	18.4	18.3	18.1			
ottom	19.7	17.1	18.4	19.8	18.9	avg	18.6	
	19.5	19.9	17.5	20.5	17.7	st dev	0.5	
avg	19.7	18.6	18.1	19.5	18.2			
st dev	0.2	1.4	0.5	1.1	0.6			
								-
3mm								
	1	2	3	4	5			
	20.3	21.3	21.4	20.5	20.5		20.0	
top	21.4	20.5	21.0	20.2	20.2	avg et dou	20.6	
200	21.3	20.1	22.2	19.7 20.1	19.7 20.1	st dev	0.1	
avg st dev	0.6	0.6	0.6	0.4	0.4			% bottom/top
n uev	0.0	0.0	0.0	0.4	0.4			94.18
								54.10
	1	2	3	4	5			
	19.7	19.3	19.6	19.1	19.1			
ottom	19.3	19.6	19.5	19.3	19.3	avg	19.4	
	19.6	19.2	20.1	19.4	19.4	st dev	0.1	,
avg	19.5	19.4	19.7	19.3	19.3			
t dev	0.2	0.2	0.3	0.2	0.2			
4mm								
	1	2	3	4	5			
	20.3	20.6	22.8	20.0	20.2			
top	21.9	21.8	20.5	20.5	21.4	avg	21.3	
	22.2	22.0	23.3	21.9	20.9	st dev	0.3	
avg	21.5	21.5	22.2	20.8	20.8			% bottom/to
st dev	1.0	0.8	1.5	1.0	0.0			61.43
	1	2	3	4	5			
	13.3	15.9	14.4	13.4	12.2			
ottom	13.8	13.8	15.4	11.1	11.8	avg	13.1	
	15.2	13.0	15.5	10.1	10.6	st dev	0.5	
avg	14.1	14.2	15.1	11.5	11.5			
st dev	1.0	1.5	0.6	1.7	0.8			
5mm								
	1	2	3	4	5			-
	22.0	20.8	20.4	20.0	22.4		21.0	
top	21.6	20.7	21.0	20.4	21.8	avg et dou	0.2	
	21.0	22.1	20.2	20.4	21.2	st dev	0.2	
avg	21.5	0.8	0.4	0.2	0.6			% bottom/to
st dev	0.5	0.0	0.4	0.2	0.0			50.62
	1	2	3	4	5			
	11.9	11.3	10.9	9.8	9.1		40.0	
oottom	10.3	9.4	11.0	12.8	10.5	avg	10.6	
	12.8	10.4	10.0	11.3	10.8	st dev	0.4	
	ARREST THE PERSON	ALC: UNKNOWN THE REAL PROPERTY.	40.0					
avg st dev	11.7	0.9	10.6	11.3	0.9			

Table 102. Photac-Fil Quick 24-Hour Hardness (KHN)

No. of Lot	47 100 (A) (A)	are a section		The second second	C/20, \$15, 35 (E)	Section in		TA S SEASON
m		2	3	4	5			
-	1 11 0	11.3	12.8	11.0	11.9			
_	11.8 11.6	10.1	12.3	12.7	11.6	avg	11.5	
p								
	12.2	11.0	11.0	10.6	11.5 11.7	st dev	0.4	
g	11.9	10.8	12.0	11.4				V
lev [0.3	0.6	0.9	1.1	0.2			% bottom/to
								81.7
-					<u>_</u>			
	1	2	3	4	5			
	9.8	9.0	10.3	8.7	7.8			
om	10.8	8.9	11.2	10.3	9.3	avg	9.4	
	10.4	9.5	10.1	9.1	8.5	st dev	0.2	
g	10.3	9.1	10.5	9.4	8.5			
lev	0.5	0.3	0.6	0.8	0.7			
						-		
m	965							
1	1	2	3	4	5			
	11.6	13.1	12.6	11.9	10.4			
p	13.3	12.4	11.7	10.9	10.1	avg	11.6	
	12.2	10.9	12.7	12.5	10.2	st dev	0.4	
g	12.4	12.1	12.3	11.8	10.2			
lev	0.9	1.1	0.6	0.8	0.2			% bottom/to
164	0.9	1,1	0.0	0.0	0.2			79.0
						-		10.0
	1	2	3	4	5			
-	1 0.7		8.7	9.0	8.0			
	9.7	9.2			9.4	~~~	9.2	
om	10.0	8.5	9.9	8.3		avg	The second second second second	
	8.7	9.8	10.0	9.6	9.8	st dev	0.1	
g	9.5	9.2	9.6	8.9	9.1			-
lev [0.7	0.7	0.7	0.7	0.9			
						and out to see	CONTROL OF THE PARTY OF	
m		10 may 194	STREET,				HEREN	E-NEW SE
	1	2	3	4	5			
	13.4	12.7	11.6	11.1	11.8			
p	12.4	12.9	10.7	12.7	12.8	avg	12.0	
	12.4	13.5	10.6	11.1	12.1	st dev	0.2	
g	12.7	13.0	11.0	11.6	12.2			
lev	0.6	0.4	0.6	0.9	0.5			% bottom/to
								70.0
-								
	1	2	3	4	5	-		
	7.6	8.2	7.4	8.1	8.9			
- H	10.1	9.3	8.1	8.0	8.4	avg	8.4	
om			7.7	8.6	9.3	st dev	0.5	
	10.4	8.6	ARREST TO THE REAL PROPERTY.			21 GBA	0.0	
g	9.4	8.7	7.7	8.2	8.9	-		
dev	1.5	0.6	0.4	0.3	0.4			
		THE REAL PROPERTY.				ALC: SOME OF THE PERSON NAMED IN	State of Manager	The Southern St.
ım		PARTIES NA	SE STATE OF	The same of	235		AND LOSS.	PERSONAL PROPERTY.
	1	2	3	4	5			
	11.7	12.5	12.6	10.7	12.0		-	
р	11.4	11.2	12.7	10.7	12.4	avg	11.6	
	13.2	11.4	11.9	9.2	11.6	st dev	0.3	
/g	12.1	11.7	12.4	10.2	12.0			
dev	1.0	0.7	0.4	0.9	0.4			% bottom/to
								75.5
	1	2	3	4	5			
	9.5	8.8	8.2	7.8	7.6			
					7.7		8.7	
om F		10.4	91	96	1 //	awn		35
tom	8.0	10.4	9.1	9.6	7.7	avg st dev		
tom /g		10.4 10.4 9.9	9.1 8.2 8.5	9.6 9.6 9.0	7.4	st dev	0.4	

Table 103. Photac-Fil Quick One-Week Hardness (KHN)

			day or manager			15,000,000	NAME OF TAXABLE PARTY.	Maria Special Control
2mm								FARMANET
	1	2	3	4	5			
	14.1	14.2	14.6	13.2	13.9		12.7	
top	13.8	12.5	14.4	13.5	14.2	avg	13,7	
0110	13.3	12.9	13.4	13.8	14.1	st dev	0.5	
avg	13.7	13.2	14.1	13.5	0.2			% hattam#-
st dev	0.4	0.9	0.6	0.3	0.2			% bottom/to 87.3
								07.3
	1	2	2		5			-
-	12.5	2	3	11.4	11.8			
	12.5	12.7	12.1				12.0	
bottom	12.2	13.1	11.9	11.6	11.5	avg	12.0	-
	12.3	11.9	12.0	11.9	11.9	st dev	0.2	-
avg	12.3	12.6	12.0	11.6	11.7			
st dev [0.2	0.6	0.1	0.3	0.2			-
							NOT RESIDENCE	
3mm		PROPERTY.					ESSER S	Contract Section 1
-	1	2	3	4	5			
	13.9	13.7	14.2	14.5	13.4			
top	13.3	12.6	14.1	13.4	13.3	avg	13.6	
	12.6	12.6	13.6	13.3	13.9	st dev	0.2	
avg	13.3	13.0	14.0	13.7	13.5			
st dev	0.7	0.6	0.3	0.7	0.3			% bottom/to
								83.3
	1	2	3	4	5			
	11.0	11.8	11.3	10.6	10.6			
ottom	10.4	11.9	11.9	11.0	11.2	avg	11.3	
	12.1	11.9	11.0	11.3	11.0	st dev	0.3	
avg	11.2	11.9	11.4	11.0	10.9			
st dev	0.9	0.1	0.5	0.4	0.3			
4mm								
	1	2	3	4	5			
	14.5	12.4	12.5	13.6	13.3			
top	13.9	12.8	11.7	13.3	13.1	avg	12.8	
	12.8	12.9	12.3	12.5	12.8	st dev	0.3	-
avg	13.7	12.7	12.2	13.1	13.1			
st dev	0.9	0.3	0.4	0.6	0.3			% bottom/to
	0.0	0.0	0.7	5.0	0.0			82.3
								02.0
1.0	1	2	3	4	5			+
			10.5	10.6	11.2			
	10.0	0.0		IU.D	11.2			
	10.6	9.6				0117	10.5	
ottom	10.9	10.9	10.4	10.9	11.0	avg	10.5	
	10.9 10.3	10.9 10.2	10.4 10.2	10.9 10.5	11.0 10.1	avg st dev	10.5 0.2	
avg	10.9 10.3 10.6	10.9 10.2 10.2	10.4 10.2 10.4	10.9 10.5 10.7	11.0 10.1 10.8			
avg	10.9 10.3	10.9 10.2	10.4 10.2	10.9 10.5	11.0 10.1			
avg	10.9 10.3 10.6	10.9 10.2 10.2	10.4 10.2 10.4	10.9 10.5 10.7	11.0 10.1 10.8			
avg st dev	10.9 10.3 10.6	10.9 10.2 10.2	10.4 10.2 10.4	10.9 10.5 10.7	11.0 10.1 10.8			
avg st dev	10.9 10.3 10.6 0.3	10.9 10.2 10.2 0.7	10.4 10.2 10.4 0.2	10.9 10.5 10.7 0.2	11.0 10.1 10.8 0.6			
avg st dev	10.9 10.3 10.6 0.3	10.9 10.2 10.2 0.7	10.4 10.2 10.4 0.2	10.9 10.5 10.7 0.2	11.0 10.1 10.8 0.6			
avg st dev	10.9 10.3 10.6 0.3	10.9 10.2 10.2 0.7	10.4 10.2 10.4 0.2	10.9 10.5 10.7 0.2	11.0 10.1 10.8 0.6		0.2	
avg st dev	10.9 10.3 10.6 0.3	10.9 10.2 10.2 0.7	10.4 10.2 10.4 0.2	10.9 10.5 10.7 0.2 4 12.3 12.0	11.0 10.1 10.8 0.6 5 12.5 13.8		12.7	
avg st dev	10.9 10.3 10.6 0.3	10.9 10.2 10.2 0.7	10.4 10.2 10.4 0.2	10.9 10.5 10.7 0.2	11.0 10.1 10.8 0.6	st dev	0.2	
avg st dev 5mm	10.9 10.3 10.6 0.3	10.9 10.2 10.2 0.7 2 13.0 13.2	10.4 10.2 10.4 0.2 3 13.0 12.4	10.9 10.5 10.7 0.2 4 12.3 12.0	11.0 10.1 10.8 0.6 5 12.5 13.8	st dev	12.7	
avg st dev 5mm top avg	10.9 10.3 10.6 0.3 1 1 12.9 12.7 13.4 13.0	10.9 10.2 10.2 0.7 2 13.0 13.2 12.9	10.4 10.2 10.4 0.2 3 13.0 12.4 12.1	10.9 10.5 10.7 0.2 4 12.3 12.0 12.5 12.3	11.0 10.1 10.8 0.6 5 12.5 13.8 12.8	st dev	12.7	% bottom/to
avg st dev 5mm top avg	10.9 10.3 10.6 0.3 1 1 12.9 12.7 13.4	10.9 10.2 10.2 0.7 2 13.0 13.2 12.9	10.4 10.2 10.4 0.2 3 13.0 12.4 12.1	10.9 10.5 10.7 0.2 4 12.3 12.0 12.5	11.0 10.1 10.8 0.6 5 12.5 13.8 12.8 13.0	st dev	12.7	% bottom/to 80.4
avg st dev 5mm top	10.9 10.3 10.6 0.3 1 1 12.9 12.7 13.4 13.0	10.9 10.2 10.2 0.7 2 13.0 13.2 12.9	10.4 10.2 10.4 0.2 3 13.0 12.4 12.1	10.9 10.5 10.7 0.2 4 12.3 12.0 12.5 12.3	11.0 10.1 10.8 0.6 5 12.5 13.8 12.8 13.0	st dev	12.7	
avg st dev 5mm top avg	10.9 10.3 10.6 0.3 12.9 12.7 13.0 0.4	10.9 10.2 10.2 0.7 2 13.0 13.2 12.9 13.0 0.2	10.4 10.2 10.4 0.2 3 13.0 12.4 12.1 12.5 0.5	10.9 10.5 10.7 0.2 4 12.3 12.0 12.5 12.3 0.3	11.0 10.1 10.8 0.6 5 12.5 13.8 12.8 13.0	st dev	12.7	
avg st dev 5mm top avg	10.9 10.3 10.6 0.3 12.9 12.7 13.4 13.0 0.4	10.9 10.2 10.2 0.7 2 13.0 13.2 12.9 13.0 0.2	3 13.0 12.4 12.1 12.5 0.5	10.9 10.5 10.7 0.2 4 12.3 12.0 12.5 12.3 0.3	11.0 10.1 10.8 0.6 5 12.5 13.8 12.8 13.0 0.7	st dev	12.7	
5mm top avg st dev	10.9 10.3 10.6 0.3 12.9 12.7 13.4 13.0 0.4	10.9 10.2 10.2 0.7 2 13.0 13.2 12.9 13.0 0.2	3 13.0 12.4 12.1 12.5 0.5	10.9 10.5 10.7 0.2 4 12.3 12.0 12.5 12.3 0.3	11.0 10.1 10.8 0.6 5 12.5 13.8 12.8 13.0 0.7	avg st dev	12.7	
avg st dev 5mm top avg	10.9 10.3 10.6 0.3 1 12.9 12.7 13.4 13.0 0.4	10.9 10.2 10.2 0.7 2 13.0 13.2 12.9 13.0 0.2	3 13.0 12.4 12.1 12.5 0.5	10.9 10.5 10.7 0.2 4 12.3 12.0 12.5 12.3 0.3	11.0 10.1 10.8 0.6 5 12.5 13.8 12.8 13.0 0.7	avg st dev	12.7 0.2	
avg st dev 5mm top avg st dev	10.9 10.3 10.6 0.3 12.9 12.7 13.4 13.0 0.4	10.9 10.2 10.2 0.7 2 13.0 13.2 12.9 13.0 0.2	3 13.0 12.4 12.1 12.5 0.5 3 10.5 9.4 10.0	10.9 10.5 10.7 0.2 4 12.3 12.0 12.5 12.3 0.3	11.0 10.1 10.8 0.6 5 12.5 13.8 12.8 13.0 0.7 5 9.8 10.6	avg st dev	12.7	
avg Smm top avg st dev	10.9 10.3 10.6 0.3 1 12.9 12.7 13.4 13.0 0.4	10.9 10.2 10.2 0.7 2 13.0 13.2 12.9 13.0 0.2	3 13.0 12.4 12.1 12.5 0.5	10.9 10.5 10.7 0.2 4 12.3 12.0 12.5 12.3 0.3	11.0 10.1 10.8 0.6 5 12.5 13.8 12.8 13.0 0.7	avg st dev	12.7 0.2	% bottom/to 80.4

Table 104. Photac-Fil Quick One-Month Hardness (KHN)

mm								
	1	2	3	4	5			
	20.1	17.1	18.2	19.2	17.8			
top	17.0	17.4	18.0	17.7	16.4	avg	17.5	
	17.0	16.4	17.2	17.5	17.6	st dev	0.5	
avg	18.0	17.0	17.8	18.1	17.3			
t dev	1.8	0.5	0.5	0.9	0.8			% bottom/top
I dev	1.0	0.5	0.5	0.3	0.0			92.9
								02.0
	1	2	3	4	5			
	16.1	17.2	16.9	16.3	16.2			
	16.1	16.5	16.4	15.8	15.6	avg	16.3	
ottom						st dev	0.1	
	15.6	16.5	16.5	15.9	15.8	Stuev	0.1	
avg	15.9	16.7	16.6	16.0	15.9			
t dev	0.3	0.4	0.3	0.3	0.3			-
		STREET, STREET, STREET, STREET, STREET, STREET, STREET, STREET, STREET, STREET, STREET, STREET, STREET, STREET,		SERVICE COM	TO A CONTRACT OF			A CONTRACTOR
mm				STATE OF THE PARTY.	SAN SAN SAN	V. 1963461		ACCRECATE A
	1	2	3	4	5			
	17.7	19.4	18.7	19.8	17.6			
top	16.7	17.3	19.0	17.2	17.3	avg	18.1	
	16.9	17.6	18.5	17.5	17.4	st dev	0.6	
avg	17.1	18.1	18.7	18.2	17.4			
t dev	0.5	1.1	0.3	1.4	0.2			% bottom/top
								86.1
	1	2	3	4	5			
1	15.4	16.7	15.8	15.7	15.9			
ttom	15.1	15.8	15.0	14.9	15.2	avg	15.6	
	15.2	16.8	14.9	14.8	15.6	st dev	0.2	
0110	15.2	16.4	15.2	15.1	15.6	51 054	0.2	
avg			0.5	0.5	0.4			
t dev	0.2	0.6	0.5	0.5	0.4			
					and the same		2001.00	17579535
mm	1	2	3	4	5			A CONTRACTOR OF THE PARTY OF TH
	18.7	19.1	18.6	18.9	19.0		10.2	
top	18.0	17.9	18.8	18.8	17.2	avg	18.2	
	17.6	18.3	17.1	17.5	17.5	st dev	0.2	
avg	18.1	18.4	18.2	18.4	17.9			
t dev	0.6	0.6	0.9	0.8	1.0			% bottom/to
								79.7
	1	2	3	4	5			
	15.3	15.5	14.9	12.2	14.5			
ottom	15.2	15.1	15.1	13.8	15.2	avg	14.5	
	15.0	14.9	15.0	14.0	14.2	st dev	0.4	
ava	15.2	15.2	15.0	13.3	14.6			
avg	0.2	0.3	0.1	1.0	0.5			-1
t dev	0.2	0.3	0.1	1.0	0.0			
imm	A. Walte	10 mg m		(a) (b) (c) (c)	STATE OF THE PARTY.	VALUE OF STREET	KEL 18	
A CONTRACTOR OF THE PARTY OF TH	1	2	3	4	5			
			17.6	19.8	17.8	1		
	17.4	16.7			17.0		17.9	
top	14.8	17.1	17.3	19.7	17.2	avg		
	17.6	18.2	17.4	18.8	17.3	st dev	0.6	
avg	16.6	17.3	17.4	19.4	17.4			
t dev	1.6	0.8	0.2	0.6	0.3			% bottom/to
								62.4
	1	2	3	4	5			
	10.9	9.7	10.4	13.5	10.4			
ottom	9.0	10.3	10.9	13.7	10.8	avg	11.2	
	9.7	10.4	11.3	12.7	10.0	st dev	0.2	
					ACCRECATION AND ADDRESS OF THE PARTY OF THE			
avg	9.9	10.1	10.9	13.3	10.4			

Table 105. Photac-Fil Quick Three-Month Hardness (KHN)

and the last of th	and the latest section in the latest section	-	Name of Street, Street, or other Designation of the Street, or oth	CONTRACTOR OF THE PARTY OF THE	Market Street, Square	NAME OF TAXABLE PARTY.	COLUMN TO LA COLUMN TO LA	CONTRACTOR OF THE PARTY.
2mm		SERVICE SERVICE		04/4/4/4/4				NAME OF STREET
	1	2	3	4	5			
	18.4	17.9	17.5	16.8	19.2		CHARLES AND ADVANCED	
top	15.3	16.5	16.5	17.2	16.9	avg	17.0	
	16.9	15.4	16.3	17.8	16.5	st dev	0.5	
avg	16.9	16.6	16.8	17.3	17.5			
t dev	1.6	1.3	0.6	0.5	1.5		9	% bottom/to
								90.56
	4	2	3	4	5			
-	1 40.4			15.6	15.8			
	16.4	16.2	16.2				15.4	
ottom	15.9	14.9	16.1	15.2	15.2	avg	CONTRACTOR AND ADDRESS OF THE PARTY OF THE P	
-	16.0	15.0	14.2	15.0	15.8	st dev	0.4	
avg	16.1	15.4	15.5	15.3	15.6			
t dev	0.3	0.7	1.1	0.3	0.3			
mm		2						CASA CONTRACTOR
	1	2	3	4	5			
	16.8	17.4	16.4	16.9	17.6		40.0	
top	17.0	16.0	17.0	16.6	16.9	avg	16.8	
	16.2	16.4	16.1	17.0	17.0	st dev	0.2	
avg	16.7	16.6	16.5	16.8	17.2			
t dev	0.4	0.7	0.5	0.2	0.4		9	% bottom/to
								80.73
	1	2	3	4	5			
	13.4	14.2	13.3	11.2	13.9			
ottors	13.4	13.6	13.8	13.1	13.3	avg	13.5	
ottom	14.3	13.4	14.0	14.5	14.2	st dev	0.5	
						St dev	0.5	
avg	13.6	13.7	13.7	12.9	13.8			
t dev	0.6	0.4	0.4	1.7	0.5			
lmm				4	5			MARKA NAMES
	1 47.0	2	3		17.6			
	17.2	16.9	16.0	17.0			40.0	
top	16.8	16.1	17.0	17.2	17.0	avg	16.9	
	16.4	16.6	16.8	16.9	17.4	st dev	0.1	
avg	16.8	16.5	16.6	17.0	17.3			
st dev	0.4	0.4	0.5	0.2	0.3		9	% bottom/to
								66.37
	1	2	3	4	5			
1	10.6	11.8	11.5	10.9	10.7			
ottom	10.8	11.8	11.4	11.5	10.1	avg	11.2	
ottoili	9.9	11.6	10.7	11.4	11.0	st dev	0.1	
	ACCOUNT WATER	ARREST VICTORIAN	Annual Control of the	11.3	10.6	5. 001		
avg	10.4	0.1	0.4	0.3	0.5			
t dev	0.5	0.1	0.4	0.3	0.3			
						AND CONTRACTOR		Carly Incide
mm	1	2	3	4	5	100000000000000000000000000000000000000	O SOUTH OF	AND DESCRIPTION OF
	17.7	16.3	17.1	17.8	18.0			-
				18.0	16.2	avg	17.3	
top	17.8	16.7	17.7				0.4	
	17.7	18.0	17.5	17.7	16.8	st dev	0.4	
avg	17.7	17.0	17.4	17.8	17.0			0/
st dev	0.1	0.9	0.3	0.2	0.9			% bottom/to 61.85
	1	2	3	4	5			
	11.2	11.3	10.9	9.9	12.2			
ottom	12.6	11.4	9.4	9.2	11.9	avg	10.7	
	11.0	10.9	10.1	10.2	11.2	st dev	0.2	
-	11.6	11.2	10.1	9.7	11.8			
avg	SECTION IN COLUMN							

APPENDIX C

SPECIFIC HEAT THERMAL ANALYSIS DATA

Table 106. Vitremer Specific Heat Results (J/gC°)

Vitremer	Specific I	leat Capac	city
MMEDIATE			
Sample	2mm	3mm	4mm
1	1.2	1.09	0.97
2	1.31	1.08	0.95
3	1.19	1.09	0.94
4	1.2	1.03	0.91
5	1.17	1.02	1.01
Mean	1.21	1.06	0.96
SD	0.06	0.03	0.04
	0.00		
24 Hours			
Sample	2mm	3mm	4mm
1	1.26	0.99	1
2	1.32	1.01	0.96
	1.32	1.12	0.96
3	1.26		0.77
4		1.07	0.95
5 Maan	1.17 1.22	1.05 1.05	0.93
Mean	0.09	0.05	0.92
SD	0.09	0.03	0.03
1 WEEK			
	20000	3mm	Amm
Sample 1	2mm 1.32	1.15	4mm 1.02
1	1.32		The second secon
2		1.04	0.99
3	1.28	1.1	0.95
4	1.28	1.07	0.94
5	1.27	1.12	1
Mean	1.28	1.10	0.98
SD	0.02	0.04	0.03
4 MONTH			
1 MONTH	2		Annu
Sample 1	2mm	3mm	4mm 1.09
2	1.25 1.24	1.11 1.09	1.09
		1.15	1.07
4	1.24	1.12	0.99
5	1.29	1.12	0.33
	1.26	1.12	1.04
Mean SD	1.26 0.02	0.02	0.05
30	0.02	0.02	0.03
3 MONTHS			
	2mm	3mm	4mm
Sample 1	1.31	1.37	1.18
2	1.31	1.18	1.10
			1.16
3	1.33	1.19	1.19
4 5	1.44 1.26		1.19
D	1.20	1.32	1.1
			1 15
Mean SD	1.36 0.08	1.23 0.11	1.15 0.04

Table 107. Fuji II LC Specific Heat Results (J/gC°)

FUJI IMMEDIAT		ific Heat (J/gC°)
	Marie Company of the Company of	Marie Park Rocal Lands	
<u>Sample</u>	<u>2mm</u>	3mm	4mm
1	1.03	0.94	0.84
2	1.02	0.96	0.78
3	1.02	0.94	0.83
4	1.06	0.89	0.86
5	1.05	0.92	0.82
Mean	1.04	0.93	0.83
SD	0.02	0.03	0.03
24 Hours		Maria Maria Santi	
<u>Sample</u>	<u>2mm</u>	3mm	4mm
1.00	1.21	0.98	0.89
2.00	1.20	0.95	0.93
3.00	1.19	1.00	0.92
4.00	1.23	0.97	0.92
5.00	1.21	0.97	0.91
Mean	1.21	0.97	0.91
SD	0.01	0.02	0.02
1 WEEK			
Sample	2mm	3mm	4mm
1.00	1.16	1.05	0.99
	1.16	1.10	1.09
2.00			
3.00	1.19	1.12	1.01
4.00	1.19	1.07	0.99
5.00	1.17	0.99	1.00
Mean	1.17	1.07	1.02
SD	0.02	0.05	0.04
1 MONTH			
Sample	2mm	3mm	4mm
1.00	1.24	0.98	0.92
	1.13	0.98	0.99
2.00			The second secon
3.00	1.22	0.95	0.97
4.00	1.03	0.92	0.97
5.00	1.18	1.00	0.95
Mean	4 40	0.97	0.96
SD	1.16		
DESCRIPTION OF THE PARTY OF THE	0.08	0.03	0.03
3 MONTHS	0.08		
3 MONTHS	0.08	0.03	0.03
Sample	0.08 S 2mm	0.03 3mm	0.03 4mm
<u>Sample</u> 1.00	0.08 <u>2mm</u> 1.16	0.03 3mm 0.96	0.03 4mm 0.98
1.00 2.00	0.08 2mm 1.16 0.77	0.03 3mm 0.96 1.02	0.03 4mm 0.98 1.00
1.00 2.00 3.00	0.08 2mm 1.16 0.77 1.01	0.03 3mm 0.96 1.02 1.01	0.03 4mm 0.98 1.00 0.87
1.00 2.00 3.00 4.00	0.08 2mm 1.16 0.77 1.01 1.16	0.03 3mm 0.96 1.02 1.01 0.96	0.03 4mm 0.98 1.00 0.87 0.99
3.00 4.00 5.00	0.08 2mm 1.16 0.77 1.01 1.16 1.14	0.03 3mm 0.96 1.02 1.01 0.96 0.99	0.03 4mm 0.98 1.00 0.87 0.99 0.99
1.00 2.00 3.00 4.00	0.08 2mm 1.16 0.77 1.01 1.16	0.03 3mm 0.96 1.02 1.01 0.96	0.03 4mm 0.98 1.00 0.87 0.99

Table 108. Photac-Fil Quick Specific Heat Results (J/gC°)

Photac-Fil Quick Specific Heat (J/gC°)

Photac-Fi		ecific Hea	t (J/gC°)
	MARKET STORY OF STREET		
Sample	2mm	3mm	4mm
1	1.07	0.86	0.82
2	1.05	0.92	0.82
3	1.05	0.97	0.83
4	1.13	0.98	0.87
5	1.08	0.90	0.81
Mean	1.08	0.93	0.83
SD	0.03	0.05	0.02
30			NAME OF THE PARTY
24 Hour			
	2mm	3mm	4mm
Sample	1.22	1.06	0.84
1			
2	1.13	1.16	0.89
3	1.21	1.10	0.94
4	1.15	1.11	0.96
5	1.22	1.12	0.81
Mean	1.19	1.11	0.89
SD	0.04	0.04	0.06
1 WEEK			
Sample	2mm	3mm	4mm
1	1.05	0.98	0.93
2	1.11	0.94	0.90
3	0.84	0.92	0.99
4	1.13	0.92	0.96
5	1.13	0.90	0.96
Mean	1.05	0.93	0.95
SD	0.12	0.03	0.03
1 MONTH			_
Sample	2mm	3mm	4mm
1	1.02	1.03	1.00
2	0.91	1.00	0.92
3	1.15	1.04	0.85
4	1.16	0.97	0.89
5	1.15	0.94	0.93
Mean	1.08	1.00	0.92
SD	0.11	0.04	0.06
3 MONTH	S	2000	
Sample	2mm	3mm	4mm
1	1.09	1.04	0.91
2	1.16	1.01	0.96
3	1.15	1.02	0.99
4		1.03	0.99
	1.17		
5	1.22	1.02	0.99
Mean	1.16	1.02	0.97
SD	0.05	0.01	0.03

APPENDIX D

THERMAL SCAN THERMAL ANALYSIS DATA

	AND DESCRIPTION OF THE PARTY OF			esults (J/	gm)
DSC SCA	N THERMA	L RESULT	S (J/gm)		
				W. 14 504	
N. BERTALL		distance of the	46 K. 2 196	ESSIE, In	
	D L. OC	2	Deak °C		Peak °C
and the first of the street, but	A SECURITION OF THE PROPERTY.	Account to the same of the sam			
					161.99 154.25
			The same of the sa		
	A CONTRACTOR OF THE PARTY OF TH		And the state of t		152.61 150.67
	经企业。1000年的企业企业企业企业				134.07
	STATE OF THE PARTY				150.72
	STATE OF THE PARTY	white the second	Name and Address of the Owner, where the Owner, which is the O		9.17
0.75	4.00	0.02	1.33	0.32	3.11
		Albert Albert			1212.350.3
2mm	Peak °C	3mm	Peak °C	4mm	Peak °C
-37.07	109.86	-46.69	131.93	-55.45	136.78
-32.23	111.53	-40.52	126.43	-50.05	136.31
-29.1	107.43	-51.13	122.41	-65.67	130.87
-16	103.84	-51.89	125.85	-51.9	132.45
-22.54	107.43	-49.53	147.43	-51.86	131.2
-27.39	108.02	47.95	130.81	-54.99	133.52
8.27	2.91	4.12	8.85	5.62	2.53
A THE SECTION				All or exc	A PERM
2mm	The second secon				Peak °C
			Annual State of the State of th		129.16
		Annual Control of the	Charles and Advanced to Control of the Control of t		133.28
					134.45
		AND THE PERSON NAMED IN COLUMN	A CONTRACTOR OF THE PARTY OF TH		133.62
	国际自己的				138.04
					133.71
6.51	3.25	6.42	5.79	8.51	2.84
	and the second	Pine Live	Marie Carrier	A CONTRACTOR OF THE PARTY OF TH	Hartingles of
2mm	Peak °C	3mm	Peak °C	4mm	Peak °C
		And the second of the second o			132.95
	ACCUMULATION OF THE PARTY OF TH		133.48	-129.96	131.92
	PROPERTY AND PROPE	-93.49	132	-92.74	132.16
	BUT TO RECEIVE THE PARTY OF THE	-101.74	132.7	-94.83	131
-52.34	114.24	-93.08	131.04	-85.27	129.4
49.72	114.63	-88.69	132.05	-103.45	131.49
4.26	2.65	10.39	0.96	16.39	1.21
all model is a	MARKET TO				
S	D 1.00	2	Dealesc		Dook °C
2mm	Peak °C	3mm	Peak °C	4mm	
2mm -102.52	122.7	-136.94	126.66	-130.69	126.83
2mm -102.52 -68.54	122.7 119.11	-136.94 -112.35	126.66 131.37	-130.69 -122.02	126.83 128.91
2mm -102.52 -68.54 -98.17	122.7 119.11 126.63	-136.94 -112.35 -104.86	126.66 131.37 131.02	-130.69 -122.02 -143.73	126.83 128.91 126.06
2mm -102.52 -68.54 -98.17 -96.78	122.7 119.11 126.63 119.91	-136.94 -112.35 -104.86 -79.61	126.66 131.37 131.02 129.56	-130.69 -122.02 -143.73 -57.48	126.83 128.91 126.06 133.62
2mm -102.52 -68.54 -98.17	122.7 119.11 126.63	-136.94 -112.35 -104.86	126.66 131.37 131.02	-130.69 -122.02 -143.73	128.91 126.06
	E 2mm -2.39 -4.14 -2.28 -3.22 -2.84 -2.97 0.75 2mm -37.07 -32.23 -29.1 -16 -22.54 -27.39 8.27 2mm -31.76 -45.33 -46.23 -45.97 -47.36 43.33 6.51 2mm -51.94 -45.7 -54.02 -44.59 -52.34 49.72	E 2mm	E 2mm	E	E 2mm

Table 110. Fuji II LC DSC Thermal Scan Results (J/gm)

FUJI II LC	DSC SCA	N THERMA	L RESULT	S		
IMMEDIAT	E	Ordina Ciferati	New Williams		ALINA-NA V	
Sample	2mm	Peak °C	3mm	Peak °C	4mm	Peak °C
1	-24.21	183.01	-33.75	174.06	-32.96	164.06
2	-25.49	183.77	-37.29	175.15	-38.49	161.60
3	-24.06	180.83	-33.73	172.40	-31.34	163.00
4	-2.25	82.65	-26.34	180.42	-35.74	164.23
5	-0.85	81.89	-36.56	171.47	-33.86	165.17
Mean	-15.37	142.43	-33.53	174.70	-34.48	163.61
SD	12.64	54.93	3.88	3.13	2.46	1.22
and the later				Service Labor		
24 hours						
Sample	2mm	Peak °C	3mm	Peak °C	4mm	Peak °C
1	-37.01	119.84	-23.58	128.37	-37.73	137.25
2	-38.71	117.39	-36.44	146.94	-40.58	143.85
3	-28.96	106.32	-40.21	151.56	-40.27	138.89
4	-41.75	142.42	-34.22	157.47	-32.74	143.35
5	-37.93	146.69	-38.80	149.04	-42.75	142.61
Mean	-36.87	126.53	-34.65	146.68	-38.81	141.19
SD	4.77	17.29	5.90	9.81	3.43	2.63
				HELLINE		
I WEEK						
Sample	2mm	Peak °C	3mm	Peak °C	4mm	Peak °C
1	-49.64	145.63	-81.90	145.78	-74.48	135.89
2	-51.97	115.01	-78.47	148.74	-74.86	136.86
3	-57.90	122.44	-78.15	147.26	-69.43	138.49
4	-52.77	131.67	-80.85	148.70	-72.22	139.33
5	-49.79	115.99	-80.76	145.95	-73.50	139.19
Mean	-52.41	126.15	-80.03	147.29	-72.90	137.95
SD	3.35	12.76	1.46	1.28	1.96	1.35
	TO THE STATE OF		NAME OF THE PARTY OF			
1 MONTH						
Sample	2mm	Peak °C	3mm	Peak °C	4mm	Peak °C
1	-30.90	128.24	-82.05	142.53	-70.14	136.68
2	-57.91	126.47	-74.32	144.23	-82.08	131.80
3	-41.79	117.13	-67.23	138.71	-84.02	142.05
4	-46.71	136.68	-77.15	144.90	-78.53	131.02
5	-53.28	134.36	-82.14	144.09	-76.62	136.58
Mean	46.12	128.58	-76.58	142.89	-78.28	135.63
SD	10.50	7.66	5.54	2.23	4.83	3.98
		SCHOOL SERVICE				A TOTAL IN
3 MONTH					MINISTER CENTER STATE OF	D 100
Sample	2mm	Peak °C	3mm	Peak °C	4mm	Peak °C
1	-54.69	136.13	-73.36	141.26	-86.27	141.12
2	-38.18	136.78	-72.26	140.06	-80.76	141.23
3	-42.03	143.60	-73.70	141.82	-73.57	136.44
4	-53.41	130.42	-79.98	142.50	-77.92	132.18
5	-49.97	117.44	-75.38	142.90	-89.20	143.32
Mean	47.66	132.87	-74.94	141.71	-81.54	138.86
SD	7.24	9.81	2.71	1.00	5.63	4.03

Table 111. Photac-Fil Quick Thermal Scan Results (J/gm)

Sample	E 2mm	Peak °C	3mm	Peak °C	4mm	Peak °C
1	-4.61	134.36	-29.49	193.12	-42.99	174.07
2	-3.44	115.69	-33.12	197.81	-42.66	176.14
3	-5.39	109.72	-30.45	196.44	-42.26	173.59
4	-3.97	59.76	-31.86	197.58	-48.59	173.60
5	-21.41	196.81	-45.49	186.72	-41.20	175.55
Mean	-7.76	123.27	-34.08	194.33	-43.54	174.59
SD	7.66	49.51	5.84	4.16	2.60	1.06
24 hours		Section 1		the section.	AL PROPERTY.	
Sample	2mm	Peak °C	3mm	Peak °C	4mm	Peak °C
1	-35.75	144.05	-65.98	158.38	-46.43	160.57
2	-32.19	162.00	-72.35	144.94	-53.12	153.53
3	-36.49	147.91	-69.73	152.55	-66.80	150.44
4	-37.70	111.14	-75.59	142.74	-69.45	153.07
5	-32.73	108.36	-70.84	148.22	-55.55	153.92
Mean	-34.97	134.69	-70.90	149.37	-58.27	154.31
SD	2.40	23.75	3.15	5.59	8.62	3.36
						Levil Do
WEEK						
Sample	2mm	Peak °C	3mm	Peak °C	4mm	Peak °C
1	-34.64	139.38	-72.18	145.72	-79.04	144.54
2	-36.54	134.54	-58.99	148.69	-72.21	143.77
3	-39.82	145.91	-69.66	145.94	-73.56	139.91
4	-64.64	143.77	-67.79	143.59	-67.58	144.99
. 5	-58.33	142.52	-69.98	142.74	-62.55	143.25
Mean	-46.79	141.22	-67.72	145.34	-70.99	143.29
SD	13.72	4.42	4.58	2.08	5.58	1.79
1 MONTH			1400	dillion tide	Lib big year	CONTRACTOR AND ADDRESS.
Sample	2mm	Peak °C	3mm	Peak °C	4mm	Peak °C
1	-48.94	139.90	-75.85	140.55	-81.85	141.00
	THE RESERVE OF THE PARTY OF THE				-81.60	142.65
/	-47 36	132 92	-/5 82	133 00		
2	-41.36 -41.01	132.92 138.24	-75.82 -74.89	139.66 139.66		
3	-41.01	138.24	-74.89	139.66	-74.78	140.31
3 4	-41.01 -53.79	138.24 138.00	-74.89 -70.63	139.66 140.19	-74.78 -74.29	140.31 140.05
3 4 5	-41.01 -53.79 -53.55	138.24 138.00 134.39	-74.89 -70.63 -72.97	139.66 140.19 140.62	-74.78 -74.29 -73.29	140.31 140.05 141.08
3	-41.01 -53.79	138.24 138.00	-74.89 -70.63	139.66 140.19	-74.78 -74.29	140.31 140.05
3 4 5 Mean	-41.01 -53.79 -53.55 -47.73	138.24 138.00 134.39 136.69	-74.89 -70.63 -72.97 - 74.03	139.66 140.19 140.62 140.14	-74.78 -74.29 -73.29 -77.16	140.31 140.05 141.08 141.02
3 4 5 Mean SD	-41.01 -53.79 -53.55 -47.73 6.28	138.24 138.00 134.39 136.69 2.91	-74.89 -70.63 -72.97 -74.03 2.00	139.66 140.19 140.62 140.14 0.42	-74.78 -74.29 -73.29 -77.16 3.76	140.31 140.05 141.08 141.02 0.91
3 4 5 Mean SD 3 MONTH Sample	-41.01 -53.79 -53.55 -47.73 6.28	138.24 138.00 134.39 136.69 2.91	-74.89 -70.63 -72.97 -74.03 2.00	139.66 140.19 140.62 140.14 0.42	-74.78 -74.29 -73.29 -77.16 3.76	140.31 140.05 141.08 141.02 0.91 Peak °C
3 4 5 Mean SD 3 MONTH Sample	-41.01 -53.79 -53.55 -47.73 6.28 S 2mm -62.47	138.24 138.00 134.39 136.69 2.91 Peak °C 125.41	-74.89 -70.63 -72.97 -74.03 2.00 3mm -89.83	139.66 140.19 140.62 140.14 0.42 Peak °C 137.58	-74.78 -74.29 -73.29 -77.16 3.76 4mm -99.04	140.31 140.05 141.08 141.02 0.91 Peak °C 139.36
3 4 5 Mean SD 3 MONTH Sample 1 2	-41.01 -53.79 -53.55 -47.73 6.28 S 2mm -62.47 -47.00	138.24 138.00 134.39 136.69 2.91 Peak °C 125.41 128.96	-74.89 -70.63 -72.97 -74.03 2.00 3mm -89.83 -90.57	139.66 140.19 140.62 140.14 0.42 Peak °C 137.58 136.56	-74.78 -74.29 -73.29 -77.16 3.76 4mm -99.04 -96.78	140.31 140.05 141.08 141.02 0.91 Peak °C 139.36 139.50
3 4 5 Mean SD 3 MONTH Sample 1 2 3	-41.01 -53.79 -53.55 -47.73 6.28 S 2mm -62.47 -47.00 -57.73	138.24 138.00 134.39 136.69 2.91 Peak °C 125.41 128.96 133.54	-74.89 -70.63 -72.97 -74.03 2.00 3mm -89.83 -90.57 -94.67	139.66 140.19 140.62 140.14 0.42 Peak °C 137.58 136.56 137.76	-74.78 -74.29 -73.29 -77.16 3.76 4mm -99.04 -96.78 -104.29	140.31 140.05 141.08 141.02 0.91 Peak °C 139.36 139.50 134.31
3 4 5 Mean SD 3 MONTH Sample 1 2 3 4	-41.01 -53.79 -53.55 -47.73 6.28 S 2mm -62.47 -47.00 -57.73 -66.28	138.24 138.00 134.39 136.69 2.91 Peak °C 125.41 128.96 133.54 133.65	-74.89 -70.63 -72.97 -74.03 2.00 3mm -89.83 -90.57 -94.67 -95.60	139.66 140.19 140.62 140.14 0.42 Peak °C 137.58 136.56 137.76 138.07	-74.78 -74.29 -73.29 -77.16 3.76 4mm -99.04 -96.78 -104.29 -100.01	140.31 140.05 141.08 141.02 0.91 Peak °C 139.36 139.50 134.31 136.68
3 4 5 Mean SD 3 MONTH Sample 1 2 3	-41.01 -53.79 -53.55 -47.73 6.28 S 2mm -62.47 -47.00 -57.73	138.24 138.00 134.39 136.69 2.91 Peak °C 125.41 128.96 133.54	-74.89 -70.63 -72.97 -74.03 2.00 3mm -89.83 -90.57 -94.67	139.66 140.19 140.62 140.14 0.42 Peak °C 137.58 136.56 137.76	-74.78 -74.29 -73.29 -77.16 3.76 4mm -99.04 -96.78 -104.29	140.31 140.05 141.08 141.02 0.91 Peak °C 139.36 139.50 134.31

APPENDIX E

THERMAL WEIGHT CHANGE DATA

Table 112. Vitremer Immediate Thermal Weight Data

im 40	µl pan										
		y illed + cure	Pre HC		Post HC 2		Post Scan	Total%loss			
wei	1 31.71 ght - pan	95.59 63.88	95.59 63.88	95.05 63.34 -0.54	94.67 62.96 -0.92	94.42 62.71 -1.17	91.4 59.69 -4.19			(mm)	stand dev
diff from p	original weight revious weight		0	-0.54	-0.92	-0.25	-3.02		Mean start weight	(mg) 72.30	9.16
	H2O loss ss during scan	%H2O loss				-1.17	4.73	6.56			
	2 31.56	95.45	95.45	94.64	94.2	93.91	90.39		Mean H20 gain	0.00	0.00
diff from (ght - pan original weight	63.89	63.89	63.08 -0.81	62.64 -1.25	62.35 -1.54	58.83 -5.06		Mean H2O loss	-1.08	0.19
diff from p	revious weight	%H2O loss	0	-0.81	-0.44	-0.29 -1.54	-3.52 2.41		Mean change HC1 Mean change HC2	-0.58 -0.40	0.14
% mat los	ss during scan						5.51	7.92	Mean change HC3	-0.26	0.02
wei	3 31.7 ght - pan	109.6 77.9	109.6 77.9	109.03 77.33	108.67 76.97	108.4 76.7	104.97 73.27		Mean change scan Mean % loss scan	-3.37 4.70	0.25 0.48
diff from	original weight revious weight		0	-0.57 -0.57	-0.93 -0.36	-1.2 -0.27	-4.63 -3.43		Mean % total loss	6.46	0.91
Est	h2O loss ss during scan	%H2O loss				-1.2	1.54	5.94			
26 mat 10	4 30.86	101.7	101.7	101.12	100.7	100.43	97.2				
wei	ght - pan	70.84	70.84	70.26	69.84	69.57	66.34				
diff from p	original weight revious weight		0	-0.58 -0.58	-0.42	-0.27	-3.23				
Est	h20 loss ss during scan	MHZU loss				-1-27	4.56	6.35			
	5 31.51	116.49	116.49	116.08	115.69	115.45	111.8				
diff from (ght - pan original weight	84.98	84.98	84.57 -0.41	84.18 -0.8	83.94 -1.04	80.29 -4.69				
diff from p	revious weight	%H2O loss	0	-0.41	-0.39	-0.24 -1.04	-3.66 1.22	1			
%	mat loss during s	can					4.30	5.52			
mm 100	D μl pan reference ref pan w		3.76	weight 7 N	lov	63.74					
		y illed + cure	Pre HC	Post HC 1	Post HC 2	Post HC 3	Post Scan				
	1 61.56 ght - pan	221.97 160.41	221.97 160.41	221.17 159.61	220.74 159.18	220.33 158.77	209.06 147.5				
diff from	original weight	100.41	0	-0.8 -0.8	-1.23 -0.43	-1.64 -0.41	-12.91 -11.27		Mean start weight	(mg) 154.07	stand dev 11.08
est	revious weight	%H2O loss		0.0	0.43	-1.54	1.02	0.0*	mdan start weight	134.07	17.00
% mat los	ss during scan			00:	-	900	7.03	8.05	Managaran	6.00	0.00
wei	2 61.39 ght - pan	221.9 160.51	221.9 160.51	221.08 159.69	220.55 159.16	220.14 158.75	208.89 147.5		Mean H20 gain Mean H2O loss	0.00 -1.26	0.00
diff from	original weight revious weight		0	-0.82 -0.82	-1.35 -0.53	-1.76 -0.41	-13.01 -11.25		Mean change HC1	-0.65	0.21
Est	h2O loss ss during scan	%H2O lose				-1.76	7.01	8.11	Mean change HC2 Mean change HC3	-0.49 -0.33	0.09 0.14
	3 61.38	199.2	199.2	198.5	198.08	197.63	188.81		Mean change scan	-9.66	1.50
diff from	ght - pan original weight	137.82	137.82	137.12	136.7	136.25 -1.57	127.43		Mean % loss scan Mean % total loss	6.26 7.22	0.79 0.95
diff from p	revious weight	%H2O loss	0	-0.7	-0.42	-0.45	-8.82				
% mat lo	ss during scan						6.40	7.54			
	4 62.05	209.5	209.5	209.19	208.57	208.29	200.28				
diff from	ght - pan original weight	147.45	147.45	147.14 -0.31	146.52 -0.93	146.24 -1.21	138.23 -9.22				
E 91	h2O lass	%H2O loss	0	-0.31	-0.62	-0.28 -1.21	-8.01 0.82				
% mat lo	ss during scan						5.43	6.25			
wei	5 61.72 ght - pan	225.86 164.14	225.86 164.14	225.23 163.51	224.8 163.08	224.68 162.96	215.73 154.01				
diff from	original weight previous weight	CARLES SANS	0	-0.63 -0.63	-1.06 -0.43	-1.18 -0.12	-10.13 -8.95			-	
Es	h20 loss mat loss during s	%H2O loss				-1.18	0.72 5.45	6.17			
mm %	acee duining 8	1									
	pan empt	y illed + cure	Pre HC 288.33	Post HC 1	Post HC 2 286.99	Post HC 3 286.31	Post Scan 272.03	1			
wei	ght - pan	288.33 227.1	227.1	226.27	225.76	225.08	210.8			(ma)	stand dev
diff from	original weight previous weight		0	-0.83 -0.83	-1.34 -0.51	-2.02 -0.68	-16.3 -14.20		Mean start weight	(mg) 239.91	10.96
est	H2O loss ss during scan	%H2O loss				-2.02	6.29	7.18			
	2 61.71	296.14	296.14	295.52	295.11	294.78	281.11		Mean H20 gain	0.00	0.00
diff from	ght - pan	234.43	234.43	233.81 -0.62	233.4 -1.03	233.07 -1.36	219.4 -15.03		Mean H2O loss	-1.38	0.26
diff from p	original weight previous weight h2O loss	%H2O loss	0	-0.62	-0.41	-0.33 -1.36	-13.67 0.58		Mean change HC1 Mean change HC2	-0.72 -0.46	0.08 0.05
	ss during scan						5.83	6.41	Mean change HC3	-0.42	0.15
	3 61.35 ght - pan	298.48 237.13	298.48 237.13	297.7 236.35	297.18 235.83	296.78 235.43	281.97 220.62		Mean change scan Mean % loss scan	6.39	1.59 0.41
diff from	original weight	201.113	0	-0.78 -0.78	-1.3 -0.52	-1.7	-16.51 -14.81		Mean % total loss	7.07	0.41
E 81	h20 loss	%H2O loss		0.70		-1.7	0.72 6.25	6.96			
% mat lo	ss during scan		617.6	210.00	310 ~	316.02		6.30			
we	4 61.76 ight - pan	317.52 255.76	317.52 255.76	316.83 255.07	316.39 254.63	254.26	298.7 236.94				
diff from	original weight previous weight		0	-0.69 -0.69	-1.13	-1.5 -0.37	-18.82 -17.32				
₹ st	h2O loss ss during scan	%H2O loss	-			-1.5	0.59 6.77	7.36			
% mat in			306.41	305.73	305.3	304.97	288.23				
% mat lo	5 613	306.41									
we	5 61.3 ight - pan original weight	306.41 245.11	245.11	244.43 -0.68	-1.11	243.67 -1.44 -0.33	226.93 -18.18				

Table 113. Vitremer 24-Hour Thermal Weight Data

nm 40 µl pan												
	empty i	lled + cure	Contract to the	Pre HC	Post HC 1	Post HC 2	Post HC 3	Post Scan	Total%los	s Total/gain%	Scan/gain%	Scan/RetH2
1 31 weight - pan	.72	105.76 74.04		106.67 74.95	105.49 73.77	104.96 73.24	104.56 72.84	100.56 68.84				
diff from original weigh	ıt	74.04		0.91	-0.27	-0.8	-1.2	-5.20				
diff from previous weigi	ht	4H2O loss		0.91	-1.18	-0.53	-0.4 -2.11	-4.00 2.82				
% mat loss during sca	in						-	5.34	8.15	671.43	439.56	-333.33
2 31 weight - pan	.35	103.34 71.99	PROPERTY OF THE PARTY OF	104.39 73.04	103.01 71.66	102.39 71.04	101.98 70.63	97.90 66.55				
diff from original weigh		7 1.00	-	1.05	-0.33 -1.38	-0.95 -0.62	-1.36 -0.41	-5.44 -4.08				
diff from previous weig		6H2O loss	net many	1.05	-1.30	-0.02	-2.41	3.30			388.57	200 00
% mat loss during sca	in							5.59	8.89	618.10	388.57	300.00
3 31 weight - pan	.16	106.38 75.22	MANUFACTURE	107.16 76	105.88 74.72	105.25 74.09	104.93 73.77	101.06 69.90				
diff from original weigh		10.41		0.78	-0.5 -1.28	-1.13 -0.63	-1.45 -0.32	-5.32 -3.87				
diff from previous weight		6H2O loss		0.70	-1,20	-0.65	-2.23	2.93		702.05	496.15	-266.90
% mat loss during sca	in			-				5.09	8.03	782.05	496.13	200.90
4 31 weight - pan	.29	110.46 79.17		111.31 60.02	110.34 79.05	109.76 78.47	109.37 78.08	105.76 74.47				
diff from original weigh	it			0.85 0.85	-0.12 -0.97	-0.7 -0.58	-1.09 -0.39	-4.70 -3.61				
diff from previous weight		6H2O loss		0.05		0.00	-1.94	2.42	6.94	652.94	424.71	-331.19
% mat loss during sca	in							4.51	0.54	032.54	424.71	331.13
5 31 weight - pan	.47	107.65 76.18	1 1200	108.66 77.19	107.56 76.09	107.22 75.75	106.8 75.33	103.46 71.99				
diff from original weigh diff from previous weig	ht ht		Name and Address of the Owner, where the Owner, which is the Ow	1.01	-0.09 -1.1	-0.43 -0.34	-0.85 -0.42	-4.19 -3.34				
Est h20 loss		6H2O loss					-1.86	2.41 4.33	6.74	514.85	330.69	392.94
% mat loss dur	-							9.33	0.74	314.03	330.03	352.54
nm 100 µl pan refer	ence w		63.76 63.74									
		lled + cure		Pre HC	Post HC 1	Post HC 2	Post HC 3	Post Scan				
1 61	.55	215.54		216.7	215.47	214.8	214.22	204.75 143.2	and the scare			
weight - pan diff from original weigh	nt	153.99		155.15	153.92 -0.07	153.25 -0.74	152.67 -1.32	-10.79				
diff from previous weig	ht	6H2O loss	SERVICE STATE	1.16	-1.23	-0.67	-0.58 -2.48	-9.47 1.60				
% mat loss during sca	in							6.10	7.70	1030.17	816.38	717.42
	1.14	228.66		229.81	228.82	228.3	227.89	218.68				
weight - pan diff from original weigh	nt	167.52		168.67	167.68 0.16	167.16 -0.36	166.75 -0.77	157.54 -9.98				
diff from previous weig	ht	6HZO Inss		1.15	-0.99	-0.52	-0.41 -1.92	-9.21 1 14				
% mat loss during sca	in							5.46	6.60	967.83	800.87	-1196.1
	1.54	213.32	ener make a	214.48 152.94	213.32 151.78	212.73 151.19	211.94 150.4	201.8 140.26				
weight - pan diff from original weigh	nt	151.78	10.00	1.16	0	-0.59	-1.38	-11.52				
diff from previous weig	ht	6H2O loss	CONTRACTOR OF THE PARTY OF THE	1.16	-1.16	-0.59	-0.79 -2.54	-10.14 1.66				
% mat loss during sca	an							6.63	8.29	1093.10	874.14	-734.78
	1.47	224.31 162.84	2010/10/00/00	225.73 164.26	224.72 163.25	224.11 162.64	223.63 162.16	212.99 151.52				
weight - pan diff fro <mark>m original weigh</mark>	nt	102.04		1.42	0.41	-0.2 -0.61	-0.68	-11.32 -10.64		1		
diff from previous weig	ht	%H2O loss		1.42	-1.01	-0.61	-0.48 -2.1	1.28			240.00	4504.7
% mat loss during sca	an							6.48	7.76	897.18	749.30	-1564.7
5 6 weight - pan	1.44	221.84 160.4		223.28 161.84	222.33 160.89	221.82 160.38	221.38 159.94	211.47 150.03				
diff from original weigh		160.4		1.44	0.49	-0.02	-0.46 -0.44	-10.37 -9.91				
diff from previous weig	ht	6H2O loss		1.44	-0.95	-0.51	-0.44	1.17				
% mat loss dur	ing sca	n						6.12	7.30	820.14	688.19	-2154.3
nm												
pan	empty	lled + cure						Post Scan				
1 6 weight - pan	1.56	310.07 248.51		311.68 250.12	249.22	310.13 248.57	248.06	292.82 231.26				
diff from original weigh	nt Int		and the same of	1.61	0.71	-0.65	-0.45 -0.51	-17.25 -16.8				
diff from previous weig est H2O loss		%H2O loss					-2.06	0.82	7.54	1171.43	1043.48	3733.3
% mat loss during sca				240.00	245.00	316.00	214.72	200 20				
2 6 weight - pan	1.21	314.71 253.5	na different	316.92 255.71	315.82 254.61	315.23 254.02	314.72 253.51	298.28 237.07 -16.43				
diff from original weigh	ht			2.21	1.11	0.52 -0.59	0.01 -0.51	-16.43 -16.44				
diff from previous weig Est h20 loss % mat loss during sca		%H2O loss					-2.2	0.86	7.29	843.44	743.89	164400.
		207.5		200.01	308.22	307.67	307.2	289.72		ļ		
weight - pan	2.03	307.54 245.51		309.21 247.18	246.19	245.64	245.17	227.69				
diff from original weigh	ht		William or comme	1.67	0.68	0.13 -0.55	-0.34 -0.47	-17.82 -17.48				
diff from previous weig		%H2O loss					-2.01	7.07	7.88	1167.07	1046.71	-5141.1
% mat loss during sc					220 70	330 13	220 4	321.96				
weight - pan	1.57	338.92 277.35		341.01 279.44	339.78 278.21	339.13 277.56	338.4 276.83	260.39				
diff from original weigh	ht	-	MARKET COMPA	2.09	0.86	0.21 -0.65	-0.52 -0.73	-16.96 -16.44				
diff from previous weig		%H2O loss					-2.61	0.93 5.88	6.82	911.48	786.60	3161.5
% mat loss during sc				1			205.01			4		
5 6 weight - pan	1.76	335.77 274.01	2	338.05 276.29	336.62 274.86	336.07 274.31	335.61 273.85	319.28 257.52				
	ht		MOST	2.28 2.28	0.85	0.3 -0.55	-0.16 -0.46	-16.49 -16.33				
diff from original weightiff from previous weigh	Name and Address of the Owner, where											

Table 114. Vitremer One-Week Thermal Weight Data

nm	40 µl pan			32.99									
		ref pan we	ight	33.09	29-Nov-05								
	1	pan empty 32.83	filled + cure 106.27		Pre HC 107.85	Post HC 1 107.51	Post HC 2 106.56	Post HC 3 106.12	Post Scan	Total%loss	Total/Gain%	Scan/Gain%	Scan/RetH2
	weight - par		73.44		75.02	74.68	73.73	73.29	68.38				
	om previous		NOTE STATE		1.58	1.24 -0.34	-0.95	-0.15 -0.44	-5.06 -4.91				
	est H2O los	В	%H2O loss					-1.73	2.31 6.54	8.85	420.25	310.76	3273.33
76 ma	at loss durin									0.00	TEURO	0.00.0	
	weight - par	32.92	116.45 83.53	district.	118.36 85.44	117.38 84.46	116.85 83.93	116.35 83.43	110.67 77.75				
diff fr	om original	veight		and the same	1.91	0.93 -0.98	-0.53	-0.1 -0.5	-5.78 -5.68				
	om previous Est h20 los	В	%H2O loss			5.55		-2.01	2.35	9.00	402.62	297.38	5680.00
% ma	at loss durin	scan							6.65	9.00	402.62	297.30	-3600.00
	3 weight - par	32.44	110.06 77.62	THE RESIDENCE OF	111.85 79.41	110.86 78.42	110.37 77.93	109.88 77.44	104.96 72.52				
diff fr	om original	veight	77.02		1.79	0.8	0.31	-0.18	-5.1 -4.92				
diff fro	m previous Est h20 los	weight	%H2O loss		1.79	-0.99	-0.49	-0.49 -1.97	2.48				
% ma	at loss durin	scan							6.20	8.68	384.92	274.86	-2733.33
	4	32.71	99.12		100.73	99.57	99.24	99.11	95.09				
diff fr	weight - par om original	veight	66.41		68.02 1.61	66.86 0.45	66.53 0.12	66.4 -0.01	62.38 -4.03				
diff fro	om previous	weight	SELECT INCO	- Des	1.61	-1.16	-0.33	-0.13	-4.02 2.38				
% ma	at loss durin	scan		,					5.91	8.29	350.31	249.69	40200.0
	5	32.68	110.21		112.31	111.04	110.3	110.04	105.89				
diff 6-	weight - par om original		77.53		79.63 2.1	78.36 0.83	77.62 0.09	77.36 -0.17	73.21 -4.32				
diff fro	om previous	weight			2.1	-1.27	-0.74	-0.26	-4.15				
	% mat loss	ss during sc	an					-2.21	5.21	8.06	305.71	197.62	-2441.1E
nm	100 µl pan			63.88									
	ιω μι pan	ref pan we		63.89		29-Nov-05							
		pan empty	filled + cure		Pre HC				Post Scan				
	1	63.53	198.8	and the same of	202.56 139.03	200.91	199.95 136.42	199.16 135.63	189.13 125.6				
diff fr	weight - par om original	veight	135.27		3.76	137.38 2.11 -1.65	1.15	0.36	-9.67				
diff fro	om previous	weight	%H2O loss		3.76	-1.65	-0.96	-0.79	-10.03 2.45				
% ma	at loss durin	scan							7.21	9.66	357.18	266.76	2786.11
	2	61.17	220.32		224.38	222.6	221.67	221.24	209.19				
diff fr	weight - par om original	veight	159.15		163.21 4.06	161.43 2.28	160.5 1.35	160.07	148.02 -11.13				
	om previous			1/2013	4.06	-1.78	-0.93	-0.43	-12.05				
% ma	at loss durin	scan	WH2U 1088						7.38	9.31	374.14	296.80	1309.78
	3	61.9	225.2		229.57	227.83	227.04	226.18	213.79				
	weight - par		163.3	10000	167.67 4.37	165.93	165.14 1.84	164.28 0.98	151.89 -11.41				
diff fro	om original om previous	weight			4.37	2.63	-0.79	-0.86	-12.39				
	et h20 los at loss durin	5	%H20 loss					-3.39	7.39	9.41	361.10	283.52	1264.29
70 1111			220.0		233.2	231.81	230.98	230.18	218.21				
	weight - par	61.55	229.2 167.65		171.65	170.26	169.43	168.63	156.66				
diff fr	om previous	weight		-	4	2.61	1.78	0.98	-10.99 -11.97				
	Est h20 los	8	%H2O loss					-3.02	6.97	8.73	374.75	299.25	1221.43
% ma	at loss durin	gscan								0.75	3,4.73	EUUIEU	
	5	61.38	213.55 152.17	ALC: MINISTER	217.2 155.82	215.85 154.47	214.78 153.4	214.09 152.71	203.83 142.45				
diff fr	weight - par rom original	weight	102.11		3.65	2.3	1.23	0.54	-9.72 -10.26				
diff fro	Est h20 lo	weight ss	%1-20 loss	I SALE	3.65		-1.07	-0.69	2.00		1		and the same of th
	% mat los		an						6.58	8.58	366.30	281.10	1900.00
mm													
		pan ampty	filled + cure		Pre HC	Post HC 1	Post HC 2	Post HC 3	Post Scan				
			313.08	NAME OF TAXABLE PARTY.	315.37	314.23 253.15	313.6	313.03 251.95	295.15 234.07				
diff fr	weight - par rom original	weight	252		254.29 2.29	1.15	252.52 0.52	-0.05	-17.93				
diff fr	om previous	weight	%H2O foss		2.29	-1,14	-0.63	-0.57 -2.34	-17.88 0.92				
% m	at loss durin	g scan							7.03	7.95	882.97	780.79	-35760.0
	2	61.34	292.25		294.69	293.12	292.44	291.94	272.8				
die c	weight - par	1	230.91	and the second	233.35 2.44	231.78 0.87	231.1 0.19	230.6 -0.31	211.46 -19.45				
diff fr	om previous	weight	WH2O loss		2.44	-1.57	-0.68	-0.5 -2.75	-19.14 1,18				
	at loss durin		MITZO 1098			, , , , ,			8.20	9.38	897.13	784.43	-6174.19
	3	61.67	290.07		292.32	291.33	290.77	290.26	273.34				
	weight - pa	1	228.4		230.65	229.66	229.1	228.59 0.19	211.67 -16.73				
diff fr	rom original om previous	weight weight			2.25	1.26	-0.56	-0.51	-16.73		-		
	Est h20 las	8	%H2O loss					-2.06	7.34	8.23	843.56	752.00	8905.26
76 m	at loss durin				246.51	242.55	217.07	217.45	300.36				-
	4 weight - pa	61.69	316.89 255.2		319.54 257.85	318.55 256.86	317.87 256.18	317.43 255.74	238.67				
diff fi	rom original	weight	and the latest designation of the latest des	-	2.65 2.65	1.66 -0.99	0.98	0.54	-16.53 -17.07				
	om previous Est h2O los	S	%H2O loss	W. Baltiman	2.05	0.55		-2.11	0.82	7	722 77	644.15	3161.11
% m	at loss durin	g scan							6.62	7.44	723.77	044.13	3101.11
	5	61.57	291.23		293.72	292.62 231.06	291.98 230.41	291.45 229.88	275.22 213.65				
diff fi	weight - pa	weight	229.66	Maria I	232.15 2.49	1.39	0.75	0.22	-16.01				
	om previous	weight	%H2O loss		2.49	-1.1	-0.64	-0.53	-16.23 0.98	Se			

Table 115. Vitremer One Month Thermal Weight Data

nm	40 µl pan	ref pan we ref pan we		32.98 33.03	21-Nov-05								
			filled + cure	Scientific Control	Pre HC					Total%loss	Total/Gain%	Total/Scan%	Scan/RetH20
	1 weight - pa	32.58	109.54 76.96	CONTRACTOR OF THE PARTY OF THE	110.52 77.94	109.78 77.2	109.4 76.82	109.11 76.53	104.73 72.15		*		
diff fr	rom original	weight	70.50		0.98	0.24	-0.14	-0.43	-4.81				
diff fro	om previous	weight	%HZQ loss	-	0.98	-0.74	-0.38	-0.29	-4.38 1.81				
% ma	at loss durin	g scan							1.81 5.62	7.43	590.82	446.94	-1018.60
	2	32.63	107.79		113.84	113.11	112.68	112.44	107.48				
	weight - pa	1	75.16		81.21	80.48	80.05	79.81	74.85				
diff fr	ro <mark>m original</mark> om previous	weight	SHEET COMMEN	-	6.05 6.05	6.32 -0.73	4.89	4.65 -0.24	-0.31 -4.96				
	Est h20 los	S	%HZO loss					-1.4	1.72	7.03	405 43	04.00	106.67
% ma	at loss durin	g scan							6.11	7.83	105.12	81.98	100.07
	3	33.02	111.3		112.31	111.6	111.1	110.95	106.18				
diff fr	weight - pa	n weight	78.28		79.29 1.01	78.58 0.3	78.08 -0.2	77.93 -0.35	73.16 -5.12				
diff fro	rom original om previous	weight			1.01	-0.71	-0.5	-0.15	-4.77				
	et loss durin	8	%H2O loss					-1.36	6.02	7.73	606.93	472.28	-1362.86
70 1111													
	4 weight - pa	32.8	100.38 67.58		101.2 68.4	100.38 67.58	99.98 67.18	99.77 66.97	95.54 62.74				
diff fr	rom original	weight			0.82	0	-0.4	-0.61	-4.84				
diff fro	om previous	weight	%H2O loss		0.82	-0.82	-0.4	-0.21 -1.43	-4.23 2.09				
% ma	at loss durin	g scan							6.18	8.27	690.24	515.85	-693.44
	5	32.9	97.15		97.86	97.16	97.16	96.54	92.29				
	weight . na	THE RESERVE OF THE PERSON NAMED IN	64.25	A MADE	64.95	64.26	64.26	63.64	59.39				
diff fr	om previous Est h20 k	weight	Name and Address of the Owner, where	-	0.7	0.01 -0.69	0.01	-0.61 -0.62	-4.86 -4.25				
Gill if	Est h20 k	ss	%H2O loss					-1.31	2.02				
	% mat los	s during sc	an						6.54	8.56	794.29	607.14	-696.72
nm	100 µl pan	reference	M	63.67									
		ref pan we		63.73		07-Dec-05							
		pan empty	filled + cure		Pre HC	Post HC 1			Post Scan				
	1	63.79	205.2	war transfer	209.29	207.5	206.28 142.49	205.69 141.9	195.36 131.57				
diff fr	weight - pa rom original		141.41		145.5 4.09	143.71	1.08	0.49	-9.84				
	om previous			-	4.09	-1.79	-1.22	-0.59	-10.33				
% ms	est H2O los at loss durin	g scan	%H2O loss			,		-3.6	7.10	9.57	340.59	252.57	2108.16
					000 77	201.7	204.42	200.66	189.94				
	2 weight - pa	63.61	199.35 135.74	AND VALUE OF	202.77 139.16	201.7 138.09	201.13 137.52	137.05	126.33				
diff fr	rom original	weight			3.42	2.35	1.78	1.31	-9.41				
diff fro	om previous	weight	%HZO loss	to the second	3.42	-1.07	-0.57	-0.47 -2.11	-10.72 1.52				
% ma	at loss durin	g scan							7.70	9.22	375.15	313.45	818.32
	3	64.04	210.5		214.49	213.27	212.71	212.04	199.19				
-	weight - pa	n	146.46	170,000	150.45	149.23	148.67	148	135.15				
diff fr	ro <mark>m original</mark> om previous	weight	NAME OF TAXABLE PARTY.	THE REAL PROPERTY.	3,99	2.77	2.21 -0.56	-0.67	-11.31 -12.85				
	Est h20 los	S	%H20 loss					-2.45	1.63	10.17	383.46	322.06	834.42
% ma	at loss durin	g scan							8.54	10.17	303.40	322.00	034.42
	4	63.75	218.63		222.94	221.9	221.04	220.5	206.8				
diff fr	weight - pa		154.88		159.19 4.31	158.15 3.27	157.29 2.41	156.75 1.87	143.05 -11.83				
diff fro	om previous		NAME OF STREET	1000	4.31	-1.04	-0.86	-0.54	-13.7				
	at loss durin	S	%H2O loss					-2.44	8.61	10.14	374.48	317.87	732.62
70 1110													
	5 weight - pa	63.76	206.35 142.59		210.67 146.91	209.5 145.74	208.76 145	207.95 144.19	195.58 131.82				
diff fr	rom original	weight	142.00		4.32	3.15	2.41	1.6	-10.77				
diff fro	em previous Est h20 li	weight	%HQO loss	OF SEALING	4.32	-1.17	-0.74	-0.81 -2.72	-12.37 1.85				
	% mat los	s during so							8.42	10.27	349.31	286.34	773.13
nm													
					_								
	1	pan empt	7 filled + cure 280.74		Pre HC 292.7	Post HC 1 287.37	286.33		Post Scan 263.16				
	weight - pa	n	216.65	a the star	228.61	223,28	222.24	221.51	199.07				
diff fr	rom original	weight		-	11.96	6.63 -5.33	5.59	4.86 -0.73	-17.58 -22.44				
	est H2O los	6	%H20 loss					-7.1	3.11			407.00	604 77
% m	at loss durin	g scan							9.82	12.92	246.99	187.63	461.73
	2	63.9	293.53		299.2								
4147.4	weight - pa	n	229.63		235.3 5.67	-63.9 -293.53	-63.9 -293.53	-63.9 -293.53	-63.9 -293.53				
om fr	rom original om previous	weight	THE REAL PROPERTY.		5.67	-299.2	0	0	0				
din iti	Est h20 les	5	%H2O loss					-299.2	#DIV/DI				
	at ioss durin												
	3	63.04	290.6	The state of	295.46 232.42	294.4 231.36	293.62 230.58	293.06 230.02	275.01 211.97				
			227.56		4.86	3.8	3.02	2.46	-15.59				
% m	weight - pa	weight	WHOO loss	-	4.86	-1.06	-0.78	-0.56	-18.05				
% m	weight - pa rom original rom previous		merrary loss			-			7.77	8.80	420.78	371.40	733.74
% m	weight - pa rom original rom previous Est h2O lo	8				298.97	298.22	297.46	278.77	-			
% m	weight - pa rom original rom previous Est N2O lo at loss durin	g scan	2015				298.22	233.41	214.72				
% m	weight - pa rom original rom previous Est A2O lo at loss durin	g scan 64.05	294.97 230.92		302.92 238.87	234.92		2.49	-16.2				
% maddiff from the made of the	weight - pa rom original rom previous Est N2O lo sat loss durir 4 weight - pa rom original	64.05 n weight	294.97 230.92		238.87 7.95	234.92	3.25						
% maddiff from the made of the	weight - pa rom original rom previous Est N2O lo sat loss durir 4 weight - pa rom original	64.05 n weight			238.87			-0.76 -5.46	-18.69 2.29				
% made diff from the diff from	weight - pa rom original om previous Est 620 lo at loss durin 4 weight - pa	64.05 n weight weight	230.92		238.87 7.95	4	3.25			10.11	303.77	235.09	750.60
% made diff from the diff from	weight - pa rom original com previous Est N20 in at loss durin 4 weight - pa rom original com previous Est N20 io at loss durin	64.05 n weight weight s	230.92 %H2O loss		238.87 7.95 7.95	-3.95	3.25 -0.75		-18.69 2.29	10.11	303.77	235.09	750.60
% made of the diff for diff for % made of the	weight - parom original or previous Est h20 lo at loss durin 4 weight - parom original om previous Est h20 lo at loss durin 5 weight - pa	64.05 n weight weight sg scan 64.33	230.92		238.87 7.95 7.95 315.38 261.05	313.38 249.05	3.25 -0.75 312.42 248.09	-0.76 -5.46 311.7 247.37	-18 69 2 29 7 82 292 03 227.7	10.11	303.77	235.09	750.60
% maddiff from the diff from t	weight - parom original com previous Est N20 lo at loss durin 4 weight - parom original com previous Est N20 lo at loss durin 4 loss durin 15	64.05 n weight weight g scan 64.33 n	230.92 5-H2O loss 309.24		238.87 7.95 7.95 315.38	313.38	3.25 -0.75 312.42	-0.76 -5.46	-18.69 2.29 7.82 292.03	10.11	303.77	235.09	750.60

Table 116. Vitremer Three-Month Thermal Weight Data

mm		ref pan we	ight	31.28									
		ref pan we	ght	31.3									NAME OF TAXABLE PARTY.
	1	pan empty 31.67	filled + cure 153.44		Pre HC 161.5	Post HC 1 158.15	157.1	156.11	143.83	Total%loss	Total/Gain%	Scan/Gain%	Scan/RetH2
eliff for	weight - pa om original	n	121.77		129.83 8.06	126.48 4.71	125.43 3.66	124.44 2.67	112.16 -9.61				
diff fro	m previous	weight		DEED SA	8.06	-3.35	-1.05	-0.99	-12.28				
% ma	t loss durin	g scan	%H2O loss					-5.33	9.46	13.61	219.23	152.36	459.93
	2	31.52	136.95		141.74	138.6	137.75	137.21	129.19				
diff for	weight - pa	n	105.43	di ser	110.22 4.79	107.08	106.23 0.8	105.69 0.26	97.67 -7.76				
diff fro	m original m previous	weight	SUPPO loss		4.79	-3.14	-0.85	-0.54	-8.02				
	t loss durin	g scan	20/12 Cilusa						7.28	11.39	262.00	167.43	3084.62
	3	31.9	136.82		140.44	139.01	138.27	137.62	128.41				
diff fro	weight - pa om original	n weight	104.92		108.54 3.62	107.11 2.19	106.37 1.45 -0.74	105.72 0.8	96.51 -8.41				
diff fro	m previous	weight	SCHOOL IONS		3.62	-1.43	-0.74	-0.65 -2.82	-9.21 2.60				
	t loss durin		20112 CHUSS						8.49	11.08	332.32	254.42	1151.25
	4	31.65	129.98		134.8	132.46	131.6	130.79	120.94				
diff fro	weight - pa m original	n weight	96.33	Mich.	103.15 4.82	100.81 2.48	99.95 1.62	99.14 0.81	89.29 -9.04				
diff fro	m previous	weight	WHO love		4.82	-2.34	-0.86	-0.81	-9.85 3.89				
	t loss durin		ACLE CO IOSS						9.55	13.44	287.55	204.36	1216.05
	5	32.09	143.63		145.82	144.67	143.93	143.38	133.78				
diff fro	weight - pa om original	weight	111.54		113.73 2.19	112.58	111.84 0.3	111.29 -0.25	101.69 -9.85				
diff fro	m previous Est h20 la	weight	%H2O loss		2.19	-1.15	-0.74	-0.55 -2.44	-9.6 2.15				
	% mat los	s during sc							8.44	10.59	549.77	438.36	-3840.00
mm	100 µl par	reference		63.76		20.0							
		ref pan we		63.78		28-Dec-05							
	1	pan empty 64.16	filled + cure 204.94		Pre HC 214.03	Post HC 1 211.1	Post HC 2 210.25	Post HC 3 209.37	Post Scan 191.65				
1140.6	weight - pa	n	140.78	Milen	149.87	146.94	146.09 5.31	145.21 4.43	127.49 -13.29				
diff fro	m original m previous	weight			9.09	6.16 -2.93	-0.85	-0.88	-17.72				
	t loss durin	8	%H2O loss					-4.66	11.82	14.93	246.20	194.94	400.00
	2	63.55	215.83		226.34	221	220.09	219.23	205.18				
.1147.6	weight - pa	n	152.28	St. M.	162.79	157.45	156.54	155.68 3.4	141.63 -10.65				
	om original on previous				10.51	5.17 -5.34	4.26 -0.91	-0.86	-14.05				
	t loss durin	8	%H2O loss					-7.11	8.63	13.00	201.33	133.68	413.24
	3	63.43	217.2		226.12	222.63	221.67	220.67	206.42				
	weight - pa	n	153.77	Sec. 1	162.69 8.92	159.2 5.43	158.24 4.47	157.24 3.47	142,99 -10.78				
	om original om previous		- 69		8.92	-3.49	-0.96		-14.25				
% ma	st h2O los it loss durin	s g scan	%H2O less					-5.45	8.76	12.11	220.85	159.75	410.66
	4	64.33	226.45		229.7	228.42	227.87	227.15	215.29				
	weight - pa	n	162.12		165.37 3.26	164.09 1.97	163.54 1.42	162.82 0.7	150.96				
diff fro	om original om previous	weight			3.25	-1.28	-0.55	-0.72	-11.86				
% ma	t loss durir	g scan	%H2O loss					-2.55	7.17	8.71	443.38	364.92	1694.29
	5	63.36	209.71		215.95	211.8	211.24	210.56	199.47				
diff for	weight - pa	n	146.35	William.	152.59 6.24	148.44	147.88	147.2 0.85	136.11 -10.24				
diff fro	m previous	weight			6.24	-4.15	-0.56	-0.68	-11.09				
	% mat los	s during sc	%H2O loss an					-5.39	7.27	10.80	264.10	177.72	1304.71
l mm													
		nan amet	filled + cure		Pre HC	Post HC 1	Post HC 2	Post HC 3	Post Scan				
	1	63.32	264.83		274.47	272.53	271.84	271.31	249.8				
diff fr	weight - pa om original	n weight	201.51		211.15 9.64	209.21 7.7	208.52 7.01	207.99 6.48	186.48 -15.03				
diff fro	om original om previous est HZO los	weight	%H2O loss	Calabara	9.64	-1.94	-0.69	-0.63 -3.16	-21.51 1.50				
% ma	at loss durin	ng scan							10.19	11.68	255.91	223.13	331.94
	2	63.15	265.23		270.75	269.31	268.74	268.31	248.02				
diff fr	weight - pa om original		202.08		207.6 5.52	206.16 4.08	205.59 3.51	205.16 3.08	184.87 -17.21				
	m previous		%HQQ loss	Santal and	5.52	-1.44	-0.57	-0.43 -2.44	-20.29 1.18				
% ma	at loss duri	ng scan							9.77	10.95	411.78	367.57	658.77
	3	64.36	253.87		270.2	261.45	260.63	259.94	238.98				
diff fr	weight - pa	weight	189.51		205.84 16.33	197.09 7.58	196.27 6.76	195.58 6.07	174.62 -14.89				
diff fro	m previous	weight	%H2O loss	A STREET	16.33	-8.75	-0.82	-0.69 -10.26	-20.96 4.98				
% ma	at loss duri	ng scan			-				10.18	15.17	191.18	128.35	345.30
	4	64.35	236.56		246.5	243.85	243.13	242.52	223.1				
diff fr	weight - pa om original		172.21		182.15 9.94	179.5 7.29 -2.65	178.78 6.57	178.17 5.96	158.75 -13.46				
diff fro	om previous	weight	%H2O loss		9.94	-2.65	-0.72	-0.61 -3.98	-19.42 2.19				
% ma	at loss duri	ng scan	ALL CONTRACTOR						10.66	12.85	235.41	195.37	325.84
		64.32	256.01		267.1	264.64	264.29	263.58	243.51				
	5				202.78	200.32	199.97	199.26	179.19				
	5 weight - pa om original	in	191.69		11.09	8.63 -2.46	8.28 -0.35	7.57	-12.5 -20.07				

Table 117. Fuji II LC Immediate Thermal Weight Data

ili II LC	Immediate	Thermal	Weight Data	Weights in	mg							
nm	40 µl pan											
		pan empty	filled + cure	Pre HC	Post HC 1	Post HC 2	Post HC 3	Post Scan	Total%loss			
	1	32.82	144.53	144.53	144.21	143.98	143.8	139.11				
diff fro	weight - par om original v	veight	111.71	111.71	-0.32	111.16 -0.55	110.98 -0.73	106.29 -5.42				and dev
diff fro	m previous	weight	% HOD love	0	-0.32	-0.23	-0.18	-4 69 0 65		Mean start weight	101.26	20.41
% ma	t loss durin	scan	ACIAL INS					0.65 4.20	4.85			
	2	33.21	160.97	160.97	160.54	160.38	160.15	154.78		Mean H20 gain	0.00	0.00
	weight - par		127.76	127.76	127.33	127.17	126.94	121.57		Mean H2O loss	-0.60	0.09
diff fro	om original v m pravious	weight		0	-0.43 -0.43	-0.59 -0.16	-0.82 -0.23	-6.19 -5.37		Mean change HC1	-0.35	0.08
E	Est h2O los	s	%H20 loss				-0.82	0.64	4 05	Mean change HC2 Mean change HC3	-0.17 -0.17	0.08
% ma	t loss durin	scan			-			4.20	4.00			
	3	32.83	136.21	136.21	135.98	135.71	135.62	131.24		Mean change scan Mean % loss scan	-3.97 3.86	1.21 0.50
diff fro	weight - par om original v	veight	103.38	103.38	103.15 -0.23	102.88 -0.5	102.79 -0.59	98.41		Mean % total loss	4.56	0.40
diff fro	m previous	weight		0	-0.23	-0.27	-0.09	-4.38				
% ma	t loss durin	scan	Delicas Iusa					4.24	4.81			
	4	32.91	108	108	107.6	107.53	107.35	104.74				
	weight - par		75.09	75.09	74.69	74.62	74.44	71.83				
diff fro	o <mark>m original v</mark> m previous	veight		0	-0.4	-0.47	-0.65 -0.18	-3.26 -2.61				-
E	Est h20 los	S	%HOU loss				-0.65	0.87				
% ma	t loss durin	scan						3.48 0.00	4.34			
	5	32.73	121.11	121.11	120.76 88.03	120.62 87.89	120.43 87.7	117.62			-	
diff fro	weight - par om original v	veight	88.38	88.38	-0.35	-0.49	-0.68	84.89 -3.49				
diff fro	m previous	weight	Carlo Mark	0	-0.35	-0.14	-0.19	-2.81				
	% mat los	s during so	an				J 68	3.18	3.95			
mm	100 µl pan			3.76								
mm	тоо рг рап	ref pan we		3.76								
			filled + cure	Pre HC	Post HC 1	Post HC 2	Post HC 3	Post Scan				
	1	60.8	231.59	231.59	231.28	231.03	230.79	223.45				
diff fro	weight - par om original v	n voight	170,79	170.79	170.48 -0.31	170.23 -0.56	169.99	162.65 -8.14			(mg) s	tand dev
diff fro	m previous	weight		Ö	-0.31	-0.25	-0.24	-7.34		Mean start weight	170.87	17.95
	esi H2O los it loss durin	s	%FQ0 lass				-0.8	0.47 4.30	4.77		-	
, ma					000 00	000.10	250.0			M U26	0.00	
-	2 weight - par	61.54	251.15 189.61	251.15 189.61	250.79 189.25	250.43 188.89	250.3 188.76	241.60 180.06		Mean H2O gain Mean H2O loss	0.00 -0.90	0.00
diff fro	om original v	weight		0	-0.36	-0.72	-0.85	-9.55				0.46
diff fro	m previous Est h20 los	weight	%H20 loss	0	-0.36	-0.36	-0.13 -0.86	-8.70 0.45		Mean change HC1 Mean change HC2	-0.54 -0.26	0.06
% ma	t loss durin	gscan						4,59	5.04	Mean change HC3	-0.21	0.06
	3	61.82	236.14	236.14	235.76	235.56	235.28	227.63		Mean change scan	-7.58	1.03
	weight - par	1	174.32	174.32	173.94	173.74	173.46	165.81		Mean % loss scan Mean % total loss	4.43 5.05	0.17
diff fro	om original v om previous	weight weight		0	-0.38 -0.38	-0.58 -0.2	-0.86 -0.28	-8.51 -7.65		mean % total loss	5.05	0.20
	Est h20 los	8	%H20 loss				-0.86	0.49	4.88			
% ma	t loss durin								4.08			
	4	61.56	202.91 141.36	202.91 141.35	201.54 139.98	201.29 139.73	201.12 139.56	195.12 133.56			-	
diff fro	weight - par om original	weight	141,00	0	-1.37	-1.62	-1.79	-7.79				
diff fro	m previous	weight	%H20 let	0	-1.37	-0.25	-0.17	-6.00 1.27				
% ma	t loss durin	gscan						4.24	5.51			
	5	61.68	239.95	239.95	239.65	239.42	239.2	230.98	-			+
-	weight - par	1	178.27	178.27	177.97	177.74	177.52	169.30				
diff fro	o <mark>m original v</mark> om previous	weight		0	-0.3 -0.3	-0.53 -0.23	-0.76 -0.22	-8.97 -8.22			+	
Sill HO	Est h20 to	SS	%HQB loss				-0.76	0.42				
	% mat los	s during so	an					4.61	5.03			
mm												
		pan empty	filled + cure	Pre HC	Post HC 1	Post HC 2	Post HC 3	Post Scan				
	1	61.83	314.96	314.96	314.57	314.35	314.14	302.5				
diff fro	weight - pa	weight	253.13	253.13 0	252.74 -0.39	252.52 -0.61	252.31 -0.82	240.67 -12.46				tand dev
diff fro	m previous	weight		0	-0.39	-0.22	-0.21	-11.64		Mean start weight	270.19	15.74
	e <mark>st H2O los</mark> at loss durin	S	%H2U loss				0.62	0.32 4.60	4.92			
			315.36	215.20	314.75	314.42	314.24	301.57		Mean H20 gain	0.00	0.00
	2 weight - pa	61.24	315.28 254.04	315.28 254.04	253.51	253.18	253	240.33		Mean H2O loss	-0.73	0.11
diff fro	om original orn previous	weight		0	-0.53 -0.53	-0.86 -0.33	-1.04 -0.18	-13.71 -12.67		Mean change HC1	-0.40	0.09
	Est N2O los	8	%H20 loss		-0.53	-0.33	-1.04	0.41		Mean change HC2	-0.25	0.05
	at loss durin	g scan						4.99	5.40	Mean change HC3	-0.20	0.02
	3	61.79	349.42	349.42	349.13	348.91	348.69	335.37		Mean change scan		0.69
,,,,,,,,	weight - pa	n	287.63	287.63	287.34 -0.29	287.12 -0.51	266.9 -0.73	273.58 -14.05		Mean % loss scan Mean % total loss	4.74 5.06	0.16
diff fro	om original om previous Est h2O los	weight	POTENTIAL PROPERTY.	Ŏ	-0.29	-0.51	-0.22	-13 32			0.00	
	Est h20 los		%H20 loss			-	-0.73	0.25 4.63	4.88		-	
% ma												
	4	61.77	342.1 280.33	342.1 280.33	341.72 279.95	341.45 279.68	341.26 279.49	328.12 266.35				
diff fr	weight - pa	weight	200.33	0	-0.38	-0.65	-0.84	-13.98				
diff fro	om previous Est h2O los	weight	%HOO loss	0	-0.38	-0.27	-0.19	-13.14 0.30	-			
	at loss durin		2012U 1088					4.69	4.99			
	5	61.09	336.9	336.9	336.51	336.28	336.06	322.87			1 - 1	
	weight - pa	n	275.81	275.81	275.42	275.19	274.97	261.78				
	THE RESERVE THE PARTY NAMED IN	weight		0	-0.39	-0.62 -0.23	-0.84 -0.22	-14.03 -13.19				
diff fr	om original	Wainbt		0	-0.39							

Table 118. Fuji II LC 24-Hour Thermal Weight Data

nm	40 ut nac	ref pan we	ight	33 03 0	6-Oct-05								
nm	40 µі рап	ref pan we		33.11 0	7-Nov-05								
		nan amnty	filled + cure		Pre HC	Post HC 1	Post HC 2	Post HC 3	Post Scan To	tal%loss To	otal/Gain% Se	can/Gain% S	can/RetH
	1	32.95	148.09		152.9	150.37	149.7	149.07	142.55				
diff fro	weight - pa m original	n weight	115.14		119.95 4.81	117.42 2.28	116.75 1.61	116.12 0.98	109.60 -5.54				
diff from	n previous	weight		OF COLUMN 2	4.81	-2.53	-0.67	-0.63	-6.52				
е	loss durin	8	%H2O loss					-3.83	5.44	8.63	215.18	135.55	665.31
70 IIIai													
	2 weight - pa	32.7	149.4 116.7		153.58 120.88	151.46 118.76	150.66 117.96	150.05 117.35	143.99 111.29				
diff fro	m original	weight			4.18	2.06	1.26	0.65	-5.41				
diff from	n previous	weight	%HDO loss	The same	4.18	-2.12	-0.8	-0.61 3.53	-6.06 2.92				
	loss durin	g scan							5.01	7.93	229.43	144.98	932.31
	3	33.12	130.99		134.58	132.32	131.75	131.29	126.52				
	weight - pa	n	97.87		101.46	99.2	98.63	98.17	93.40				
diff fro	m original	weight	NO SECURITOR DE LA COMPA	-	3.59	1.33	0.76 -0.57	0.3	-4.47 -4.77				
E	n previous st h20 los	8	%H2O loss					-3.29	3.24		224.54	400.07	4500.0
% mat	loss durin	g scan							4.70	7.94	224.51	132.87	1590.0
	4	32.94	122.61		125.75	123.61	123.06	122.6 89.66	117.93 84.99				
diff fro	weight - pa m original	weight	89.67	200	92.81	90.67	90.12 0.45	-0.01	-4.68				
diff from	n previous	weight	SEAT WATE		3.14	-2.14	-0.55	-0.46	-4.67				
	loss durin		MITEO IGSS						5.03	8.43	249.04	148.73	46700.0
	5	32.59	117.92		119.98	118.84	118.3	117.84	113.58				
	weight - pa	n	85.33		87.39	86.25	85.71	85.25	80.99				
diff fro	m original m previous	weight		NAME OF TAXABLE PARTY.	2.06	0.92	0.38 -0.54	-0.08 -0.46	-4.34 -4.26				
	Est h20 k	SS	%H2O loss		-		The second second	-214	2.45				
	% mat los	s during sc	an						4.87	7.32	310.68	206.80	-5325.0
nm	100 µl par	reference	wt	63.76									
	F. F. T.	ref pan we		63.74		07-Nov-05							
		pan empty	illed + cure		Pre HC			Post HC 3	Post Scan				
	1	63.75	235.42		236.82	235.91	235.52	235.21 171.46	226.67		-		
	weight - pa m original		171.67	-	173.07	172.16 0.49	171.77 0.1	171.46 -0.21	162.92 -8.75				
diff from	n previous	weight	WING TO SERVICE STATE OF THE PARTY OF THE PA		1.4	-0.91	-0.39	-0.31	-8.54				
	st H2O los loss durin	S	%H2O loss					-1,61	4.93	5.86	725.00	610.00	4066.6
			250.00		259.78	258.96	258.5	258.17	248.78				
	2 weight - pa	63.79	258.08 194.29	B1 77827	195.99	195.17	194.71	194.38	184.99				
diff fro	m original m previous	weight			1.7	0.88	0.42 -0.46	0.09	-9.30				
diff from	n previous	weight	%H2O loss		1.7	-0.82	-0.46	-0.33 -1.61	-9.39 0.82				
	loss durin	g scan							4.79	5.61	647.06	552.35	10433.3
	3	61.56	228.4		230.26	229.31	228.98	228.5	219.97				
	weight - pa	n	166.84		168.7	167.75	167.42	166.94 0.1	158.41				
diff fro	m original	weight		MANAGED IN	1.86	-0.95	-0.33	-0.48	-8.53				
	n previous		%H2O loss					-1.76	1.04	6.10	553.23	458.60	8530.0
76 mat	loss durin	g scan							3.00	0.10	GGGEG	140.00	-
	4	61.71	228.86	-	230.86	229.84 168.13	229.39 167.68	228.96 167.25	220.80 159.09				
diff fro	weight - pa m original	n weight	167,15	200000000000000000000000000000000000000	169.15	0.98	0.53	0.1	-8.06				
diff from	n previous	weight	%H2O loss		2	-1.02	-0.45	-0.43	-8.16				
	loss durin		76H2O 1088						4.82	5.95	503.00	408.00	8160.0
	5		258.52		260.63	259.71	259.28	258.85	248.94				
	weight - pa	61.42 n	197.1	MATERIAL SERVICE	199.21	198.29	197.86	197.43	187.52				
diff fro	m original	weight	-		2.11	1.19	0.76	0.33	-9.58 -9.91				
diff from	m previous Est h20 l	ss	%H2O loss	The second second		5/5/5		-1 78	0.89				
		s during so	an					-	4.97	5.87	554.03	469.67	3003.0
nm													
			70-4		Den 110	Best HO 4	Boot HC 3	Post HC 3	Post Sear				
	1	61.84	334.34		336.51	335.58	335.03	334.47	Post Scan 320.90				
	weight - pa	n	272.5		274.67	273.74	273.19	272.63	259.06				
diff fro	m original m previous	weight	DESCRIPTION OF THE PARTY OF THE	and the same of	2.17	1.24 -0.93	0.69 -0.55	0.13 -0.56	+13.44 +13.57				
6	est H2O to	s	%H2O less					-2.04	0.74	E CO.	719.35	625.35	10438.
% mat	loss durin	g scan						****		5.68	7 19.33	623.33	10436
	2	61.48	322.86		325.56	324.58	323.95	323.49	309.77				
	weight - pa		261.38	MA TO SERVICE	264.07	263.1 1.72	1.09	262.01 0.63	248.29 -13.09				
diff fro	m original m previous	weight		INCOME.	2.69	-0.97	-0.63	-0.46	-13.72				
E	t loss durir	S	%H2O loss					-2,00	5.20	5.98	586.62	510.04	2177.7
70 mai							0.00	0.17					
	3 weight - pa	61.31	346.65 285.34	distributed	350.04 288.73	348.83 287.52	348.01 286.7	347.5 286.19	332.58 271.27				
diff fro	m original	weight	200.04		3.39	2.18	1.36	0.85	-14.07				
diff fro	m original m previous st h20 lo	weight	%H2O loss		3.39	-1.21	-0.82	-0.51 -2.54	-14.92 0.88				
	t loss durin		-						0.88 5.17	6.05	515.04	440.12	1755.2
	4	62.07	348.96		352.01	350.62	349.87	349.26	335.74				
	weight - pa	in	348.96 266.89	APPENDED	289.94	288.55	287.8	287.19	273.67				
diff fro	m original	weight	-		3.05	1.66	0.91 -0.75	0.3 -0.61	-13.22 -13.52				
diff fro	m previous	weight	%H2O loss		9,00	1,00	0.75	-2.75	0.95				
	t loss duri								4.66	5.61	533.44	443.28	4506.6
	5	61.45	314.63		317.53	316.36	315.68	315.26	302.29				
	weight - pa	an	253.18		256.08	254.91	254.23	253.81	240.84				
		weight		-	2.9	1.73	1.05	-0.42	-12.34 -12.97		-		
diff fro	m previous	SACO LOS LA											

Table 119. Fuji II LC One-Week Thermal Weight Data

ım	40 µl pan ref pan w			03-Jan-06								
	ref pan v			11-Jan-06								
	pan emp 1 32.65	tyilled + cure 112.82		116.94	114.24	113.85	113.53	108.67	otal%loss To	otal/Gain% Si	can/Gain%	Scan/H2Ore
diff fr	weight - pan om original weight	80.17		84.29 4.12	81.59 1.42	81.2 1.03	80.88 0.71	76.02 -4.15				
diff fro	om previous weight	WINDO INC.		4.12	-2.7	-0.39	-0.32	-4.86 4.05				
	at loss during scan	MIZU IOSS						4.28	9.81	200.73	117.96	684.51
	2 33.12	125.34		128.28	127.07	126.68	126.36	121.00				
diff fr	weight - pan rom original weight	92.22		95.16 2.94	93.95 1.73	93.56 1.34	93.24 1.02	67.66 -4.34				
diff fro	om previous weight	%H2O loss		2.94	-1.21	-0.39	-0.32	-5.36 2.02				
	at loss during scan	MILEO IUSS						4.24	7.65	247.62	182.31	525.49
	3 32.82	140.87		146.3	143.53	142.98	142.43	135.60				
diff fr	weight - pan rom original weight	108.05		113.48 5.43	110.71 2.66	110.16 2.11	109.61 1.56	102.78 -5.27				
diff fro	om original weight om previous weight Est h2O loss	%H2O loss		5.43	-2.77	-0.55	-0.55 -3.07	-6.83 3.41				
	at loss during scan							4.80	9.43	197.05	125.78	437.82
	4 33.03	136.92		141.27	139.4	138.78	138.26	131.94				
diff fr	weight - pan rom original weight	103.89		108.24 4.35	106.37 2.48	105.75 1.86	105.23	98.91 -4.98				
diff fro	om previous weight	%H20 loss	1254	4.35	-1.87	-0.62	-0.52 -3.01	-6.32 2.78				
	at loss during scan							4.57	8.62	214.48	145.29	471.64
	5 33.31	123.05		129.48	124.47	124.46	124.11	118.74				
diff fr	weight - pan rom original weight	89.74		96.17 6.43	91.16 1.42	91.15	90.8 1.06	95.43 -4.31				
diff fro	om previous weight	WH20 less		6.43	-5.01	-0.01	-0.35 -5.37	-5.37 6.58				
	% mat loss during	can						4.33	11.17	167.03	83.51	506.60
ım	100 µl pan reference		61.76		10 :- 0-							
	ref pan v		61.78		12-Jan-06							
	pan emp 1 61.1	ty illed + cure 240.63		Pre HC 248.1	Post HC 1 246.1	Post HC 2 245.06	Post HC 3 244.49	Post Scan 231.79				
	weight - pan	179.53		187 7.47	185	183.96 4.43	183.39 3.86	170.69 -8.84				
diff fro	rom original weight om previous weight est H2O loss			7.47	5.47	-1.04	-0.57	-12.70				
	est H2O loss at loss during scan	%H2O loss					-3.61	5.19	8.72	218.34	170.01	329.02
	2 61.55	222.16		229.6	226.54	225.92	225.4	214.31				
	weight - pan	160.61		168.05	164.99	164.37	163.85	152.76 -7.85				
diff fro	rom original weight om previous weight Est h2O loss		SECTION S	7.44	4.38 -3.06	3.76 -0.62	3.24 -0.52	-11.09				
	Est h2O loss at loss during scan	%H2O less					-4.2	4.92	9.10	205.51	149.06	342.28
	3 61.34	229.68		237.8	234.65	233.86	233.24	221.40				
	weight - pan	168.34		176.46 8.12	173.31 4.97	172.52 4.18	171.9 3.56	160.06 -8.28				
diff fro	rom original weight om previous weight Est h2O loss			8.12	-3.15	-0.79	-0.62	-11.84				
	et h20 loss at loss during scan	%H2O loss	-				-4.56	5.00	9.29	201.97	145.81	332.58
	4 61.22	249.4		258.64	254.2	253.62	253.09	240.06				
4000	weight - pan	188 18		197.42	192.98 4.8	192.4 4.22	191.87 3.69	178.84 -9.34				
diff fro	rom original weight om previous weight			9.24	-4.44	-0.58	-0.53	-13.03				
	at loss during scan	%HZO loss					-9.00	5.15	9.41	201.08	141.02	353.12
	5 61.33	254.61		260.3	258.8	258.3	257.83	244.63				
	weight - pan	193.28		198.97 5.69	197.47 4.19	196.97 3.69	196.5 3.22	183.30 -9.98				
diff fr	rom original weight om previous weight			5.69	-1.5	-0.5	-0.47	-13.20				
	Est h20 loss % mat loss during	%H2O toss scan					-2.47	5.12	7.88	275.40	231.99	409.94
nm												
		tuilled t		Pre HO	Post HC	Post HC	Post HC 3	Post Scan				
	1 61.37		AND DESCRIPTION OF	334.8	333.36	332.57	331.83	313.36				
diff fr	weight - pan rom original weight	265.72		273.43 7.71	271.99 6.27	271.2 5.48	270.46 4.74	251.99 -13.73				
diff fre	om previous weight	%H2O loss		7.71	-1.44	-0.79	-0.74	-18.47 1.09				
% ma	at loss during scan							5.57	7.84	278.08	239.56	389.66
	2 61.78	317.62		328.22	324.43	323.58	323.01	305.01				
diff fe	weight - pan rom original weight	255.84		266.44 10.6	262.65 6.81	261.8 5.96	261.23 5.39	243.23 -12.61				
diff fr	rom previous weight	%H2O loss	ALUANIA	10.6	-3.79	-0.05	-0.57 -5.21	-18.00 1.96				
% m	at loss during scan							5.57	8.71	218.96	169.81	333.95
	3 61.69	305.1		317.93	310.4	309.77	309.27	293.38				
diff 6	weight - pan from original weight	243.41		256.24 12.83	248.71 5.3	248.08 4.67	247.58 4.17	231.69 -11.72				
diff fr	rom previous weight	%H2O foss		12.83	-7.53	-0.63	-0.5 -8.66	-15.89 3.38				
	nat loss during scan	The state of the s						5.14	9.58	191.35	123.85	381.06
	4 61.64			334.85	333.08	332.47	331.93	314.00				
diff 6	weight - pan from original weight	265.57	10000	273.21 7.64	271.44 5.87	270.83 5.26	270.29 4.72	252.36 -13.21				
um II	rom previous weight		-	7.64	-1.77	-0.61	-0.54	-17.93				
diff fr	Est h20 loss nat loss during scan	%H2O loss						5.40	7.63	272.91	234.69	379.87
	The second second			322.46	319.88	319.2	318.62	302.12				
	5 61 76	314.65										
% m	5 61.76 weight - pan from original weight	314.65 252.89		260.7 7.81	258.12 5.23	257.44 4.55	256.86 3.97	240.36 -12.53				

Table 120. Fuji II LC One-Month Thermal Weight Data

nm	40 µl pan	ref pan we	eight		06-Oct-05								
		ref pan we	ight		11-Nov-05								
			filled + cure	the Armes	Pre HC					otal%loss To	tal/Gain% S	can/Gain%S	can/H2O
	1 weight - pa	33.01	96.51 63.5	or some man	100.14 67.13	97.84 64.83	97.39 64.38	96.95 63.94	93.41 60.40				
diff fro	m original	weight			3.63	1.33	0.88	0.44 -0.44	-3.10 -3.54				
е	m previous	S	%H2O loss	A SHALL PLOTE	3.03	-2.3	-0.45	-3.19	4.75				
% mat	loss duri	ng scan							5.27	10.03	185.40	97.52	804.55
	2	32.93	145.72		149.35	148.11	147.6	147.2	140.29				
diff from	weight - pa m original	weight	112.79		116.42 3.63	115.18 2.39	114.67 1.88	114.27 1.48	107.36 -5.43				
diff from	n previous	weight	MICO GEO	100	3.63	-1.24	-0.51	-0.4	-6.91				
	loss duri		%H2O loss	B 6 7 89				-2.10	5.94	7.78	249.59	190.36	466.89
	3	32.98	112.32		115.04	113.96	113.62	113.09	108.38				
	weight - pr	in	79.34		82.06	80.98	80.64	80.11	75.40				
diff from	m original	weight	NAME OF TAXABLE PARTY.	and the same of	2.72	1.64	1.3 -0.34	0.77	-3.94 -4.71				
	n previous st h2O lo		%H2O loss					-1.95	2.38			470.40	****
% mat	loss duri	ng scan	-						5.74	8.12	244.85	173.16	611.69
	4	32.96	136.21		139.92 106.96	138.88 105.92	138.11 105.15	137.61 104.65	130.89 97.93				
diff from	weight - pr m original	weight	103.25		3.71	2.67	1.9	1.4	-5.32				
diff from	n previous	weight	WHO have		3.71	-1.04	-0.77	-0.5	-6.72				
	loss duri								6.28	8.44	243.40	181.13	480.00
	5	32.75	126.68		131.12	129.25	128.4	127.98	122.55				
	weight - pa	an	93.93		98.37	96.5	95.65	95.23	89.80				
diff from	m original n previous	weight	and the latest the lat	NAME OF TAXABLE PARTY.	4.44	2.57	1.72 -0.85	1.3	-4.13 -5.43				
	Est h20 I	ess	%H2O less					-3.14	3.19 5.52	8.71	193.02	122.30	417.69
	% mat lo	ss during sc	an						5.02	8.71	193.02	122.30	417.65
nm	100 µl pa	n reference		63.76		06-Oct-05							
		ref pan we	eight	63.74		07-Nov-05							
			filled + cure	225 100	Pre HC		Post HC 2		Post Scan 248.04				
	1 weight - pa	63.26 an	258.29 195.03	E-17 (12.2)	264.8 201.54	263.65 200.39	263.1 199.84	262.34 199.08	184.78				
diff from	m original	weight		THE PARTY	6.51	5.36 -1.15	4.81 -0.55	4.05	-10.25				
е	n previous	S	%H2O loss	hi Nedati	6.51	-1.15	-0.55	-0.76 -2.46	1 22				
% mat	loss duri	ng scan							7.10	8.32	257.45	219.66	353.09
	2	62.62	235.92		243.35	240.75	240.21	239.58	227.04				
dist 4	voight - no	en .	173.3		180.73 7.43	178.13 4.83	177.59 4.29	176.96 3.66	164.42 -8.88				
diff from	m original m previous st h2O lo	weight		200000	7.43	-2.6	-0.54	-0.63	-12.54				
94 mar	st h2O lo	s scar	%H2O loss					377	6.94	9.02	219.52	168.78	342.62
∞ mát					001	000 00	251.12	200 20					
	3 weight - pa	64 an	247.68 183.68		254.37 190.37	252.03 188.03	251.48 187.48	250.86 186.86	238.07 174.07			-	
diff from	m original	weight		THE S	6.69	4.35	3.8	3.18	-9.61				
diff from	n previous	weight	%H2O loss	a modern	6.69	-2.34	-0.55	-0.62 -3.51	-12.79 1.84				
	loss duri								6.72	8.56	243.65	191.18	402.20
	4	63.25	257.45		264.43	263	262.18	261.77	247.67				
	weight - pa	an	194.2	B	201.18	199.75	198.93	198.52 4.32	184.42 -9.78				
diff from	m original m previous	weight		No. of Concession,	6.98	5.55 -1.43	4.73 -0.82	-0.41	-14.10				
	SHE WALL TO	8 5	%H2O loss					-2.56	1.32 7.01	8.33	240.11	202.01	326.39
% mat	loss duri	ng scan	-							0.00			
	5	64.26	261.87 197.61	e subsectored	268.33 204.07	267.2 202.94	266.64 202.38	266.11 201.85	251.31 187.05				
diff from	weight - pa m original	weight	197.01	SIL I	6.46	5.33	4.77	4.24	-10.56				
diff from	n previous	weight oss	W100 H		6.46	-1.13	-0.56	-0.53	-14.80				
	% mat lo	ess ss during sc	an						7.25	8.34	263.47	229.10	349.06
mm		1											
					_				5-15				
	1	pan empty 64.27	322.79		330.88	Post HC 1 329.6	Post HC 2 329.04	328.48	Post Scan 309.52	-			
	weight - p	an	258.52	1000	266.61	265.33	264.77	264.21	245.25				
diff fro	m original m previous	weight			8.09	6.81	6.25 -0.56	5.69 -0.56	-13.27 -18.96				
6	est H2O to	8 S	%H2O loss					-2.4	0.90 7.11	8.01	264.03	234.36	333.2
% mat	loss duri	ng scan	-							0.01	204.03	234.30	333.2
	2	64.67	314.62	-	322.61	321.09	320.55 255.88	320.1 255.43	301.61 236.94				
diff fro	weight - pa m origina	weight	249.96	2	257.94 7.99	256.42 6.47	5.93	5.48	-13.01				
diff from	m previous	weight	%HQC loss		7.99	-1.52	-0.54	-0.45 -2.51	0.97				
	loss duri		THE STREET						7.17	8.14	262.83	231.41	337.4
	3	64.05	300.28		308.14	306.37	306.63	305.11	267.66				
	weight - p	an	236.23		244.09	242.32	241.58	241.06	223.61				
diff fro	m origina m previou	weight		The state of the s	7.86 7.86	6.09	5.35 -0.74	4.83 -0.52	-12.62 -17.45				
E	st h2O lo	es s	%H2O loss		4			-3.03	7.15	8.39	260.56	222.01	361.2
% mat	t loss duri	ng scan								6.39	200.36	222.01	301.2
	4	63.21	314.03		322.12	320.89	320.22	319.54	301.62				
	weight - p om origina		250.82	NEW YORK	258.91 8.09	257.68 6.86	257.01 6.19	256.33 5.51	238.41 -12.41				
diff from	m previou	weight			8.09	-1.23	-0.67	-0.68	-17.92				
	t loss duri		%H2O loss					The second second	6.92	7.92	253.40	221.51	325.2
		I	210.0		220 4	300.04	325.84	325.37	306.91				
		62.89	319.9 257.01	14 3 4 KM	328.1 265.21	326.61 263.72	262.95	262.48	244.02				
	5 weight - p	an							-12.99				
diff fro	weight - p om origina m previou	weight	257,01		8.2 8.2	6.71	-0.77	5.47	-12.99				

Table 121. Fuji II LC Three-Month Thermal Weight Data

nm	40 μl pan	ref pan we ref pan we			04-Oct-05 04-Jan-06								
	1	32.83	filled + cure 113.6	Edward III	115.66	114.73	114.39	114.06	109.48	tal%loss To	otal/Gain% S	can/Gain% Sc	an/H2Oref
diff fr	weight - pa	weight	90.77		82.83 2.06	81.9 1.13	81.56 0.79	81.23 0.46	76.65 -4.12				
diff fro	om previous	weight	%H2O loss	On the last	2.06	-0.93	-0.34	-0.33 -1.6	-4.58 1.93				
% m	at loss durin	g scan							4.02	7.46	300.00	222.33	995.
	2	32.91	108.7		110.52	109.85	109.47	109.19 76.28	105.16				
diff fr	weight - pa	weight	75.79		77.61 1.82	76.94 1.15	76.56 0.77	0.49	72.25 -3.54				
diff fro	om previous Est h20 los	weight	%H2O loss		1.82	-0.67	-0.38	-0.28 -1.33	1.71				
	at loss durin								3.69	6.91	294.51	221.43	822.
	3	32.95	108.77	-	110.12	109.32	108.99	108.62 75.67	104.77 71.82				
diff fr	weight - pa rom original	weight	75.82		77.17 1.35	76.37 0.55	76.04 0.22	-0.15	-4.00				
diff fro	om previous	weight	%H2O loss	100000	1.35	-0.8	-0.33	-0.37 -1.5	-3.85 1.94				
% ma	at loss durin	g scan							3.54	6.93	396.30	285.19	-2566
	4	33.42	115.58		117.6	116.87	116.68	116.23	111.47				
diff fr	weight - pa		82.16		84.18 2.02	83.45 1.29	83.26 1.1	82.81 0.65	78.05 -4.11				
	om previous		%H2O love		2.02	-0.73	-0.19	-0.45 -1.37	-4.76 1.63				
% m	at loss durin	g scan							4.10	7.28	303.47	235.64	732
	5	32.81	126.35		128.7	127.95	127.55	127.29	121.77				
diff fr	weight - pa		93.54		95.89 2.35	95.14	94.74	94.48 0.94	88.96 -4.58				
diff fro	om previous Est h2O l	weight	%H20 loss		2.35	-0.75	-0.4	-0.26 -1.41	-5.52				
	% mat los	s during sc							4.34	7.23	294.89	234.89	587
nm	100 µl par	reference	wt		04-Oct-05								
		ref pan we	ight	63.47	04-Jan-06								
			filled + cure		Pre HC		Post HC 2	Post HC 3	Post Scan 208.24				
	1 weight - pa	64.19 n	216.01 151.82	Section 1	219.63 155.44	218.96 154.77	154.31	153.95	144.05				
diff fr	rom original om previous	weight	SERVICE STATE		3.62	2.95	2.49	2.13 -0.36	-7.77 -9.90				
	est H2O tot	S	%H2O loss					-1.49	0.96 4.54	7.33	314.64	273,48	464
% ma	at loss durin	g scan								7.33	314.04	273.40	404
	2 weight - pa	63.95	233.64 169.69	erstrete total	238.22 174.27	237.2 173.25	236.54 172.59	236.2 172.25	225.00 161.05				
diff fr	rom original	weight		77 150	4.58 4.58	3.56 -1.02	2.9	2.56 -0.34	-8.64 -11.20				
	om previous Est h20 los	s	%H2O loss		4.00	-1.02	-0.00	-2.02	1.16		200 00		407
% m	at loss durin	g scan							4.74	7.59	288.65	244.54	437
	3	64.06	208.48 144.42	LA TORONISMO	212.03 147.97	211.29 147.23	211.01 146.95	210.56 146.5	201.01 136.95				
diff fr	weight - pa rom original	weight	144.42	1000	3.55	2.81	2.53	2.08	-7.47				
diff fro	om previous Est h20 las	weight	%H2O loss	A PAGE	3.55	-0.74	-0.26	-0.45 -1.47	-9.55 0.99				
% m	at loss durin	g scan							4.54	7.45	310.42	269.01	459
	4	64.24	245.48	-	250.2 185.96	249.4 185.16	248.91 184.67	248.4 184.16	236.26 172.02				
diff fr	weight - pa		181.24		4.72	3.92	3.43	2.92	-9.22				
diff fr	Fat h20 los	weight	%H2O loss		4.72	-0.8	-0.49	-0.51 -1.8	-12.14 0.97				
% m	at loss durin	g scan							4.89	7.50	295.34	257.20	415
	5	64.64	237.05		241.35	240.56	240.03	239.7	228.39				
diff fr	weight - pa	weight	172.41		176.71	175.92 3.51	175.39 2.98	175.06 2.65	163.75 -8.66				
diff fr	rom original rom previous Est h20 l	weight	%H20 (gss		4.3	-0.79	-0.53	-0.33	-11.31				
		s during so	REAL PROPERTY.						4.72	7.33	301.40	263.02	426
nm													
		nan emnts	filled + cure		Pre HC	Post HC 1	Post HC 2	Post HC 3	Post Scan				
	1	63.61	301.86		308.5	307.71	307.24 243.63	306.8 243.19	289.98 226.37				
diff fi	weight - pa fro <mark>m original</mark>	weight	238.25		244.89 6.64	244.1 5.85	5.38	4.94	-11.88				
diff fr	est H2O le	weight	%H2O loss	ALC: N	6.64	-0.79	-0.47	-0.44 -1.7	-16.82 0.69				
% m	at loss durin	ig scan				-			5.48	7.56	278.92	253.31	340
	2	63.82	305.74		311.67	310.92	310.59	310.22 246.4	293.84 230.02				
diff fi	weight - pa from original	weight	241.92		247.85 5.93	247.1 5.18	246.77 4.85	4.48	-11.90				
diff fr	rom previous	weight	%H2O loss	220	5.93	-0.75	-0.33	-0.37 -1.45	-16.38 0.59				
% m	at loss durin	ng scan				-	-		5.28	7.19	300.67	276.22	36
	3	64.01	312.58		318.12	317.4	317.04	316.72	301.40				
diff 6	weight - pa		248.57		254.11 5.54	253.39 4.82	253.03 4.46	252.71 4.14	237.39 -11.18				
	rom previous		%H2O loss		5.54	-0.72	-0.36	-0.32 -1.4	-15.32 0.66				
	nat loss duri	ng scan							0:56 4.84	6.58	301.81	276.53	37
diff fr		63.59	304.05		310.39	309.36	308.74	308.25	292.23				
diff fr	4	in	240.46	ALC: N	246.8 6.34	245.77 5.31	245.15 4.69	244.66 4.2	228.64 -11.82				
diff fr % m	weight - pa				6.34	-1.03	-0.62	-0.49	-16.02 0.87				
diff fr % m					-	1			5.20	7.36	286.44	252.68	38
diff fr % m diff fr	weight - pa	weight s	%H2U loss										
diff fr % m diff fr	weight - pa from original rom previous Est h20 lo nat loss durin	weight s ng scan	289 38		295.56	294.68	294.16	293.74	277.28				
diff fr % m diff f diff fr % m	weight - pa from original from previous Est h20 la nat loss durin 5 weight - pa	s weight s ng scan 63.56	289.38 225.82		232	231.12	230.6	230.18	213.72				
diff fr % m diff f diff fr % m	weight - pa from original rom previous Est h20 le nat loss duri weight - pa from original	s weight s ng scan 63.56 an weight					294.16 230.6 4.78 -0.52		277.28 213.72 -12.10 -16.46				

Table 122. Photac-Fil Quick Immediate Thermal Weight Data

nm 40	μl pan ref pan we	ight	33.1								
	pan empty	filled + cure 106.04	Pre HC 106.04	Post HC 1 105.71	Post HC 2 105.45	Post HC 3 105.19	Post Scan 102.05	Total%los	8		
	ght - pan	75.06	75.06	74.73	74.47	74.21 -0.85	71.07			(mg)	stand de
diff from p	original weight previous weight		ŏ	-0.33	-0.26	-0.26	-3.14		Mean start weight	79.37	8.2
% mat lo	ss during scan	mrizOjuss					1.13 4.18	5.32			
	2 31.3	113.59	113.59 82.29	113.28 81.98	113.02 81.72	112.85 81.55	109.48 78.18		Mean H20 gain Mean H2O loss	-0.54 -0.66	1.2
diff from	ght - pan original weight	82.29	0	-0.31	-0.57 -0.26	-0.74	-4.11 -3.37			-0.36	0.0
E st	h2O loss	%H2O loss	0	-0.31	-U.20	-0.17 -0.74	0.90	4.99	Mean change HC1 Mean change HC2 Mean change HC3	-0.22 -0.19	0.0
% mat lo	ss during scan	101.10	101.46	404.0	400.07	100.70		4.33			0.0
	3 32.85 ght - pan	104.16 71.31	68.61	101.2 68.35	100.97 68.12	100.78 67.93	97.82 64.97		Mean change scan Mean % loss scan	-3.40 4.29	0.6
diff from p	original weight previous weight		-2.7	-2.96 -0.26	-3.19 -0.23	-3.38 -0.19	-6.34 -2.96		Mean % total loss	5.27	0.3
% mat lo	ss during scan	%H2O loss				U.60	4,31	5.31			
	4 33.44	109.23	109.23	108.74	108.53	108.38	105.44				
	ght - pan original weight	75.79	75.79 0	75.3 -0.49	75.09 -0.7	74.94 -0.85	72.00 -3.79				
diff from p	revious weight	%H2O loss	0	-0.49	-0.21	-0.15 -0.85	-2.94 1.12				
% mat lo	ss during scan						3.88	5.00			
wai	5 31.54 ght - pan	123.95 92.41	123.95 92.41	123.58 92.04	123.42 91.88	123.26 91.72	118.67 87.13				
diff frame	original weight		0	-0.37 -0.37	-0.53 -0.16	-0.69 -0.16	-5.28 -4.59				
Es	revious weight t h20 loss mat loss during sc	%H2O loss				-0.69	0.75	5.71			
			0.74				4.01	3.71			
mm 10	0 μl pan reference		3.74								
	pan empty 1 61.26	illed + cure 259.12	Pre HC 269.12	258.61	Post HC 2 258.36	258.1	248.24				
	ght - pan original weight	197.86	197.86	197.35 -0.51	197.1 -0.76	196.84 -1.02	186.98 -10.88			(mg)	stand de
diff from p	previous weight	%H2O loss	0	-0.51	-0.25	-0.26 -1.02	-9.86 0.52		Mean start weight	202.13	16.5
	ss during scan						4.98	5.50			
	2 61.25	268.09 206.84	268.09 206.84	267.53 206.28	267.22 205.97	267.04 205.79	256.50 195.25		Mean H20 gain Mean H2O loss	0.00	0.00
diff from	ght - pan original weight	200.04	0	-0.56 -0.56	-0.87 -0.31	-1.05 -0.18	-11.59 -10.54		Mean change HC1	-0.51	0.0
E st	h2O loss	%H2O loss	·	-0.56	-0.31	-1.05	0.51 5.10		Mean change HC2	-0.28	0.0
% mat lo	ss during scan							5.60	Mean change HC3	-0.21	0.0
wei	3 61.3 ght - pan	254.46 193.16	264.46 193.16	253.9 192.6	253.68 192.38	253.44 192.14	247.64 186.34		Mean change scan Mean % loss scan	-8.78 4.29	2.8
diff from	original weight previous weight		0	-0.56 -0.56	-0.78 -0.22	-1.02 -0.24	-6.82 -5.80		Mean % total loss	4.79	1.1
Est	h2O loss ss during scan	%H2O loss				-1.02	0.53 3.00	3.53			
-	4 61.31	246.11	246.11	245.62	245.33	245.15	239.40				
	ght - pan original weight	184.8	184.8	184.31 -0.49	184.02 -0.78	183.84	178.09 -6.71				
	previous weight	WHICHDAR	0	-0.49	-0.29	-0.18 -0.96	-5.75				
% mat lo	ss during scan	menzo uss					3.11	3.63			
	5 61.33	289.3	289.3	288.85	288.5	288.31	276.34				
diff from	ght - pan original weight	227.97	227.97	227.52 -0.45	227.17 -0.8	226.98 -0.99	215.01 -12.96				
diff from p	revious weight	%H2O loss	0	-0.45	-0.35	-0.19 -0.99	-11.97 0.43				
%	mat loss during so	an					5.25	5.68			
MM	nan emnty	filled + cure	Pre HC	Post HC 1	Post HC 2	Post HC 3	Post Scan				
	1 61.73 ight - pan	364.49	364.49 302.76	364.03 302.3	363.84 302.11	363.64 301.91	346.67 284.94				
diff from	original weight	302.70	0	-0.46 -0.46	-0.65 -0.19	-0.85 -0.2	-17.82 -16.97		Mean start weight	(mg) 300.34	stand de
	H2O loss	%H2O loss		-0.40	-0.15	-0.85	0.28	5.89	moun start trong	300.01	
% mat lo	ss during scan						3.61	3.09		0.00	0.0
wei	2 61.15 ight - pan	346.35 285.2	346.35 285.2	345.98 284.83	345.68 284.53	345.48 284.33	329.29 268.14		Mean H20 gain Mean H2O loss	-0.83	
diff from	original weight	AND DESCRIPTION OF THE PERSON	0	-0.37 -0.37	-0.67 -0.3	-0.87 -0.2	-17.06 -16.19		Mean change HC1	-0.41	
	h2O loss s during scan	%H20 loss				-0.87	D.31 5.68	5.98	Mean change HC2 Mean change HC3	-0.30 -0.23	
	3 61.8	363.51	363.51	363.16	362.82	362.58	346.10		Mean change scan	-16.90	0.8
Wei	ight - pan	301.71	301.71	301.36	301.02	300.78	284.30		Mean % loss scan Mean % total loss	5.63 5.94	0.1
diff from	original weight previous weight	WINOT CO.	0	-0.35	-0.34	-0.93 -0.24 -0.93	-16.48				
E 81	h2O loss ess during scan	%H2O loss				0.90	5.46	5.77			
	4 61.11	376.56	376.56	376.08	375.58	375.33	356.99				
diff from	ight - pan original weight	315.45	315.45 0	314.97 -0.48	314.47 -0.98	314.22 -1.23	295.88 -19.57				
diff from	h20 loss	%H2O loss	0	-0.48	-0.5	-0.25 -1.23	-18.34 0.39				
% mat lo	ss during scan						5.81	6.20			
	5 62.15	358.75	358.75 296.6	358.37 296.22	358.18 296.03	357.93 295.78	341.41 279.26				
We we	ight - pan original weight	296.6	0	-0.38 -0.38	-0.57 -0.19	-0.82	-17.34				
din irom	previous weight				1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	-0.25	-16.52				

Table 123. Photac-Fil Quick 24-Hour Thermal Weight Data

	, 1	pan empty 31.71	illed + cure 106.75	Pre HC 107.31	Post HC 1 106.29	Post HC 2 105.9	Post HC 3 105.55	Post Scan	Total%loss	Total/Gain%	Scan/Gain%	Scan/RetH2
diff from	ight - par original v		75.04	75.6 0.56	74.58 -0.46	74.19 -0.85	73.84	69.51 -6.53				
diff from	previous v	weight		0.56	-1.02	-0.39	-1.2 -0.35	-4.33				
	oss during		96H2U loss				-1.76	5.73	8.06	1087.50	773.21	360.83
	2	30.92	107.73	109.02	108.23	107.74	107.48	103.25				
	ight - pan original v		76.81	78.1 1.29	77.31 0.5	76.82 0.01	76.56 -0.25	72.33 -4.48				
diff from	previous v	weight	%H2O loss	1.29	-0.79	-0.49	-0.26	-4.23				
	oss during		MH2U IOSS			-	-1.04	5.42	7.39	447.29	327.91	-1692.00
	3	32.13	104.39	105.7	104.8	104.42	104.14	100.04				
diff from	ight - pan original v	veight	72.26	73.57 1.31	72.67 0.41	72.29 0.03	72.01 -0.25	67.91 -4.35				
diff from	previous v	weight	WHOO loss	1.31	-0.9	-0.38	-0.28	-4.10				
	oss during		201120 1055					5.57	7.69	432.06	312.98	-1640.00
	4	31.98	127.82	129.23	128.28	127.87	127.5	122.27				
diff from	ight - pan original v	veight	95.84	97.25 1.41	96.3 0.46	95.89 0.05	95.52 -0.32	90.29 -5.55				
diff from	previous	weight		1.41	-0.95	-0.41	-0.37	-5.23				
% mat lo	n20 los		BERZO IOSS					5.38	7.16	493.62	370.92	-1634.38
	5	31.45	110.66	112	111.23	110.8	110.51	106.91				
diff from	ight - pan	No. of Concession, Name of Street, or other transferred to the Concession of the Con	79.21	80.55	79.78 0.57	79.35 0.14	79.06 -0.15	75.46 -3.75				
diff from	original v previous v st h20 to	weight	WIND IS	1.34	-0.77	-0.43	-0.29	-3.60				
%		ss during sc	an				1.45	4.47	6.32	379.85	268.66	-2400.00
nm 10	00 µl pan	reference v	vt 6	3.74								
			filled + cure	Pre HC	Post HC 1	Post HC 2	Post HC 3	Post Scan				
	1	61.69	252.68	255.13	254.03	253.43	252.89	239.91				
diff from	ight - pan original v	veight	190.99	193.44 2.45	192.34 1.35	191.74 0.75	191.2 0.21	178.22 -12.77				
diff from	original v previous v H2O los	weight	%H2O lass	2.45	-1.1	-0.6	-0.54 -2 24	-12.98 1.16				
	oss during							6.71	7.87	621.22	529.80	6180.95
	2	61.45	242.27	247.6	244.46	243.75	243.14	229.54				
diff from	ight - par original v	veight	180.82	186.15 5.33	183.01	182.3 1.48	181.69 0.87	168.09 -12.73				
diff from	previous v	weight	WHOCH lower	5.33	-3.14	-0.71	-0.61	-13.60 2.40				
	oss during				****			7.31	9.70	338.84	255.16	1563.22
	3	62.08	247.55	251.48	249.78	249.09	248.64	235.16				
diff from	eight - par original v	veight	185.47	189.4 3.93	187.7 2.23	187.01 1.54	186.56	173.08 -12.39				
diff from	previous	weight	%H2O loss	3.93	-1.7	-0.69	-0.45 -2.84	-13.48 1.50				
	oss during							7.12	8.62	415.27	343.00	1236.70
	4	61.85	248.78	253.62	251.67	250.89	250.32	235.99				
diff from	ight - par original v	veight	186.93	191.77 4.84	189.82 2.89	169.04 2.11 -0.78	188.47 1.54	174.14 -12.79				
diff from	previous t	weight	%H2O loss	4.84	-1.96	-0.78	-0.57 -3 3	-14.33 1.72				
	oss during							7.47	9.19	364.26	296.07	930.52
	5	61.12	262.42	266.1	264.34	263.71	263.23	248.82				
diff from	eight - par original v	veight	201.3	204.98 3.68	203.22 1.92	202.59 1.29	202.11 0.81	187.70 -13.60				
diff from	previous st h2O lo	weight	%H2O loss	3.68	-1.76	-0.63	-0.48 -2.87	1.40				
%		during sc						7.03	8.43	469.57	391.58	1779.01
IM								D				
	1	pan empty 61.63	filled + cure 386.52	989.91	388.14	387.34	386.68	Post Scan 367.86				
diff from	eight - par original v		324.89	328.28 3.39	326.51 1.62	325.71 0.82	325.05 0.16	306.23 -18.66				
diff from	previous	weight	%H2O loss	3.39	-1.77	-0.8	-0.66	-18.82				
e sa	oss during	5	20 1055				3.23	5.73	6.72	650.44	555.16	11762.5
	2	60.99	375.78	379.55	377.83	376.73	376.12	356.88				
diff fro	eight - par original v		314.79	318.56 3.77	316.84 2.05	315.74 0.95	315.13 0.34	295.89 -18.90				
diff from	previous t h20 los	weight	THE RESIDENCE	3.77	-1.72	-1.1	-0.61	-19.24				
	t h20 los oss during		%H2O loss				-3.43	6.04	7.12	601.33	510.34	5658.82
	3	61.31	334.76	341.2	337.99	336.94	336.24	316.96				
WE	eight - par		273.45	279.89 6.44	276.68 3.23	275.63 2.18	274.93 1.48	255.65 -17.80				
diff from	original v	weight		6.44	-3.21	-1.05	-0.7	-19.28				
Es	h20 los	S	%H2O loss				-4.96	6.89	8.66	376.40	299.38	1302.70
	4	62.03	335.89	340.01	338.55	337.66	336.96	317.13				
	eight - par	1	273.86	277.98	276.52	275.62	274.93	255.10 -18.76				
diff from	original v	weight weight		4.12	2.66 -1.46	1.76	1.07 -0.69	-18.76				
diff from	t b20 los	S	%H2O loss				-3.05	7.13	8.23	555.34	481.31	1853.27
diff from	oss durin	gscan					1					
diff from Es % mat le	oss durin		366.21	368 87	357.2	356 32	355.77	337 26		+		
% mat l	5 eight - par	61.39	365.21 293.82	358.87 297.48 3.66	357.2 295.81 1.99	356.32 294.93 1.11	355.77 294.38 0.56	337.26 275.87 -17.95	-			

Table 124. Photac-Fil Quick One-Week Thermal Weight Data

m 40 µl pan	ref pan we		33.1	13-Oct-05								
	ref pan we	ight	33.07	20-Oct-05								
	pan empty	illed + cure	- profiler						Total%loss	Total/Gain% S	Scan/Gain%	Scan/RetH
1 weight - pa	n	125.49 92.87	97.00	126.76 94.14	126.14 93.52	125.88 93.26	126.15 93.53	120.12 87.5				
diff from original	weight	MINE STATE	and the last	1.27	-0.62	0.39	0.66	-5.37 -6.03				
diff from previous est H2O los		%H2O less					-0.61	0.65	7.05	522.83	474.80	91:
% mat loss durir	ig scan							6.41	7.03	322.03	474.00	31.
2 weight - pa	32.96	140.04 107.08		141.82 108.86	140.89 107.93	140.55 107.59	140.2 107.24	133.57 100.61				
diff from original	weight	107.00	19.00	1.78	0.85	0.51	0.16	-6.47				
diff from previous	8	%H2O loss		1.78	-0.93	-0.34	-0.35 -1.62	-6.63 1.49				
% mat loss durin	g scan							6.09	7.58	463.48	372.47	414
3	32.8	127.3		128.93	127.93	127.62	127.27	121.92				
weight - pa diff from original		94.5	Edition.	96.13 1.63	95.13 0.63	94.82 0.32	94.47	69.12 -5.38				
diff from previous	weight		di della	1.63	-1	-0.31	-0.35	-5.35				
% mat loss durin	g scan	%H2O loss					-1.86	5.57	7.29	430.06	328.22	-1783
4	32.93	127.61		129.12	128.36	128.02	127.73	121.91				
weight - pa	n	94.68		96.19	95.43	95.09	94.8 0.12	88.98				
diff from original diff from previous	weight	NISTI HIND	NAME OF TAXABLE PARTY.	1.51	0.75 -0.76	-0.34	-0.29	-5.70 -5.82				
Est h2O for	s	%H2O less					-1.39	1.45 6.05	7.60	477.48	385.43	485
% mat loss durin						Total Control			7.50	477.40	303.43	463
5 weight - pa	32.64 n	122.31 89.67	contract.	123.84 91.2	122.91 90.27	122.51 89.87	122.12 89.48	116.78 84.14				
diff from original	weight	05.31	1700	1.53	0.6	0.2	-0.19	84.14 -5.53				
diff from previous	weight	%H2O loss	14/4/8	1.53	-0.93	-0.4	-0.39 -1.72	-5.34 1.89				
	s during sc	an						5.86	7.74	461.44	349.02	-281
nm 100 µl par	reference			13-Oct-05								
	ref pan we		63.88	20-Oct-05								
	pan empty	illed + cure				Post HC 2		Post Scan				
1 weight - pa	63.26	236.38 173.12		239.36 176.1	238.48 175.22	237.98 174.72	237.51 174.25	226.03 162.77				
diff from original	weight	110.112		2.98	2.1	1.6	1.13	-10.35				
diff from previous	weight	%H2O loss	Sec. 1	2.98	-0.88	-0.5	-0.47	1.06				
% mat loss durin	g scan							6.52	7.57	447.32	385.23	101
2	63.71	261.84		264.76	264.12	263.55	263.19	250.46				
weight - pa	n	198.13		201.05	200.41	199.84	199.48	186.75 -11.38				
diff from previous	weight		THE REAL PROPERTY.	2.92	-0.64	-0.57	-0.36	-12.73				
% mat loss durin	8	%H2O less					-1,57	6.33	7.11	489.73	435.96	94
		250.03		252.82	252.04	251.46	251.11	238.60				
3 weight - pa	64.58 n	250.03 185.45		188.24	187.46	186.88	186.53	174.02				
diff from original	weight			2.79	2.01	1.43	1.08	-11.43 -12.51				
Est h20 lo	S	%H2O loss					-1.71	0.91			***	
% mat loss durin	g scan							6.65	7.55	509.68	448.39	115
4	63.52	254.51	and the same	257.79	256.61	256.28	255.77 192.25	243.53 180.01				
weight - pa diff from original		190.99		194.27 3.28	193.09	192.76	1.26	-10.98				
diff from previous	weight	%H2O less	100	3.28	-1.18	-0.33	-0.51	1.04				
% mat loss durin	ng scan							6.30	7.34	434.76	373.17	97
5	62.96	248.61		252.09	251.05	250.6	250.14	237.90				
weight - pa	in	185.65		189.13	188.09	187.64	187.18	174.94				
diff from original diff from previous	weight		NAME OF	3.48 3.48	-1.04	1.99	1.53	-10.71 -12.24				
Est (20)	SS	%H2O loss					-1.95	1.03	7.50	407.76	351.72	80
	s during sc	an						6.47	7.50	407.76	331.72	80
М	nan carri	filled + cure		Pre HC	Post HC	Poet HC	Post HC 3	Post Scan				
1	63.49	301.34		305.51	304.52	304.03	303.58	287.01				
weight - pa	weight	237.85		242.02 4.17	241.03 3.18	240.54 2.69	240.09	223.52 -14.33				
diff from previous	weight		No. of Lot	4.17	-0.99	-0.49	-0.45	-16.57				
% mat loss during	S	%H2O loss					-1.93	6.85	7.64	443.65	397.36	73
		200.24		300.1	299.21	298.66	298.25	282.70				
2 weight - pa	64.36 in	296.24 231.88	The same	235.74	234.85	234.3	233.89	218.34				
diff from original	weight	-		3.86 3.86	2.97	2.42 -0.55	2.01	-13.54 -15.55				
diff from previous Est h2O to	s	%H2O less		3.00	0.05	1	-1.85	8.78			402.85	77
% mat loss during	ng scan							6.60	7.38	450.78	402.85	11
3	63.94	294.26	No. of the last of	298.04	297.13 233.19	296.59 232.65	296.12 232.18	279.93 215.99				
weight - pa diff fro <mark>m original</mark>	weight	230.32		234.1 3.78	2.87	2.33	1.86	-14.33				
diff from previous	weight	KHOO III	A STATE OF	3.78	-0.91	-0.54	-0.47	-16.19 0.82	-			
% mat loss duri		MITE O TORS						6.92	7.74	479.10	428.31	87
4	63.2	302.01		305.78	304.8	304.18	303.79	288.36	-			
weight - pa	in	238.81		242.58	241.6	240.98	240.59	225.16				
diff from original diff from previous	weight			3.77	2.79 -0.98	2.17 -0.62	-0.39	-13.65 -15.43				
Est h2O lo	s	%H2Q loss					-1.99	0.82	7.18	462.07	409.28	86
% mat loss duri	ng scan	h						6.36	7.18	402.07	409.20	00
5	63.45	300.19		303.72	302.73 239.28	302.35 238.9	301.92 238.47	286.22 222.77				
weight - pa	weight	236.74	THE PARTY NAMED IN	240.27 3.53	2.54	2.16	1.73	-13.97 -15.70				
diff from original diff from previous				3.53		-0.38	-0.43					

Table 125. Photac-Fil Quick One-Month Thermal Weight Data

	ef pan wei ef pan wei		3.13 13-Nov-05								
					Post HC 2	Boot MC 3	Poet Scan T	otal%loss T	ntal/Gain% S	can/Gain% So	an/PatH20
1	33.21	130.15	Pre HC 131.28	130.59	130.23	129.97	123.87	Utal 761088 II	otal/Gain % Si	Larv Gain 76 St	andredirize
weight - pan diff from original w	eight	96.94	98.07 1.13	97.38 0.44	97.02 0.08	96.76 -0.18	90.66 -6.28				
diff from previous v	reight	WH2O loss	1.13	-0.69	-0.36	-0.26	-6.10				
% mat loss during		MATZO IOSS				-1.51	6.22	7.56	655.75	539.82	-338
2	32.56	123.81	125.48	124.54	124.02	123.78	118.54				
weight - pan		91.25	92.92	91.98	91.46	91.22	85.98				
diff from original w	eight reight		1.67 1.67	0.73 -0.94	0.21 -0.52	-0.03 -0.24	-5.27 -5.24				
% mat loss during		%H2O loss				-1.7	1.83 5.64	7.47	415.57	313.77	-1746
											10-1-1
weight - pan	32.74	123.73 90.99	125.4 92.66	124.57 91.83	124.14 91.4	123.77 91.03	118.23 85.49				
diff from original w	eight		1.67	0.84 -0.83	0.41	0.04 -0.37	-5.50 -5.54				
diff from previous v		%H2O loss	1.07	-0.03	-0.40	-1.63	1.76				
% mat loss during	scan						5.98	7.74	429.34	331.74	1385
4	33.4	141.69	143.81	142.82	142.38	142.1	134.90 101.50				
weight - pan diff from original w	eight	108.29	110.41 2.12	109.42	108.98 0.69	108.7 0.41	-6.79				
diff from previous v	veight	%H2O loss	2 12	-0.99	-0.44	-0.28 -1.71	-7.20 1.55				
% mat loss during							6.52	8.07	420.28	339.62	175
5	33.3	144.2	146.56	145.59	145.15	144.82	137.62				
weight - pan diff from original w		110.9	113.26 2.36	112.29	111.85 0.95	111.52 0.62	104.32 -6.58				
diff from previous v	veight		2.36	-0.97	-0.44	-0.33	-7.20				
Est h20 los	S	%H2O loss				-1.74	6.36	7.89	378.81	305.08	116
			63 AF 15 0								
m 100 µl pan	eference v ef pan wei		63.45 15-Oct-05 63.47 15-Nov-05								
			Pre HC		Post HC 2	Post HC 3	Post Scan				
1	64.98	illed + cure 242.47	246.06	244.99	244.41	243.96	231.42				
weight - pan	MARKET SER	177.49	181.08 3.59	180.01 2.52	179.43 1.94	178.98	166.44 -11.05				
diff from previous v	veight		3.59	-1.07	-0.58	-0.45	-12.54				
% mat loss during		%H2O loss				-41	6.93	8.08	407.80	349.30	84
		242.25	246.73	245.8	245.24	244.91	232.29				
2 weight - pan	62.37	243.25 180.88	184.36	183.43	182.87	182.54	169.92				
diff from original w	eight	ne or trap assessment	3.48 3.48	2.55	1.99 -0.56	1.66	-10.96 -12.62				
Est h20 loss		%H2O loss				-1.82	0.99	7.83	414.94	362.64	76
% mat loss during	scan						6.85	7.83	414.94	302.04	/6
3	64.07	246.88	250.38 186.31	249.43 185.36	249 184.93	248.75 184.68	235.97 171.90				
weight - pan diff from original w	eight	182.81	3.5	2.55	2.12	1.87	-10.91				
diff from previous v	veight	%H2O loss	3.5	-0.96	-0.43	-0.25 -1.63	-12.78 0.87				
% mat loss during	scan						6.86	7.73	411.71	365.14	68
. 4	63.4	259.55	263.13	262.26	261.76	261.41	248.12				
weight - pan		196.15	199.73 3.58	198.86	198.36 2.21	198.01 1.86	184.72				
diff from original w diff from previous v Est h2O loss	veight		3.58	-0.87	-0.5	-0.35	-13.29				
% mat loss during	scar	%H2O inss		-		-1.72	0.86 6.65	7.52	419.27	371.23	71
		240.55	252.55	252.42	251.75	251.42	238.52				
5 weight - pan	62.93	249.33 186.4	253.05 190.12	252.12 189.19	251.75 188.82	188.49	175.59				
diff from original w	reight	No. of Concession, Name of Street, or other Designation, Name of Street, or other Designation, Name of Street,	3.72 3.72	2.79 -0.93	2.42 -0.37	2.09 -0.33	-10.81 -12.90				
diff from previous y	ss	%H20 loss	THE PARTY OF THE P	0.00	- CICITAL	163	0.86				Market Activities
% mat loss	during sc	an					6.79	7.64	390.59	346.77	61
м							D 0				
1	64 32	307.44	913.54	Post HC 1	911.55	311.11	Post Scan 292.97				
weight - pan		243.12	249.22	247.68	247.23	246.79 3.67	228.65 -14.47				
diff from original w	veight	STATE OF THE PARTY.	6.1 6.1	4.56	4.11 -0.45	-0.44	-18.14				
est H2O loss	3	%H2O luss				-2.43	7.28	8.25	337.21	297.38	49
				207 7	207.25	200.00					
weight - pan	64.12	304.29 240.17	309.23 245.11	307.79 243.67	307.33 243.21	306.93 242.81	289.03 224.91				
diff from original w	reight	-	4.94	3.5	3.04	2.64	-15.26 -17.90				
diff from previous y	3	%H2O loss	4.34	1,44	0.40	-2.3	0.94			200 20	
% mat loss during	scan						7.30	8.24	408.91	362.35	67
3	64.09	307.15	311.87	310.97	310.43	310.03	292.72				
weight - pan diff from original v	veight	243.06	247.78 4.72	246.88 3.82	246.34 3.28	245.94 2.88	228.63 -14.43				
diff from previous v	weight		4.72	-0.9	-0.54	-0.4	-17.31				
% mat loss during		%H2O loss					6.99	7.73	405.72	366.74	60
		307.53	312.42	311.45	310.89	310.44	292.76				
weight - par	63.53	307.53 244	248.89	247.92	247.36	246.91	229.23				
diff from original v	veight		4.89	3.92	3.36 -0.56	2.91 -0.45	-14.77 -17.68				
diff from previous	s and the same	%H2O loss		0.0		-1.98	0.60	7.00	402.04	361.55	60
% mat loss during	scan						7.10	7.90	402.04	301.33	bl
5	63.75	303.05	306.44	305.99	305.63	305.23	287.95				
weight - par		239.3	242.69 3.39	242.24	241.88 2.58	241.48 2.18	224.20 -15.10 -17.28				
diff from original v		THE RESERVE AND ADDRESS OF THE PARTY OF THE	3.39	-0.45	-0.36	-0.4	THE RESERVE OF THE PERSON NAMED IN				

Table 126. Photac-Fil Quick Three-Month Thermal Weight Data

n 40 µl pan	ref pan we	ight	33.1								
	pan empty	illed + cure	Pre HC	Post HC 1 106.29	Post HC 2 105.9	Post HC 3 105.55	Post Scan 101.22	Total%loss	sTotal/Gain%	Scan/Gain%	Scan/RetH20
1 weight - pa	1	106.75 75.04	107.31 75.6	74.58	74.19	73.84	69.51				
diff from original diff from previous	weight weight		0.56 0.56	-0.46 -1.02	-0.85 -0.39	-1.2 -0.35	-5.53 -4.33				
est H20 los % mat loss durin	S	%H2O loss				-1 76	2:33 5.73	8.06	1087.50	773.21	-360.83
2	30.92	107.73	109.02	108.23	107.74	107.48	103.25				
weight - pa	n	76.81	78.1	77.31	76.82	76.56	72.33				
diff from original diff from previous	weight weight		1.29	0.5 -0.79	0.01 -0.49	-0.26 -0.26	-4.48 -4.23				
% mat loss durin	s n scan	%H2O loss				-1.54	5.42	7.39	447.29	327.91	-1692.00
3		104.39	105.7	104.8	104.42	104.14	100.04				
weight - pa		72.26	73.57	72.67	72.29	72.01	67.91				
diff from previous	weight weight		1.31	-0.9	0.03 -0.38	-0.25 -0.28	-4.35 -4.10				
% mat loss durin	s	%H2O loss				-1.56	5.57	7.69	432.06	312.98	-1640.00
4	31.98	127.82	129.23	128.28	127.87	127.5	122.27				
weight - pa	1	95.84	97.25	96.3	95.89	95.52	90.29				
diff from original diff from previous	weight weight		1.41	0.46	0.05 -0.41	-0.32	-5.55 -5.23				
% mat loss durin	S	%H2O loss	-			-1.73	1.78 5.38	7.16	493.62	370.92	-1634.38
	31.45	110.66	112	111.23	110.8	110.51	106.91				
5 weight - pa	n	79.21	80.55	79.78	79.35	79.06	75.46				
diff from original diff from previous	weight weight		1.34	-0.77	0.14 -0.43	-0.15 -0.29	-3.75 -3.60				
Est 1/20 la	ss s during sc	%H2O loss	-		-	-1.49	1.85	6.32	379.85	268.66	-2400.00
			~~~		-						
m 100 µl pan	reference		3.74								
1	pan empty 61.69	filled + cure 252.68	Pre HC 255.13	Post HC 1 254.03	Post HC 2 253.43	Post HC 3 252.89	Post Scan 239.91				
weight - pa	n	190.99	193.44 2.45	192.34	191.74 0.75	191.2 0.21	178.22 -12.77				
diff from original diff from previous	weight weight		2.45	-1.1	-0.6	-0.54	-12.98				
% mat loss durin	S	%H2O lass				-2.24	6.71	7.87	621.22	529.80	6180.95
2	61.45	242.27	247.6	244.46	243.75	243.14	229.54				
weight - pa	n	180.82	186.15 5.33	183.01	182.3	181.69 0.87	168.09 -12.73				
diff from original	weight weight		5.33	-3.14	-0.71	-0.61	-13.60				+
% mat loss durin	S	%H2C loss				-4.46	7.31	9.70	338.84	255.16	1563.22
3	62.08	247.55	251.48	249.78	249.09	248.64	235.16				
weight - pa	n	185.47	189.4 3.93	187.7	187.01	186.56 1.09	173.08				
diff from original diff from previous	weight		3.93	-1.7	-0.69	-0.45	-13.48				
% mat loss durin	S	%H2O loss				-2.84	7.12	8.62	415.27	343.00	1236.70
4	61.85	248.78	253.62	251.67	250.89	250.32	235.99				
weight - pa	n	186.93	191.77	189.82	189.04	188.47	174.14 -12.79				
diff from previous	weight weight		4.84	2.89 -1.96	-0.78	1.54	-14.33				
% mat loss durin	S	%H2O lass			-	-3.3	7.47	9.19	364.26	296.07	930.52
		262.42	266.1	264.34	263.71	263.23	248.82				
5 weight - pa	61.12 n	201.3	204.98	203.22	202.59	202.11	187.70				
diff from previous	weight		3.68	1.92	1.29 -0.63	0.81 -0.48	-14.41				
diff from previous	ss s during sc	%H2O loss				-2.87	7.03	8.43	469.57	391.58	1779.01
	s during sc										
u	pan empty	filled + cure		Post HC 1	Post HC 2	Post HC 3	Post Scan				
1	61.63	386.52 324.89	389.91 328.28	388.14 326.51	387.34	386.68 325.05	367.86 306.23				
diff from original	weight		3.39	1.62	0.82	0.16 -0.66	-18.66 -18.82				
diff from previous	S	%H2O loss	3.33	- Loft	3.0	-3.23	0.98	6.72	650.44	555.16	11762.50
% mat loss durin	g scan						5.73	0.12	030.44	333.10	11702.30
2	60.99	375.78 314.79	379.55 318.56	377.83 316.84	376.73 315.74	376.12 315.13	356.88 295.89				
weight - pa	weight		3.77	2.05	0.95	0.34	-18.90 -19.24				
diff from previous	S	%H2O loss	3.77			-3.43	1 08	7 42	604 22	510.34	5658.82
% mat loss durin	g scan							7.12	601.33	310.34	3636.62
3	61.31	334.76 273.45	341.2 279.89	337.99 276.68	336.94 275.63	336.24 274.93	316.96 255.65				
weight - pa diff from original	weight	2/3.45	6.44	3.23	2.18	1.48	255.65 -17.80 -19.28				
diff from previous	weight	%H2O loss	6.44	-5.21	-1.05	-4.96	1.77			200 00	4999
% mat loss durin	g scan			1			6.89	8.66	376.40	299.38	1302.70
4	62.03	335.89	340.01	338.65	337.65	336.96 274.93	317.13 255.10				
weight - pa	n weight	273.86	277.98 4.12	276.52 2.66	275.62 1.76	1.07	-18.76				
diff from previous Est h2O los	weight	MH2O loss	4.12	-1.46	-0.9	-0.69 -3.06	1 10				
% mat loss durir		- III C IUS			1		7.13	8.23	555.34	481.31	1853.27
5	61.39	365.21	358.87	357.2	356.32	355.77	337.26 275.87				
weight - pa	n weight	293.82	297.48 3.66	295.81 1.99	294.93	294.38 0.56	-17.95				
Con month original	weight		3.66	-1.67	-0.88	-0.55	-18.51				